

## folding instructions

1. Press spring button (B) on right side of handlebar and lift handlebar.
2. Hold handlebar in one hand and with the other hand, swing left upright outward and squeeze into Z position until fully folded.
3. Drop handlebar over folded Gait Trainer.
4. Secure both uprights by putting the ring (A) over the left upright and under the spring button.

## accessories

- CE 1035 FP - adjustable forearm platforms  
CE 1080 S - trunk support, small, for pediatric and tyke gait trainers  
CE 1080 L - trunk support, large, for adult gait trainer  
CE 1260 - extended uprights for pediatric and tyke gait trainers  
CE 1270 - extended uprights for adult gait trainer

## lifetime limited warranty

Your Drive brand product is warranted to be free of defects in materials and workmanship for the lifetime of the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This Lifetime Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your Drive device or this warranty, please contact an authorized Drive dealer.



© 2012 Medical Depot, Inc. All rights reserved.  
Drive is a trademark of Medical Depot, Inc.  
Port Washington N.Y. 11050 USA Assembled in USA

## comet anterior gait trainer

item # CO 2000, CO 2100, CO 2200

WenzelRe/hab<sup>TM</sup>  
drive



EU Authorized Representative  
Drive Medical LTD  
Ainley's Industrial Estate  
Elland, West Yorkshire,  
United Kingdom HX5 9JP

The Comet Gait Trainer can be customized to suit user's growth and needs as their motor skills change by adjusting or removing the positioning features.

### set up instructions

1. Remove Comet Gait Trainer from carton.
2. Unfold Comet Gait Trainer and open frame. Lift handlebar and position over right upright. Make sure spring button (B) snaps into hole on the handlebar underside.




### position gait trainer before transfer

1. Position handlebar frame at user's waist by opening both height adjustment knobs (C), push in spring buttons and raise or lower uprights to desired height. Tighten knobs.
2. Unfasten both front Fastex clips (E) on seat harness straps.
3. Open Fastex clip (G) of pelvic stabilizer strap.
4. Adjust pelvic stabilizer width by loosening knobs (H) on back of pelvic stabilizer to desired width. Tighten knobs.
5. Position pelvic stabilizer depth by loosening the two screws(I)on the mounting brackets. Move positioning bar (J) forward or backward. Tighten screws.

### transfer into gait trainer

1. Place user into gait trainer and close Fastex clip (G) of pelvic stabilizer strap. Adjust strap so the user is comfortable and adequately supported.
2. Bring up the seat harness between the legs and close Fastex clip (E). Adjust length of front and back straps.
3. Adjust variable resistance tabs on rear wheels by turning the knobs (N) for desired speed and gait alignment.
4. To change wheels from swivel to non-swivel, loosen knob (O), slide bracket backwards and tighten knob when not positioned over the wheel.



Ankle prompts (D):  
Loosen Velcro™ straps and place around users ankles. Length of belt can be adjusted by pulling at the buckle Length of stride can be adjusted by loosening collars (F) and moving it further apart or closer.