



## STANDARD ADJUSTABLE FOAM HANDLE

### EXERCISE SUGGESTIONS USES AND FEATURES

- ▲ Accessory attachment for Cando® exercise bands and tubing
- ▲ Soft, comfortable, yet firm foam grip handle
- ▲ Adds challenge to resistive exercise
- ▲ Heavy-duty webbing for durability
- ▲ Unique fastening system secures all bands or tubing
- ▲ Easy to use! Thread band/tube through adjustable webbing strap and pull down fastener to secure band/tube in place.



### EXERCISE ROUTINE

The exercise routine should be determined by a doctor, therapist, coach, or trainer.

### DETERMINING APPROPRIATE CAN-DO® BAND AND TUBE RESISTANCES

Have your doctor or therapist select a Can-Do® band/tubing or Tube for your present stage of rehabilitation. The Can-Do® band/tubing selected should allow you to adhere to all established exercise guidelines and execute each exercise proficiently.


### WARNING

- Avoid using Can-Do® band/tubing in areas where product could accidentally strike others.
- Read directions and carefully follow exercise instructions at all times.
- Inspect the band/tubing before each exercise session for any nicks or tears that may occur from continued use.
- Always check to see that the band/tubing is securely attached to Foam Handle before any exercise is performed.
- Never tie multiple strength bands/tubings together!
- Avoid looking directly at the Can-Do® band/tubing while performing exercises.
- Immediately consult your doctor or therapist if any discomfort is experienced while performing exercise(s).

# ADJUSTABLE FOAM HANDLE EXERCISE SUGGESTIONS

Attach Cando® exercise band or tubing by loosening straps to form 2 loops.

Thread band or tubing through one loop and out the other as shown in diagram. Pull straps taut to secure band or tubing.



## WRIST FLEXION

### Start

- Sit in chair facing away from and to the side of attachment site.
- Grasp Foam Handle with exercise hand, placing band/tubing under same side foot, bend arm and place lower arm on top of same upper leg.
- Position elbow directly below shoulder and allow hand to hang unsupported in front of upper leg with palm facing up and thumb pointing outward.
- Keep foot firmly positioned on band/tubing and allow wrist to bend back and down slightly.

### Finish

- Close hand, bend wrist up and back toward upper arm with thumb continuing to face outward.
- Keep upper torso stationary with buttock and lower back in contact with chair.
- Avoid bending elbow and raising lower arm from top of upper leg.
- Slowly return to starting position.
- Modify placement of band/tubing under foot if increased resistance is desired.



and modify placement of foot on band/tubing if increased resistance is desired.

## WRIST EXTENSION

### Start

- Sit in chair facing away from and to the side of attachment site.
- Grasp Foam Handle with exercise hand, placing band/tubing under same side foot, bend arm and place lower arm on top of same upper leg.
- Position elbow directly below shoulder and allow hand to hang unsupported in front of upper leg with palm facing up and thumb pointing inward.
- Keep foot firmly positioned on band/tubing and allow wrist to bend down and forward in front of same side of knee.

### Finish

- Close hand, bend wrist up and back toward upper arm with thumb continuing to face inward.
- Keep upper torso stationary with buttock and lower back in contact with chair.
- Avoid bending elbow and raising lower arm from top of upper leg.
- Slowly return to starting position.
- Modify placement of band/tubing under foot if increased resistance is desired.



## SHOULDER FLEXION

### Start

- Stand, grasp Foams Handle with hand of exercise arm and face away from attachment site.
- Step on band/tubing with same side foot. Straighten arm and allow it to hang comfortably along side of body with palm of hand facing backward.
- Keep knees slightly bent and wrist firm.

### Finish

- Lift arm up and directly forward in front of body until hand reaches shoulder height and end with palm of hand facing floor.
- Keep head and upper torso stationary and wrist firm.
- Avoid arching lower back and (or) lifting upper shoulder toward ear as arm is lifted upward.
- Slowly return to starting position.
- Move body further away from attachment site



## SHOULDER ABDUCTION

### Start

- Stand with exercise arm positioned furthest from attachment site. Grasp Foam Handle and step on band/tubing with foot furthest from attachment site.
- Straighten exercise arm and allow it to hang comfortably along side of body with palm of hand facing body.
- Keep knees slightly bent and wrist firm.

### Finish

- Lift arm up and directly out from side of body until hand reaches shoulder height and end with palm of hand facing floor.



- Keep head stationary and wrist firm.
- Avoid bending upper torso toward attachment site as arm is lifted away from side of body .
- Slowly return to starting position.
- Move body further away from attachment site and modify placement of foot on band/tubing if increased resistance is desired.

## SCAPULA ELEVATION

### Start

- Stand with exercise arm positioned closest to attachment site. Grasp Foam Handle and step on band/tubing with foot closest to attachment site.
- Straighten exercise arm and allow it to hang comfortably at side of body with palm of hand facing backward.
- Keep knees slightly bent, wrist firm and elbow locked.



### Finish

- Lift and elevate arm and shoulder up toward ear.
- Keep head stationary, arm straight and avoid bending upper torso away from attachment site as shoulder is elevated.
- Slowly return to starting position.
- Modify placement of foot on band/tubing if increased resistance is desired.
  - OPTION: Above exercise may be upgraded by allowing elbow to bend and perform an “up-right row” movement following the “Elevation” exercise.

## ELBOW FLEXION

### Start

- Stand with exercise arm positioned furthest from attachment site and grasp Foam Handle.
- Step on band/tubing with same side foot straighten arm and allow it to hang comfortably along side of body with palm of hand facing side of leg.
- Keep knees slightly bent and wrist firm.




### Finish

- Bend elbow and bring hand up and out in front of same side shoulder with thumb facing outward.
- Keep head stationary and wrist firm.
- Avoid moving upper arm away from side of body and (or) arching lower back as arm is bent.
- Slowly return to starting position.
- Move body further away from attachment site and modify placement of foot on band/tubing if increased resistance is required.

## SHOULDER HORIZONTAL EXTENSION

### Start

- Stand with exercise arm positioned furthest from attachment site and grasp Foam Handle.
- Straighten exercise arm across front of body at shoulder height with hand positioned directly out in front of opposite side shoulder and palm facing backward.
- Keep knees slightly bent, arm straight and wrist firm.



### Finish

- Pull arm back and across upper body at shoulders height and end with palm of hand facing forward.
- Keep head stationary and wrist firm.
- Avoid rotating upper body back and away from attachment site as arm is pulled across front of body.
- Slowly return to starting position.
- Move body further away from attachment site if increased resistance is desired.

## SHOULDER D-2 PATTERN

### Start

- Stand with exercise arm positioned furthest from, and body slightly behind attachment site.
- Grasp Foam Handle, straighten exercise arm and position hand out in front of opposite side hip with thumb facing down and back.
- Bend non-exercise arm and place hand on same hip.
- Keep knees slightly bent and wrist firm.



### Finish

- Pull arm up, across and behind body in a diagonal pattern while rotating shoulder backward.
- End with band/tubing across upper chest, arm straight and thumb pointing backward.
- Keep head stationary and wrist firm.
- Avoid rotating upper torso in direction which exercise arm is moving.
- Slowly return to starting position.
- Move body further away from and behind attachment site if increased resistance is desired.

## OTHER USES

### HAND AND WRIST REHABILITATION




Pronation  
Supination  
Ulnar Deviation  
Radial Deviation

### FITNESS & REHABILITATION


Biceps Curl  
Abdominal Crunch  
Obliques  
Lunge Press



Simple, Cost Effective, Portable Exercise Products



1. Cando® exercise band with adjustable handle
2. AccuForce™ band
3. Tubing with Handles™
4. TBS™ bar
5. Exercise web
6. Digi-Extend™
7. Digi-Flex™
8. Theraputty™ material
9. Wate™ bar
10. Dumbbell
11. Cuff® weight
12. Handy ball
13. Medicine ball
14. Balance pad
15. Closed cell fitness mat
16. Open cell fitness mat
17. Jumbo™ foam roller
18. MVP™ balance board
19. Vestibular disc
20. Inflatable ball



Visit us at

[www.CandoProducts.net](http://www.CandoProducts.net)

for more exercise products

Authorized CE representative:

RMS UK Ltd.

28 Trinity Road

Nailsea, Somerset BS48 4NU (UK)

