
Instructions for Use

COLD RUSH[®] COMPACT

3

EN | Instructions for Use




5

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16

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29

EN - Caution: This product has been designed and tested based on single patient usage and is not recommended for multiple patient use. If any problems occur with the use of this product, immediately contact your medical professional.

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ATTENTION - See operating instructions

Meets IEC type BF leakage current requirements

IP21

Protected against solid foreign objects of 12.5mmØ and greater. Protected against vertically falling water drops.

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ENGLISH

INSTRUCTIONS FOR PATIENTS AND MEDICAL PROFESSIONALS

Patients and Medical Professionals should carefully read all of the following instructions and warnings before operating the Cold Rush Compact Cold Therapy System.

INDICATION FOR USE

- The Cold Rush Compact Cold Therapy System is designed to reduce localized pain and swelling following surgical procedures or trauma.
- The Cold Rush Compact Cold Therapy System is intended for single patient use only.

WARNING: RX ONLY

This device should only be used with a prescription from a physician that includes the following treatment information: (1) the number of days or weeks that the treatment should last and (2) the length and frequency of product use (and breaks) during treatment. **DO NOT USE THIS DEVICE IF YOU DO NOT HAVE A PRESCRIPTION FROM YOUR PHYSICIAN THAT INCLUDES ALL SUCH TREATMENT INFORMATION.**

A blank treatment protocol to be filled out by your healthcare provider is provided on top of the Cold Rush Compact Cold Therapy Device.

WARNING: POSSIBLE COLD-RELATED INJURIES

This device reduces the temperature of the skin and tissue. If not used properly and in accordance with the prescribed instructions from your physician, this device can cause serious injury, including frostbite, tissue necrosis and nerve damage.

WARNING: SKIN INSPECTION AND ADVERSE REACTIONS

You or your healthcare provider must check your skin at least every 1–2 hours during use of this device, regardless of whether you are asleep or awake, for any change in skin condition, including for any increase in pain, burning, blistering, itching, increased swelling, discoloration of the skin, increased redness, welts or any other changes in the appearance of the skin. If you experience any of these reactions, you should immediately discontinue use of your Cold Rush Compact Cold Therapy Device and contact your healthcare provider. Any dressing, casting, bracing or wrapping must not interfere with your ability to check your skin condition. If you are unable to check your skin condition for any reason, do not use the Cold Rush Compact Cold Therapy Device.

WARNING: NO DIRECT CONTACT WITH SKIN; INSULATION BARRIER REQUIRED

The Cold Rush Pad should never make direct contact with the skin. The Cold Rush pad operating temperature is too cold to be placed directly on the skin. An insulation barrier between the pad and the skin should be used at all times without exception so as to prevent the skin from becoming too cold. Failure to use an insulation barrier between the pad and the skin can lead to a cold-related injury.

WARNING: USE ONLY ÖSSUR COLD RUSH PADS

Competitive brand cold therapy pads should not be used with the Cold Rush Compact Cold Therapy Device as the temperature may become too cold when using competitive pad and the pads may not otherwise

function properly, causing injury. Only Össur Cold Rush Pads should be used with the Cold Rush Compact Cold Therapy Device.

WARNING: NO MOISTURE IN INSULATION BARRIER

If the moisture is present in the insulation barrier, the skin may become colder than intended. The insulation barrier should be inspected regularly to check for moisture caused by bleeding, sweating or condensation. If moisture is found on the insulation barrier between the pad and the skin, immediately discontinue use of the Cold Rush until the moisture is removed and a fresh barrier is placed between the pad and the skin.

CONTRAINDICATIONS

The Cold Rush Compact Cold Therapy Device should not be used on patients with the following conditions:

- History of cold-related injury or adverse reactions to cold therapy;
- Unconsciousness or incoherence due to general anesthesia, sedation or coma;
- Impairment of circulation or potential impairment of healing in the treatment area, including infection;
- Circulatory disease, including Reynaud's disease, Burger's disease, peripheral vascular disease, vasospastic disorders, sickle cell anemia, or hypercoagulable clotting disorders; or
- Polyneuropathy or other nerve damage.

RISK FACTORS

Physicians should carefully consider the following conditions or factors before before prescribing the Cold Rush Compact Cold Therapy Device:

- Sensitivity to cold;
- General impaired circulation;
- Use of medications that may negatively affect peripheral vascular circulation, including beta adrenergic blockers and local epinephrine use (such as local anesthetics);
- Poor nutrition, smoking, tobacco use, excessive use of caffeine or alcohol and any other behavior negatively affecting circulation;
- Desensitization of the treatment area due to local anesthesia or regional nerve blockers;
- Cognitive impairment and use of medications that have a negative effect upon mental capacity or judgment;
- Moisture at the application site due to excessive bleeding, sweating or condensation;
- Diabetes;
- Use of product on hand, wrist, foot or ankle; and
- Young children and elderly.

DEVICE PARTS IDENTIFICATION (Figure 1)

1. Cold Rush Compact Cold Therapy System
2. Hose
3. Couplings
4. DC Connector 12V Power Supply
5. 12V Power Supply
6. Cold Rush Pad (Sold Separately)

APPLICATION/OPERATION

1. Do NOT plug unit into wall socket until steps 2–7 are completed.
2. Fill reservoir with water to the indicated fill line (**Figure 2**).

Warning: water level should not exceed the reference mark on the

sticker inside the bucket to avoid risk of overflow when closing the cover.

3. Add chopped or cubed ice to the reservoir up to the indicated fill line **(Figure 3)**.

Warning: Make sure the ice level does not exceed the reference mark on the sticker. Failure to comply will cause device overflow when closing the lid.

4. Close the reservoir. With handle in the open position, press lid firmly onto unit making certain that the lid is in full contact with the vessel and the seal is engaged, then pull handle into the locked position **(Figure 4)**
5. Attach the provided 12V power supply to DC power port, but do not plug into wall **(Figure 5)**.
6. Apply insulation barrier over the patient's skin in the area to be treated **(Figure 6)**.
7. Connect Cold Rush pad to hose extending from the device. To ensure proper connection, push couplings together until they snap into place **(Figure 7)**.
8. Apply Cold Rush Pad over insulation barrier. The Cold Rush Pad should never be applied directly over the patient's skin **(Figure 8)**.
9. Plug the provided 12V power supply to wall socket **(Figure 9)**.
10. To turn off the cold therapy device, unplug the power supply unit from wall socket.
11. Always turn the cold therapy device off before disconnecting the Cold Rush Pad or any hoses. To disconnect the hoses, press the release tabs on the coupling and pull the connectors apart.
12. To open lid, press handle down to disengage lid seal, then lift to remove **(Figure 10)**
13. Drain all of the water out of the cold therapy device after each use.
14. Follow directions 2–8 before each use.

TROUBLESHOOTING

Pump will not turn on

Ensure that the power supply connector properly inserted in to port on unit and that the wall plug has power.

Pad will not get cold

Visually inspect pad to see if water is flowing.

If water is flowing: check pad for pockets of trapped air. Gently massage the pad to help remove the trapped air. Once water flow appears to moving freely through the pad, resume usage

If no water is flowing, make certain the unit is properly filled. Refer to filling instructions.

Check the hose connection to ensure that the coupling is fully engaged. Check to ensure that the hose is not kinked or pinched by any object.

Water leaking at connector

Unplug device, disconnect hose, visually inspect to ensure coupling connectors are in the open position and then reconnect hose. Snap or clicking sound indicates proper connection **(Figure 4 & 7)**.

Leaking

Note: certain conditions such as hot or humid climates may cause condensation to develop on unit and or pads.

Leaking Pad -Remove pad, towel dry and then visually inspect pad. If a leak is detected, then the pad will need to be replaced. See warranty regarding pad replacement.

Leaking Unit – visually inspect to determine the location of the leak. Refer to limited warranty for damaged replacement requirements. If no leak is detected, it may be condensation.

HEALTHCARE PROVIDER RESPONSIBILITIES

- The prescribing physician must determine (1) the number of days or weeks that the treatment should last and (2) the length and frequency of product use (and breaks) during treatment. All such information should be included in the patient's prescription.
- Appropriate training must be provided by the healthcare provider, in the proper application, use and operation and care of the Cold Rush Compact Cold Therapy System.
- The healthcare provider must monitor the patient's use of the device to assure compliance with the prescribed protocol, appropriate use, proper application and operation of the Cold Rush Compact Cold Therapy System, including but not limited to application and maintenance of an insulation barrier between the patient's skin and the Cold Rush Pad.

Treatment Time Period	Frequency and Duration of Use	Inspect Skin Every
Day	Through	
Day		
Day	Through	
Day		
Day	Through	
Day		

ADDITIONAL WARNINGS

- The Cold Rush Compact Cold Therapy Device should never be left unattended when plugged in.
- Do not place the tubes or the power cord where your or other's feet can get entangled, causing a fall.
- Never drop or insert any object into any opening or hose.
- Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
- Do not attempt to drag or reposition the Cold Rush Compact Cold Therapy Device by pulling or tugging on the insulated hoses, as this may cause damage to the device and/or cause the hoses to disconnect. Use handle when moving the Cold Rush Compact Cold Therapy Device.
- Do not place the tubes or the power cord by the top of your bed where they could twist around your neck while you are sleeping.
- The pump in Cold Rush Compact Cold Therapy Device is designed to run with water. Running the device without water will cause

permanent damage to the pump.

- A wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect the Cold Rush Compact Cold Therapy Device and should be kept at least a distance 3.3 m away from the Cold Rush Compact Device.

SAFETY STANDARDS

Cold Rush Compact Cold Therapy System is tested and certified to comply with the IEC60601-1 standard of electrical safety of medical devices and IEC/EN60601-1-2, electromagnetic compatibility for medical electrical devices. IEC 60601-1-11, requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment. The company fulfills the requirements of ISO13485.

ELECTRIC SHOCK HAZARD

To create the safest operating conditions possible, Össur has designed the Cold Rush Compact Cold Therapy System to operate on a 12 volt power supply, significantly reducing the chance of dangerous electric shock. Avoid use near any water source.

- Only use the power supply provided with the device.
- Unplug device before filling it with ice and water.
- To avoid danger of electric shock be certain that hands are dry before inserting or removing the power supply from wall socket.
- Avoid getting water or ice on the power adapter.
- Always unplug this product immediately after each use.
- Never use Cold Rush Compact Cold Therapy device while bathing.
- Do not place or store product where it can fall or be pulled into a tub or sink.
- Do not place in or drop into water or other liquid.
- Do not reach for a product that has fallen into water. Unplug immediately.
- Keep cord & device away from heated surfaces.
- Never operate this product if:
 - The device has a damaged electrical cord or plug;
 - The device has been dropped or damaged;
 - The has been dropped into water; or
 - The device is not working properly.
- If any of the above have occurred, immediately discontinue use of the Cold Rush Compact Cold Therapy Device and return the product to your healthcare provider for examination and repair.

TECHNICAL SPECIFICATION: ELECTRICAL - POWER SUPPLY

- Model no: UES06WU-120050SPA
- Input: 120V US
- Input frequency: 60 Hz US
- Output power: 6W
- Output Voltage: 12V
- Output Current: 0.5A

MAINTENANCE

Do not attempt to repair your Cold Rush Compact Device as this may create a user hazard and will automatically void the warranty.

CLEANING AND PROLONGED STORAGE

- When you have finished all of your cold therapy treatments, prior to prolonged storage, fill the device so the water level is at the ice fill line.
- Add 2 tablespoons of liquid bleach or 3 ounces of 3% hydrogen peroxide to the water and mix.
- Connect hoses to the Cold Rush Pad, plug in the device to start.
- Allow water mixture to circulate for 5 minutes.
- To turn device off, unplug power supply.
- Carefully drain the water mixture out of the device and store in a dry, cool, and dark place leaving the lid slightly open. (Exposure to sunlight and extreme heat may damage hoses and device.)
- When reusing the device after storage, fill the reservoir with clean water and turn the pump on to circulate the water through the hoses and flush the remaining water mixture out of the system (use old pad for this procedure, discard the pad after cleaning is completed).
- Dispose water and refill the reservoir with fresh ice and water combination to begin cold therapy treatment. Follow operating instructions 2–9 before each use.
- Before you use the device, inspect hoses and new pad for any tears or breaks. Contact your healthcare provider to purchase new Cold Rush Pad.
- Clean or replace insulation barrier as recommended by your healthcare practitioner.

ENVIRONMENTAL

- Operating Temperature: +10 °C (+50 °F) to 40 °C (104 °F)
- Operating Humidity: 30% - 75% Relative Humidity
- Operating Atmospheric Pressure: 700–1060 h Pa
- Shipping and Storage Temperature: -20 °C (-4 °F) to 60 °C (140 °F)
- Shipping and Storage Humidity: 10% - 90% Relative Humidity, non- condensing
- Shipping and Storage Atmospheric Pressure: 700–1060 h Pa

DISPOSAL

Follow local ordinance or regulations for proper disposal of device, accessories, and packaging.

ACCESSORIES

The following accessories are available for your Cold Rush Compact Cold Therapy System at an additional charge. To order, contact your healthcare provider.

- Cold Rush Ankle Pad
 - Cold Rush Hip Pad Left
 - Cold Rush Hip Pad Right
 - Cold Rush Knee Pad
 - Cold Rush Lumbar Pad
 - Cold Rush Rectangular Pad Small
 - Cold Rush Shoulder Pad Small
 - Cold Rush Shoulder Pad Large
 - Cold Rush Universal Pad
-
- Insulation barrier used for testing thermal properties of the system is a single layer of 3-ply Seal-In Sock, Össur product numbers (K-531310, K-531313, K-531316, K-531320).

LIMITED WARRANTY


Össur warrants, to the original purchaser, that this product is free from defects in materials and workmanship for 6 months from initial purchase. This warranty does not apply if this product: 1 was not purchased from Össur or its authorized dealer 2 has been altered in any way; or 3 has not been used in accordance with Össur's Instructions For Use (IFU). This Limited Warranty does not cover damage due to accidents, neglect, misuse, or operation beyond capacity, parts damaged by improper installation, substitution of parts not approved by Össur, or any other alteration or repair by others that, in Össur's judgment, materially or adversely affects the product or part. This product expected service life is 6 months from the date of first use.
For complete warranty guidelines please visit www.ossur.com

ELECTROMAGNETIC COMPATIBILITY PRECAUTIONS

The Cold Rush Compact Cold Therapy System needs special precautions regarding electromagnetic compatibility (EMC). Specifically it needs to be installed and put into service according to the EMC information provided as follows:

- The Cold Rush Compact Cold Therapy System should not be used adjacent to or stacked with other equipment. In case adjacent or stacked use is necessary, the Cold Rush Compact Cold Therapy System should be observed to verify normal operation in the configuration in which it will be used.
- The Cold Rush Compact Cold Therapy System may be susceptible to electromagnetic interference from portable and mobile RF communications devices such as mobile (cellular) telephones.
- The Cold Rush Compact Cold Therapy System may be interfered with by other equipment, even if that other equipment complies with cisp emission requirements.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or user of device should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	
Harmonics Emissions IEC 61000-3-2	Class B	
Voltage Fluctuations/ Flicker Emissions vIEC 61000-3-3	Complies	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions			
The device is intended for use in the electromagnetic environment specified below. The customer or user of device should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, Short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle. 40 % U_T (60 % dip in U_T) for 5 cycles. 70 % U_T (30 % dip in U_T) for 25 cycles. <5 % U_T (>95% dip in U_T) for 5 s.	<5 % U_T (>95 % dip in U_T) for 0.5 cycle. 40 % U_T (60 % dip in U_T) for 5 cycles. 70 % U_T (30 % dip in U_T) for 25 cycles. <5 % U_T (>95% dip in U_T) for 5 s.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U_T is the a.c. mains voltage prior to application of the test level.			

Recommended separation distances between portable and mobile RF communications equipment and the device.			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz d=1,2√P	80 MHz to 800 MHz d=1,2√P	800 MHz to 2,5 GHz d=2,3√P
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1	At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.		
NOTE 2	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or user of the device should assure that it is used in such an environment.		
Immunity Test	IEC 60601 Test Level	Compliance Level
Conducted RF IEC 61000-4-6	3 Vrms, 150 kHz to 80 MHz	3 Vrms
Radiated RF IEC 61000-4-3	3 V/m, 80 MHz to 2.5 GHz	3 V/m
		Note 1
		Note 2
		*
		**

SAVE THESE INSTRUCTION

Electromagnetic Environment – Guidance
Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
$d=1.2\sqrt{P}$
$d=1.2\sqrt{P}$ 80 MHz to 800 MHz $d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz
where “P” is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey*, should be less than the compliance level in each frequency range**.
Interference may occur in the vicinity of equipment marked with the following symbol: 
At 80 MHz and 800 MHz, the higher frequency range applies.
These guidelines may not apply in situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.
Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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