

# Push® MetaGrip®

Provides CMC joint stability for pain-free thumb mobility

## *The Sophisticated Solution to Thumb CMC Osteoarthritis*

- Can be worn during work and household activities as well as golf, tennis, knitting, swimming, etc.
- Imbedded, contoured metal insert offers custom fit
- Provides CMC stability during pinch & grip without impeding other joint movement
- Provides long term durability, resisting abrasion
- Can be worn under a glove
- Can be machine washed on gentle cycle at low temperature
- Won't deform if left in hot environment

## Measuring / Sizing / Fitting

The Push MetaGrip must be correctly SIZED and SQUEEZED to work! Visit **BraceLab.com** for instructional videos on sizing, fitting, and care or consult your therapist.

- To determine the appropriate size, measure circumference of the hand just below the knuckles, excluding the thumb (see below). This indicates hand size, but not necessarily thumb size. Please take this into consideration when choosing your MetaGrip size. If measurement is near the size break point, try on both sizes to determine best fit.
- Ideal fitting around the thumb is accomplished by trying on the Push MetaGrip. Professional measurement and fitting by a therapist / physician is suggested for optimal results.
- Be sure to follow instructions for squeezing metal insert to fit snugly around your thumb.

### Sizing: Right or Left Hand

SIZE	INCHES	CENTIMETERS
1	6 <sup>3</sup> / <sub>8</sub> - 7 <sup>3</sup> / <sub>4</sub>	16 - 19.5
2	7 <sup>3</sup> / <sub>4</sub> - 9	19.5 - 22.5
3	9 - 10 <sup>1</sup> / <sub>4</sub>	22.5 - 26



*The Push MetaGrip must be correctly SIZED and SQUEEZED to work!*

*Please see next page for insurance filing information.*

---

## Insurance Information

*PLEASE NOTE: This information is provided by BraceLab as a courtesy. We are not a medical provider and we do not file insurance. We are unable to provide any information beyond what is included below. Please contact your insurer with any questions about claims/filing.*

**Medicare** - You must use the CMS form for DME reimbursement, which can be found with instructions at <https://www.cms.gov/Medicare/CMS-Forms/CMS-Forms/Downloads/CMS1490S-ENGLISH.pdf>. BraceLab's receipt (order confirmation) must be attached which has the correct code for filing (L3924-CG: Hand finger orthosis, without joints, may include soft interface, straps, prefabricated, off-the-shelf).

Choose the option: ***"The provider or supplier is not enrolled with Medicare."*** BraceLab does not have a National Provider Identifier Number (NPI number).

**Medicare Advantage, Medicare Supplemental Plans and Other Insurers** - Contact your insurer to obtain the correct form to request reimbursement. Whether your insurer requires the "CG" modifier will vary, so you will need to ask them which code is appropriate: L3924-CG or L3924.

*Following these guidelines does not guarantee reimbursement; that is the decision of each insurer.*