Quick Start Guide

Digital Heart Rate And Pedometer Watch

REF MDSP3044

To watch a video walkthrough of this item, scan the QR code (right) or type in the URL (below) on your computer:
<https://www.medline.com/sku/item/MDPMDSP3044>

Turning on the Watch

- Press & hold all 4 buttons. Watch screen will show "HoLd"
- Continue holding all 4 buttons
- Release buttons when watch displays numbers

Using the Pedometer Mode

NOTE: When you turn the watch ON, the pedometer will default to ON

- Start from time screen (the Home screen)
- Press and hold **UP** until an icon of a runner appears in the bottom right corner of the watch screen

Using the Heart Rate Mode

NOTE: You can measure heart rate at any time/from any screen while using the watch

- Wear the watch for at least five minutes before measuring
- Use two fingers to cover both silver metal pulse sensors
- Cover the sensors fully and exert light & even pressure. Keep your fingers on the sensors
- The heart rate icon will flash when a heartbeat is detected
- The heart rate value (bottom) will update in real-time. The percentage (top) is a reflection of how hard your heart is working
- Press **MODE** to return to time screen

Setting the Date & Time

NOTE: When setting the Date & Time, **UP** moves the value forward and **SET** moves the value backward.

- Starting from time screen, press **SET**
- Press and hold **MODE** until you see "12H" or "24H" flash at the top of the watch screen. Then, release **MODE**
- Press **SET** to switch "12H" (AM/PM time) and "24H" (Military time)
- Press **MODE** to switch to hour selection
- Set the hour by pressing **UP** or **SET**

- Press **MODE** to switch to minute selection
- Set the minute by pressing **UP** or **SET**
- Press **MODE** twice to switch to date format selection
- Press **SET** to switch "M-D" (Month-Day) and "D-M" (Day-Month)
- Press **MODE** to switch to year selection
- Set the year by pressing **UP** or **SET**

- Press **MODE** to switch to month selection
- Set the month by pressing **UP** or **SET**
- Press **MODE** to switch to day selection
- Set the day by pressing **UP** or **SET**
- Press **MODE** to switch to hourly chime selection
- Press **SET** to switch "SIG OFF" (no chime) and "SIG ON" (hourly chime)
- Press and hold **MODE** until you return to the time screen

Please see User Manual for more details. If you have any questions or are experiencing difficulties, please call the Medline Service Center at **1-800-MEDLINE (1-800-633-5463)**. V3 RH22POA

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