

WEARING GUIDE RING POSITIONING



Position ring over 'belly' of finger

Palm side of hand

- Slide your finger into the CIRCUL+ with the raised part of the ring facing away from your body. Fit the CIRCUL+ onto any finger so that the ring is comfortable, but not easy to rotate (i.e. not too loose).
- Be sure the square sensor is positioned over 'belly' of finger.

HEARST GLOBAL LICENSING & BRAND DEVELOPMENT 300 West 57 Street, New York, NY 1003 P 212.492.1309 F 646.280.1109

WEARING GUIDE
RING POSITIONING



SPACER RING
The 'Spacer Ring' is worn at the finger root

- If your 'CIRCL+' ring does not remain snug in position on your finger, use the 'Spacer Ring' as pictured.
- This will keep the 'CIRCL+' ring in position and will ensure your 'Sleep' recording is the most accurate.

HEARST
global licensing &
brand development

Partner:	PREVENTION BIOMEDICS	Project:	Prevention_CIRCUL+
Date:	08.27.21	File Name:	PREV_CIRCUL+_InsertCardB_01_OL
Program:	Illustrator CC	Fonts:	Outlined
Ink/Effects:	<div><div>PMS 7458</div><div>PMS COOL GRAY 10C</div><div>PMS 473</div></div>		