AIRZONE™ Peak Flow Meter

How the Meter works

The AIRZONE™ Peak Flow Meter measures your “Peak Expiratory Flow Rate” (PEFR) which is the fastest speed a person can blow air out of the lungs after taking as big a breath as possible. When you blow into the mouthpiece of the meter the indicator will move along the scale. The number next to the final position is your peak flow measurement. The faster you blow, the higher the reading. The peak flow number tells you how well air is moving through the airways in your lungs. When you use the AIRZONE™ regularly, you will be able to detect changes in your peak flow. Changes in your peak flow will tell you and your healthcare professional what is happening with your lungs. These changes may require special treatment of your condition according to the treatment plan given to you by your healthcare professional. Your healthcare professional will tell you when and how often to use your Peak Flow Meter. It is generally recommended to use it two or three times a day, in the morning when you wake up and before you go to bed. You should always use your meter when you are feeling symptoms of breathing problems.

⚠ Warnings

- When the AIRZONE™ is used for home monitoring of conditions such as asthma, the user should be under the care of a physician or other licensed healthcare professional. A healthcare professional's advice is recommended to understand the meaning and importance of the measurements and to decide on an appropriate treatment plan.
- The treatment plan given to you by your healthcare professional will tell you what action to take when there are changes in your peak flow readings.
- No matter what your peak flow readings are, if you have symptoms such as chest tightness, shortness of breath, coughing or wheezing you should contact your healthcare professional.
- Instructions for using the Peak Flow Meter and the patient diary must be followed carefully. If you are unable to obtain a reading contact your healthcare professional.

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- The Peak Flow Meter is intended for single patient use. The Peak Flow Meter can be used by more than one person if it is cleaned thoroughly after each use.
- If you have any questions about your Peak Flow Meter and its use, contact your healthcare professional.

Troubleshooting

Should you take your peak flow measurement and get an unusually low reading, it could mean that your Peak Flow Meter is broken, or it could mean that the reading is accurate and your asthma is getting worse. Check to make sure that the meter is not broken. Slide the indicator along the slot to make sure it moves freely. If you feel comfortable that you are following the instructions correctly and that your meter is not broken, follow the instructions in your treatment plan for low readings and contact your healthcare professional. If you have any questions about whether your meter is working correctly, call 1-800 848 8923.

Instructions for Use

Step 1

Remove the AIRZONE™ from its packaging. The AIRZONE™ must be cleaned as shown in the care and cleaning section before your initial PEF trial and periodically thereafter.

Step 2

Gently reset the pointer to the end of the slot nearest the mouthpiece.


Step 3

Hold the Meter lightly, so as not to obstruct the slot in which the pointer slides, the holes on the end of the Meter or the vents either side of the mouthpiece. Take as deep a breath as possible and place the mouthpiece end into your mouth, making sure to form a tight seal with your lips.

Step 4

Blow into the AIRZONE™ as quickly and as hard as you can. The correct method is a quick, hard blow, rather than a long slow exhalation.


Step 5

The pointer will slide up the slot, stopping at a number corresponding to your Peak Expiratory Flow (PEF) for that test.

Step 6

Since each test session should consist of three PEF trials, repeat STEPS 1-5 twice more. Remember to reset the pointer to the end of the slot nearest the mouthpiece before each trial. The highest of the three values reached should be recorded, along with the time and date of the test session. DO NOT AVERAGE THE THREE TRIAL VALUES. Some patients may wish to rest for a few minutes between trials.

Determining your Individual Best Reading

In order to understand the importance of the changes in your peak flow reading, it is necessary for your healthcare professional to determine your best peak flow reading. The best reading is the highest peak flow reading that you can reach when you are in your healthiest physical condition. The importance of any changes in your daily reading are dependent upon how much they are different from your best reading.

Healthcare professionals will typically use one of two methods to identify an individual's best reading, upon which to set up a patient's treatment plan:

- Your healthcare professional may set up your treatment plan based upon changes from your own individual best. In this case your healthcare professional will work with you to determine what your individual best is by determining the reading you can reach when you are in your best physical condition, as for example after a period of peak flow monitoring while you are under effective treatment.
- Your healthcare professional may set up your treatment plan based upon changes from a predicted normal peak flow value for people of your same sex, age and height. In this method the predicted normal value becomes your individual best for your treatment plan. Should your healthcare professional prefer this method, we have included a set of Normal Predicted Average Peak Expiratory Flow charts for both adults and children. If you have trouble determining your predicted normal value from the charts, ask your healthcare professional for assistance.

It is important to know that these predicted normal peak flow values are average numbers for large groups of people. You may have a higher peak flow number than the average number and you may not be healthy. Or you may have a lower peak flow number than the average and be healthy. The best way to determine what is a healthy number for you is to discuss this with your physician or other licensed healthcare professional.

No matter which method your healthcare professional prefers to use, it is important that you clearly understand the meaning of your individual best and how it relates to the treatment plan he/she will provide you.

Once your best reading has been determined, record it at the beginning of this pamphlet in the space provided. Also, record your name in the space provided.

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Normal Adult Predicted Average Peak Expiratory Flow (LPM)

Age (yrs)	Men Height					Women Height				
	60"	65"	70"	75"	80"	55"	60"	65"	70"	75"
20	554	602	649	693	740	390	423	460	496	529
25	543	590	636	679	725	385	418	454	490	523
30	532	577	622	664	710	380	413	448	483	516
35	521	565	609	651	695	375	408	442	476	509
40	509	552	596	636	680	370	402	436	470	502
45	498	540	583	622	665	365	397	430	464	495
50	486	527	569	607	649	360	391	424	457	488
55	475	515	556	593	634	355	386	418	451	482
60	463	502	542	578	618	350	380	412	445	475
65	452	490	529	564	603	345	375	406	439	468
70	440	477	515	550	587	340	369	400	432	461

Data from: Leiner GC, et al.:Expiratory Peak Flow Rate. Standard values for normal subjects. Use as a clinical test of ventilatory function. Am. Rev. Resp. Dis. 88:644, 1963.

Normal Child and Adolescent Predicted Average Peak Expiratory Flow (LPM)

Data from: Polger, G, Promedhat V: Pulmonary function testing in children: Techniques and standards. Philadelphia, W.B. Saunders, 1971.

Note: These tables are averages and are based on tests with a large number of people. An individual's Peak Expiratory Flow Rate (PEFR) may vary widely. Further, many individuals' PEFR values are consistently higher or lower than the average values. It is recommended that PEFR objectives for therapy be based upon each "individual's best", which is established after a period of PEFR monitoring while the individual is under effective treatment.

The above note, normal adult and normal child and adolescent tables from: "Guidelines for the Diagnosis and Management of Asthma", U.S. Department of Health and Human Services, Public Health Service, National Institute of Health, Publication No. 97 - 4053, October 1997.

It is important to watch for changes in your peak flow readings from one measuring to the next, and to follow the actions you are to take according to the treatment plan provided to you by your healthcare professional.

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