

# PRODUCT DESCRIPTION BULLETIN

claire

## C-915 GENERAL PURPOSE DISINFECTANT WIPE

### PRODUCT DESCRIPTION:

Wipe formulated to protect hard surfaces from the spread of germs. General Purpose Disinfectant Wipe effectively cleans, disinfects and deodorizes any hard, non-porous, inanimate environmental surface. Ideal for use in hospitals & healthcare settings, clinics, doctor's offices, home health care, emergency rooms, nursing homes, dentist offices, correctional facilities, health clubs, schools, animal kennels, veterinary clinics, gyms, operating rooms, etc.

### EFFECTIVE AGAINST:

Influenza A2/Japan - Escherichia coli - Herpes Simplex Virus 1 & 2 - Campylobacter jejuni  
Klebsiella pneumoniae - Listeria monocytogenes - Pseudomonas aeruginosa  
Salmonella enterica - Staphylococcus aureus - Streptococcus pyogenes



Effective against:  
**H1N1,  
CA-MRSA,  
MRSA,  
VRE**

### SPECS:

Wipe Size:	6" x 10"
Count:	60 wipes per tub
Case Count:	6 tubs per case
Net Wt.:	15.5 oz.
UPC Number:	7 13014 11915 5
Product Number:	C-915
Other:	EPA registered

### REGULATORY:

Recyclable Package:	Yes
VOC Compliant CARB:	Yes
VOC Compliant OTC:	Yes

### PHYSICAL CHARACTERISTICS:

Color:	White cloth
Fragrance:	Lemon

### DIRECTIONS FOR USE:

**CLEANING:** Use one wipe to completely remove surface of all dirt from surfaces to be treated. **DISINFECTING:** To disinfectant hard, non-porous, inanimate environmental surface, thoroughly wet the surface to be treated. Treated surface must remain visibly wet for disinfecting.

### FEATURES + BENEFITS:



- Kills 99.9% of bacteria\* in 30 seconds.
- Kills Flu Virus<sup>††</sup>
- Phenol-free, alcohol-free, bleach-free and non-staining

\*Staphylococcus aureus (Staph), Klebsiella pneumoniae, Listeria monocytogenes, Escherichia coli (E. coli)

††Influenza A2/Japan - H1N1 Strain

CLAIRE MANUFACTURING COMPANY • Addison, IL 60101  
800-252-4731 • www.clairemfg.com

C915\_120109