




TRADITIONAL WOUND CARE

Reference Guide



VALUE AND QUALITY
FOR HEALTH CARE



COTTON WOVEN GAUZE SPONGES



Q: When should I use cotton gauze?

- ✓ Cleaning
- ✓ Scrubbing
- ✓ Packing
- ✓ Debriding
- ✓ Prepping

NON-STERILE

The most cost-effective choice for cleaning, debriding and covering wounds.

- Packaged in sleeves of 200 sponges
- The preferred choice of most nursing homes

Item Number	Item Description	Case Pkg
PRM21208C	2" x 2", 8-ply (5.08 cm x 5.08 cm)	200/sleeve, 25 pk/cs, 5,000/cs
PRM21208CZ	2" x 2", 8-ply (5.08 cm x 5.08 cm)	200/sleeve
PRM21212C	2" x 2", 12-ply (5.08 cm x 5.08 cm)	200/sleeve, 40 pk/cs, 8,000/cs
PRM21312C	3" x 3", 12-ply (7.62 cm x 7.62 cm)	200/sleeve, 20 pk/cs, 4,000/cs
PRM21312CZ	3" x 3", 12-ply (7.62 cm x 7.62 cm)	200/sleeve
PRM21408C	4" x 4", 8-ply (10.16 cm x 10.16 cm)	200/sleeve, 20 pk/cs, 4,000/cs
PRM21408CZ	4" x 4", 8-ply (10.16 cm x 10.16 cm)	200/sleeve
PRM91408C	4" x 4", 8-ply, economy (10.16 cm x 10.16 cm)	200/sleeve, 20 pk/cs, 4,000/cs
PRM21412C	4" x 4", 12-ply (10.16 cm x 10.16 cm)	200/sleeve, 10 pk/cs, 2,000/cs
PRM21416C	4" x 4", 16-ply (10.16 cm x 10.16 cm)	200/sleeve, 10 pk/cs, 2,000/cs

Cotton Filled Dental Sponges

PRM22DS	2" x 2" (5.08 cm x 5.08 cm)	200/sleeve, 25 pk/cs, 5,000/cs
---------	--------------------------------	--------------------------------

STERILE

Ideal for wound debriding, prepping, packing, dressing and general wound care.

- 100% cotton sponges constructed of fine mesh gauze
- Low lint manufacturing process assures debris-free packing
- "C fold" design minimizes loose threads and lint

Sterile 2's

Item Number	Item Desc.	Case Pkg
PRM2208	2" x 2", 8-ply (5.08 cm x 5.08 cm)	50 pk/bx, 30 bx/cs, 1,500 pks/cs
PRM2208Z	2" x 2", 8-ply (5.08 cm x 5.08 cm)	50 pk/bx
PRM3312	3" x 3", 12-ply (7.62 cm x 7.62 cm)	40 pk/bx, 30 bx/cs, 1,200 pks/cs
PRM3312Z	3" x 3", 12-ply (7.62 cm x 7.62 cm)	40 pk/bx
PRM4408	4" x 4", 8-ply (10.16 cm x 10.16 cm)	50 pk/bx, 12 bx/cs, 600 pks/cs
PRM4408Z	4" x 4", 8-ply (10.16 cm x 10.16 cm)	50 pk/bx
PRM4412	4" x 4", 12-ply (10.16 cm x 10.16 cm)	25 pk/bx, 24 bx/cs, 600 pks/cs
PRM4412H	4" x 4", 12-ply (10.16 cm x 10.16 cm)	25 pk/bx
PRM9412/Econ	4" x 4", 12-ply (10.16 cm x 10.16 cm)	25 pk/bx, 24 bx/cs, 600 pks/cs

Sterile 1's

Item Number	Item Desc.	Case Pkg
PRM21419	2" x 2", 12-ply (5.08 cm x 5.08 cm)	100 pk/bx, 24 bx/cs, 2,400 pks/cs
PRM21423	4" x 4", 12-ply (10.16 cm x 10.16 cm)	100 pk/bx, 12 bx/cs, 1,200 pks/cs



CARING® SELF-ADHERENT WRAP



- Self-adherent, (cohesive) wrap, sticks to itself but not to skin, hair or dressings. No pins or clips required.
- Secondary dressing that functions like a tape to retain dressings without tension around digits or limbs.
- Can be applied to provide light compression if desired.
- Promotes pain-free removal that will not damage fragile skin.
- Comfortable, porous and lightweight construction encourages patient compliance.
- Ideal for wrapping splints, gauze pads and dressings securement, compression wrap post blood draw, I.V. securement
- Each bandage is individually wrapped.

Tan, LATEX

Item Number	Item Desc.	Case Pkg
PRM086001	1" X 5 yds stretched, Latex (5.08 cm x 4.57 m)	30 rolls/cs
PRM086002	2" X 5 yds stretched, Latex (5.08 cm x 4.57 m)	36 rolls/cs
PRM086003	3" X 5 yds stretched, Latex (5.08 cm x 4.57 m)	24 rolls/cs
PRM086004	4" X 5 yds stretched, Latex (5.08 cm x 4.57 m)	18 rolls/cs
PRM086006	6" X 5 yds stretched, Latex (5.08 cm x 4.57 m)	12 rolls/cs

Tan, LATEX-FREE

Item Number	Item Desc.	Case Pkg
PRM088001	1" X 5 yds stretched, Latex-free (2.54 x 4.57 m)	30 rolls/cs
PRM088002	2" X 5 yds stretched, Latex-free (5.08 cm x 4.57 m)	36 rolls/cs
PRM088003	3" X 5 yds stretched, Latex-free (7.62 cm x 4.57 m)	24 rolls/cs
PRM088004	4" X 5 yds stretched, Latex-free (10.16 cm x 4.57 m)	18 rolls/cs
PRM088006	6" X 5 yds stretched, Latex-free (15.24 cm x 4.57 m)	12 rolls/cs

COTTON GAUZE SPONGES

USP TYPE VII, Sterile 2's

Item Number	Item Desc.	Case Pkg
PRM21420	2" X 2", 8-Ply (5.08 x 5.08 m)	50 pk/box, 30 bx/cs, 1500 pks/cs
PRM21422	3" x 3", 12-Ply (7.62 cm x 7.62 m)	25 pk/bx, 48 bx/cs, 1200 pks/cs
PRM21424	4" X 4", 12-Ply (10.16 cm x 10.16 m)	25 pk/box, 24 bx/cs, 600 pks/cs
PRM21448	4" X 4", 8-Ply (10.16 cm x 10.16 m)	25 pk/box, 24 bx/cs, 600 pks/cs

USP TYPE VII X-Ray Detectable, Non-Sterile

Item Number	Item Desc.	Case Pkg
PRM25430	4" x 4", 16-Ply (10.16 cm x 10.16 m)	10 sponges/band, 100/sleeve 2000/cs
PRM25432	4" x 8", 12-Ply (10.16 cm x 20.3 m)	10 sponges/band, 100/sleeve 2000/cs
PRM25433	4" x 8", 16-Ply (10.16 cm x 20.3 m)	10 sponges/band, 100/sleeve 1000/cs
PRM25434	4" x 8", 24-Ply (10.16 cm x 20.3 m)	10 sponges/band, 50/sleeve 1000/cs

USP TYPE VII, Sterile 10's

Item Number	Item Desc.	Case Pkg
PRM21426	4" x 4", 12-Ply (10.16 cm x 10.16 m)	10 sponges/pk, 128 pks/cs, 1280/cs
PRM21428	4" x 4", 16-Ply (10.16 cm x 10.16 m)	10 sponges/pk, 128 pks/cs, 1280/cs

USP TYPE VII X-Ray Detectable, Sterile

Item Number	Item Desc.	Case Pkg
PRM21430LF	4" x 4", 16-Ply (10.16 cm x 10.16 m)	10 sponges/pk, 135 pks/cs, 1350/cs
PRM21432LF	4" x 8", 12-Ply (10.16 cm x 20.3 m)	10 sponges/pk, 80 pks/cs, 800 pks/cs
PRM21436LF	4" x 8", 16-Ply (10.16 cm x 20.3 m)	10 sponges/pk, 64 pks/cs, 640/cs
PRM21438LF	4" x 8", 24-Ply (10.16 cm x 20.3 m)	10 sponges/pk, 48 pks/cs, 480/cs



To order, call 1-800-633-5463

NON-WOVEN GAUZE SPONGES



More absorbency and less linting than 100% cotton gauze.

- Unique rayon/polyester formed fabric
- Significantly less linting than cotton gauze
- Virtually non-adherent
- More cost-effective than cotton gauze

Sterile 2's

Item Number	Item Desc.	Case Pkg
PRM21224	2" x 2", 4-ply (5.08 cm x 5.08 cm)	50pk/box, 30 bx/cs, 1500 pks/cs
PRM21444	4" x 4", 4-ply (10.16 cm x 10.16 cm)	50pk/box, 12 bx/cs, 600 pks/cs

NON-STERILE Bulk

Item Number	Item Desc.	Case Pkg
PRM25224	2" x 2", 4-ply (5.08 cm x 5.08 cm)	200/sleeve, 40 sleeves/cs, 8000/cs
PRM25444	4" x 4", 4-ply (10.16 cm x 10.16 cm)	200/sleeve, 10 sleeves/cs, 2000/cs

Sterile 1's

Item Number	Item Desc.	Case Pkg
DG21334	3" x 3", 4-ply (7.62 x 7.62 cm)	1ea/pouch, 10 pouches/box, 24ea/cs

Q: When should I use non-woven?

- ✓ Absorption
- ✓ Low-linting
- ✓ Fast wicking
- ✓ Wound release

ABD/COMBINE PADS

For heavily draining wounds.

- Soft, non-woven rayon/polyester outer layer is virtually non-adherent
- Compressed cellulose center quickly absorbs heavy exudate
- Hydrophobic back retards strike through
- Sealed, finely finished edges prevent linting

Sterile 1's

Item Number	Item Desc.	Case Pkg
PRM21450	5" x 9" (12.7 cm x 22.9 cm)	25/box, 16 bx/cs, 400/cs
PRM21450Z	5" x 9" (12.7 cm x 22.9 cm)	25/box
PRM21453	8" x 7.5" (20.3 x 19.1 cm)	20/box, 12 bx/cs, 240/cs
PRM21454	8" x 10" (20.3 x 25.4 cm)	16/box, 20 bx/cs, 320/cs

NON-STERILE Bulk

Item Number	Item Desc.	Case Pkg
PRM21451	5" x 9" (12.7 cm x 22.9 cm)	576/cs
PRM21456	8" x 10" (20.3 cm x 25.4 cm)	432/cs



Most Caring Traditional Wound Care products are latex-free

To order, call 1-800-633-5463

COTTON GAUZE BANDAGE ROLLS

Ideal for bandaging heads, limbs and difficult to dress wounds.

- 6-ply, 100% cotton
- Unique crimped, open-weave design means maximum absorbency and cushioning
- Gently stretches and conforms to body contours
- Used as either a primary or secondary dressing

Item Number	Item Desc.	Case Pkg
PRM25865	Sterile, 4.5" x 4.1 yds (11.4 cm x 3.75 m)	100 ea/cs
PRM25865H	Sterile, 4.5" x 4.1 yds (11.4 cm x 3.75 m)	1 ea
PRM25855	Non-sterile, 4.5" x 4.1 yds (11.4 cm x 3.75 m)	100 ea/cs
DG25865	Sterile, 4.5" x 4 yds (11.4 cm x 3.66 m)	1ea/pouch, 24ea/cs



SUPRA FORM™ CONFORMING BANDAGES



Gently secures primary dressings without restriction of movement.

- One-ply, rayon/polyester weave conforms to hard to wrap areas
- One-way stretch with memory accommodates swelling
- Finished edges for low linting

STERILE - Unstretched

Item Number	Item Desc.	Case Pkg
PRM25496	2" x 75" (5.08 cm x 190.5 m)	12/bx, 8 bx/cs, 96/cs
PRM25497	3" x 75" (7.62 cm x 190.5 m)	12/bx, 8 bx/cs, 96/cs
PRM25497Z	3" x 75" (7.62 cm x 190.5 m)	12/bx
PRM25498	4" x 75" (10.16 cm x 190.5 m)	12/bx, 8 bx/cs, 96/cs
PRM25498Z	4" x 75" (10.16 cm x 190.5 m)	12/bx
PRM25499	6" x 80" (15.2 cm x 203.2 m)	8/bx, 6 bx/cs, 48/cs

NON-STERILE - Unstretched

Item Number	Item Desc.	Case Pkg
PRM25492	2" x 75" (5.08 cm x 190.5 m)	12/bag, 8 bag/cs, 96/cs
PRM25493	3" x 75" (7.62 cm x 190.5 m)	12/bag, 8 bags/cs, 96/cs
PRM25493Z	3" x 75" (7.62 cm x 190.5 m)	12/bag
PRM25494	4" x 75" (10.16 cm x 190.5 m)	12/bag, 8 bags/cs, 96/cs
PRM25494Z	4" x 75" (10.16 cm x 190.5 m)	12/bag
PRM25495	6" x 80" (15.2 cm x 203.2 m)	12/bag, 4 bags/cs, 48/cs


DRAIN SPONGES

Unique opening helps to wick away more fluid from tube site.

- Excellent fenestration more closely hugs drain tubes, keeping tube sites drier
- 6-ply, rayon/polyester fabric provides maximum absorbency and low linting
- Pre-cut sponge eliminates the need for sterile scissors and conserves nursing time

Sterile 2's

Item Number	Item Desc.	Case Pkg
PRM256001	2" x 2", 6-ply (5.08 cm x 5.08 cm)	35 pk/bx, 20 bx/cs, 700 pks/cs
PRM256000	4" x 4", 6-ply (10.16 cm x 10.16 cm)	25 pk/bx, 12 bx/cs, 300 pks/cs
PRM256000H	4" x 4", 6-ply (10.16 x 10.16 cm)	25 pk/bx



CARING® SURGICAL TAPES



PAPER TAPE

Gentle, breathable tape for sensitive skin.

- Gentle, air-permeable, non-woven material
- Especially recommended for geriatric patients or others with fragile skin
- Conforms and adheres even to damp skin
- Secures lightweight dressing and tubing, ostomy pouches and facial dressings

Item Number	Item Desc.	Case Pkg
PRM260001	1" x 10 yds (2.54 cm x 9.14 m)	12 rolls/bx, 12 bx/cs, 144 rolls/cs
PRM260001H	1" x 10 yds (2.54 cm x 9.14 m)	12 rolls/bx
PRM260002	2" x 10 yds (5.08 cm x 9.14 m)	6 rolls/bx, 12 bx/cs, 72 rolls/cs
PRM260002H	2" x 10 yds (5.08 cm x 9.14 m)	6 rolls/bx
PRM260003	3" x 10 yds (7.62 cm x 9.14 m)	4 rolls/bx, 12 bx/cs, 48 rolls/cs
PRM260003H	3" x 10 yds (7.62 cm x 9.14 m)	4 rolls/bx

TRANSPARENT TAPE

Permits skin examination without tape removal.

- Clear, breathable plastic material
- Micro-perforation allows for bi-directional tearing
- Excellent for IV sites and tubing retention

Item Number	Item Desc.	Case Pkg
PRM260201	1" x 10 yds (2.54 cm x 9.14 m)	12 rolls/bx, 12 bx/cs, 144 rolls/cs
PRM260201H	1" x 10 yds (2.54 cm x 9.14 cm)	12 rolls/bx
PRM260202	2" x 10 yds (5.08 cm x 9.14 m)	6 rolls/bx, 12 bx/cs, 72 rolls/cs
PRM260202H	2" x 10 yds (5.08 cm x 9.14 m)	6 rolls/bx
PRM260203	3" x 10 yds (7.62 cm x 9.14 m)	4 rolls/bx, 12 bx/cs, 48 rolls/cs
PRM260203H	3" x 10 yds (7.62 cm x 9.14 m)	4 rolls/bx

CLOTH "SILK" TAPE

A strong tape for affixing dressings.

- Sturdy rayon/acetate fabric helps secure dressings in place
- Serrated edges facilitate tearing
- Durable adhesion with minimal residue on skin
- Water repellent

Item Number	Item Desc.	Case Pkg
PRM260101	1" x 10 yds (2.54 cm x 9.14 m)	12 rolls/bx, 12 bx/cs, 144 rolls/cs
PRM260101H	1" x 10 yds (2.54 cm x 9.14 cm)	12 rolls/bx
PRM260102	2" x 10 yds (5.08 cm x 9.14 m)	6 rolls/bx, 12 bx/cs, 72 rolls/cs
PRM260102H	2" x 10 yds (5.08 cm x 9.14 m)	6 rolls/bx
PRM260103	3" x 10 yds (7.62 cm x 9.14 m)	4 rolls/bx, 12 bx/cs, 48 rolls/cs

PLASTIC ADHESIVE BANDAGES



Ideal for bandaging minor wounds and scrapes.

- Tan plastic backing is water resistant
- Long-lasting adhesive with minimal residue
- Absorbent non-adherent continuous pad
- Sterile

Item Number	Item Desc.	Case Pkg
PRM25500	3/4" x 3" (1.91cm x 7 cm)	100/bx, 12 bx/cs, 1200/cs
PRM25500H	3/4" x 3" (1.91cm x 7 cm)	100/bx
PRM25600	1" x 3" (2.54 cm x 7 cm)	100/bx, 12 bx/cs, 1200/cs
PRM25600H	1" x 3" (2.54 cm x 7.62 cm)	100/bx

Need Traditional Wound Care items you can't find in this brochure?

If you can't find something you need, please see Medline's Traditional Wound Care brochure or ask your Medline sales representative for assistance.

Call 1-800-MEDLINE (633-5463), or visit our on-line catalog at www.medline.com

1-800-MEDLINE (633-5463) www.medline.com ©2009 Medline Industries, Inc., Medline is a registered trademark of Medline Industries, Inc., One Medline place, Mundelein, IL 60060.

We reserve the right to correct any errors that may occur within this brochure. Avant Gauze®, Caring® and Sheer-Gard® are registered trademarks of Medline Industries, Inc.

Supra Form™ is a trademark of Medline Industries, Inc.

Medline Canada
1-800-396-6996
www.medline.ca
canada@medline.com

Medline México
01-800-831-0898
www.medlinemexico.com
mexico@medline.com