



# DIAGNOSTIC PRODUCTS

Simple Solutions, Significant Results™





# Medline Diagnostics

*Medline is proud to present our latest line of diagnostics. We have built a reputation on providing diagnostics that are highly reliable, extremely accurate, and built to last.*

Blood Pressure Monitors



## Index

Blood Pressure Monitors	3-5
Sphygmomanometers	6-7
Accessories	8-9
Stethoscopes	10-13
Compass Diabetes Program	14-15
Diagnostic Testing	16
Lancets	17
Sharps Containers	18-21
Patient Monitoring	22
Thermometers	21-26
Defibrillator	27
Bladder Scanner	28

If you have questions about any of the products in this brochure, contact your Medline representative or call us at: 1-800-MEDLINE.

**New!**



MDS2007

## Desktop Blood Pressure Monitor Features Average Mode Technology!

These premium monitors detect the appearance of an irregular heartbeat during measurement and gives a warning signal with your reading.

- Measurement includes patented Irregular Heartbeat Detector.
- Monitor also automatically stores 99 readings with date and time.
- Average mode feature takes three consecutive measurements and averages them for more accurate results.

The MDS93AC1PC comes with easy-to-use blood pressure analyzer software with point and click help feature.

Item #	Description	Pkg.
<b>MDS93AC1PC</b>	Premium Advanced Monitor. Includes blood pressure software, USB cable for PC download, 4-AA batteries, medium arm cuff and deluxe storage case.	ea
<b>MDS2007</b>	Premium Monitor Includes 4-AA batteries, medium arm cuff and deluxe storage case.	ea

Pressure Monitor Accessories Replacement cuffs and adapters.

Item #	Description	Pkg.
MDS9872	Adult Cuff, (fits arms 22-30 cm, 8.58"-11.70")	ea
MDS9873	Large Adult Cuff, (fits arms 30-38 cm, 11.70"-14.82")	ea
M-IADE01	AC Adapter	ea



Real-time graphic display of blood pressure being measured, on-screen instructions.



MDS2006

## Average Mode Technology Wrist Blood Pressure Monitor

This advanced blood pressure monitor uses a unique, patented Average Mode Technology to take three consecutive measurements then display a single, averaged result — a procedure practiced by many physicians for a more accurate result.

- Real-time graphic display of blood pressure being measured, on-screen instructions.
- 5 languages.
- 30 memories for 2 separate users with date and time.
- Average mode feature takes three consecutive measurements and averages them for more accurate results.

Item #	Description	Pkg.
MDS2006	Wrist Monitor	ea



MDS2003

## Wrist Blood Pressure Monitor

- One-button operation
- Automatically inflates and deflates
- Digital readout of systolic and diastolic pressure
- Digital readout of pulse
- Memory with date and time
- Operates on two AAA batteries

Item #	Description	Pkg.
MDS2003	Digital BP Wrist Monitor	ea

## Automatic Digital Blood Pressure Unit

Our revolutionary automatic digital blood pressure units have automatic inflation and error detection thanks to the proprietary “fuzzy logic” programming. The easy one-button design provides accurate and reliable readings.

- “Fuzzy logic” circuitry allows this unit to customize inflation levels for individual patients
- Eight-bit processor provides unsurpassed accuracy
- Auto pulse readings
- Memory
- Large digital readout
- Operates on four-AA batteries



MDS2001

Item #	Description	Pkg.
MDS2001	Automatic Inflation w/Std Adult Cuff (fits arms 22-30 cm, 8-12 inches)	ea
MDS2001LA	Automatic Inflation w/Lrg Adult Cuff (fits arms 30-38 cm, 12-15 inches)	ea
MDS9873	Large Adult Cuff (fits arms 30-38 cm, 12-15 inches)	ea
MDS9872	Adult Cuff	ea
M-IADE01	AC Adapter	ea

## Digital Blood Pressure Unit - Manual Inflate

Manually inflated digital blood pressure units have the convenience of a large, digital readout. Also provides accurate and reliable readings.

- Manual inflation
- Auto pulse readings
- Memory
- Large digital readout
- Operates on four-AA batteries



MDS2002

Item #	Description	Pkg.
MDS2002	Manual Inflation w/Std Adult Cuff (fits arms 22-30 cm, 8-12 inches)	ea
MDS2002LA	Manual Inflation w/Lrg Adult Cuff (fits arms 30-38 cm, 12-15 inches)	ea
MDS9873	Large Adult Cuff (fits arms 30-38 cm, 12-15 inches)	ea
MDS9872	Adult Cuff	ea
M-IADE01	AC Adapter	ea

## Digital Blood Pressure Monitors

Large, raised numbers and a luminescent (glow in the dark) coating make the Nite Shift Wall-Mount Aneroid easy to read in almost any type of lighting.

- Eight-foot coiled tubing and pivoting bracket extend service area for more convenient use.
- Inflation system and bulb assembly store easily behind gauge.
- Features Medline's nylon range finder cuff.
- Available in choice of Latex or Latex-Free

Item #	Description	Pkg.
MDS9400	Adult, Latex	ea
MDS9400LF	Adult, Latex-Free Neoprene	ea
MDS9472	Wall mount basket (not shown)	ea
MDS9418	Mobile Stand	ea

The Mobile Aneroid maintains all of the features of its Wall-Mount counterpart, the 9400, while adding the convenience of portability and height adjustment.

- Adjusts from 29" to 48" (74 cm to 122 cm) to accommodate a wide range of viewing levels.
- Includes inflation system and bulb assembly storage basket.

Item #	Description	Pkg.
MDS9407	Adult, Latex	ea
MDS9407LF	Adult, Latex-Free Neoprene	ea



MDS9400LF



MDS9407LF

Blood Pressure Monitors



### Blood Pressure Cuff Barriers

Single-use wrap for added protection from cross-contamination.

- Protection against bacteria pass through
- Will not alter blood pressure readings
- Latex-Free
- Available in 3 sizes to fit any size blood pressure cuff

Item #	Description	Pkg.
MDS916113	Small	300/cs
MDS916114	Adult	300/cs
MDS916115	Large	150/cs



### Blood Pressure Cuff-Guard

Help reduce the potential of cross-contamination by using our protective, disposable blood pressure cuff covers.

- Latex-free
- Made of Tyvek® material
- Quick and easy to replace or remove the cover!

Item #	Description	Pkg.
MDS9157	Child Cuff-Guard	20/box
MDS9158	Adult Cuff-Guard	20/box
MDS9159	Large Adult Cuff-Guard	20/box

#### Soft Fabric for Single Time Use

MDS9167	Child Soft Cuff-Guard	20/bx
MDS9168	Adult Soft Cuff-Guard	20/bx
MDS9169	Lg. Adult Cuff Guard	20/bx



**MDS9114**



**MDS9123**



**MDS9117**

The Compli-Mate Aneroid Sphygmomanometers are designed for many years of reliable service.

- Designer-colored cuffs offer a choice of vibrant colors.
- Comes complete with oversized matching nylon carrying case which allows the user to carry additional items.

Item #	Description	Pkg.
MDS9104	Red	ea
MDS9106	Navy Blue	ea
MDS9107	Lavender	ea
MDS9108	Hunter Green	ea
MDS9109	Orchid	ea
MDS9112	Teal	ea
MDS9113	Purple	ea
MDS9114	Royal Blue	ea

Combining the Compli-Mates Aneroid Sphygmomanometers with the very popular Sprague Rappaport Stethoscope makes this Combination Kit popular with nurses and EMT's.

- Kit comes with an aneroid sphygmomanometer and a 22" (56 cm) Sprague Rappaport stethoscope.
- Packed together in a matching nylon carrying case.

Item #	Description	Pkg.
MDS9121	Magenta	ea
MDS9123	Purple	ea
MDS9124	Royal Blue	ea
MDS9125	Black	ea
MDS9127	Lavender	ea
MDS9128	Hunter Green	ea

The Compli-Mates Dual-Head Combination Kit combines an aneroid sphygmomanometer and a dual-head stethoscope.

- Ideal for nurses and students.
- Kit comes complete with an aneroid sphygmomanometer and matching dual-head stethoscope.
- Packed together in a matching nylon carrying case.

Item #	Description	Pkg.
MDS9116	Magenta	ea
MDS9117	Teal	ea
MDS9118	Purple	ea
MDS9119	Royal Blue	ea
MDS9133	Hunter Green	ea
MDS9135	Burgundy	ea
MDS9136	Red	ea

The Latex-Free Handheld Aneroid's lightweight design and portability make it a popular handheld unit.

- No-stop pin action allows a visual calibration check, ensuring accuracy.
- Nylon rangefinder cuff provides quick, convenient patient fitting.
- Latex-Free construction.
- Field adjustable for long-term accuracy.

PVC			Neoprene		
Item #	Description	Pkg.	Item #	Description	Pkg.
MDS9380	Adult	ea	MDS9380LF	Adult	ea
MDS9386	Infant	ea	MDS9386LF	Infant	ea
MDS9387	Child	ea	MDS9387LF	Child	ea
MDS9388	Large Adult	ea	MDS9388LF	Large Adult	ea



### Our Best Selling Sphygmomanometer.

Our finest handheld aneroid, the Nite-Shift Premier, has a luminescent dial that facilitates reading over a wide range of lighting conditions.

- Durable zinc-alloy shell protects precision aneroid mechanism.
- Shatter resistant acrylic faceplate surrounds the luminescent dial.
- Field adjustable for long-term accuracy.

Item #	Description	Pkg.
MDS9410	Adult	ea
MDS9411	Child	ea
MDS9413	Large Adult	ea



Sphygmomanometers

## FOR HOME USE


Medline's Home Blood Pressure Kits feature D-ring cuffs to simplify home use

- Both units feature an easy-to-read aneroid gauge, complete instructions, recording chart, and carrying case.
- MDS9301 features an attached stethoscope while the two-person MDS9300 features a separate single-head stethoscope.
- Ideal for home use.

Item #	Description	Pkg.
MDS9300	With detached stethoscope	ea
MDS9301	With attached stethoscope	ea

Simple and convenient







MDS91420LF



MDS91410



MDS9505LF



MDS9205



MDS9420

**Two-Tube Inflation Bag & Nylon Range finder Cuff**

Item #	Description	Size	Pkg.
<b>Latex-Free (PVC)</b>			
MDS91420	Adult	5¾" x 22", (15 cm x 56 cm)	ea
MDS91421	Child	4¼" x 15", (11 cm x 38 cm)	ea
MDS91422	Infant	2⅞" x 10¼", (7 cm x 26 cm)	ea
MDS91423	Large Adult	6⅞" x 26", (17 cm x 66 cm)	ea
MDS91424	Thigh	8½" x 30½", (22 cm x 77 cm)	ea
<b>Neoprene Latex-Free</b>			
MDS91420LF	Adult	5¾" x 22", (15 cm x 56 cm)	ea
MDS91421LF	Child	4¼" x 15", (11 cm x 38 cm)	ea
MDS91422LF	Infant	2⅞" x 10¼", (7 cm x 26 cm)	ea
MDS91423LF	Large Adult	6⅞" x 26", (17 cm x 66 cm)	ea

**One-Tube Inflation Bag & Nylon Range finder Cuff**

Item #	Description	Size	Pkg.
<b>Latex-Free (PVC)</b>			
MDS91410	Adult	5¾" x 22", (15 cm x 56 cm)	ea
MDS91411	Child	4¼" x 15", (11 cm x 38 cm)	ea
MDS91412	Infant	2⅞" x 10¼", (7 cm x 26 cm)	ea
MDS91413	Large Adult	6⅞" x 26", (17 cm x 66 cm)	ea
MDS91414	Thigh	8½" x 30½", (22 cm x 77 cm)	ea
<b>Neoprene Latex-Free</b>			
MDS91410LF	Adult	5¾" x 22", (15 cm x 56 cm)	ea
MDS91411LF	Child	4¼" x 15", (11 cm x 38 cm)	ea
MDS91412LF	Infant	2⅞" x 10¼", (7 cm x 26 cm)	ea
MDS91413LF	Large Adult	6⅞" x 26", (17 cm x 66 cm)	ea
MDS91414LF	Thigh	8½" x 30½", (22 cm x 77 cm)	ea

**Two-Tube Seamless Inflation Bag**

Item #	Description	Size	Pkg.
<b>Latex</b>			
MDS9505	Adult	4¾" x 8⅞", (12 cm x 23 cm)	ea
MDS9506	Large Adult	6" x 12¾", (15 cm x 32 cm)	ea
<b>Neoprene Latex-Free</b>			
MDS9505LF	Adult	4¾" x 8⅞", (12 cm x 23 cm)	ea
MDS9506LF	Large Adult	6" x 12¾", (15 cm x 32 cm)	ea

**One-Tube Seamless Inflation Bag**

Item #	Description	Size	Pkg.
<b>Latex</b>			
MDS9205	Adult	4¾" x 8⅞", (12 cm x 23 cm)	ea
MDS9206	Large Adult	6" x 12¾", (15 cm x 32 cm)	ea
<b>Neoprene Latex-Free</b>			
MDS9205LF	Adult	4¾" x 8⅞", (12 cm x 23 cm)	ea
MDS9206LF	Large Adult	6" x 12¾", (15 cm x 32 cm)	ea

**Nylon Range finder Velcro Cuffs**

Item #	Description	Size	Pkg.
MDS9420	Adult	5¾" x 22", (15 cm x 56 cm)	ea
MDS9421	Child	4¼" x 15", (11 cm x 38 cm)	ea
MDS9423	Large Adult	6⅞" x 26", (17 cm x 66 cm)	ea



**Gauges and gauge parts**

Item #	Description	Pkg.
MDS9414	Premier Gauge Only	ea
MDS9415	Face Glass for Premier Gauge	ea
MDS9381	Standard Gauge Only	ea
MDS9507	Face Glass for Standard and Economy Gauge	ea



MDS9414



MDS9381


**Bulbs and Valves**

Item #	Description	Pkg.
MDS9508	Valve Only	ea
MDS91481	Bulb Only, Latex-Free (PVC)	ea
MDS91481LF	Bulb Only, Latex-Free (Neoprene)	ea
MDS91482	Bulb and Valve Assembly, Latex-Free (PVC)	ea
MDS9630	Bulb and Valve Assembly, Premium Latex	ea
MDS9560	Bulb and Valve, Large, Latex	ea

MDS9630



MDS9508



MDS91481LF

**Tubing**

Item #	Description	Pkg.
MDS9473	4-Foot Coiled Tubing w/Connector, Latex	ea
MDS9475	8-foot Coiled Tubing w/Connectors, Latex	ea
MDS9475LF	8-foot Coiled Tubing w/Connectors, Latex-Free (Neoprene)	ea



MDS9475LF

**Other parts**

Item #	Description	Pkg.
MDS9472	Inflation Systems Holder for Wall Mounted or Mobile Sphygmomanometers	ea
MDS9481	Metal Connector	ea



MDS9472



MDS9481



## Our Best Value Stethoscopes!

Quality and economy makes our single and dual-head stethoscopes excellent values.

- Both stethoscopes feature 22" (56 cm) tubing and lightweight design.
- Thick-walled tubing minimizes background interference.
- Dual-heads feature a non-chill ring on the bell for increased patient comfort.

Item #	Description	Pkg.
<b>Single-Head Stethoscopes</b>		
MDS926101	Black	ea
MDS926102	Gray	ea
MDS926103	Blue	ea
MDS926105	Lavender	ea
MDS926106	Red	ea
MDS926107	Hunter Green	ea
MDS926108	Yellow	ea
<b>Dual-Head Stethoscopes</b>		
MDS926201	Black	ea
MDS926202	Gray	ea
MDS926203	Blue	ea
MDS926205	Lavender	ea
MDS926206	Red	ea
MDS926207	Hunter Green	ea



## This Stethoscope has EXCELLENT Sound!

The Sprague Rappaport Stethoscope incorporates versatility into an excellent sound-transmitting device.


- Each stethoscope comes complete with a component kit that includes two sets of eartips, large and small diaphragm, adult bell, pediatric bell, and infant bell.
- The dual 22" (56 cm) tubing design transmits even the faintest sounds.

Item #	Description	Pkg.
MDS926301	Black	ea
MDS926302	Gray	ea
MDS926303	Blue	ea
MDS926305	Lavender	ea
MDS926306	Red	ea
MDS926307	Hunter Green	ea



### Single or Dual-Head Parts

Item #	Description	Pkg.
MDS9544	White Plastic Ear Tips	ea
MDS9547	Soft Rubber Ear Tips	ea
MDS9552	Diaphragm	ea



**Medline's Stainless Steel Stethoscopes offer unsurpassed quality and value**

- Each stethoscope is handcrafted from high tensile strength German stainless steel. Features 22" (56 cm) tubing.
- The color coordinated non-chill rings provide excellent patient comfort.
- All models come in a foam-lined box with one extra set of ear tips.
- Available in adult, pediatric, neonatal, and cardiology styles.
- Cardiology Stethoscope has a deep, cone-shaped bell — a unique feature which provides excellent sound transmission and acoustic sensitivity. Features 17" (44 cm) tubing.

Item #	Description	Pkg.
<b>Adult</b>		
MDS92250	Gray	ea
MDS92260	Black	ea
MDS92270	Blue	ea
MDS92280	Burgundy	ea
<b>Pediatric</b>		
MDS92200	Gray	ea
<b>Neonatal</b>		
MDS92150	Gray	ea
MDS92160	Pink	ea
<b>Cardiology</b>		
MDS92500	Black	ea

Our Bowles Stethoscope features thick-walled 22" (56 cm) tubing, a chrome-plated binaural and a chrome-plated, raised, metal-stem chest piece.

Item #	Description	Color	Pkg.
MDS9515	Bowles Stethoscope	Black	ea

The Pediatric and Neonatal Stethoscopes have all the features of the dual-head stethoscopes, but are sized to fit the needs of the individual practitioner and patient.

- Both stethoscopes feature 22" (56 cm) tubing and lightweight design.
- Thick-walled tubing minimizes background interference.
- Aluminum binaurals.

Item #	Description	Color	Pkg.
MDS9557	Pediatrics	Gray	ea
MDS9565	Neonatal	Gray	ea





### MDF® 2-in-1 Deluxe Sprague Rappaport Stethoscope

The advanced X-configuration acoustic tubing on the MDF® 767X eliminates all tube rubbing noises associated with the traditional Sprague-type stethoscopes and provides superior acoustic performance making it the most advanced Sprague Rappaport available today.

The chestpiece enables users to switch between adult and pediatric-sided bowels/bells with ease and interchange both sides of the chest piece with adult, pediatric and infant fittings (2 diaphragms and 3 bells). The 2-in-1 tubes in the X configuration converge at the middle of the tube and retain a dual lumen and dual input from the chestpiece. This innovation provides superior acoustic quality never available before on a Sprague Rappaport type stethoscope.

Item #	Color	Pkg.
MIF767X04	Navy Blue	ea
MIF767X08	Purple	ea
MIF767X10	Royal Blue	ea
MIF767X17	Burgundy	ea
MIF767XBO	BlackOut (All Black)	ea



MIF767XBO - BlackOut

MIF767X04 - Navy Blue

### MDF® Pulse Time™ Stethoscope

The original, patented MDF® 740 Pulse Time Stethoscope, integrated with a digital timing device, enables medical professionals to maintain eye contact with patients during diagnosis while monitoring elapsed time. U.S.A. PAT. NO. 5252787

This innovative and patented Pulse Time stethoscope has a liquid-crystal analog (LCA) quartz timepiece equipped on top of the integrated bowels-type chestpiece.

The advanced chestpiece is made of lightweight anodized aluminum and the diaphragm is fitted with the unique, raised, ultra-thin fiber diaphragm for greater amplification and transmission of heart and lung sounds.

Item #	Color	Pkg.
MIF74008	Purple	ea
MIF74010	Royal Blue	ea
MIF74011	Black	ea
MIF74017	Burgundy	ea
MIF740BO	BlackOut (All Black)	ea



MIF74010 - Royal Blue




MIF740BO - BlackOut

The Teaching/Training Stethoscope is designed for multi-function teaching applications.

- The dual head design offers the necessary flexibility to use this scope in all clinical and educational conditions.
- Comes complete with two binaurals, two pieces of Y-tubing, and an extension tube.

Item #	Description	Pkg.
MDS9750	Teaching/Training Stethoscope	ea



The Fetal Stethoscope is designed for superior sound sensitivity.

- Chrome-plated inner spring binaural.
- Vinyl coated headrest.
- 22" (56 cm) tubing.


Item #	Description	Pkg.
MDS9775	Fetal Stethoscope	ea



Our Disposable Stethoscope features plastic binaural and a sensitive chest piece that provides excellent sound quality.

- Designed for single patient use.
- Helps reduce cross-contamination in infectious areas.
- High visibility yellow.

Item #	Description	Pkg.
MDS9543	Disposable Stethoscope	ea





# Improving Diabetes Management in Long-Term Care

## Medline's Compass Program A Diabetes Resource for Long-Term Care

Diabetes is a devastating disease affecting approximately 23.6 million people in the U.S. 5.7 million cases continue to go undiagnosed<sup>1</sup> for many years without treatment until severe complications occur. As many as 1 in 4 nursing home residents have diabetes.<sup>2</sup> The time and cost associated with caring for these individuals greatly impact the long term care facilities.

Medline has developed a comprehensive tool to help long-term care facilities provide better care for residents with diabetes by compiling and organizing data acquired during the course of working with healthcare professionals for more than 40 years. "Compass- A Diabetes Resource for Long-Term Care" emphasizes the importance of proper diabetes management to prevent complications and also illustrates how to achieve it.

Considering that heart disease, stroke, high blood pressure, blindness, kidney disease, dental disease, nervous system disease, and even amputations are all complications of diabetes<sup>1</sup>, it is easy to see why healthcare providers strive to provide the best possible care for residents. Additionally, knowing that the Centers for Medicare & Medicaid Services is reviewing Pay for Performance measures and evaluating avoidable hospital admissions, like those due to complications from diabetes, it is vital for nursing home staff to understand current standards of practice.

**With 26 percent of the nursing home population battling diabetes and its complications on a daily basis, Medline is providing educational tools for residents, their families and their caregivers.**





### What's in the Compass Diabetes box?

- One **Director of Nursing (DON) Instruction Manual**—  
A teacher's guide to the program, which includes 4 CE contact hours.
- Three additional **Glucose Monitoring Training** resource books—  
Put them on treatment carts or in nurses' stations for easy access.
- Five additional **Family Education** books—Give them to residents' families.
- One **Educational Compact Disc and** and one **Educational DVD**  
To educate staff on using the OptiumEZ™ Glucose Monitoring System.

A look at a few things from the box:

Resource Books



Instructional Manual



Educational CD and DVD



Family Educational Books

## Medline OptiumEZ Glucose Monitoring System

The Medline OptiumEZ Glucose Monitoring System requires no coding and provides easy, accurate, and reliable results! The OptiumEZ test strips are designed to prevent inaccurate tests and test strip waste caused by interferences from many common substances\* and/or application of an inadequate blood sample volume. The OptiumEZ Glucose Monitoring System is designed to handle all of your blood glucose needs.



No Coding Required



### Easy to Test

- No Coding Required
- Simple 2 Step Testing
- Results in 5 seconds
- Small blood sample size- 0.6 uL

### Easy, Accurate and Reliable Results

- Unique chemistry minimizes the potential interfering effects of many common medications.
- Allows for the ability to reapply more blood on the same test strip.
- Individually wrapped test strips preserve the integrity of the test strips, which ensures an accurate glucose reading.

Item #	Description	Pkg.
MIS71039z	OptiumEZ Blood Glucose Meter	ea
MIS71041	OptiumEZ Glucose Test Strips	12 bx/cs, 50 ea/bx
MIS71041z	OptiumEZ Glucose Test Strips	50 ea/bx
MIS71042	OptiumEZ Glucose Test Strips	6 bx/cs, 100 ea/bx
MIS71042z	OptiumEZ Glucose Test Strips	100 ea/bx
MIS70550z	Precision Control Solution 1 vial high, 1 vial low per box	1 bx

Point of Care Testing

## INRatio2 PT/INR Monitor

#1 selling brand in point of care PT/INR Monitoring for Three Years in a Row\*\*  
The INRatio2 PT/INR testing system combines reliable results with practical convenience making it an excellent solution for anticoagulation patient management.



### 1 Minute Test Time

- 1 button operation with simple icon interface.
- 1 step auto-on with test strip insertion.
- 1 drop of blood from a fingerstick.

### 2 Level On Board Quality Control

Rigorous clot-based and electronic quality control checks are performed automatically with every test, ensuring a reliable result

### Individually-wrapped Strips/No Refrigeration Required

Individually foil-wrapped strips maintain freshness and stability, and can be stored at room temperature.

### Small Sample Size

1 unmeasured drop of blood applied directly to the strip performs both the test and quality controls.

Item #	Description	Pkg.
P-T200431	INRatio2 Professional Kit (includes meter)	ea.
P-T100139Z	INRatio Test Strips	48/bx
P-T0200235A	Blood Collection Tubes	100/bg
MPHSAFETY21Z	Lancet, Normal Flow, 21 Gauge 1.8mm	200/bx

Optium EZ \*  
Data on file, Abbott Diabetes Care Inc.  
INRatio2 \*\*

HPIS comparative sales data, point of care PT/INR monitors sold Q4 2004-Q4 2008.






DIB0310510

## Hemoglobin A1C Test

The A1cNow® InView System offers incredibly fast, quantitative A1C results (in just 5 minutes) utilizing micro-optical technology that is integrated into a proprietary reusable monitor with single use cartridges. Accuracy is assessed by the same standards used by laboratories.

Item #	Description	Pkg.
DIB0310510	A1cNow InView Monitor	10/kit



MSC0002

## Monofilament

This monofilament can be used as a routine test to assess a resident's foot sensitivity.

Item #	Description	Pkg.
MSC0002	Foot Monofilament, Diabetic Neuropathy Test	25/bx




MPHSAFETY28

## Medline Safety Lancets

Medline safety lancets are simple and safe for both the healthcare professional and the patient. The needle retracts automatically and locks out for extra safety, protecting patients and healthcare professionals from needlestick injuries and reducing cross-contamination.

Item #	Description	Pkg.
MPHSAFETY21	Lancet, Normal Flow, 21 Gauge 1.8 mm	18 bx/cs; 200 ea/bx
MPHSAFETY21Z	Lancet, Normal Flow, 21 Gauge 1.8 mm	200 ea/bx
MPHSAFETY28	Lancet, Micro Flow, 28 Gauge 1.6 mm	18 bx/cs; 200 ea/bx
MPHSAFETY28Z	Lancet, Micro Flow, 28 Gauge 1.6 mm	200 ea/bx



OMMAT07 series

OMMAT1004Z

## Unistik® Safety Lancets

Offers maximum comfort for your patients. The Unistik needlepoint is hidden before use and automatically retracts after use to avoid needle stick injuries and help reduce cross-contamination.

Item #	Description	Pkg.
OMMAT0747Z	Unistik 2 Gentle	200/bx
OMMAT0704Z	Unistik 2 Normal Depth	200/bx
OMMAT0704	Unistik 2 Normal Depth	6000/cs
OMMAT0714Z	Unistik 2 Extra Depth	200/bx
OMMAT1004Z	Unistik 3 Normal Depth	200/bx

## Lancets

Ultra Sharp Lancets have fine-point tip for smooth and fast procedures. Choose from our general lancets in 100 or 200 count. These lancets can be used with all general lancing devices.

Item #	Description	Pkg.
MPH150515	Lancet, 23 Gauge Needle	2400/cs
MPH150515Z	Lancet, 23 Gauge Needle	200/bx
MPH150514	Lancet, 26 Gauge Needle	2400/cs
MPH150514Z	Lancet, 26 Gauge Needle	100/bx
MPH281001DME	Lancet, 28 Gauge Needle	20,000/cs
MPH281001DMEZ	Lancet, 28 Gauge Needle	100/bx



MPH281001DMEZ



MDS705153 Shown with SWD85161H cabinet

## Medline has a variety of sharps disposal systems to fit your needs.

And if you're looking to cut costs, you'll find savings, too, without having to sacrifice quality. Medline Sharps Containers are among the safest containers on the market. Easy-to-see fill levels reduce overfilling and accidental needle sticks. Containers lock easily for final disposal and safe transportation. Specially designed lids limit unwanted access to container contents.

### Patient Room Sharps Containers

These sharps disposal containers are designed for use in patient rooms, on med carts and in high traffic areas.

- Lid design facilitates one-handed operation
- Horizontal drop design for maximum capacity to fill line
- Easy-to-see fill levels
- Large openings accommodate larger syringes
- Multiple cabinet and bracket options

Item #	Description
MDS705125	2.5 Quart (2.4 L), 7.5" H x 4.25" L x 11" W, Wall Mount, 30/cs
MDS705153	5 Quart (4.7 L), 11" H x 4.75" L x 10.75" W Wall Mount, 20/cs
PLA141020	5 Quart (4.7 L), 10.75" H x 4.75" L x 10.5" W Wall Mount/Counter-Balanced Lid, 20/cs
MDS705151	5 Quart (4.7 L), 14" H x 6" L x 13" W Wall Mount/Counter-Balanced Lid, 14/cs
MDS705203	3 Gallon (11.4 L), 20.5" H x 6" L x 14" W Wall Mount/Counter-Balanced Lid, 12/cs

### Phlebotomy Sharps Containers

Portable sharps containers are designed to fit into blood drawing trays for point-of-use disposal. Square base for added stability.

- Translucent top allows easy identification of contents
- Choice of single or dual openings to accommodate a variety of sharps
- Containers lock for final disposal and safe transport
- Brackets and holders available

MDS705110	1 Quart Phlebotomy Container (0.9 L), 6.25" H x 4.25" L x 4.5" W (159 x 108 x 114 mm)
MDS705115	1.5 Quart (1.4 L) Transportable Container 10.5" H x 3.5" L x 3.5" W (267 x 89 x 89 mm)



MDS705151

MDS705203



MDS705115

MDS705110 shown with Bracket MDS707912



**MDS705210**



**MDS705201**



**MDS705202**

## Multipurpose Sharps Containers

These containers are designed for use in restricted access areas.

- Containers lock for final disposal
- Nestable containers save space and reduce shipping costs
- Non-tortuous lid design accommodates a variety of sharps sizes
- Horizontal drop maximizes use of container volume
- Brackets available

### Sharps / Biohazard

MDS705201	1 Gallon (3.8 L) 6.25" H x 7" L x 10.5" W, 32/cs
MDS705202	2 Gallon (7.6 L) 10.25" H x 7" L x 10.5" W, 20/cs
MDS705210	10 Gallon (37.9 L), 21.5" H x 12" L x 15.5" W, 6/cs

### Chemotherapy

MDS706202	2 Gallon (7.6 L) 10" H x 7.25" L x 10.5" W, 20/cs
-----------	--

## Mail-Back Sharps Containers

Each kit includes all components required for collection and transportation of sharps and regulated medical waste.

- Compact design uses a minimum of counter space
- Pre-paid, pre-addressed shipping labels and pre-paid disposal
- Lids provide temporary and final locking closure
- Easy-to-use enclosed outer shipping box
- Step-by-step instructions
- Tracking documents

Item #	Description
SCC1QV2	1 Quart (0.9 L), 1/ea
MDS705000	1 Gallon (3.8 L), 1/ea
SCC2GV2	2 Gallon (7.6 L), 1/ea
SCC3GV2	3 Gallon (11.4 L), 1/ea
SCC5GV2	5 Gallon (18.9 L), 1/ea
SCC18GV2	18 Gallon (68.1 L), 1/ea



## Sharps Containers and Accessories at-a-Glance

Item #	Description	Locking Bracket	Non-Locking Bracket	Cabinet	Glove Box for Cabinet	Holder	Cart
MDS705110	1 Quart Phlebotomy, Red		MDS707912			MDS707910	
MDS705115	1.5 Quart Transportable, Red						
MDS705125	2.5 Quart, Wall Mount, Red	MDS707960		MDS707549	SWD8550LG and SWD8554B		
MDS705153	5 Quart, Wall Mount, Red	SWD8518X	SWD8519C	SWD85165H or SWD85161H	SWD8555SA		
PLA141020	5 Quart (4.7 L), Wall Mount/ Counter-Balanced Lid, Red			SWD85165H or SWD85161H	SWD8555SA		
MDS705151	5 Quart, Wall Mount, Counter-balanced Lid, Red	MDS707054	MDS707062	MDS707553	SWD307013 or SWD320032	MDS707955	
MDS705203	3 Gallon, Wall Mount, Counter-balanced Lid, Red	MDS707054	MDS707062	MDS707005	SWD307013 or SWD320032	MDS707016	
MDS705201	1 Gallon, Sliding Lid, Red	MDS707963	MDS707975			MDS707974	
MDS705202	2 Gallon, Sliding Lid, Red	MDS707963	MDS707975			MDS707974	
MDS705210	10 Gallon, Hinged Lid, Red						MDS707109
MDS706202	2 Gallon Chemotherapy, Hinged Lid, Yellow	MDS707963	MDS707975			MDS707974	
MDS706318	18 Gallon, Chemotherapy, Port Lid, Yellow						MDS705218CART


Sharps Containers

### Medication Disposal System

This is a complete disposal solution for unused medications. The medication-disposal system includes the cost of the medication-disposal container, return freight and total destruction costs. The unused medications can be over-the-counter (OTCs) drugs or prescription drugs, so long as the prescription drugs are not controlled substances.

Item #	Description/Pkg.
MDS707001	Medication Disposal, Envelope 8" x 11", Ea
MDS707010	Medication Disposal, Container 10 Gallon, Ea
MDS707020	Medication Disposal, Container 20 Gallon, Ea





**SWD511110Z**

## Safety Insulin Syringes

These safety syringes have a secure safety shield that covers the entire needle length.

Item #	Description/Pkg.
SWD511110Z	1 mL syringe, 29 gauge x 1/2" needle, 100/bx
SWD5111110	1 mL syringe, 29 gauge x 1/2" needle, 500/cs
SWD511136Z	1/2 mL syringe, 29 gauge x 1/2" needle, 100/bx
SWD511136	1/2 mL syringe, 29 gauge x 1/2" needle, 500/cs
SWD511144Z	3/10 mL syringe, 29 gauge x 1/2" needle, 100/bx
SWD511144	3/10 mL syringe, 29 gauge x 1/2" needle, 500/cs

## Medline's Tympanic Thermometer with Quick Probe Release

Get fast, accurate and easy-to-read results with Medline's **Tympanic Thermometer with Quick Probe Release**. This ear thermometer is designed as a comfortable and convenient method for taking a person's body temperature. No embarrassing or invasive procedures; you simply place the thermometer in the ear to get a reading within seconds. It's perfectly suited for both home and clinical use on patients of all ages. Plus, its unique mechanical probe release helps reduce cross-contamination, allowing you to easily change probe covers with every use!

### Why an ear thermometer?

It turns out that the eardrum is an extremely accurate point to measure body temperature from because it is recessed inside the head (just like your tongue). Instead of touching the eardrum with a thermometer, the tympanic thermometer uses an infrared sensor to measure the temperature of energy radiating from the eardrum.

### Innovative, quick-response ear thermometer features:


- ✓ One-touch operation
- ✓ Mechanical probe cover quick-release allows easy probe cover change
- ✓ Includes 100 probe covers
- ✓ Safe to use on patients of any age
- ✓ Readings in Fahrenheit and Celsius
- ✓ Takes an accurate, non-invasive temperature from the tympanic membrane

Item #	Description	Pkg.
MDS9700	Thermometer	ea
MDS9701	Probe covers	100/bx



*Apply disposable cover onto the tympanic thermometer.*

*Turn device "on" by pressing button on front.*



*Place thermometer in ear and depress top button to take reading.*

*After reading, flip the blue tab to release the disposable probe cover without contact.*



SunTech 247 Patient Monitor  
**An Expandable Diagnostic Station with BP, Temperature and SpO2**

This patient monitor was designed to measure clinical grade vital signs, while allowing the unit to be customized to the facility's needs. The 247 Patient Monitor is completely modular, meaning that you can purchase the functionalities you need today and then add additional parameters, if needed. In addition, it combines the efficiency of an automated monitor with the convenience of a manual mode that allows the clinician to take a BP reading with a stethoscope whenever needed!

Item #	Description	Pkg.
SNT99006200	BP with Mobile Stand	1/ea
SNT99006300	BP & Temperature with Mobile Stand	1/ea
SNT99006400	BP & SpO2 with Mobile Stand	1/ea
SNT99006100	BP, SpO2, Temperature with Mobile Stand	1/ea
<b>Wall Mount System</b>		
SNT99006600	BP	1/ea
SNT99006700	BP & Temperature	1/ea
SNT99006800	BP & SpO2	1/ea
SNT99006500	BP, Temperature & SpO2	1/ea
<b>Table Top System with Battery</b>		
SNT99007000	BP	1/ea
SNT99007100	BP & Temperature	1/ea
SNT99007200	BP & SpO2	1/ea
SNT99006900	BP, Temperature & SpO2	1/ea
<b>Modules</b>		
SNT98012800	Temperature Module	1/ea
SNT98012900	SpO2 Module	1/ea

**Accessories Included**

**BP:** 8' (2.5 meters) blood pressure hose, adult & large adult all purpose cuffs, mains/power supply, mains/power cable, resource cd, quick start guide, reference card

**BP + Temperature:** oral/axillary probe (blue), 1 box of Filac Fastemp™ probe covers

**BP + SpO2:** reusable adult finger sensor, 6' (1.8 meters) extension cable

**Additional Accessories (not included)**

- SWD202020Z Probe Covers for the FastTemp Thermometer
- STCUS1320DM Child Replacement Cuff, 5/bx
- STCUS2635DM Adult Replacement Cuff, 5/bx
- STCUS3242DM Large Adult Replacement Cuff, 5/bx
- STCUS4250DM Thigh Replacement Cuff, 5/bx

## Welch Allyn SureTemp® Model 690

The SureTemp Plus 690 Thermometer delivers oral results in just four seconds! This electronic thermometer is easy to use, compact and dependable. It delivers oral results in four seconds, axillary temperatures (for children under four) in 10 seconds and rectal temperatures in 10 seconds with the use of an optional, interchangeable probe. Large LCD display features Fahrenheit or Celsius readouts and last temperature recall. This affordable thermometer is ideal for alternate care environments. Three 1.5V AA batteries provide 6000 readings. Automatic probe cover ejection; Two-year limited warranty.

Item #	Description	Pkg.
W-A01690200	Oral/Axillary Thermometer	ea
W-A02892000	Rectal Probe	ea
W-A5031101	Disposable Probe Covers	1000/cs
W-A5031101Z	Disposable Probe Covers	250/bx



W-A01690200

## Thermoscan® Thermometers

The NEW Braun® Thermoscan Pro 4000 Thermometer patented technology is an ergonomically designed ear thermometer engineered to take temperatures easily and accurately with the click of a button.

- Exclusive ExacTemp® technology provides an active user feedback light that signals that the probe has been positioned in a stable way throughout the measurement, avoiding inaccurate measurements by handling errors
- Preheated probe tip minimizes the cool down effect experienced when taking repeat temperatures, improving accuracy and repeatability
- Automatic probe cover eject button for quick and easy removal of used probe covers
- Built-in cross-contamination safeguards – probe cover sensor reminds operator to use a new probe cover
- Wide range of patient temperature – 68°F-108°F
- Wide range of ambient temperature – 50°F-104°F ambient range allows for use in most environments
- Durable – tested for damage resistance
- Convenient wall-holder for Pro 4000 unit w/probe cover storage is available
- Three-year limited warranty coverage

Item #	Description	Pkg.
TMS40002	Thermoscan PRO 4000	ea
TMSPC800	Probe Covers	40 bx/cs



TMS40002

## Exergen® Temporal Scanner™ Thermometer

Precision infrared scanner provides a truly non-invasive method for measuring temperature. Simply scan across the temporal artery area of the forehead, touch behind the earlobe, and an accurate reproducible temperature is measured in about three seconds.

- Convenient one-handed operation
- Completely non-invasive
- Stainless steel probe for highest durability
- Easy snap-on probe caps
- Three-year warranty

Item #	Description	Pkg.
EXGTAT5000	Infrared Temporal Artery Thermometer	ea
EXGTAT2000	Infrared Temporal Artery Thermometer	ea
EXG134203	Disposable Probe Caps	1000 ea/cs



EXGTAT5000



### Dual Scan Modes

Forehead and object/liquid

### Easy-To-Read

Nite-Glo® display

### Quick. Accurate. Reliable

#### Infrared Thermometer with Digital Readout

- Professional accuracy
  - Fast 1-second readout
  - Wide range readout
    - 50°F-122°F person
    - 32°F-221°F object/liquid
  - Nite-Glo® display
  - Fever alarm
  - 10-reading memory recall
  - °F or °C switchable
  - Automatic shut-off
- Includes:
- Convenient storage bag
  - Detailed guidebook in English and Spanish
  - 2 AAA batteries
  - 5-year limited warranty

MBH18535000 Infrared Thermometer, each



### Get Both Speed And Accuracy with The FILAC® 3000 Features Both!

The FILAC 3000 Electronic Thermometer utilizes advanced technology to ensure the most accurate and fastest readings are provided in easy-to-use formats. These cutting edge thermometers address the need for fast, accurate readings in under 4 seconds or an advanced accuracy reading in under 6 seconds. Compact, user friendly design will stand upright on a flat surface; intuitive icon-based interface system with backlight LCD display; unique progress bar indicates status of temperature measurement; last temperature recall function. Patented isolation chamber design and color schemes along with the durable probe covers help minimize the chance of cross-contamination and enhance infection control measures. Operates on 4 AA batteries for 6000 measurements. Equipped with electronic anti-theft system, self-calibrating electronics and auto on-off function. Includes oral probe and push button probe cover ejection.

Item #	Description Pkg.
SWD504000	Filac 3000 Thermometer, each
SWD202020	Probe Covers for the Filac Fastemp Thermometer & 3000, 5,000/cs
SWD202020Z	Probe Covers for the Filac, Fastemp Thermometer & 3000, 500/box



### Genius2 Thermometer

Measures temperature from the tympanic membrane, which shares the same blood flow as the hypothalamus, the body's natural thermostat.

- Utilizes infrared technology to derive the most accurate body temperature from the tympanic membrane.
- Records and analyzes up to 100 temperature measurements, and selects the highest and most accurate in less than 2 seconds, as compared to 32 in Genius and 1-8 measurements in competitive thermometers.
- Reduces the chance of an inaccurate reading due to improper technique or patient movement.
- Durable tip-based sensor technology Improves durability and accuracy, minimizing the need for recalibration.
- Tapered probe design
- Rapid response time 1 - 2 second
- Real time temperatures
- Programmable display
- Reduces cross-contamination
  - No contact with mucus membranes
  - No-touch probe cover application and release

Item #	Description/Pkg.
SWD303000	Genius2 Thermometer each
SWD303030Z	Probe Covers 96 ea/bx
SWD303030	Probe Covers 2,112/case



### Six Seconds\*, Temple Digital Thermometer

The NEW Medline Instant Read\* Digital Temple Thermometer accurately measures body heat using Patented R.A.T.E.™ technology (Rapid, Accurate, Temperature, Establishment).

- Sensors rapidly track heat flow generated from blood vessels to the skins surface and convert measurements to body temperature.
- Easy-to-use while asleep or awake.
- Meets accuracy standards.
- Large digital display.
- Switchable from Fahrenheit to Celsius.
- Beeps when ready.
- Memory stores previous reading.

Item #	Description	Pkg.
MDS9698	Temple Thermometer	ea
MDS9699	Probe Covers	250ea/bx



Gentle, non-invasive sensor accurately reads in F° & C°

Power/memory control button

Large Digital Display

2-AAA Batteries included



Features Patented R.A.T.E. Technology

\*Readings in six to ten seconds depending on individual.

### FeverScan™ Disposable Forehead Thermometers

FeverScan™ products utilize thermo chromic liquid crystal (TLS) technology for the measurement of clinical temperature trends in the healthcare environment. Liquid crystals are chemical compounds and mixtures that exhibit the mechanical properties of liquids and the optical characteristics of solids.

FeverScan™ thermometers are composed of a series of discrete temperature events, each event shows a green color as the indicated temperature is attained.

#### Anesthesia & PACU (with adhesive backing)

Item #	Description	Pkg.
ATC430	°F 94°F-104°F	100/dispenser
ATC431	°C 35°C-40°C	100/dispenser

#### Hypo/Hyperthermic Monitoring (with adhesive backing)

ATC434	Extended Range °F, 84°F-106°F	100/dispenser
ATC435	Extended Range °C, 30°C-41°C	100/dispenser

#### Reusable (no adhesive)

ATC432	°F 94°F-104°F	100/dispenser
ATC433	°C 35°C-40°C	100/dispenser



- Continuous Display
- Forehead Temperature Monitor
- Anesthesia through PACU
- Watch for Hypo/Hyperthermia
- Pre and Post-Op Surveillance
- Disposable, Latex-Free
- Core Adjusted Temperature



New!

MDS96801



MDS9650



MDS9628



MDS9654B



MDS9629

### Digital Thermometers


Our digital thermometers provide accurate temperatures in about 60 sec.

- Large LCD screen for easy readings.
- Simple push-button start.
- Disposable sheaths allow for sanitary oral, rectal, or axillary use.

Item #	Description	Pkg.
MDS9650	Premier Oral °F & °C	ea
MDS9652	Premier Rectal °F	ea
MDS9654B	Standard Oral °F	200/cs
MDS9628	Oral Kit w/20 Sheaths °F	ea
MDS9629	Rectal Kit w/20 Sheaths °F	ea
MDS9607	Oral Sheaths	100/bx
MDS9608	Rectal Sheaths	100/bx
MDS96801	Flex-tip, waterproof, °C or °F	ea
MDS9655	Standard, Oral, Bulk, °F & °C	200/cs
MDS9655Z	Standard, Oral, °F & °C	20/bx
MDS9650B	Premier Oral, Bulk, °F & °C	200/cs
MDS9656	Standard, Rectal, Bulk, °F & °C	200/cs
MDS9656Z	Standard, Rectal, °F & °C	20/bx



Patient Monitoring



### Disposable Penlights

Lightweight and disposable, these high-intensity exam lights are perfect for the physician, nurse, and EMT. They're constructed of high-impact molded styrene and can be clipped on pockets.

Item #	Description	Pkg.
MDS131040	Disposable Penlight, White	6/pk

## AED External Defibrillator

Cardiac Science® PowerHeart® AED G3 Automated External Defibrillator with voice prompts and text screen displays guide you through the rescue process. When cardiac arrest occurs, every second counts. The advantage of the PowerHeart AED G3 is that it may be used by anyone who has received basic training on its operation. Electrodes are pre-connected and disposable. Semi-automatic defibrillator with biphasic variable waveform delivers an effective shock in less than 55 seconds. Daily, weekly, and monthly automated self-checks ensure that PowerHeart is always ready to go. Internal memory holds 60 minutes of ECG data. Includes one pair of defibrillation electrodes, instructional CD-ROM and training video, RescueLink® and MDLink® software and serial communication cable. Seven-year warranty on device & four-year full replacement on battery.

Item #	Description	Pkg.
AED9300SEMICC	Semi Automatic AED Kit	ea
AED9300AUTOCC	Fully automatic AED Kit	ea
Includes:	1 Replacement set of adult defibrillator electrode pads	
	1 Soft carrying case	
	1 Ready Kit	



Powerheart G3 AED

### References:

1. National Diabetes Fact Sheet, 2007. Centers for Disease Control and Prevention website. Available at [www.cdc.gov/diabetes/pubs/pdf/ndfs\\_2007.pdf](http://www.cdc.gov/diabetes/pubs/pdf/ndfs_2007.pdf). Accessed on April 24, 2009
2. Johnson EL, Brosseau JD, Soule M, et al. Treatment of diabetes in long-term care facilities: a primary care approach. *Clinical Diabetes*. 2008;26(4):152-156

AlcNow is a registered trademark of Metrika, Inc.  
 Braun is a registered trademark of Braun GmbH Corporation  
 Cardiac Science is a registered trademark of Cardiac Science Corporation  
 ExacTemp is a registered trademark of Braun GmbH Corporation  
 Exergen is a registered trademark and Temporal Scanner is a trademark of Exergen Corp.  
 Fastemp is a trademark of Westek Corporation  
 Feverscan is a trademark of LCR Hallcrest  
 FILAC is a registered trademark of Covidien AG Corp.  
 INRatio is a registered trademark of HemoSense, Inc.  
 MDLink is a registered trademark of Cardiac Science Corporation  
 MDF is a registered trademark and PulseTime is a trademark of MDF Instruments USA, Inc.  
 Medline is a registered trademark of Medline Industries, Inc.  
 Nite-Glo is a registered trademark of MABIS Healthcare, Inc.  
 OptiumEZ is a trademark of Abbott Laboratories  
 PowerHeart is a registered trademark of Cardiac Science Corp.  
 R.A.T.E. is a trademark of BestMed  
 RescueLink is a registered trademark of Cardiac Science Corporation  
 Simple Solutions. Significant Results. is a trademark of Medline Industries, Inc.  
 Suntech is a registered trademark of Suntech Medical  
 Thermoscan is a registered trademark of Braun GmbH Corp.  
 Tyvek is a registered trademark of E. I. Du Pont De Nemours & Co.  
 Unistik is a registered trademark of Owen Mumford Ltd.  
 Welch Allyn and SureTemp are registered trademarks of Welch Allyn, Inc.

# BioCon-500 Scanner safely measures bladder volume.

## The BioCon-500 Scanner is brought to you exclusively by Medline

In our commitment to bring you high quality products at a good value, Medline has become the exclusive distributor of the BioCon-500—the only portable, 3-D ultrasound device with a pre-scan feature. The BioCon-500 is a quick, safe and non-invasive way to measure bladder volume and post-void urinary retention.

## Help Save time, Increase Accuracy with Our Bladder Scanner

Incontinence is widespread among the elderly population, especially those who live in nursing homes, where 50 percent or more of residents are incontinent.<sup>1</sup> The BioCon-500 can be helpful in developing individualized urinary continence programs.

As part of a urinary continence program, the BioCon-500 Bladder Scanner allows caregivers to:

- Scan the resident's bladder to determine if catheterization or toileting is necessary.
- Show residents the results so they can actually see how much fluid is in their bladder. This helps residents re-learn the sensation of needing to urinate.
- Develop individualized toileting schedules based on each resident's bladder function pattern rather than at fixed times, such as every two hours. This will save staff valuable time.
- More accurately determine if a resident's bladder volume and voided amounts are consistently low, prompting an increase in the resident's fluid intake.



Item #	Description	Pkg.
MDSB10500	BioCon-500 Bladder Scanner	each
MDS3110100032	Mobile Stand	each
MDS092005	Ultrasound Gel	12 ea/cs

Medline United States  
1-800-MEDLINE (633-5463)

www.medline.com | info@medline.com

Medline Canada  
1-800-396-6996

www.medline.ca | canada@medline.com

Medline México  
01-800-831-0898

www.medlinemexico.com | mexico@medline.com



©2010 Medline Industries, Inc. One Medline Place, Mundelein, IL 60060  
Medline is a registered trademark of Medline Industries, Inc.  
MKT210058/LIT007A/10M/SG/2