



Protective Apparel

Disposable All-Around Protection





Why are more hospitals and purchasing groups choosing Medline for their all-around protection?

Our products provide the highest protection with the greatest comfort at the lowest cost.

Medline is a leader in protective apparel with the broadest product line and low acquisition costs. If high costs are plaguing your facility, Medline can offer you solutions to bring those costs under control without sacrificing the protection of your caregivers.

Our protective apparel line includes generously sized footwear, protective headwear, gowns that meet AAMI PB70 requirements, protective multi-purpose gowns, a full range of protective coveralls and comfortable, colorful lab coats and lab jackets.

We know you rely on disposable apparel to keep your staff protected and comfortable. For this reason, we use the finest materials available for protection, offering the greatest level of comfort. We apply the highest quality standards to our products and in our manufacturing facilities.


 **Medline is ISO 13485-2003 certified. Medline has been manufacturing protective apparel for over 40 years.**

Table of Contents

2

- LEVELS OF BARRIER PROTECTION 4
- AAMI ISOLATION GOWNS 5-6
- ISOLATION GOWNS 6-8
- PROTECTIVE GOWNS 9-12
- PROCEDURE & PROTECTIVE GOWNS 12
- COVERALLS 13
- FOOTWEAR 14-15
- HEADWEAR 16-17
- LAB JACKETS AND COATS 18
- SCRUBS, APRONS AND SLEEVES 19

We created easy to use icons to help identify location of use based on our recommendations.

- Lab 
- Visitors 
- Intensive Care Unit 
- Food Services 
- Decontamination 
- Operating Room 
- Radiology 
- Maternity 
- Emergency 
- Nursing 
- Airborne Precautions 
- Maintenance 
- Contact Precautions 
- Fluid Resistance 

Guide to Protective Materials

To eliminate confusion, here's a brief description of the main materials we use to produce Medline protective apparel. Refer back to this page as you look through the brochure.

Spunbond Polypropylene

The standard fabric for basic protective apparel; formed by bonding fibers together to form a single layer of breathable, woven-like material. Its main advantages are economy and comfort, but it also provides light fluid protection.

Polypropylene Protection and Comfort Levels			
	Good	Better	Best
Protection			
Comfort			

Multi-Ply Material

This fabric consists of inner layers of meltblown polypropylene sandwiched between outer layers of spunbond polypropylene. Strong, fluid-resistant and breathable, this multi-layer material provides an optimal mix of protection and comfort.

Appropriate for extended wear situations involving low to moderate fluid contact. Available in lightweight and medium weight options.

Medium Weight Multi-Ply Protection and Comfort Levels			
	Good	Better	Best
Protection			
Comfort			

Lightweight Multi-Ply Protection and Comfort Levels			
	Good	Better	Best
Protection			
Comfort			

Coated Polypropylene

This fabric is created by coating soft, spunbond polypropylene with a layer of impervious polyethylene. This material is ideal for any situation demanding a high level of fluid protection. Available in medium weight and heavy weight materials.

Coated Protection and Comfort Levels			
	Good	Better	Best
Protection			
Comfort			

Breathable Laminate

Microporous film bonded to nonwoven material provides fluid-proof protection and great breathability. This material is ideal for extended wear conditions where fluid contact is expected.

Breathable Microporous Laminate Protection and Comfort Levels			
	Good	Better	Best
Protection			
Comfort			

Polyethylene

This single layer of polyethylene film provides tremendous barrier properties. This material is ideal for situations when high levels of fluids are anticipated.

Polyethylene Protection and Comfort Levels			
	Good	Better	Best
Protection			
Comfort			

Prevention Plus

This impervious and breathable material consists of a hybrid film laminated between two layers of bi-component spunbond material. This fabric is ideal when impervious protection is needed.

Prevention Plus Protection and Comfort Levels			
	Good	Better	Best
Protection			
Comfort			



Levels of Barrier Protection

(from AAMI PB70:2003)

Medline isolation gowns with color-coded neck bindings meet the test requirements set forth by the Association for the Advancement of Medical Instrumentation (AAMI) PB70.

The chart below outlines AAMI PB70 guidelines. The color-coding has been developed by Medline to differentiate each level and corresponds with the color-coded neck bindings on Medline's AAMI Compliant Isolation Gowns.

ANSI / AAMI PB70 Barrier Performance	Test Method Water Resistance	Test Definition AATCC	Requirement
Level 1	AATCC 42 Impact Penetration	AATCC 42 <i>Measures the resistance of fabrics to the liquid penetration of water by impact.</i>	Water Impact ≤ 4.5 g
Level 2	AATCC 42 Impact Penetration AATCC 127 Hydrostatic Pressure	AATCC 42 <i>Measures the resistance of fabrics to the liquid penetration of water by impact.</i> AATCC 127 <i>Measures the resistance of fabrics to the liquid penetration of water by impact under constant and increasing hydrostatic pressure.</i>	Spray Impact ≤ 1.0 g Hydrostatic Pressure ≥ 20 cm
Level 3	AATCC 42 Impact Penetration AATCC 127 Hydrostatic Pressure	AATCC 42 <i>Measures the resistance of fabrics to the liquid penetration of water by impact.</i> AATCC 127 <i>Measures the resistance of fabrics to the liquid penetration of water by impact under constant and increasing hydrostatic pressure.</i>	Spray Impact ≤ 1.0 g Hydrostatic Pressure ≥ 50 cm

Classic Yellow Isolation Gowns

Breathable lightweight multi-ply material maximizes comfort



Item No.	Size	Pkg
Elastic Wrists, Tape Tab Neck, Flat Pack		
NONLV100	Regular/Large	10/bg, 50/cs
NONLV100XL	X-Large	10/bg, 50/cs

Classic Yellow Medium Weight Multi-Ply Gowns

Comfortable, flexible and strong enough for tough situations



Item No.	Size	Pkg
Elastic Wrists, Tape Tab Neck, Flat Pack		
NONLV200M	Medium	10/bg, 100/cs
NONLV200	Regular/Large	10/bg, 100/cs
NONLV200XL	X-Large	10/bg, 100/cs
NONLV200XXL	XX-Large	10/bg, 100/cs
Elastic Wrists, Tape Tab Neck, Individually Folded-Dispensing box		
NONLV200F	Regular/Large	8/bx, 40/cs
NONLV200FXL	X-Large	8/bx, 40/cs
Thumb Loop Wrists, Tape Tab Neck, Individually Folded-Dispensing Box		
NONLV200TL	Regular/Large	8/bx, 40/cs
NONLV200TLXL	X-Large	8/bx, 40/cs

Thumb Loop Wrists, Over-the-Head Neck, Side Ties, Flat Pack		
NONLV240*	Regular/Large	10/bg, 100/cs
NONLV240XL**	X-Large	10/bg, 100/cs



* NONLV240 features plain neckline

** Size X-Large Color Identification (NONLV240 series only)



NONLV200



NONLV240



NONLV350



Premium Yellow Multi-Ply Gowns

Comfortable and breathable gown provides a high level of protection against fluids

Level 3 Protection

Item No.	Size	Pkg
Elastic Wrists, Tape Tab Neck, Flat Pack		
NONLV315	Regular/Large	10/bg, 100/cs
NONLV315XL	XL	10/bg, 100/cs



Premium Blue Barrier Protection

Level 3 Protection

Item No.	Size	Pkg
Knit Cuffs, Hook and Loop Neck, Flat Pack		
NONLV325	Regular/Large	10/bg, 50/cs
NONLV325XL	XL	10/bg, 50/cs



Microporous Breathable White Gowns

Ideal for heavy fluids where comfort is needed and protection is expected

Level 3 Protection

Item No.	Size	Pkg
Knit Cuffs, Tape Tab Neck, Flat Pack		
NONLV350	Regular/Large	10/bg, 50/cs
NONLV350XL	X-Large	10/bg, 50/cs



NON27SMS2

Our most popular isolation gown

Medium Weight, Classic Style Neck and Waist Ties with Elastic Wrists

- Made from medium weight fluid-resistant multi-ply material
- Generously sized for full coverage and flexibility
- Waist and neck ties ensure a secure fit
- Available in yellow or blue

Item No.	Size	Pkg
Yellow, Neck and Waist Ties, Elastic Wrists, Flat Pack		
NON27SMS2	Regular/Large	10/bg, 100/cs
NON27SMS2XL	X-Large	10/bg, 100/cs
Blue, Neck and Waist Ties, Elastic Wrists, Flat Pack		
NON27SMS3	Regular/Large	10/bg, 100/cs
NON27SMS3XL	X-Large	10/bg, 100/cs



Medline heard what you had to say about the importance of easy donning and simple removal in facilities today. We are pleased to introduce the newest additions to our isolation gown line. Over-the-head neck styles and thumb loop wrists have become a staple in effortless application; we have added new features based on your feedback. These new styles will help promote simple and safe compliance among the staff in your facilities.

Medium Weight, Over-the-Head Style Gown with Thumb-Loop Wrists

- Made from medium weight, AAMI level 2 barrier material
- Features over-the-head neck style, open back
- Thumb loop wrists feature elastic for a good gown and glove fit



Item No.	Size	Pkg
Yellow, Over-the-Head Neck, Thumb-Loop Wrists, Flat Pack		
NON27SMS5	Regular/Large	10/bg, 100/cs
NON27SMS5XL	X-Large	10/bg, 100/cs



NON27SMS5

Medium Weight, Over-the-Head Style Gown, Full Coverage with Stretch Sizing

- Made from medium weight AAMI Level 2 barrier material
- Features over-the-head neck, thumb loop wrists, and full back coverage
- Sides feature elastic gathering to provide increased sizing



Item No.	Size	Pkg
Yellow, Over-the-Head Neck, Thumb Loop Wrists, Side Gatherings, Flat Pack		
NON27550	Regular/Large	10/bg, 100/cs
NON27550XL	X-Large	10/bg, 100/cs



NON27550



NON27236

Lightweight Fluid-Resistant Multi-Ply Gowns

- Made from lightweight fluid-resistant multi-ply material
- Available with elastic wrists or knit cuffs
- Generously sized for full coverage and flexibility
- Waist and neck ties ensure a secure fit



Item No.	Size	Pkg
Yellow, Neck and Waist Ties, Elastic Wrists, Flat Pack		
NON27236	Regular/Large	10/bg, 50/cs
NON27236XL	X-Large	10/bg, 50/cs
Yellow, Neck and Waist Ties, Knit Cuffs, Flat Pack		
NON27239Y	Regular/Large	10/bg, 50/cs
NON27239YXL	X-Large	10/bg, 50/cs



CRI4000

Basic Polypropylene Gowns

- Made from lightweight, fluid-repelling polypropylene material
- Generously sized for full coverage and flexibility
- Styles available in classic neck and waist ties with elastic wrists or over-the-head with thumb loop wrists



Item No.	Size	Color	Pkg
Classic Neck and Waist Ties, Elastic Wrists, Flat Pack			
CRI4000	Regular/Large	Yellow	10/bg, 50/cs
CRI4001	X-Large	Yellow	10/bg, 50/cs
CRI4000B	Regular/Large	Blue	10/bg, 50/cs
CRI4001G	X-Large	Green	10/bg, 50/cs
Over-The-Head Neck and Thumb-Loop Wrists, Flat Pack			
CRI4200	Regular/Large	Yellow	10/bg, 50/cs
CRI4201	X-Large	Yellow	10/bg, 50/cs

Standard Yellow Polyethylene-Coated Polypropylene Gowns

- Made from medium weight yellow polyethylene-coated polypropylene material
- High level of fluid protection
- Waist and neck ties ensure gown remains securely in place
- Closed back gowns are available with elastic wrists or knit cuffs



Item No.	Size	Pkg
Elastic Cuffs, Flat Pack		
CRI4010	Regular/Large	10/bg, 50/cs
CRI4011	X-Large	10/bg, 50/cs
Knit Cuffs, Flat Pack		
CRI4010CF	Regular/Large	10/bg, 50/cs
CRI4011CF	X-Large	10/bg, 50/cs
CRI4012CF	XX-Large	10/bg, 50/cs

CRI4010



CRI5000

Stay Protected with Standard Polyethylene Thumbs Up® gowns

- Made from blue impervious polyethylene material
- Material passed ASTM1670 and ASTM1671 test methods
- Apron-style neck enables gown to be placed on and off quickly
- Thumb loop style wrists keep sleeves from sliding
- Individually folded and wrapped gowns in dispenser box



Item No.	Size	Pkg
CRI5000	Regular/Large	15/bx, 75/cs
CRI5001	X-Large	15/bx, 75/cs
CRI5003	XXX-Large	15/bx, 75/cs



NON27114



NON27279OHTL



Back View

Classic Polyethylene-Coated Polypropylene Gowns feature 'Closed Back'

Offers excellent fluid protection for worry-free job performance

- Made from blue polyethylene-coated polypropylene material
- Generously sized for full coverage and flexibility
- Designed with classic neck and waist ties
- Available with elastic wrists or knit cuffs



Item No.	Size	Pkg
Elastic Cuffs, Flat Pack		
NON27114	Regular/Large	10/bg, 50/cs
NON27114XL	X-Large	10/bg, 50/cs
Knit Cuffs, Flat Pack		
NON27116	Regular/Large	10/bg, 50/cs
NON27116XL	X-Large	10/bg, 50/cs
NON27116XXL	XX-Large	10/bg, 50/cs

Classic Polyethylene-Coated Polypropylene 'Open Back' Gowns

High fluid protection that keeps you comfortable longer

- Made from blue polyethylene-coated polypropylene material
- Generously sized for full coverage and flexibility
- Designed with classic neck ties or over-the-head style neck
- Available with knit cuffs or thumb loop wrists



Item No.	Size	Pkg
Knit Cuffs, Flat Pack		
NON27279	Regular/Large	10/bg, 50/cs
NON27279XL	X-Large	10/bg, 50/cs
Knit Cuffs, Over-The-Head, Flat Pack		
NON27279OH	Regular/Large	10/bg, 50/cs
NON27279OHXL	X-Large	10/bg, 50/cs
Thumb Loop Wrists, Over-The-Head, Flat Pack		
NON27279OHTL	Regular/Large	10/bg, 50/cs
NON27279OHTLX	X-Large	10/bg, 50/cs

Protective Over-the-Head Gown in White Microporous Breathable Material

- High level fluid protection
- Apron-style neck enables gown to be placed on and off quickly
- Thumb loop style wrists keep sleeves from sliding up when putting on gloves
- Individually folded and packaged in a dispenser box



Item No.	Size	Pkg
NONTH350	Regular/Large	8/bx, 40/cs
NONTH400	X-Large	8/bx, 40/cs

NONTH350



Comfortable and Protective White Microporous Breathable Material in Our Classic Gown Style

- Created for high level fluid protection where comfort is needed
- Available with neck ties, waist ties and elastic wrists for a secure fit
- Generously sized for full coverage and flexibility
- Ten gowns per bag, flat pack



Item No.	Size	Pkg
NON27314	Regular/Large	10/bg, 50/cs
NON27314XL	X-Large	10/bg, 50/cs



NON27314

Prevention Plus Blue Procedure Gowns

- Made from blue impervious breathable material
- Material passed ASTM 1671 testing methods
- Hook-and-loop neck closure for easy on and off
- Available in blue with knit cuffs

Chemotherapy Drug Tested

Prevention Plus gowns are tested against the top chemotherapy drugs.

The Prevention Plus Material has been tested against the following chemotherapy drugs per ASTM F739-99a:

- Dacarbazine
- Carmustine (BiCNU)
- Cyclophosphamide (Cytoxan)
- Doxorubicin Hydrochloride
- 5-Fluorouracil
- Cisplatin
- Etoposide (Toposar)
- Paclitaxel (Taxol)
- Thio-Tepa
- Mitomycin C
- Vincristine Sulfate

(Test data available upon request)



Item No.	Size	Pkg
NON27457M	Medium	10/bg, 30/cs
NON27457	Regular/Large	10/bg, 30/cs
NON27457XL	X-Large	10/bg, 30/cs
NON27457XXL	XX-Large	10/bg, 30/cs

Thumbs Up[®] Polyethylene Gowns

- Made from blue impervious polyethylene material
- Material passed ASTM1670 and ASTM1671 test methods
- Apron-style neck enables gown to be placed on and off quickly
- Thumb loop style prevents sleeves from sliding up when putting on gloves
- Available in bulk and individually folded and wrapped in dispenser box



Item No.	Size	Color	Pkg
NONTH150	Regular/Large	Blue	15/bx, 75/cs
NONTH150Y	Regular/Large	Yellow	15/bx, 75/cs
NONTH200	X-Large	Blue	15/bx, 75/cs
NONTH100B	Regular/Large	Blue	(Bulk) 100/cs
NONTH200B	X-Large	Blue	(Bulk) 75/cs



NON27457

NONTH150

A wide selection of coveralls for every situation

Heavy fluids and tough jobs require a strong and protective material such as our *microporous breathable material*.

Need coverage for light jobs or to cover clothing?
Try our *Spunbond Coveralls*.

Breathable Coveralls

- Made from microporous breathable laminate material
- Maximum protection - ideal for maternity or the OR
- Comfortable for extended periods
- Available in white



Item No.	Description	Pkg
NONCV300L-XXL*	Straight Wrist and Ankle	5/bg, 25/cs
NONCV700L-XXXL*	Elastic Wrist and Ankle	5/bg, 25/cs

Multi-Layer Coveralls

- Made from fluid-resistant material
- Good combination of comfort and protection
- Available in white



Item No.	Description	Pkg
NONCV250L-XXL*	Elastic Wrist and Straight Ankle	5/bg, 25/cs

Spunbond Coveralls

- Made from heavyweight polypropylene material
- Breathable fluid-repellent material
- Available in white



Item No.	Description	Pkg
NONCV150L-XXXXL*	Elastic Wrist and Straight Ankle	5/bg, 25/cs

Standard Spunbond Coveralls

- Made from medium weight polypropylene material
- Highly breathable, fluid repellent material
- Available in white



Item No.	Size	Description	Pkg
CRI3000	X-Large	Straight Wrist and Straight Ankle	5/bg, 25/cs
CRI3001	XX-Large	Straight Wrist and Straight Ankle	5/bg, 25/cs



CRI3001

* When ordering, add the correct size letter code to the product number.

For example:
NONCV150XL = X-Large



NON28858

Non-Skid Multi-Layer Shoe Covers

- Made from breathable, fluid-resistant material
- Sport size available for larger, athletic-type shoes

Item No.	Size	Pkg
NON28858	Regular Size (up to men's size 12)	100/bx, 300/cs
NON28859	X-Large (up to men's size 15)	100/bx, 200/cs
NON28852	Sports Size (fits athletic shoes)	100/bx, 200/cs



NON28759

New! Non-Skid, Multi-Layer Shoe Covers

- Made from breathable, fluid-resistant material
- New, stripe pattern non-skid for increased traction

Item No.	Size	Pkg
NON29858	Regular Size (up to men's size 12)	100/bx, 300/cs
NON29859	XL (up to men's size 15)	100/bx, 200/cs
NON29852	Sports Size (fits athletic shoes)	100/bx, 200/cs



CRI2003

Our most popular line

Comfortable, Protective Shoe Covers

- Made from breathable, fluid-repellent spunbond material
- Sport size available for larger, athletic-type shoes

Item No.	Size	Pkg
NON28758	Regular Size (up to men's size 12)	100/bx, 300/cs
NON28759	X-Large (up to men's size 15)	100/bx, 200/cs
NON28752	Sports Size (fits athletic shoes)	100/bx, 200/cs

Boundary[®] Non-Skid Shoe Covers

- Made from breathable, fluid-repellent spunbond material
- Non-skid bottom

Item No.	Size	Pkg
CRI2002	Regular Size (up to men's size 12)	100/bx, 300/cs
CRI2003	X-Large (up to men's size 15)	100/bx, 200/cs

Boundary[®] Smooth Bottom Shoe Covers

- Made from breathable, fluid-repellent spunbond material
- Smooth bottom for situations when traction is not essential

Item No.	Size	Pkg
CRI2000	Regular Size (up to men's size 12)	100/bx, 300/cs
CRI2001	X-Large (up to men's size 15)	100/bx, 200/cs

Boundary[®] Polyethylene Shoe Covers

- Made from CPE plastic material that resists high levels of fluid
- Smooth bottom

Item No.	Size	Pkg
CRI2010	One Size Fits Most	1000/cs



CRI2010

Prevention Plus Impervious Breathable Knee-High Boot Covers

- Rated number one for breathability and complete protection
- Made from Blue Prevention Plus material
- Hook-and-loop closure strap allows for a snug fit
- Foam strip non-skid bottom

Item No.	Size	Pkg
NON27348P	Regular Size (up to men's size 12)	50/bx, 150/cs
NON27348PXL	X-Large (up to men's size 15)	50/bx, 150/cs



NON27348P

Non-Skid Fluid-Proof Boot Covers

- Made from fluid-proof coated polypropylene material
- High fluid protection on the bottom
- Convenient hook-and-loop strap allows wearer to adjust the fit
- Strap is heat-sealed to boot
- Available with non-skid pattern or foam strip



NON27150

Item No.	Size	Pkg
Non Skid Pattern		
NON27144	Regular (up to men's size 12)	50/bx, 150/cs
NON27144XL	X-Large (up to men's size 15)	50/bx, 150/cs
Foam Strip		
NON27150	Regular (up to men's size 12)	50/bx, 150/cs
NON27150XL	X-Large (up to men's size 15)	50/bx, 150/cs



NON27144

Non-Skid Multi-Layer/Poly Knee-High Boot Covers

- Breathable, fluid-resistant material above the ankle
- High fluid protection on the bottom
- Heat-sealed seams on foot reduce the chance of fluid penetration

Item No.	Size	Pkg
NON27143	Regular (up to men's size 12)	50/bx, 150/cs
NON27143XL	X-Large (up to men's size 15)	50/bx, 150/cs



NON27143XL

Non-Skid Multi-Layer/Poly Ankle Covers

- Designed for users who want full foot coverage without wearing a boot cover
- Impervious bottom with breathable material above the ankle
- Same heat-sealed seams and non-skid tread as boot cover design

Item No.	Size	Pkg
NON27145	Regular (up to men's size 12)	100/bx, 300/cs
NON27145XL	X-Large (up to men's size 15)	100/bx, 300/cs



NON27145

Pro Series Bouffant Caps

For hours of comfort and protection

- Lint-free, latex-free caps made from lightweight multi-ply and heavyweight polypropylene
- Stitched band keeps the cap flat and comfortable against the forehead
- Sold in solid blue and a variety of prints
- Available in 21", 24" and 28"
- All caps comply with NFPA standards

Item No.	Size	Description	Pkg
Lightweight Multi-Ply			
NON28232B	21"	Blue	100/bx, 500/cs
NON28233B	24"	Blue	100/bx, 500/cs
NON28223B	28"	Blue	100/bx, 500/cs
NON28226	24"	U.S. Flag Print	100/bx, 500/cs
NON28228	21"	Brush Stroke Print	100/bx, 500/cs
NON28229	24"	Brush Stroke Print	100/bx, 500/cs
NON28230R	24"	Breast Cancer Ribbon Print	100/bx, 500/cs
NON28233F	24"	Butterfly Print	100/bx, 500/cs
Heavyweight Polypropylene			
NON28234P	24"	Multi-Color Print	100/bx, 500/cs




NON28229

Add
Color & Style
to your
facility!

Medline's Boundary Bouffant Caps

- Choose from lightweight single-ply polypropylene or heavyweight spunbond material
- Available in solid blue, green, red and white
- Available in 21" and 24"
- All caps comply with NFPA standards

Item No.	Size	Description	Pkg
Lightweight Polypropylene			
CRI1000	21"	Green	100/bx, 500/cs
CRI1001	21"	Blue	100/bx, 500/cs
CRI1002	21"	White	100/bx, 500/cs
CRI1003	24"	White	100/bx, 500/cs
CRI1004	24"	Blue	100/bx, 500/cs
CRI1005	24"	Green	100/bx, 500/cs
CRI1009	24"	Red	100/bx, 500/cs
Heavyweight Polypropylene			
CRI1203	24"	White	100/bx, 500/cs



CRI1004

Sheer-Guard Surgeons' Caps

You're sure to find the style you need in our expanded line of surgeons' caps.

Choose from three different materials:

- SRM** - Scrim reinforced material, a paper-like absorbent material
- Multi-Layer Polypropylene** - Highly fluid-resistant, yet breathable
- Spunlace** - Extremely soft, comfortable and most absorbent

Item No.	Side Material	Crown Material	Color	Pkg
NON28625	SRM	Polypropylene	Light Blue	100/bx, 500/cs
NON28626	Multi-Layer	Polypropylene	Dark Blue	100/bx, 500/cs
NON28628	Spunlace	Polypropylene	Medium Blue	100/bx, 500/cs

NON28625



All caps are one size fits most and tie in back. All caps are packaged 100/bx, 500/cs.

Medline Surgeons' Hood

- Complete head coverage
- Available in multi-layer material
- Ties in back of head

Item No.	Description	Pkg
NONSH100C	Multi-Layer Material	100/bx, 300/cs

NONSH100C



Medline Surgeons' Head and Beard Covers

- Complete head coverage
- Made of multi-layer material with an elastic neck
- Beard cover is made of spunbond polypropylene, stays in place with an elastic band
- One size fits most

Item No.	Description	Pkg
NONSH600	Head and Beard Cover	100/bx, 300/cs
NONSH400	Beard Cover	100/bx, 1000/cs

NONSH600





Medline's classic lab coats and lab jackets are made from a fluid resistant multi-layer material to provide comfort and protection for everyday needs.

Medline's disposable products are available with a five snap closure or zip front, reinforced pockets, back vent for extended wear and a convenient loop inside for hanging. Enjoy our variety to fit your facility's needs.

Classic Lab Jackets with Knit Collar and Knit Cuffs

Item No.	Description	Pkg
NONRP500S-XXL*	White, Snap Front	10/bg, 30/cs
NONRP600S-XXL*	Blue, Snap Front	10/bg, 30/cs

Classic Lab Coats with Knit Collar and Knit Cuffs

Item No.	Description	Pkg
NONSW500S-XXXL*	White, Snap Front	10/bg, 30/cs
NONSW600S-XXL*	Blue, Snap Front	10/bg, 30/cs

Classic Lab Coats with Traditional Collar and Knit Cuffs

Item No.	Description	Pkg
NONSW100S-XXXL*	White, Snap Front	10/bg, 30/cs
NONSW400S-XXL*	Blue, Snap Front	10/bg, 30/cs
NONSW480S-XXXXL*	Blue, Zipper Front	10/bg, 30/cs

* When ordering, add the correct size letter code (*S-XXXL for Small through Triple Extra Large*) to the product number. **For example:** NONSW100L = Large

Classic Scrub Wear

An alternative to reusable and laundered scrubs.

- Made from breathable, fluid-resistant multi-layer material
- Generously cut shirts have one chest pocket and are available in either round neck or v-neck
- Pants have one rear pocket, reinforced crotch area, and are available with drawstring or elastic waist
- Available in blue

Item No.	Description	Pkg
Shirts		
NON27212S-XXXL*	Scrub Shirt, Round Neck	10/bg, 30/cs
NON27202S-XXL*	Scrub Shirt, V-Neck	10/bg, 30/cs
Pants		
NON27203S-XXL*	Scrub Pants, Drawstring Waist	10/bg, 30/cs
NON27213S-XXXL*	Scrub Pants, Elastic Waist	10/bg, 30/cs

NON27202XL
NON27203XL

NON27212M
NON27213M



Heavyweight Blue Apron

Provides highest level of protection against fluids.

- Made from heavyweight 4 mil plastic
- Complete protection for the OR, ER, decontamination, morgue or any area where heavy fluids are present
- Styled to easily pull over the head and tie in back

Item No.	Size	Pkg
NON24280	41" x 62"	10/bg, 100/cs



Poly Disposable White Aprons

Effective protection against liquids and waste.

- Made from polyethylene material
- Available in two lengths to protect the torso
- Convenient pullover styling for quick and easy removal

Item No.	Description	Pkg
NON24272	Lightweight, 24" x 42"	100/bx, 1000/cs
NON24274	Medium weight, 28" x 46"	100/bx, 500/cs
NON24274W	Medium weight, 28" x 46" Individually Wrapped, Dispenser Box	50/bx, 250/cs

NON24280

NON24274

Protective Sleeve Covers

- Get complete protection for your arms
- Available in two materials to fit your needs
- Covers are 18" and are made with elastic bands at both ends for wrist to upper arm coverage and a snug fit

Item No.	Description	Pkg
NON18SLEEVE1	Microporous Breathable Material, White	200/cs
NON18P2SLEEVE	Prevention Plus Material, Blue	200/cs



NON18P2SLEEVE

All-Around Protection in Protective Apparel

Medline offers many options to help you select the right protective apparel for your facility


To see samples of the materials used in our protective apparel, particularly isolation gowns, ask for a copy of our Isolation Gown Fabric brochure.

We can help you organize your protective apparel and other personal protective equipment (PPE) with one of our organizer units.

For more information, ask for a copy of our Personal Protection Organizers brochure.

Questions? Comments?

Contact your local Medline sales representative or call 1-800-MEDLINE. You can also visit us online at www.medline.com.



LIT109



Medline Industries, Inc.

One Medline Place
Mundelein, IL 60060

Medline United States
1-800-MEDLINE (633-5463)

www.medline.com | info@medline.com

Medline Canada
1-800-396-6996

www.medline.ca | canada@medline.com

Medline México
01-800-831-0898

www.medlinemexico.com | mexico@medline.com

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure. ©2011 Medline Industries, Inc. Medline, Boundary and Thumbs Up are registered trademarks of Medline Industries, Inc.