

Lumbar Support/Neck Roll

Three layers provides support and comfort to neck or lumbar area

Latex Free

Jelly roll design:

1. Inner high density foam cylinder core
2. Non-aqueous specialized visco-elastic gel layer surrounding core, reduces tissue interface pressure at critical points
 - Proprietary gel is a flexible solid that will not leak
3. Surrounded by 1.5" of visco-elastic memory foam

Cover features

- Fluid-proof stretch top cover, reduces shear/friction
- Antimicrobial
- Flame retardant
- Removable and washable
- End handle and convenient strap provide easy handling and positioning
- Limited warranty 18 months

8200 16"L x 6"D cylinder

