



TENA® *InstaDri Air*™ Underpad

**Maximum Absorbency
Underpad**



SOFT NONWOVEN TOPSHEET
provides extra strength and prevents
ripping or tearing.

**SUPERABSORBENT
MICROBEADS**
wick fluid away to help keep
skin dry and comfortable.



SPECIAL PIN-CORE STRUCTURE
allows even distribution of fluid
which eliminates clumps and blotches.
Moisture lock controls odors and
bacteria.

STRONG NONWOVEN BACKSHEET
resists bunching to promote comfort.

SPECIAL BREATHABLE FILM
allows moisture vapor and
excess heat to escape, helping
to prevent skin breakdown.
Prevents liquids from leaking
through onto clothing or bed
linens.



All TENA® Underpads
are Latex Free



TENA® InstaDri Air™ Underpads

TENA® Underpads	Size	Order #	Packaging	SCS/OCC
InstaDri Air™	23" x 36" / 58 x 91cm	359*	10 bags of 7/70 per case	10768702003596
	30" x 36" / 76 x 91cm	360*	8 bags of 5/40 per case	10768702003602

OTHER AVAILABLE TENA® UNDERPADS

Air Flow For use on air loss therapy beds	23" x 36" / 58 x 91cm	370*	6 bags of 10/60 per case	10768702003701
Regular	17" x 24" / 43 x 61cm	350	12 bags of 25/300 per case	10768702003503
	23" x 24" / 58 x 61cm	351	8 bags of 25/200 per case	10768702003510
	23" x 36" / 58 x 91cm	352	6 bags of 25/150 per case	10768702003527
Extra	17" x 24" / 43 x 61cm	353	12 bags of 25/300 per case	10768702003534
	23" x 24" / 58 x 61cm	354	8 bags of 25/200 per case	10768702003541
	23" x 36" / 58 x 91cm	355	6 bags of 25/150 per case	10768702003558
	36" x 36" / 91 x 91cm	361*	10 bags of 10/100 per case	10768702003619
	28" x 30" / 71 x 76cm	366	15 bags of 10/150 per case	10768702003664
	28" x 36" / 71 x 91cm	368	10 bags of 10/100 per case	10768702003688
Large	29.5" x 29.5" / 75 x 75cm	61310	10 bags of 15/150 per case	10768702613108
Ultra	23" x 36" / 58 x 91cm	357*	4 bags of 25/100 per case	10768702003572
Ultra Plus	28" x 30" / 71 x 76cm	365*	10 bags of 10/100 per case	10768702003657
	28" x 36" / 71 x 91cm	367*	10 bags of 10/100 per case	10768702003671
Premium	29.5" x 29.5" / 75 x 75cm	61312*	5 bags of 15/75 per case	10768702613122

*Contains Superabsorbent microbeads

Questions? Contact Us:

1.800.992.9939 | info.usa@tena.com | www.tena-usa.com