

INVACARE® ROLLATORS



Invacare® Rollators will meet all your needs with a focus on style and affordability. These rollators feature easy-to-operate ergonomic hand brakes and easy folding for storage and transport. Invacare Rollators offer great maneuverability and stability, and are height adjustable for a custom fit. Baskets are standard on all rollators. Trays are standard on models P429/2 and 65550.

Invacare® FR300™ Rollator Features

Model # 67100

Tray and Basket included

Cross Brace

Streamlined design

Curved backrest

Curb climber

65650 / 65650R Rollator Features

Model # 65650 Blue, 65650R Red

Easy-to-operate, ergonomic hand brakes lock the rear wheels for security

Flip-up padded seat with zipper compartment for storage of valuables

Removable foam padded curved backrest enhances comfort and can be easily removed

Lightweight aluminum frame

Rollator folds easily with basket attached

Universal basket fits underneath seat or hooks easily on front seat brace

HCPCS Code: E0143 & E0156

67100



65650: Blue (pictured), 65650R: Red

65851R / 65851B Adult Rollator Features

Model # 65851R Red, 65851B Blue

Easy-to-operate, ergonomic hand brakes lock the rear wheels for security

Flip-up seat

Lightweight aluminum frame

Basket hooks easily on front seat brace

Removable foam-padded straight backrest enhances comfort and can be easily removed

65350GR Junior Rollator Features

Model # 65350GR

Easy-to-operate, ergonomic hand brakes lock the rear wheels for security

Flip-up seat

Lightweight aluminum frame

Basket hooks easily on front seat brace

Removable foam-padded straight backrest enhances comfort and can be easily removed

HCPCS Code: E0143 & E0156

65851R: Red (pictured), 65851B: Blue Adult Rollator



65350GR: Junior Rollator

NOTE: Rollators are not intended to be used as wheelchairs.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further details.

P429/2 Rollator Features

Model # P429/2

Easy-to-operate, ergonomic hand brakes

Rollator folds easily with tote bag attached

Indoor/outdoor 7½" x 2" solid rubber wheels

Basket standard, attaches easily to the frame

Tray standard, attaches easily to the basket

Tote bag standard, attaches easily to the frame

65550 Rollator Features

Model # 65550

Easy-to-operate, ergonomic hand brakes lock the rear wheels for security

Sturdy seat

Rollator folds easily with basket attached

Basket and tray standard

**HCPCS Code: E0143 Model # P429/2
E0143 & E0156 Model # 65550**



P429/2: Three-Wheel Rollator



65550: Four-Wheel Rollator

66550 Bariatric Rollator Features

Model # 66550

500 lb. weight capacity

Wide, deep frame

Flip-up padded seat

Removable padded backrest

Folds for ease of storage and transport

Large 8" wheels for improved mobility over rougher terrain

Height adjustable loop handbrakes with foam handgrips

Durable steel construction

Pouch included

HCPCS Code: E0149



1. Bariatric size

NOTE: Rollators are not intended to be used as wheelchairs.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further details.



66550: Bariatric Rollator



Model #	Overall Dimensions	Dimension Inside Hand Grips	Seat Dimensions	Seat-to-Floor-Height	Caster Dimensions	Product Weight	Weight Capacity
67100	24.4" W 26.3" D	16"	15.5" W 8" D	23.6"	8"	17.5 lb.	285 lb.
65650, 65650R	24.25" W 28.28" D	18"	13.75" W 16.375" D	23.25"	7.5"	19 lb.	300 lb.
65851R, 65851B	28.5" W 27" D	16"	14" W 13" D	22"	6"	13.4 lb.	300 lb.
65350GR	26" W 24.5" D	13.5"	14" W 13" D	19.5"	6"	13 lb.	250 lb.
P429/2	26" W 24" D	14"	14" W 13" D	19.5"	6"	17 lb.	294 lb.
65550	23.25" W 27" D	16.5"	15.25" W 12.5" D	24"	7.5"	25 lb.	285 lb.
66550	27" W 29.5" D	22"	18" W 14" D	22"	7.5"	26 lb.	500 lb.



	Patient Height Range	Units Per Carton	Limited Warranty
67100	N/A	1	5 Years Brakes 1 Year Components*
65650, 65650R	5'5"-6'3"	1	5 Years Brakes 1 Year Components*
65851R, 65851B	5'6"-6'3"	1	5 Years Brakes 1 Year Components*
65350GR	5'-5'6"	1	5 Years Brakes 1 Year Components*
P429/2	5'3"-6'	1	5 Years Brakes 1 Year Components*
65550	5' 3"-6'2"	1	5 Years Brakes 1 Year Components*
66550	N/A	1	5 Years Brakes 1 Year Components*

*Excluding tires and wheels
NOTE: Rollators are not intended to be used as wheelchairs

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.
Form No. 05-139
Rev. 03/12
111259

