

INVACARE®  
**TRANSPORT**  
WHEELCHAIRS

---



---

**Invacare® Transport Wheelchairs** are designed for safe, efficient patient transport. Their lightweight frames and fold down backs make the transport chairs easy to travel with and store.

---

The **Invacare® Transport Chair** is designed for safe and efficient patient transport. The Transport chair has the added advantage of sharing interchangeable parts with the other chairs in the transport series. The Transport chair is compact and lightweight.

### Features








- Stylish, Silver Vein, powder-coated carbon steel frame
- Footrests are included with the chair and match frame finish
- Black nylon upholstery is flame retardant and bacteria resistant
- Full-length arms are padded for comfort
- Carry straps on seat allow for ease of folding and lifting
- Fold-down back makes transport and storage easy








**HCPCS Code: E1038**



1. Seating options available



							
Model #	Overall Width	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Width	Seat Depth	Seat-to-Floor Height
TRAN17FR	19"	39"	27"	7.25"	17"	16"	19"
TRAN19FR	21"	39"	27"	7.75"	19"	16"	19"

							
	Back Height	Caster Height	Wheel Height	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty
TRAN17FR	19" Fold-down	8"	8"	250 lb.	34 lb.	25"x12"x31"	1 Year
TRAN19FR	19" Fold-down	8"	8"	250 lb.	34 lb.	25"x12"x31"	1 Year

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

The **Invacare® Lightweight Aluminum Transport Chair** is designed for safety, ease and comfortable patient transport. The lightweight frame weighs only 19 lb. and incorporates a fold-down back, making these chairs easy to transport and store.

### Features

- Black nylon upholstery is flame retardant and bacteria resistant
- Full-length arms are padded for comfort
- Carry straps on seat allow for ease of folding and lifting
- Footrest painted to match frame finish
- Fold-down back makes transport and storage easy
- Weighs 19 lb. without riggings

**HCPCS Code: E1038**



1. Seating options available



Model #	Overall Width	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Width	Seat Depth	Seat-to-Floor Height
ALR(B)19HBFR	23.5"	37.5"	24.25"	8.5"	17", 19"	15.5", 16"	19.25"
LTTR(B)19FR	21.5"	36.5"	24.25"	8.5"	19"	15.5"	19.25"
LTB17FR	23.5"	36.5"	24.25"	7.75"	17"	15.5"	19.25"
	Back Height	Caster Height	Wheel Height	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty
ALR(B)19HBFR	17.5"	8"	12"	250 lb.	30 lb.	27"x10"x 31"	13 Months
LTTR(B)19FR	17.5"	8"	8"	250 lb.	27 lb.	26"x10"x 32"	13 Months
LTB17FR	17.5"	8"	8"	250 lb.	27 lb.	26"x10"x 31"	13 Months

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

The **Invacare® Heavy-Duty Transport Chair** is designed for safe and efficient patient transport. The Transport chair has the added advantage of sharing interchangeable parts with the other chairs in the transport series. It is compact and lightweight.

## Features

Removable padded desk arms  
Carbon steel frame with dual crossbraces  
Swingaway footrest with aluminum footplates

**HCPCS Code: E1039**



1. Seating options available
2. Bariatric size



**Invacare Corporation**  
www.invacare.com

## USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900















## Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2012

Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols ™ and  
®. All trademarks are  
owned by or licensed  
to Invacare Corporation  
unless otherwise noted.  
Specifications are subject  
to change without  
notification.  
Form No. 11-223  
111237



Model #							
TRHD22FR	Overall Width 28"	Overall Height N/A	Overall Depth w/o riggings 36.5"	Folded Width N/A	Seat Width 22"	Seat Depth 18"	Seat-to-Floor Height 19.5"
							
TRHD22FR	Back Height 20"	Caster Height 8"	Wheel Height 12"	Weight Capacity 400 lb.	Shipping Weight 51 lb.	Carton Dimensions 29"x11"x 33"	Limited Warranty 13 Months

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.