



DIAGNOSTIC PRODUCTS

Simple Solutions, Significant Results™



Medline Diagnostics

Medline is proud to present our latest line of diagnostics.

We have built a reputation on providing diagnostics that are highly reliable, extremely accurate, and built to last.



Index

| | |
|---------------------------|-------|
| Blood Pressure Monitors | 3 |
| Sphygmomanometers | 6 |
| Accessories | 8-9 |
| Stethoscopes | 10-13 |
| Compass Diabetes Program | 14-15 |
| Diagnostic Testing | 16 |
| Lancets | 17 |
| Patient Monitoring | 18 |
| Thermometers | 19-22 |
| Miscellaneous Diagnostics | 23 |

If you have questions about any of the products in this brochure, contact your Medline representative or call us at: 1-800-MEDLINE.

New!



MDS2007

Desktop Blood Pressure Monitor Features Average Mode Technology!

These premium monitors detect the appearance of an irregular heartbeat during measurement and gives a warning signal with your reading.

- Measurement includes patented Irregular Heartbeat Detector.
- Monitor also automatically stores 99 readings with date and time.
- Average mode feature takes three consecutive measurements and averages them for more accurate results.

The MDS93AC1PC comes with easy-to-use blood pressure analyzer software with point and click help feature.

| Item # | Description | Pkg. |
|------------|---|------|
| MDS93AC1PC | Premium Advanced Monitor. Includes blood pressure software, USB cable for PC download, 4-AA batteries, medium arm cuff and deluxe storage case. | ea |
| MDS2007 | Premium Monitor Includes 4-AA batteries, medium arm cuff and deluxe storage case. | ea |

Pressure Monitor Accessories Replacement cuffs and adapters.

| Item # | Description | Pkg. |
|----------|---|------|
| MDS9872 | Adult Cuff, (fits arms 22-30 cm, 8.58"-11.70") | ea |
| MDS9873 | Large Adult Cuff, (fits arms 30-38 cm, 11.70"-14.82") | ea |
| M-IADE01 | AC Adapter | ea |



Real-time graphic display of blood pressure being measured, on-screen instructions.



MDS2006

Average Mode Technology Wrist Blood Pressure Monitor

This advanced blood pressure monitor uses a unique, patented Average Mode Technology to take three consecutive measurements then display a single, averaged result — a procedure practiced by many physicians for a more accurate result.

- Real-time graphic display of blood pressure being measured, on-screen instructions.
- 5 languages.
- 30 memories for 2 separate users with date and time.
- Average mode feature takes three consecutive measurements and averages them for more accurate results.

| Item # | Description | Pkg. |
|---------|---------------|------|
| MDS2006 | Wrist Monitor | ea |

latexfree



Digital readout of systolic, diastolic pressure and pulse

MDS2003

Wrist Blood Pressure Monitor

- One-button operation
- Automatically inflates and deflates
- Digital readout of systolic and diastolic pressure
- Digital readout of pulse
- Memory with date and time
- Operates on two AAA batteries

| Item # | Description | Pkg. |
|---------|--------------------------|------|
| MDS2003 | Digital BP Wrist Monitor | ea |

Automatic Digital Blood Pressure Unit

Our revolutionary automatic digital blood pressure units have automatic inflation and error detection thanks to the proprietary "fuzzy logic" programming. The easy one-button design provides accurate and reliable readings.

- "Fuzzy logic" circuitry allows this unit to customize inflation levels for individual patients
- Eight-bit processor provides unsurpassed accuracy
- Auto pulse readings
- Memory
- Large digital readout
- Operates on four-AA batteries



MDS2001

| Item # | Description | Pkg. |
|-----------|---|------|
| MDS2001 | Automatic Inflation w/Std Adult Cuff (fits arms 22-30cm, 8-12 inches) | ea |
| MDS2001LA | Automatic Inflation w/Lrg Adult Cuff (fits arms 30-38cm, 12- 15 inches) | ea |
| MDS9873 | Large Adult Cuff (fits arms 30-38cm, 12- 15 inches) | ea |
| MDS9872 | Adult Cuff | ea |
| M-IADE01 | AC Adapter | ea |

Digital Blood Pressure Unit - Manual Inflate

Manually inflated digital blood pressure units have the convenience of a large, digital readout. Also provides accurate and reliable readings.

- Manual inflation
- Auto pulse readings
- Memory
- Large digital readout
- Operates on four-AA batteries



MDS2002

| Item # | Description | Pkg. |
|-----------|--|------|
| MDS2002 | Manual Inflation w/Std Adult Cuff (fits arms 22-30cm, 8-12 inches) | ea |
| MDS2002LA | Manual Inflation w/Lrg Adult Cuff (fits arms 30-38cm, 12- 15 inches) | ea |
| MDS9873 | Large Adult Cuff (fits arms 30-38cm, 12- 15 inches) | ea |
| MDS9872 | Adult Cuff | ea |
| M-IADE01 | AC Adapter | ea |

Digital Blood Pressure Monitors

Large, raised numbers and a luminescent (glow in the dark) coating make the Nite Shift Wall-Mount Aneroid easy to read in almost any type of lighting.

- Eight-foot coiled tubing and pivoting bracket extend service area for more convenient use.
- Inflation system and bulb assembly store easily behind gauge.
- Features Medline's nylon range finder cuff.
- Available in choice of Latex or Latex-Free

| Item # | Description | Pkg. |
|-----------|----------------------------------|------|
| MDS9400 | Adult, Latex | ea |
| MDS9400LF | Adult, Latex-Free Neoprene | ea |
| MDS9472 | Wall mount basket (not shown) | ea |



MDS9400LF

The Mobile Aneroid maintains all of the features of its Wall-Mount counterpart, the 9400, while adding the convenience of portability and height adjustment.

- Adjusts from 29" to 48" (74cm to 122cm) to accommodate a wide range of viewing levels.
- Includes inflation system and bulb assembly storage basket.

| Item # | Description | Pkg. |
|-----------|-------------------------------|------|
| MDS9407 | Adult, Latex | ea |
| MDS9407LF | Adult, Latex-Free Neoprene | ea |

MDS9407LF



Blood Pressure Cuff Barriers

Single-use wrap for added protection from cross contamination.

- ∴ Protection against bacteria pass through
- ∴ Will not alter blood pressure readings
- ∴ Latex Free
- ∴ Available in 3 sizes to fit any size blood pressure cuff

| Item # | Description | Pkg. |
|-----------|-------------|--------|
| MDS916113 | Small | 300/cs |
| MDS916114 | Adult | 300/cs |
| MDS916115 | Large | 150/cs |



Protection from Cross Contamination

Blood Pressure Cuff-Guard

Help reduce the potential of cross-contamination by using our protective, disposable blood pressure cuff covers.

- ∴ Latex-free
- ∴ Made of Tyvek® material
- ∴ Quick and easy to replace or remove the cover!

| Item # | Description | Pkg. |
|--|---------------------------|--------|
| MDS9157 | Child Cuff-Guard | 20/box |
| MDS9158 | Adult Cuff-Guard | 20/box |
| MDS9159 | Large Adult Cuff-Guard | 20/box |
| Soft Fabric for Single Time Use | | |
| MDS9167 | Child Soft Cuff-Guard | 20/bx |
| MDS9168 | Adult Soft Cuff-Guard | 20/bx |
| MDS9169 | Lg. Adult Cuff Guard | 20/bx |



MDS9114

The Compli-Mate Aneroid Sphygmomanometers are designed for many years of reliable service.

- Designer-colored cuffs offer a choice of vibrant colors.
- Comes complete with oversized matching nylon carrying case which allows the user to carry additional items.

| Item # | Description | Pkg. |
|---------|--------------|------|
| MDS9104 | Red | ea |
| MDS9106 | Navy Blue | ea |
| MDS9107 | Lavender | ea |
| MDS9108 | Hunter Green | ea |
| MDS9109 | Orchid | ea |
| MDS9112 | Teal | ea |
| MDS9113 | Purple | ea |
| MDS9114 | Royal Blue | ea |



MDS9123

Combining the Compli-Mates Aneroid Sphygmomanometers with the very popular Sprague Rappaport Stethoscope makes this Combination Kit popular with nurses and EMT's.

- Kit comes with an aneroid sphygmomanometer and a 22" (56 cm) Sprague Rappaport stethoscope.
- Packed together in a matching nylon carrying case.

| Item # | Description | Pkg. |
|---------|--------------|------|
| MDS9121 | Magenta | ea |
| MDS9123 | Purple | ea |
| MDS9124 | Royal Blue | ea |
| MDS9125 | Black | ea |
| MDS9127 | Lavender | ea |
| MDS9128 | Hunter Green | ea |



MDS9117

The Compli-Mates Dual-Head Combination Kit combines an aneroid sphygmomanometer and a dual-head stethoscope.

- Ideal for nurses and students.
- Kit comes complete with an aneroid sphygmomanometer and matching dual-head stethoscope.
- Packed together in a matching nylon carrying case.

| Item # | Description | Pkg. |
|---------|--------------|------|
| MDS9116 | Magenta | ea |
| MDS9117 | Teal | ea |
| MDS9118 | Purple | ea |
| MDS9119 | Royal Blue | ea |
| MDS9133 | Hunter Green | ea |
| MDS9135 | Burgundy | ea |
| MDS9136 | Red | ea |

The Latex-Free Handheld Aneroid's lightweight design and portability make it a popular handheld unit.

- No-stop pin action allows a visual calibration check, ensuring accuracy.
- Nylon rangefinder cuff provides quick, convenient patient fitting.
- Latex-Free construction.
- Field adjustable for long-term accuracy.

| PVC | | | Neoprene | | |
|---------|-------------|------|-----------|-------------|------|
| Item # | Description | Pkg. | Item # | Description | Pkg. |
| MDS9380 | Adult | ea | MDS9380LF | Adult | ea |
| MDS9386 | Infant | ea | MDS9386LF | Infant | ea |
| MDS9387 | Child | ea | MDS9387LF | Child | ea |
| MDS9388 | Large Adult | ea | MDS9388LF | Large Adult | ea |



Our Best Selling Sphygmomanometer.

Our finest handheld aneroid, the Nite-Shift Premier, has a luminescent dial that facilitates reading over a wide range of lighting conditions.

- Durable zinc-alloy shell protects precision aneroid mechanism.
- Shatter resistant acrylic faceplate surrounds the luminescent dial.
- Field adjustable for long-term accuracy.

| Item # | Description | Pkg. |
|---------|-------------|------|
| MDS9410 | Adult | ea |
| MDS9411 | Child | ea |
| MDS9413 | Large Adult | ea |



FOR HOME USE

Medline's Home Blood Pressure Kits feature D-ring cuffs to simplify home use

- Both units feature an easy-to-read aneroid gauge, complete instructions, recording chart, and carrying case.
- MDS9301 features an attached stethoscope while the two-person MDS9300 features a separate single-head stethoscope.
- Ideal for home use.

| Item # | Description | Pkg. |
|---------|---------------------------|------|
| MDS9300 | With detached stethoscope | ea |
| MDS9301 | With attached stethoscope | ea |

Simple and convenient





MDS91420LF

Two-Tube Inflation Bag & Nylon Range finder Cuff


| Item # | Description | Size | Pkg. |
|----------------------------|-------------|---------------------------|------|
| <i>Latex-Free (PVC)</i> | | | |
| MDS91420 | Adult | 5¾" x 22", (15cm x 56cm) | ea |
| MDS91421 | Child | 4¼" x 15", (11cm x 38cm) | ea |
| MDS91422 | Infant | 2⅞" x 10¼", (7cm x 26cm) | ea |
| MDS91423 | Large Adult | 6⅞" x 26", (17cm x 66cm) | ea |
| MDS91424 | Thigh | 8½" x 30½", (22cm x 77cm) | ea |
| <i>Neoprene Latex-Free</i> | | | |
| MDS91420LF | Adult | 5¾" x 22", (15cm x 56cm) | ea |
| MDS91421LF | Child | 4¼" x 15", (11cm x 38cm) | ea |
| MDS91422LF | Infant | 2⅞" x 10¼", (7cm x 26cm) | ea |
| MDS91423LF | Large Adult | 6⅞" x 26", (17cm x 66cm) | ea |



MDS91410

One-Tube Inflation Bag & Nylon Range finder Cuff

| Item # | Description | Size | Pkg. |
|----------------------------|-------------|---------------------------|------|
| <i>Latex-Free (PVC)</i> | | | |
| MDS91410 | Adult | 5¾" x 22", (15cm x 56cm) | ea |
| MDS91411 | Child | 4¼" x 15", (11cm x 38cm) | ea |
| MDS91412 | Infant | 2⅞" x 10¼", (7cm x 26cm) | ea |
| MDS91413 | Large Adult | 6⅞" x 26", (17cm x 66cm) | ea |
| MDS91414 | Thigh | 8½" x 30½", (22cm x 77cm) | ea |
| <i>Neoprene Latex-Free</i> | | | |
| MDS91410LF | Adult | 5¾" x 22", (15cm x 56cm) | ea |
| MDS91411LF | Child | 4¼" x 15", (11cm x 38cm) | ea |
| MDS91412LF | Infant | 2⅞" x 10¼", (7cm x 26cm) | ea |
| MDS91413LF | Large Adult | 6⅞" x 26", (17cm x 66cm) | ea |
| MDS91414LF | Thigh | 8½" x 30½", (22cm x 77cm) | ea |



MDS9505LF

Two-Tube Seamless Inflation Bag

| Item # | Description | Size | Pkg. |
|----------------------------|-------------|--------------------------|------|
| <i>Latex</i> | | | |
| MDS9505 | Adult | 4¾" x 8⅞", (12cm x 23cm) | ea |
| MDS9506 | Large Adult | 6" x 12¾", (15cm x 32cm) | ea |
| <i>Neoprene Latex-Free</i> | | | |
| MDS9505LF | Adult | 4¾" x 8⅞", (12cm x 23cm) | ea |
| MDS9506LF | Large Adult | 6" x 12¾", (15cm x 32cm) | ea |



MDS9205

One-Tube Seamless Inflation Bag

| Item # | Description | Size | Pkg. |
|----------------------------|-------------|--------------------------|------|
| <i>Latex</i> | | | |
| MDS9205 | Adult | 4¾" x 8⅞", (12cm x 23cm) | ea |
| MDS9206 | Large Adult | 6" x 12¾", (15cm x 32cm) | ea |
| <i>Neoprene Latex-Free</i> | | | |
| MDS9205LF | Adult | 4¾" x 8⅞", (12cm x 23cm) | ea |
| MDS9206LF | Large Adult | 6" x 12¾", (15cm x 32cm) | ea |



MDS9420

Nylon Range finder Velcro Cuffs

| Item # | Description | Size | Pkg. |
|---------|-------------|--------------------------|------|
| MDS9420 | Adult | 5¾" x 22", (15cm x 56cm) | ea |
| MDS9421 | Child | 4¼" x 15", (11cm x 38cm) | ea |
| MDS9423 | Large Adult | 6⅞" x 26", (17cm x 66cm) | ea |

Gauges and gauge parts

| Item # | Description | Pkg. |
|---------|---|------|
| MDS9414 | Premier Gauge Only | ea |
| MDS9415 | Face Glass for Premier Gauge | ea |
| MDS9381 | Standard Gauge Only | ea |
| MDS9507 | Face Glass for Standard and Economy Gauge | ea |



MDS9414



MDS9381

Bulbs and Valves

| Item # | Description | Pkg. |
|------------|---|------|
| MDS9508 | Valve Only | ea |
| MDS91481 | Bulb Only, Latex-Free (PVC) | ea |
| MDS91481LF | Bulb Only, Latex-Free (Neoprene) | ea |
| MDS91482 | Bulb and Valve Assembly, Latex-Free (PVC) | ea |
| MDS9630 | Bulb and Valve Assembly, Premium Latex | ea |
| MDS9560 | Bulb and Valve, Large, Latex | ea |

MDS9630



MDS9508



MDS91481LF

Tubing

| Item # | Description | Pkg. |
|-----------|--|------|
| MDS9473 | 4-Foot Coiled Tubing w/Connector, Latex | ea |
| MDS9475 | 8-foot Coiled Tubing w/Connectors, Latex | ea |
| MDS9475LF | 8-foot Coiled Tubing w/Connectors, Latex-Free (Neoprene) | ea |



MDS9475LF

Other parts


| Item # | Description | Pkg. |
|---------|---|------|
| MDS9472 | Inflation Systems Holder for Wall Mounted or Mobile Sphygmomanometers | ea |
| MDS9481 | Metal Connector | ea |



MDS9472



MDS9481



Our Best Value Stethoscopes!

Quality and economy makes our single and dual-head stethoscopes excellent values.

- Both stethoscopes feature 22" (56 cm) tubing and lightweight design.
- Thick-walled tubing minimizes background interference.
- Dual-heads feature a non-chill ring on the bell for increased patient comfort.

| Item # | Description | Pkg. |
|---------------------------------|--------------|------|
| Single-Head Stethoscopes | | |
| MDS926101 | Black | ea |
| MDS926102 | Gray | ea |
| MDS926103 | Blue | ea |
| MDS926105 | Lavender | ea |
| MDS926106 | Red | ea |
| MDS926107 | Hunter Green | ea |
| MDS926108 | Yellow | ea |
| Dual-Head Stethoscopes | | |
| MDS926201 | Black | ea |
| MDS926202 | Gray | ea |
| MDS926203 | Blue | ea |
| MDS926205 | Lavender | ea |
| MDS926206 | Red | ea |
| MDS926207 | Hunter Green | ea |



This Stethoscope has EXCELLENT Sound!

The Sprague Rappaport Stethoscope incorporates versatility into an excellent sound-transmitting device.

- Each stethoscope comes complete with a component kit that includes two sets of eartips, large and small diaphragm, adult bell, pediatric bell, and infant bell.
- The dual 22" (56 cm) tubing design transmits even the faintest sounds.

| Item # | Description | Pkg. |
|-----------|--------------|------|
| MDS926301 | Black | ea |
| MDS926302 | Gray | ea |
| MDS926303 | Blue | ea |
| MDS926305 | Lavender | ea |
| MDS926306 | Red | ea |
| MDS926307 | Hunter Green | ea |



Single or Dual-Head Parts

| Item # | Description | Pkg. |
|---------|------------------------|------|
| MDS9544 | White Plastic Ear Tips | ea |
| MDS9547 | Soft Rubber Ear Tips | ea |
| MDS9552 | Diaphragm | ea |



Medline's Stainless Steel Stethoscopes offer unsurpassed quality and value

- Each stethoscope is handcrafted from high tensile strength German stainless steel. Features 22" (56 cm) tubing.
- The color coordinated non-chill rings provide excellent patient comfort.
- All models come in a foam-lined box with one extra set of ear tips.
- Available in adult, pediatric, neonatal, and cardiology styles.
- Cardiology Stethoscope has a deep, cone-shaped bell — a unique feature which provides excellent sound transmission and acoustic sensitivity. Features 17" (44 cm) tubing.

| Item # | Description | Pkg. |
|-------------------|-------------|------|
| Adult | | |
| MDS9225 | Gray | ea |
| MDS9226 | Black | ea |
| MDS9227 | Blue | ea |
| MDS9228 | Burgundy | ea |
| Pediatric | | |
| MDS9220 | Gray | ea |
| Neonatal | | |
| MDS9215 | Gray | ea |
| MDS9216 | Pink | ea |
| Cardiology | | |
| MDS9250 | Black | ea |



Our Bowles Stethoscope features thick-walled 22" (56 cm) tubing, a chrome-plated binaural and a chrome-plated, raised, metal-stem chest piece.

| Item # | Description | Color | Pkg. |
|---------|--------------------|-------|------|
| MDS9515 | Bowles Stethoscope | Black | ea |



The Pediatric and Neonatal Stethoscopes have all the features of the dual-head stethoscopes, but are sized to fit the needs of the individual practitioner and patient.

- Both stethoscopes feature 22" (56 cm) tubing and lightweight design.
- Thick-walled tubing minimizes background interference.
- Aluminum binaurals.

| Item # | Description | Color | Pkg. |
|---------|-------------|-------|------|
| MDS9557 | Pediatrics | Gray | ea |
| MDS9565 | Neonatal | Gray | ea |



MIF767X04 - Navy Blue

MIF767XBO - BlackOut

MDF® 2-in-1 Deluxe Sprague Rappaport Stethoscope

The advanced X-configuration acoustic tubing on the MDF® 767X eliminates all tube rubbing noises associated with the traditional Sprague-type stethoscopes and provides superior acoustic performance making it the most advanced Sprague Rappaport available today.

The chestpiece enables users to switch between adult and pediatric-sided bowels/bells with ease and interchange both sides of the chest piece with adult, pediatric and infant fittings (2 diaphragms and 3 bells). The 2-in-1 tubes in the X configuration converge at the middle of the tube and retain a dual lumen and dual input from the chestpiece. This innovation provides superior acoustic quality never available before on a Sprague Rappaport type stethoscope.

| Item # | Color | Pkg. |
|-----------|-----------------------|------|
| MIF767X04 | Navy Blue | ea |
| MIF767X08 | Purple | ea |
| MIF767X10 | Royal Blue | ea |
| MIF767X17 | Burgundy | ea |
| MIF767XBO | BlackOut® (All Black) | ea |



MIF74010 - Royal Blue



MIF740BO - BlackOut

MDF® Pulse Time™ Stethoscope

The original, patented MDF® 740 Pulse Time Stethoscope, integrated with a digital timing device, enables medical professionals to maintain eye contact with patients during diagnosis while monitoring elapsed time. U.S.A. PAT. NO. 5252787

This innovative and patented Pulse Time stethoscope has a liquid-crystal analog (LCA) quartz timepiece equipped on top of the integrated bowels-type chestpiece.

The advanced chestpiece is made of lightweight anodized aluminum and the diaphragm is fitted with the unique, raised, ultra-thin fiber diaphragm for greater amplification and transmission of heart and lung sounds.

| Item # | Color | Pkg. |
|----------|----------------------|------|
| MIF74008 | Purple | ea |
| MIF74010 | Royal Blue | ea |
| MIF74011 | Black | ea |
| MIF74017 | Burgundy | ea |
| MIF740BO | BlackOut (All Black) | ea |



Customize Our Silver-infused Stethoscope Cover

- Fits ALL stethoscope models
- Easy-to-use design
- For use with or without your current diaphragm without reduction in sound quality
- Customize with YOUR LOGO!

| Item # | Description | Pkg. |
|-----------|-------------|--------|
| MDS926400 | Stethocap | 25/box |

The Teaching/Training Stethoscope is designed for multi-function teaching applications.

- The dual head design offers the necessary flexibility to use this scope in all clinical and educational conditions.
- Comes complete with two binaurals, two pieces of Y-tubing, and an extension tube.

| Item # | Description | Pkg. |
|---------|-------------------------------|------|
| MDS9750 | Teaching/Training Stethoscope | ea |



MDS9750

The Fetal Stethoscope is designed for superior sound sensitivity.

- Chrome-plated inner spring binaural.
- Vinyl coated headrest.
- 22" (56 cm) tubing.

| Item # | Description | Pkg. |
|---------|-------------------|------|
| MDS9775 | Fetal Stethoscope | ea |



MDS9775

Our Disposable Stethoscope features plastic binaural and a sensitive chest piece that provides excellent sound quality.

- Designed for single patient use.
- Helps reduce cross-contamination in infectious areas.
- High visibility yellow.

| Item # | Description | Pkg. |
|---------|------------------------|------|
| MDS9543 | Disposable Stethoscope | ea |



MDS9543



Improving Diabetes Management in Long-Term Care


Medline's Compass Program A Diabetes Resource for Long-Term Care

Diabetes is a devastating disease affecting approximately 23.6 million people in the U.S. 5.7 million cases continue to go undiagnosed¹ for many years without treatment until severe complications occur. As many as 1 in 4 nursing home residents have diabetes.² The time and cost associated with caring for these individuals greatly impact the long term care facilities.

Medline has developed a comprehensive tool to help long-term care facilities provide better care for residents with diabetes by compiling and organizing data acquired during the course of working with healthcare professionals for more than 40 years. "Compass- A Diabetes Resource for Long-Term Care" emphasizes the importance of proper diabetes management to prevent complications and also illustrates how to achieve it.

Considering that heart disease, stroke, high blood pressure, blindness, kidney disease, dental disease, nervous system disease, and even amputations are all complications of diabetes¹, it is easy to see why healthcare providers strive to provide the best possible care for residents. Additionally, knowing that the Centers for Medicare & Medicaid Services is reviewing Pay for Performance measures and evaluating avoidable hospital admissions, like those due to complications from diabetes, it is vital for nursing home staff to understand current standards of practice.

With 26 percent of the nursing home population battling diabetes and its complications on a daily basis, Medline is providing educational tools for residents, their families and their caregivers.






What's in the Compass Diabetes box?

- One **Director of Nursing (DON) Instruction Manual**—
A teacher's guide to the program, which includes 4 CE contact hours.
- Three additional **Glucose Monitoring Training** resource books—
Put them on treatment carts or in nurses' stations for easy access.
- Five additional **Family Education** books—Give them to residents' families.
- One **Educational Compact Disc** and one **Educational DVD**
To educate staff on using the OptiumEZ™ Glucose Monitoring System.

A look at a few things from the box:




Instructional Manual



Educational CD and DVD

Resource Books



Family Educational Books

**New!!
No Coding
Required**



Medline OptiumEZ Glucose Monitoring System

The New Medline OptiumEZ Glucose Monitoring System requires no coding and provides easy, accurate, and reliable results! The OptiumEZ test strips are designed to prevent inaccurate tests and test strip waste caused by interferences from many common substances* and or/ application of an inadequate blood sample volume. The OptiumEZ Glucose Monitoring System is designed to handle all of your blood glucose needs.

Easy to Test

- No Coding Required
- Simple 2 Step Testing
- Results in 5 seconds
- Small blood sample size- 0.6 uL

Easy, Accurate and Reliable Results

- Unique chemistry minimizes the potential interfering effects of many common medications.
- Allows for the ability to reapply more blood on the same test strip.
- Individually wrapped test strips preserve the integrity of the test strips, which ensures an accurate glucose reading.

| Item # | Description | Pkg. |
|-----------|---|--------------------|
| MIS71039z | OptiumEZ Blood Glucose Meter | ea |
| MIS71041 | OptiumEZ Glucose Test Strips | 12 bx/cs, 50 ea/bx |
| MIS71041z | OptiumEZ Glucose Test Strips | 50 ea/bx |
| MIS71042 | OptiumEZ Glucose Test Strips | 6 bx/cs, 100 ea/bx |
| MIS71042z | OptiumEZ Glucose Test Strips | 100 ea/bx |
| MIS70550 | Precision Control Solution 1 vial high, 1 vial low per box | 6 bx/cs |
| MIS70550z | Precision Control Solution 1 vial high, 1 vial low per box | 1 bx |

Introducing the INRatio2 PT/INR Monitor

#1 selling brand in point of care PT/INR Monitoring for Three Years in a Row**
The INRatio2 PT/INR testing system combines reliable results with practical convenience making it an excellent solution for anticoagulation patient management.

1 Minute Test Time

- 1 button operation with simple icon interface.
- 1 step auto-on with test strip insertion.
- 1 drop of blood from a fingerstick.

2 Level On Board Quality Control

Rigorous clot-based and electronic quality control checks are performed automatically with every test, ensuring a reliable result

Individually wrapped Strips/No Refrigeration Required

Individually foil-wrapped strips maintain freshness and stability, and can be stored at room temperature.

Small Sample Size

1 unmeasured drop of blood applied directly to the strip performs both the test and quality controls.

| Item # | Description | Pkg. |
|--------------|--|--------|
| P-T200431 | INRatio2 Professional Kit (includes meter) | ea. |
| P-T100139Z | INRatio Test Strips | 48/bx |
| P-T0200235A | Blood Collection Tubes | 100/bg |
| MPHSAFETY21Z | Lancet, Normal Flow, 21 Gauge 1.8mm | 200/bx |



Optium EZ *
Data on file, Abbott Diabetes Care Inc.
INRatio2 **
HPIS comparative sales data, point of care PT/INR monitors sold Q4 2004-Q4 2008.



DIB0310510

Hemoglobin A1c Test

The A1cNow® InView System offers incredibly fast, quantitative A1c results (in just 5 minutes) utilizing micro-optical technology that is integrated into a proprietary reusable monitor with single use cartridges. Accuracy is assessed by the same standards used by laboratories.

| Item # | Description | Pkg. |
|------------|-----------------------|--------|
| DIB0310510 | A1cNow InView Monitor | 10/kit |



MSC0002

Monofilament

This monofilament can be used as a routine test to assess a resident's foot sensitivity.

| Item # | Description | Pkg. |
|---------|---|-------|
| MSC0002 | Foot Monofilament, Diabetic Neuropathy Test | 25/bx |



MPHSAFETY28

Medline Safety Lancets

Medline safety lancets are simple and safe for both the healthcare professional and the patient. The needle retracts automatically and locks out for extra safety, protecting patients and healthcare professionals from needlestick injuries and reducing cross-contamination.

| Item # | Description | Pkg. |
|--------------|-------------------------------------|---------------------|
| MPHSAFETY21 | Lancet, Normal Flow, 21 Gauge 1.8mm | 18 bx/cs; 200 ea/bx |
| MPHSAFETY21Z | Lancet, Normal Flow, 21 Gauge 1.8mm | 200 ea/bx |
| MPHSAFETY28 | Lancet, Micro Flow, 28 Gauge 1.6mm | 18 bx/cs; 200 ea/bx |
| MPHSAFETY28Z | Lancet, Micro Flow, 28 Gauge 1.6mm | 200 ea/bx |



OMMAT07 series

OMMAT1004Z

Unistik® Safety Lancets

Offers maximum comfort for your patients. The Unistik needlepoint is hidden before use and automatically retracts after use to avoid needle stick injuries and help reduce cross-contamination.

| Item # | Description | Pkg. |
|------------|------------------------|---------|
| OMMAT0704Z | Unistik 2 Normal Depth | 200/bx |
| OMMAT0704 | Unistik 2 Normal Depth | 6000/cs |
| OMMAT0714Z | Unistik 2 Extra Depth | 200/bx |
| OMMAT1004Z | Unistik 3 Normal Depth | 200/bx |



MPH281001DMEZ

Lancets

Ultra Sharp Lancets have fine-point tip for smooth and fast procedures. Choose from our general lancets in 100 or 200 count. These lancets can be used with all general lancing devices.

| Item # | Description | Pkg. |
|---------------|-------------------------|-----------|
| MPH1505515 | Lancet, 23 Gauge Needle | 2400/cs |
| MPH1505515Z | Lancet, 23 Gauge Needle | 200/bx |
| MPH1505514 | Lancet, 26 Gauge Needle | 2400/cs |
| MPH1505514Z | Lancet, 26 Gauge Needle | 100/bx |
| MPH281001DME | Lancet, 28 Gauge Needle | 20,000/cs |
| MPH281001DMEZ | Lancet, 28 Gauge Needle | 100/bx |



SunTech 247 Patient Monitor An Expandable Diagnostic Station with BP, Temperature and SpO2

This patient monitor was designed to measure clinical grade vital signs, while allowing the unit to be customized to the facilities needs. The 247 Patient Monitor is completely modular, meaning that you can purchase what the functionalities you need today and then add additional parameters, if needed. In addition, it combines the efficiency of an automated monitor with the convenience of a manual mode that allows the clinician to take a BP reading with a stethoscope whenever needed!

| Item # | Description | Pkg. | Accessories Included |
|--------------------------------------|---|------|--|
| SNT99006200 | BP with Mobile Stand | 1/ea | BP: 8' (2.5 meters) blood pressure hose, adult & large adult all purpose cuffs, mains/power supply, mains/power cable, resource cd, quick start guide, reference card BP + Temperature: oral/axillary probe (blue), 1 box of Filac Fastemp [™] probe covers BP + SpO2: reusable adult finger sensor, 6' (1.8 meters) extension cable |
| SNT99006300 | BP & Temperature with Mobile Stand | 1/ea | |
| SNT99006400 | BP & SpO2 with Mobile Stand | 1/ea | |
| SNT99006100 | BP, SpO2, Temperature with Mobile Stand | 1/ea | |
| Wall Mount System | | | Additional Accessories (not included) SWD202020Z Probe Covers for the FastTemp Thermometer STCUS1320DM Child Replacement Cuff, 5/bx STCUS2635DM Adult Replacement Cuff, 5/bx STCUS3242DM Large Adult Replacement, Cuff 5/bx STCUS4250DM Thigh Replacement Cuff, 5/bx |
| SNT99006600 | BP | 1/ea | |
| SNT99006700 | BP & Temperature | 1/ea | |
| SNT99006800 | BP & SpO2 | 1/ea | |
| SNT99006500 | BP, Temperature & SpO2 | 1/ea | |
| Table Top System with Battery | | | |
| SNT99007000 | BP | 1/ea | |
| SNT99007100 | BP & Temperature | 1/ea | |
| SNT99007200 | BP & SpO2 | 1/ea | |
| SNT99006900 | BP, Temperature & SpO2 | 1/ea | |
| Modules | | | |
| SNT98012800 | Temperature Module | 1/ea | |
| SNT98012900 | SpO2 Module | 1/ea | |

Welch Allyn SureTemp® Model 690

The SureTemp Plus 690 Thermometer delivers oral results in just four seconds! This electronic thermometer is easy to use, compact and dependable. It delivers oral results in four seconds, axillary temperatures (for children under four) in 10 seconds and rectal temperatures in 10 seconds with the use of an optional, interchangeable probe. Large LCD display features Fahrenheit or Celsius readouts and last temperature recall. This affordable thermometer is ideal for alternate care environments. Three 1.5V AA batteries provide 6000 readings. Automatic probe cover ejection; Two-year limited warranty.

| Item # | Description | Pkg. |
|-------------|---------------------------|---------|
| W-A01690200 | Oral/Axillary Thermometer | ea |
| W-A02892000 | Rectal Probe | ea |
| W-A5031101 | Disposable Probe Covers | 1000/cs |
| W-A5031101Z | Disposable Probe Covers | 250/bx |



W-A01690200

Thermoscan® Thermometers

The NEW Braun Thermoscan Pro 4000 Thermometer patented technology is an ergonomically designed ear thermometer engineered to take temperatures easily and accurately with the click of a button.

- Exclusive ExacTemp® technology provides an active user feedback light that signals that the probe has positioned in a stable way throughout the measurement, avoiding inaccurate measurements by handling errors
- Preheated probe tip minimizes the cool down effect experienced when taking repeat temperatures, improving accuracy and repeatability
- Automatic probe cover eject button for quick and easy removal of used probe covers
- Built-in cross-contamination safeguards – probe cover sensor reminds operator to use a new probe cover
- Wide range of patient temperature – 68°F-108°F patient range allows for rescue of intensive care
- Wide range of ambient temperature – 50°F-104°F ambient range allows for use in most environments
- Durable – tested for damage resistance
- Convenient wall-holder for Pro 4000 unit w/probe cover storage is available
- Three-year limited warranty coverage

| Item # | Description | Pkg. |
|----------|---------------------|----------|
| TMS40002 | Thermoscan PRO 4000 | ea |
| TMSPC800 | Probe Covers | 40 bx/cs |



TMS40002

Exergen® Temporal Scanner™ Thermometer

Precision infrared scanner provides a truly non-invasive method for measuring temperature. Simply scan across the temporal artery area of the forehead, touch behind the earlobe, and an accurate reproducible temperature is measured in about three seconds.

- Convenient one-handed operation
- Completely non-invasive
- Stainless steel probe for highest durability
- Easy snap-on probe caps
- Three-year warranty

| Item # | Description | Pkg. |
|------------|--------------------------------------|------------|
| EXGTAT5000 | Infrared Temporal Artery Thermometer | ea |
| EXGTAT2000 | Infrared Temporal Artery Thermometer | ea |
| EXG134203 | Disposable Probe Caps | 1000 ea/cs |



EXGTAT5000



Genius2 Thermometer

Measures temperature from the tympanic membrane, which shares the same blood flow as the hypothalamus, the body's natural thermostat.



- Utilizes infrared technology to derive the most accurate body temperature from the tympanic membrane.
- Records and analyzes up to 100 temperature measurements, and selects the highest and most accurate in less than 2 seconds, as compared to 32 in Genius and 1-8 measurements in competitive thermometers.
- Reduces the chance of an inaccurate reading due to improper technique or patient movement.
- Durable tip-based sensor technology Improves durability and accuracy, minimizing the need for recalibration.
- Tapered probe design
- Rapid response time 1 - 2 second
- Real time temperatures
- Programmable display
- Reduces cross-contamination No contact with mucus membranes No touch probe cover application and release

| Item # | Description | Pkg. |
|------------|---------------------|------------|
| SWD303000 | Genius2 Thermometer | ea |
| SWD303030Z | Probe Covers | 96 ea/bx |
| SWD303030 | Probe Covers | 2,112/case |

Get Both Speed And Accuracy with The FILAC® 3000 Features Both!



The FILAC® 3000 Electronic Thermometer utilizes advanced technology to ensure the most accurate and fastest readings are provided in easy-to- use formats. These cutting edge thermometers address the need for fast, accurate readings in under 4 seconds or an advanced accuracy reading in under 6 seconds. Compact, user friendly design will stand upright on a flat surface; intuitive icon based interface system with backlight LCD display; unique progress bar indicates status of temperature measurement; last temperature recall function. Patented isolation chamber design and color schemes along with the durable probe covers help minimize the chance of cross-contamination and enhance infection control measures. Operates on 4 AA batteries for 6000 measurements. Equipped with electronic anti-theft system, self-calibrating electronics and auto on-off function. Includes oral probe and push button probe cover ejection.

| Item # | Description | Pkg. |
|------------|---|----------|
| SWD504000 | Filac 3000 Thermometer | ea |
| SWD202020 | Probe Covers for the Filac Fastemp Thermometer & 3000 | 5,000/cs |
| SWD202020Z | Probe Covers for the Filac Fastemp Thermometer & 3000 | 500/box |

Six Seconds*, Temple Digital Thermometer

The NEW Medline Instant Read* Digital Temple Thermometer accurately measures body heat using Patented R.A.T.E.™ technology (Rapid, Accurate, Temperature, Establishment).

- Sensors rapidly track heat flow generated from blood vessels to the skins surface and convert measurements to body temperature.
- Easy-to-use while asleep or awake.
- Meets accuracy standards.
- Large digital display.
- Switchable from Fahrenheit to Celsius.
- Beeps when ready.
- Memory stores previous reading.

| Item # | Description | Pkg. |
|---------|--------------------|----------|
| MDS9698 | Temple Thermometer | ea |
| MDS9699 | Probe Covers | 250ea/bx |



Gentle, non-invasive sensor accurately reads in F° & C°

Power/memory control button

Large Digital Display

2-AAA Batteries included

Features Patented R.A.T.E. Technology

*Readings in six to ten seconds depending on individual.

FeverScan™ Disposable Forehead Thermometers

FeverScan™ products utilize thermo chromic liquid crystal (TLS) technology for the measurement of clinical temperature trends in the healthcare environment. Liquid crystals are chemical compounds and mixtures that exhibit the mechanical properties of liquids and the optical characteristics of solids.

FeverScan™ thermometers are composed of a series of discrete temperature events, each event shows a green color as the indicated temperature is attained.

Anesthesia & PACU (with adhesive backing)

| Item # | Description | Pkg. |
|--------|---------------|---------------|
| ATC430 | °F 94°F-104°F | 100/dispenser |
| ATC431 | °C 35°C-40°C | 100/dispenser |

Hypo/Hyperthermic Monitoring (with adhesive backing)

| | | |
|--------|-------------------------------|---------------|
| ATC434 | Extended Range °F, 84°F-106°F | 100/dispenser |
| ATC435 | Extended Range °C, 30°C-41°C | 100/dispenser |

Reusable (no adhesive)

| | | |
|--------|---------------|---------------|
| ATC432 | °F 94°F-104°F | 100/dispenser |
| ATC433 | °C 35°C-40°C | 100/dispenser |



- Continuous Display
- Forehead Temperature Monitor
- Anesthesia through PACU
- Watch for Hypo/Hyperthermia
- Pre and Post-Op Surveillance
- Disposable, Latex-Free
- Core Adjusted Temperature



New!


MDS96801



MDS9650



MDS9628



MDS9654B



MDS9629

Digital Thermometers

Our digital thermometers provide accurate temperatures in about 60 sec.

- Large LCD screen for easy readings.
- Simple push-button start.
- Disposable sheaths allow for sanitary oral, rectal, or axillary use.

| Item # | Description | Pkg. |
|----------|---------------------------------|--------|
| MDS9650 | Premier Oral °F & °C | ea |
| MDS9652 | Premier Rectal °F | ea |
| MDS9654B | Standard Oral °F | 200/cs |
| MDS9628 | Oral Kit w/20 Sheaths °F | ea |
| MDS9629 | Rectal Kit w/20 Sheaths °F | ea |
| MDS9607 | Oral Sheaths | 100/bx |
| MDS9608 | Rectal Sheaths | 100/bx |
| MDS96801 | Flex-tip, waterproof, °C or °F | ea |
| MDS9655 | Standard, Oral, Bulk, °F & °C | 200/cs |
| MDS9655Z | Standard, Oral, °F & °C | 20/bx |
| MDS9650B | Premier Oral, Bulk, °F & °C | 200/cs |
| MDS9656 | Standard, Rectal, Bulk, °F & °C | 200/cs |
| MDS9656Z | Standard, Rectal, °F & °C | 20/bx |



Disposable Penlights

Lightweight and disposable, these high-intensity exam lights are perfect for the physician, nurse, and EMT. They're constructed of high-impact molded styrene and can be clipped on pockets.

| Item # | Description | Pkg. |
|-----------|----------------------------|------|
| MDS131020 | Disposable Penlight, White | 6/pk |



AED External Defibrillator

Cardiac Science PowerHeart® AED G3 Automated External Defibrillator with voice prompts and text screen displays guide you through the rescue process. When cardiac arrest occurs, every second counts. The advantage of the PowerHeart AED G3 is that it may be used by anyone who has received basic training on its operation. Electrodes are pre-connected and disposable. Semi-automatic defibrillator with biphasic variable waveform delivers an effective shock in less than 55 seconds. Daily, weekly, and monthly automated self-checks ensure that PowerHeart is always ready to go. Internal memory holds 60 minutes of ECG data. Includes one pair of defibrillation electrodes, instructional CD-ROM and training video, RescueLink® and MDLink® software and serial communication cable. Seven-year warranty on device & four-year full replacement on battery.

| Item # | Description | Pkg. |
|---------------|---|------|
| AED9300SEMICC | Semi Automatic AED Kit | ea |
| AED9300AUTOCC | Fully automatic AED Kit | ea |
| Includes: | 1 Replacement set of adult defibrillator electrode pads | |
| | 1 Soft carrying case | |
| | 1 Ready Kit | |



Powerheart G3 AED

References:

1. National Diabetes Fact Sheet, 2007. Centers for Disease Control and Prevention website. Available at www.cdc.gov/diabetes/pubs/pdf/nf07_2007.pdf. Accessed on April 24, 2009
2. Johnson EL, Brosseau JD, Soule M, et al. Treatment of diabetes in long-term care facilities: a primary care approach. *Clinical Diabetes*. 2008;26(4):152-156

A1cNow is a registered trademark of Metrika, Inc.
 Thermoscan is a registered trademark of Braun GmbH Corp.
 Exergen is a registered trademark and Temporal Scanner is a trademark of Exergen Corp.
 Feverscan is a trademark of LCR Hallcrest
 FILAC is a registered trademark of Covidien AG Corp.
 INRatio is a registered trademark of HemoSense, Inc.
 MDF is a registered trademark and PulseTime is a trademark of MDF Instruments USA, Inc.
 Medline is a registered trademark of Medline Industries, Inc.
 OptiumEZ is a trademark of Abbott Laboratories
 PowerHeart is a registered trademark of Cardiac Science Corp.
 R.A.T.E. is a trademark of BestMed
 Suntech is a registered trademark of Suntech Medical
 Unistik is a registered trademark of Owen Mumford Ltd.
 Welch Allyn and SureTemp are registered trademarks of Welch Allyn, Inc.



Four easy steps. One accurate reading.

Medline's Tympanic Thermometer with Quick Probe Release

Get fast, accurate and easy-to-read results with Medline's Tympanic Thermometer with Quick Probe Release. This ear thermometer is designed as a comfortable and convenient method for taking a person's body temperature. No embarrassing or invasive procedures; you simply place the thermometer in the ear to get a reading within seconds. It's perfectly suited for both home and clinical use on patients of all ages. Plus, its unique mechanical probe release helps reduce cross contamination, allowing you to easily change probe covers with every use!

Why an ear thermometer?

It turns out that the eardrum is an extremely accurate point to measure body temperature from because it is recessed inside the head (just like your tongue). Instead of touching the eardrum with a thermometer, the tympanic thermometer uses an infrared sensor to measure the temperature of energy radiating from the eardrum.



Apply disposable cover onto the tympanic thermometer.

Turn device "on" by pressing button on front.



Place thermometer in ear and depress top button to take reading.

After reading, flip the blue tab to release the disposable probe cover without contact.

Innovative, quick-response ear thermometer features:

- ✓ One-touch operation
- ✓ Mechanical probe cover quick-release allows easy probe cover change
- ✓ Includes 100 probe covers
- ✓ Safe to use on patients of any age
- ✓ Readings in Fahrenheit and Celsius
- ✓ Takes an accurate, non-invasive temperature from the tympanic membrane

| Item # | Description | Pkg. |
|---------|--------------|--------|
| MDS9700 | Thermometer | ea |
| MDS9701 | Probe covers | 100/bx |



Medline United States
1-800-MEDLINE (633-5463)

www.medline.com | info@medline.com

Medline Canada
1-800-396-6996

www.medline.ca | canada@medline.com

Medline México
01-800-831-0893

www.medlinemexico.com | mexico@medline.com



©2009 Medline Industries, Inc. One Medline Place, Mundelein, IL 60060
Medline is a registered trademark of Medline Industries, Inc.
MKT209138/LIT007A/10M/ZPN/2