



# Alternating Pressure, Low Air Loss and Self-Adjusting Systems

Designed Specifically  
for Long Term Care  
and Homecare





# Medline Specialty Support Services

*Medline is proud to introduce our comprehensive line of Specialty Support Surfaces. We understand that every facility has a variety of clinical needs, and that is why we offer options in every major therapeutic modality including:*

- “True” Low Air Loss
- Alternating Pressure
- Alternating Pressure with Low Air Loss
- Non Powered Self Adjusting
- Self Adjusting with Alternating Pressure
- Bariatric Alternating Pressure or Low Air Loss

*Medline’s complete line of Specialty Support Surfaces are recommended to help prevent and treat Stage 1 through Stage 4 pressure ulcers. All our products are incredibly easy to use and are virtually maintenance free. Our systems are reliable, easy to clean and built to last. Let Medline help your facility break your reliance on expensive rental mattresses. Come experience what hundreds of facilities already have by owning Medline Support Surfaces:*

- Huge Costs Savings versus rentals
- Enhanced Quality
- Positive Clinical Outcomes
- Lowest Pricing Available – Guaranteed!



## Medline Supra DPS

*Advanced Alternating Pressure & Low Air Loss Therapy*

This state-of-the-art system features a digital pump, which is microprocessor controlled. Its advanced pressure sensor technology reacts to resident movement automatically readjusting air- flow for optimal performance.

- Functions include Alternating Pressure/ Low Air Loss/Static Float/Auto Firm/Seat inflation
- Choice of four alternating pressure cycle times
- Upgraded low friction, anti-shear “all-way” stretch cover

- Cell on cell mattress design prevents “bottoming out” and stays inflated during power outages
- Panel lock out feature to protect pressure settings
- Quick Connector allows for easy set up and keeps mattress inflated during transport
- Power failure & low pressure alarms with audio and visual indicators
- Comes with carry bag for easy transport or storage
- 400 lb. weight capacity

Medline Supra DPS  
MDT24SUPRADPS




Cell on Cell design prevents bottoming out



Ventilated low air loss air cells help manage moisture



User-friendly digital pump allows you to make adjustments with the press of a button



Quick connector allows for easy set up and keeps mattress inflated during transport



Medline Supra CXC  
MDT24SUPRACXC



## Medline Supra CXC

*Alternating Pressure & Low Air Loss Therapy*

This easy to use mattress replacement system is suitable for any Nursing Home or Homecare wound treatment protocol. This feature driven system is both trouble free and cost effective – and is Medline’s best selling powered mattress.

- Functions include Alternating Pressure/ Low Air Loss/ Static float
- Easy to clean, waterproof, breathable nylon cover reduces shear and friction
- Cell on cell mattress design prevents “bottoming out” and stays inflated during power outages
- Whisper quiet pump with easy to use weight dial for customized therapy
- Low pressure alarm with audio and visual indicators
- Comes with carry bag for easy transport or storage
- 300 lb. weight capacity



Ventilated low air loss air cells help manage moisture



Weight dial makes set up fast and easy



Comes with nylon carry bag for storage or easy transport

## Medline MedTech 8000®

*Alternating Pressure Therapy*

This is Medline’s most economical mattress replacement system. It features twenty horizontal eight-inch air cells which provides optimal comfort and weight distribution.

- Functions include Alternating Pressure/Static Float
- Alternating cells provide extremely low interface pressures
- Easy to clean, waterproof, breathable nylon cover reduces shear and friction
- Whisper quiet pump, features static mode designed for dressing changes and pain management
- Low pressure alarm with audio and visual indicators
- 300 lb. weight capacity

Medline MedTech 8000  
MDT248000



Control unit provides easy to use weight dial for customized therapy

## Medline MedTech Air

### "True" Low Air Loss Therapy

This blower based "True" low air loss mattress is our most advanced low air loss system. It not only effectively manages moisture, it provides exceptional protection against pressure, shearing and friction. Recommended for patients with any stage pressure ulcer, myocutaneous skin flaps and skin grafts.

- Functions include "True" low air loss to manage moisture, "firm mode" for dressing changes and repositioning, and "upright mode" for added support when the head of the bed is elevated
- "True" low air loss is distributed equally through the 16 individual air cells wicking moisture away from the patient
- User friendly microprocessor controlled digital pump provides rapid inflation (less than 60 seconds) and easy set up
- Integrated side bolsters keep patients centered in bed reducing risk of entrapment
- Foot pillow helps to prevent unnecessary skin breakdown due to patient migration
- Low shear and friction cover is easy to clean and fully encases the air cells
- 400 lb. weight capacity

Medline MedTech Air  
MDT24MEDAIR



Digital pump is easy to use, quiet and durable

Raised Edge mattress available with our MedTech Air and MedTech Air Max models. Call your Medline representative for details.

## Medline MedTech Air Max

### "True" Low Air Loss Therapy

MedTech Air Max provides all the benefits of the MedTech Air plus these special features:

- System Lock: Secures the control panel and pressure settings preventing accidental adjustments
- Auto Comfort: When activated, senses patient's weight distribution and automatically adjusts pump to optimal pressure settings
- Alternating Pressure (pulsation): Provides three cycle settings that continually readjust the system to ensure patient comfort
- 400 lb. weight capacity
- MDT24MEDAIRMX



Medline MedTech Air Max System MDT24MEDAIRMX

## Medline MedTech Bari

### Bariatric "True" Low Air Loss

MedTech Bari provides all the benefits of the MedTech Air plus these special features:

- Larger width and depth mattress to accommodate patients up to 700 lbs
- Quickset: Allows you to input the patient's actual weight so you can preset the pressure
- Digital Pump: Offers specialized software to customize therapy
- Alternating Pressure (pulsation): Provides three cycle settings that continually readjust the system to ensure patient comfort



Medline Bari System MDT24B39, MDT24B42, MDT24B48, MDT24B54

MDT24B39	39" Low Air Loss System
MDT24B42	42" Low Air Loss System
MDT24B48	48" Low Air Loss System
MDT24B54	54" Low Air Loss System



Digital pump is easy to use, quiet and durable

Bariatric APM  
MSCBARIAPM42  
MSCBARIAPM48

## Medline Bariatric APM Alternating Pressure Therapy

Our most economical bariatric mattress replacement system provides both comfort and therapy. Features including adjustable cycle times and digital display, this system is a cost effective skin protection and treatment product.

- Functions include Alternating Pressure/Static Float/Autofirm
- Adjustable cycle times and pressure settings for customized fit
- Bottom cells stay inflated during power outages to protect patient
- One hand quick connect tubing
- Low pressure alarm
- Easy to clean, waterproof, breathable nylon cover reduces shear and friction
- Available in 42 inch (MSCBARIAPM42) and 48 inch (MSCBARIAPM48) widths
- 750 lb. weight capacity

## Medline MedTech 5000® Overlay Alternating Pressure & Low Air Loss Therapy

When budgets are an issue, our MedTech 5000 Mattress Overlay is the answer. This high quality, dual therapy product features 5 inch air cells versus a 4 inch air cell mattress. That's 20% more air and better flotation than competitive models.

- Functions include Alternating Pressure/Low Air Loss
- Whisper quiet pump with easy to use weight dial for customized therapy
- Low pressure alarm with audio and visual indicators
- Easy to clean, waterproof, breathable nylon cover reduces shear and friction
- 300 lb. weight capacity

Medline MedTech 5000  
MDT245000



Ventilated low air loss air cells help manage moisture



MedTech 5000 overlay fits on any standard healthcare mattress


Control unit features easy to use weight dial for customized therapy

## Medline's EqualizeAir Convertible Mattress Systems

*A new era in self adjusting, non powered pressure redistribution mattresses*

Medline is proud to offer the versatility of the EqualizeAir mattresses. Offered in two versions, both operate as a dynamic self-adjusting non-powered mattress replacement system or an alternating pressure mattress replacement with the addition of the optional control unit. The EqualizeAir mattresses can be utilized for high- risk patients for prevention through stage IV pressure ulcers. Now you can have it all -one low cost, highly effective pressure redistribution product to meet all of your support surface needs.

- Unique valve system constantly adjusts pressures by allowing air to move in and out of the mattress to conform to the specific weight, shape and contour of each patient.
- Constantly adjusting to movement, the valves maintain optimal pressure redistribution.
- Independent zones, created by using different densities of foam in the air tubes, form an anatomically designed profile that mirrors the shape of the human body.
- The simplicity of a therapeutic foam mattress – easy to clean, easy to use, fully enclosed and maintenance free – with the therapy of an air mattress.



Entire top layer features conforming visco elastic memory foam



Heel section features gently sloping air tubes which significantly reduce interface pressures in this sensitive area

EqualizeAir 12000  
MDTE12C3576 35" x 76" x 7"  
MDTE12C3580 35" x 80" x 7"  
MDTE12C3584 35" x 84" x 7"

## Medline EqualizeAir 12000

*Self-Adjusting Therapy*

- (12) 5 inch air cells dynamically adjust to patient profile
- Six anatomically designed pressure relief zones
- Two layer foam topper features enveloping visco elastic memory foam
- Heel zone with sloping air tubes reduces pressure in critical heel area
- Firm perimeter helps keep patients centered and aids in ingress and egress
- Zippered stretch cover allows for better patient immersion and features a non skid bottom with transport handles
- Optional pump adds alternating pressure therapy
- 500 lb. weight capacity

## Medline EqualizeAir 9000

*Self-Adjusting Therapy*

- (9) 5 inch air cells dynamically adjust to patient profile
- Four anatomically designed pressure relief zones
- Two inch foam top comfort layer features super soft heel zone with built in heel slope
- Firm perimeter helps keep patients centered and aids in ingress and egress
- Zippered stretch cover allows for better patient immersion and features a non skid bottom with transport handles
- Optional pump adds alternating pressure therapy
- 500 lb. weight capacity



By adding optional pump, the EqualizeAir 9000 or 12000 can be converted to an alternating pressure system featuring a 5 minute cycle time

EqualizeAir 9000  
MDTE9C3576 35" x 76" x 7"  
MDTE9C3580 35" x 80" x 7"  
MDTE9C3584 35" x 84" x 7"



# Support Surfaces Features Guide

	Mattress Replacement	Alternating Pressure	Low Air Loss	Blower Based Low Air Loss	Digital Pump	Lock Out Feature	Non-Powered Self Adjusting	Bariatric
MedTech 5000								
MedTech 8000								
Supra CXC								
Supra DPS								
Medline Bariatric APM								
MedTech Air								
MedTech Air Max								
Medline MedTech Bari								
EquilizeAire 9000		*						
EquilizeAire 1200		*						

\*Optional.

*Support Surfaces are a critical part of any wound care protocol. In addition to the low air loss and alternating pressure mattresses shown here, Medline Industries, Inc offers pressure redistribution therapeutic foam mattresses and a full line of pressure reduction seating and positioning products.*

*Call your Medline Representative for more information.*




1-800-MEDLINE (633-5463)  
www.medline.com


©2008 Medline Industries, Inc.  
Medline is a registered trademark of Medline Industries, Inc.  
MedTech 5000® and MedTech 8000® are registered trademarks of Medline Industries, Inc.

We reserve the right to correct any errors that may occur within this brochure.

MKT208060/LIT054R/10M



Medline seating and positioning products – a comprehensive line to meet all of your pressure redistribution needs.



Medline Advantage Mattresses – Introducing Octa-zoned technology in Health Care Mattresses from Medline. The Advantage is Yours!