



# Disposable Protective Apparel

ALL-AROUND PROTECTION





# Why are more hospitals and purchasing groups choosing Medline for their all-around protection?

Medline is a leader in protective apparel with the broadest product line and low acquisition costs. If high costs are plaguing your facility, Medline can offer you solutions to bring those costs under control without sacrificing the protection of your caregiver.

Our protective apparel line includes generous-sized footwear, protective headwear, gowns that meet AAMI PB70 requirements, protective gowns for all needs, a full range of protective coveralls and comfortable, colorful lab coats and lab jackets.

We know you rely on disposable apparel to keep your staff both protected and comfortable. For this reason,

**OUR PRODUCTS PROVIDE  
the highest protection  
with the greatest comfort  
at the lowest cost**

we use the finest materials available for protection which offer the greatest level of comfort. We apply the highest quality standards to our products and in our manufacturing facilities. Medline is ISO 13485-2003 certified. And we know our business. Medline has been manufacturing protective apparel for over 40 years.

**AAMI Education** 4  
**AAMI Isolation Gowns** 4-5



**Isolation and Protective Gowns** 6-8



**Footwear** 9-10



**Headwear** 11-12



**Scrubwear** 12



**Coveralls** 13



**Lab Coats** 14



**Lab Jackets/Aprons** 15



# Guide to Protective Materials

To eliminate confusion, here's a brief description of the main materials we use to produce Medline protective apparel. Refer back to this page as you look through the brochure.

## Spunbond Polypropylene

The standard fabric for basic protective apparel; formed by bonding fibers together to form a single layer of breathable, woven-like material. Its main advantages are economy and comfort, but it also provides light fluid protection.

## Multi-Ply Material

This fabric consists of inner layers of meltblown polypropylene sandwiched between outer layers of spunbond polypropylene. Strong, fluid-resistant and breathable, this multi-layer material provides an optimal mix of protection and comfort. Appropriate for extended wear situations involving low to moderate fluid contact. Available in lightweight, medium weight and environmentally preferred weight.

## Coated Polypropylene

This fabric is created by coating soft, spunbond polypropylene with a layer of impervious polyethylene. This material is ideal for any situation demanding a high level of fluid protection. Available in medium weight and heavy weight materials.

## Breathable Laminate









Microporous film bonded to nonwoven material provides fluid-proof protection and great breathability. This material is ideal for extended wear conditions where fluid contact is expected.

## Polyethylene

This single layer of polyethylene film provides tremendous barrier properties. This material is ideal for situations when high levels of fluids are anticipated.

## Prevention Plus

This impervious and breathable material consists of a patent-pending hybrid film laminated between two layers of bi-component spunbond material. This fabric is ideal when impervious protection is needed.



# Comfort Drives Compliance/AAMI Gowns

## Levels of Barrier Protection (from AAMI PB70:2003)

Medline Isolation Gowns with color-coded neck bindings meet the test requirements set forth by the Association for the Advancement of Medical Instrumentation (AAMI) PB70.

The chart below outlines AAMI PB70 guidelines. The color-coding has been developed by Medline to differentiate each level and corresponds with the color-coded neck bindings on Medline's AAMI Compliant Isolation Gowns.

ANSI/AAMI PB70 Barrier Performance	Test Method	Test Definition	Requirement
<b>Level 1</b>	Water Resistance: Impact Penetration AATCC 42	AATCC 42 Measures the resistance of fabrics to the liquid penetration of water by impact.	<b>Water Impact <math>\leq</math> 4.5 g</b>
<b>Level 2</b>	Water Resistance: Impact Penetration AATCC 42  Water Resistance: Hydrostatic Pressure AATCC 127	AATCC 42 Measures the resistance of fabrics to the liquid penetration of water by impact.  AATCC 127 Measures the resistance of fabrics to the liquid penetration of water by impact under constant and increasing hydrostatic pressure.	<b>Spray Impact <math>\leq</math> 1.0 g</b> <b>Hydrostatic Pressure <math>\geq</math> 20 cm</b>
<b>Level 3</b>	Water Resistance: Impact Penetration AATCC 42  Water Resistance: Hydrostatic Pressure AATCC 127	AATCC 42 Measures the resistance of fabrics to the liquid penetration of water by impact.  AATCC 127 Measures the resistance of fabrics to the liquid penetration of water by impact under constant and increasing hydrostatic pressure.	<b>Spray Impact <math>\leq</math> 1.0 g</b> <b>Hydrostatic Pressure <math>\geq</math> 50 cm</b>
<b>Level 4</b>	Viral Penetration ASTM F1671	ASTM F1671 Measures the resistance of fabrics to bloodborne pathogens using viral penetration at 2 psi and ambient pressure.	<b>Totally Impervious</b>



NONLV100

## Level 1 Protection

**Breathable Polypropylene maximizes comfort in Medline's Classic Isolation Gowns**

Item No.	Description	Pkg
NONLV100	Lightweight Multi-Ply Material, Elastic Wrists, Tape Tab Neck, Yellow, Regular Size, Flat Pack	10/bg, 50/cs
NONLV100XL	Lightweight Multi-Ply Material, Elastic Wrists, Tape Tab Neck, Yellow, X-Large Size, Flat Pack	10/bg, 50/cs

## Level 2 Protection

Medline's Classic Yellow Medium Weight Multi-Ply Gowns are breathable, flexible and strong enough for the toughest situations



NONLV200

Item No.	Description	Pkg
NONLV200	Elastic Wrists, Tape Tab Neck, Regular Size, Flat Pack	10/bg, 100/cs
NONLV200XL	Elastic Wrists, Tape Tab Neck, X-Large Size, Flat Pack	10/bg, 100/cs
NONLV200XXXL	Elastic Wrists, Tape Tab Neck, XXX-Large Size, Flat Pack	10/bg, 100/cs
NONLV200F	Elastic Wrists, Tape Tab Neck, Regular Size, Individually Folded-Dispensing Box	8/bx, 40/cs
NONLV200FXL	Elastic Wrists, Tape Tab Neck, X-Large Size Individually Folded-Dispensing Box	8/bx, 40/cs
NONLV200TL	Thumb Loop Wrists, Tape Tab Neck, Regular Size Individually Folded-Dispensing Box	8/bx, 40/cs
NONLV200TLXL	Thumb Loop Wrists, Tape Tab Neck, X-Large Size Individually Folded-Dispensing Box	8/bx, 40/cs

New! Level 2 Environmentally Preferred Gowns reduce waste in landfills and provide greater comfort



Item No.	Description	Pkg
NONLV210	Medium weight Multi-Ply Material, Elastic Wrists, Tape Tab Neck, Yellow, Regular Size, Flat Pack	10/bg, 100/cs
NONLV210XL	Medium weight Multi-Ply Material, Elastic Wrists, Tape Tab Neck, Yellow, X-Large Size, Flat Pack	10/bg, 100/cs



NONLV350

## Level 3 Protection

Premium comfortable and breathable Multi-Ply Gowns provide a high level of protection against fluids

Item No.	Description	Pkg
NONLV315	Elastic Wrists, Tape Tab Neck, Yellow Regular Size, Flat Pack	10/bg, 100/cs
NONLV315XL	Elastic Wrists, Tape Tab Neck, Yellow X-Large Size, Flat Pack	10/bg, 100/cs

Medline's unique Microporous Breathable Gowns are ideal for decontamination

Item No.	Description	Pkg
NONLV350	Microporous Breathable Laminate Material, Knit Cuffs, Tape Tab Neck, White, Regular Size, Flat Pack	10/bg, 50/cs
NONLV350XL	Microporous Breathable Laminate Material, Knit Cuffs, Tape Tab Neck, White, Regular Size, Flat Pack	10/bg, 50/cs

# Isolation Gowns



NON27SMS2

## The gowns on this page are recommended as:

- Cover gowns
- Standard isolation/Nursing gowns
- Visitor gowns

## Medline's Classic, Fluid-Resistant Multi-Ply Gowns are breathable, flexible and strong enough for tough duties

- Made from lightweight or medium weight fluid-resistant multi-ply material.
- Generously sized for full coverage and flexibility.
- Waist and neck ties ensure a secure fit.
- Lightweight gowns are available with elastic wrist or knit cuff.
- Ten gowns per bag, flat pack.

Item No.	Description	Pkg
NON27SMS2	Medium weight, Yellow, Regular Size	100/cs
NON27SMS2XL	Medium weight, Yellow, X-Large Size	100/cs
NON27SMS2XXXL	Medium weight, Yellow, XXX-Large Size	100/cs
NON27236	Lightweight, Yellow, Regular Size	50/cs
NON27236XL	Lightweight, Yellow, X-Large Size	50/cs
NON27239Y	Lightweight, Yellow, Knit Cuffs, Regular Size	50/cs
NON27239YXL	Lightweight, Yellow, Knit Cuffs, X-Large Size	50/cs

## Medline's Basic Polypropylene Gowns offer cover protection

- Made from lightweight, fluid-repelling polypropylene material.
- Gowns are generously sized for full coverage and flexibility.
- Available with neck ties and elastic wrists or in an apron-style over-the-head design with thumb-loop wrists.
- Ten gowns per bag, flat pack.

Neckties and Elastic Wrists		
Item No.	Description	Pkg
CRI4000	Yellow, Regular Size	50/cs
CRI4001	Yellow, X-Large Size	50/cs
CRI4000B	Blue, Regular Size	50/cs

Over-the-Head Neck and Thumb-loop Wrists		
Item No.	Description	Pkg
CRI4200	Yellow, Regular Size	50/cs
CRI4201	Yellow, X-Large Size	50/cs



CRI4000



## Questions? Comments?

Contact your local Medline sales representative or call 1-800-MEDLINE. You can also visit us online at [www.medline.com](http://www.medline.com)



NON27114

**The gowns on this page are recommended for:**

- O.B.
- Laboratory
- Decontamination
- High fluid situations

**Medline’s Classic Polyethylene-Coated Polypropylene Gowns offer excellent fluid protection for worry-free job performance**

- Made from blue polyethylene-coated polypropylene material.
- High level of fluid protection.
- Waist and neck ties ensure gown remains securely in place.
- Available in open or closed back with knit or elastic cuffs.
- Ten gowns per bag, flat pack.

Item No.	Description	Pkg
NON27114	Closed Back, Elastic Cuffs, Regular Size	50/cs
NON27114XL	Closed Back, Elastic Cuffs, X-Large Size	50/cs
NON27116	Closed Back, Knit Cuffs, Regular Size	50/cs
NON27116XL	Closed Back, Knit Cuffs, X-Large Size	50/cs
NON27116XXL	Closed Back, Knit Cuffs, XX-Large Size	50/cs
NON27279	Open Back, Knit Cuffs, Regular Size	50/cs
NON27279XL	Open Back, Knit Cuffs, X-Large Size	50/cs
NON27279OH	Open Back, Over-the-Head Style, Knit Cuffs, Regular Size	50/cs
NON27279OHTLX	Open Back, Over-the-Head Style, Knit Cuffs, X-Large Size	50/cs

**Medline’s Standard Polyethylene-Coated Polypropylene Gowns provide protection against high levels of fluid**

- Made from medium weight yellow polyethylene-coated polypropylene material.
- High level of fluid protection.
- Waist and neck ties ensure gown remains securely in place.
- Available with closed back and knit or elastic wrist cuffs.
- Ten gowns per bag, flat pack.

Item No.	Description	Pkg
CRI4010	Elastic Cuffs, Yellow, Regular Size	50/cs
CRI4011	Elastic Cuffs, Yellow, X-Large Size	50/cs
CRI4010CF	Knit Cuffs, Yellow, Regular Size	50/cs
CRI4011CF	Knit Cuffs, Yellow, X-Large Size	50/cs

CRI4010



# Protective Gowns



NON27457

## Prevention Plus Gowns offer protection against heavy fluids

- Made from blue impervious breathable material.
- High level of fluid protection.
- Hook and loop neck closure for easy on and off.

### Chemotherapy drug tested

Prevention Plus gowns are tested against the top chemotherapy drugs. See chart below.

The Prevention Plus Material has been tested against the following chemotherapy drugs per ASTM F739-99a:

- |                              |                       |
|------------------------------|-----------------------|
| • Dacarbazine                | • Etoposide (Toposar) |
| • Carmustine (BiCNU)         | • Paclitaxel (Taxol)  |
| • Cyclophosphamide (Cytoxan) | • Thio-Tepa           |
| • Doxorubicin Hydrochloride  | • Mitomycin C         |
| • 5-Fluorouracil             | • Vincristine Sulfate |
| • Cisplatin                  |                       |

*Test data available upon request*

Item No.	Description	Pkg
NON27457	Knit Cuffs, Blue, Regular Size	30/cs
NON27457XL	Knit Cuffs, Blue, X-Large Size	30/cs
NON27457XXL	Knit Cuffs, Blue, XX-Large Size	30/cs

## Save time and prevent “sleeve-slide” with Thumbs Up® Polyethylene Gowns

- Made from blue polyethylene material.
- High-level fluid protection.
- Apron-style neck enables gown to be placed on and off quickly.
- Specially designed thumb loop style prevents sleeves from sliding up when putting on gloves.
- Specially folded and packaged in individual pouches in a box.

Item No.	Description	Pkg
NONTH150	Thumbs Up style, Regular Size	15/bx, 75/cs
NONTH200	Thumbs Up style, X-Large Size	15/bx, 75/cs
NONTH100B	Thumbs Up style, Regular Size, Bulk	100/cs
NONTH200B	Thumbs Up style, X-Large Size, Bulk	75/cs



NONTH150



CRI5000

## Stay protected with Standard Polyethylene Thumbs Up® gowns

- Made from blue coated polyethylene material.
- High-level fluid protection.
- Apron-style neck enables gown to be placed on and off quickly
- Thumb loop style wrists keep sleeves from sliding up when putting on gloves.
- Individually folded and packaged in a dispenser box.

Item No.	Description	Pkg
CRI5000	Thumbs Up style, Regular Size	15/bx, 75/cs
CRI5001	Thumbs Up style, X-Large Size	15/bx, 75/cs
CRI5002	Thumbs Up style, XX-Large Size	15/bx, 75/cs



### Non-Skid Multi-Layer Shoe Covers are ideal for extended wear situations involving low to moderate fluid content

- Made from highly breathable, fluid-resistant material.
- Sport size available for larger, athletic-type shoes.

Item No.	Description	Pkg
NON28858	Regular Size (up to men's size 12)	100/bx, 300/cs
NON28859	X-Large Size (up to men's size 15)	100/bx, 200/cs
NON28852	Sport Size (fits athletic shoes)	100/bx, 200/cs

### Comfortable, Protective Shoe Covers - our most popular line

- Made from breathable, fluid-repellent spunbond material.
- Sport size available for larger, athletic-type shoes.
- Our most popular line of shoe covers.

Item No.	Description	Pkg
NON28758	Regular Size (up to men's size 12)	100/bx, 300/cs
NON28759	X-Large Size (up to men's size 15)	100/bx, 200/cs
NON28752	Sport Size (fits athletic shoes)	100/bx, 200/cs

NON28759

### Boundary® Shoe Covers with non-skid bottom

- Made from breathable, fluid-repellent spunbond material.
- Non-skid bottom.

Item No.	Description	Pkg
CRI2002	Regular Size (up to men's size 12)	100/bx, 300/cs
CRI2003	X-Large Size (up to men's size 15)	100/bx, 200/cs

CRI2003

### Boundary® Shoe Covers with smooth bottom

- Made from breathable, fluid-repellent spunbond material.
- Smooth bottom for situations when traction is not essential.

Item No.	Description	Pkg
CRI2000	Regular Size (up to men's size 12)	100/bx, 300/cs
CRI2001	X-Large Size (up to men's size 15)	100/bx, 200/cs

### Boundary® Polyethylene Shoe Covers resist high levels of fluids

- Made from CPE plastic material.
- Smooth bottom.

Item No.	Description	Pkg
CRI2010	One Size Fits Most, Blue	1000/cs

CRI2010

# Footwear



**NON27348P**

## Prevention Plus Impervious Breathable Boots

- Rated number one for breathability and complete protection
- Made from Prevention Plus material, which is impervious yet breathable.
- Hook-and-loop closure ensures a snug fit.
- Foam strip non-skid bottom.
- Regular up to men's size 12, XL up to men's size 15.

Item No.	Description	Pkg
NON27348P	Knee High, Regular Size	50/bx, 150/cs
NON27348PXL	Knee High, X-Large Size	50/bx, 150/cs



**NON27144**

## Non-Skid Fluid-Proof Boot Covers

- Ultrasonically sealed on all sides.
- Seamless; no holes that allow leaking.
- Convenient hook-and-loop strap allows wearer to adjust the fit.
- Strap is heat-sealed to boot.
- Available with non-skid pattern or foam strip.
- Regular up to men's size 12, XL up to men's size 15.

Item No.	Description	Pkg
NON27144	Regular Size Non-Skid bottom	50/bx, 150/cs
NON27144XL	X-Large Non-Skid bottom	50/bx, 150/cs
NON27150	Regular Foam Strip Non-Skid bottom	50/bx, 150/cs
NON27150XL	X-Large Foam Strip Non-Skid bottom	50/bx, 150/cs



**NON27143XL**

## Non-Skid Multi-Layer/Poly Boot Covers

- Breathable, fluid-resistant material above the ankle.
- Impervious, coated polypropylene below the ankle for maximum fluid protection.
- Heat-sealed seams reduce the chance of fluid penetration.
- Regular up to men's size 12, XL up to men's size 15.

Item No.	Description	Pkg
NON27143	Knee High, Regular Size	50/bx, 150/cs
NON27143XL	Knee High, X-Large	50/bx, 150/cs

## Non-Skid Multi-Layer/Poly Ankle Covers

- Designed for users who want full foot coverage without wearing a boot cover.
- Impervious bottom with breathable material above the ankle.
- Same heat-sealed seams and non-skid tread as boot cover design.
- Regular up to men's size 12, XL up to men's size 15.

Item No.	Description	Pkg
NON27145	Ankle High, Regular Size	100/bx, 300/cs
NON27145XL	Ankle High, X-Large	100/bx, 300/cs



**NON27145**



NON28625

## Sheer-Guard Surgeons' Caps

You're sure to find the style you need in our expanded line of surgeons' caps. Choose from three different materials:

**SRM** - Scrim reinforced material, a paper-like absorbent material.

**Multi-Layer Polypropylene** - Highly fluid-resistant, yet breathable.

**Spunlace** - Extremely soft, comfortable and most absorbent.

Reorder No.	Side Material	Crown Material	Absorbent	Color
NON28625	SRM	Polypropylene	Yes	Medium Blue
NON28626	Multi-Layer	Polypropylene	No	Dark Blue
NON28628	Spunlace	Polypropylene	Yes	Medium Blue

*All caps tie in back. All caps are packaged 100/bx, 500/cs.*



NONSH100C

## Medline Surgeons' Hoods

- Complete head coverage.
- Available in multi-layer material or comfortable rayon
- Ties in back of head.

Item No.	Description	Pkg
NONSH100C	Multi-Layer Material	100/bx, 300/cs
NONSH500	Rayon Material	100/bx, 300/cs



NONSH600

## Medline Surgeons' Head and Beard Covers

- Complete head coverage.
- Head and Beard Cover is made of multi-layer material with an elastic neck.
- Beard cover is made of spunbond polypropylene. Stays in place with an elastic band.

Item No.	Description	Pkg
NONSH600	Head and Beard Cover	100/bx, 300/cs
NONSH400	Beard Cover	100/bx, 1000/cs



CRI1004

## Medline's Boundary Line of Bouffant Caps

- Lightweight single-ply polypropylene material.
- Available in solid blue, green and white.
- Available in 21-in. and 24-in.
- All caps comply with NFPA standards.

Item No.	Description	Pkg
CRI1000	21" Green	100/bx, 500/cs
CRI1001	21" Blue	100/bx, 500/cs
CRI1002	21" White	100/bx, 500/cs
CRI1003	24" White	100/bx, 500/cs
CRI1004	24" Blue	100/bx, 500/cs
CRI1005	24" Green	100/bx, 500/cs



NON28229

**Pro Series Bouffant Caps provide hours of comfortable protection**

- Lint-free, latex-free caps made from lightweight multi-ply and heavyweight polypropylene.
- Stitched band keeps the cap flat and comfortable against the forehead.
- Available in solid blue and a variety of prints.
- Available in 21-, 24- and 28-in.
- All caps comply with NFPA standards.




**Lightweight Multi-Ply**

Item No.	Description	Pkg
NON28232B	21" Blue	100/bx, 500/cs
NON28233B	24" Blue	100/bx, 500/cs
NON28223B	28" Blue	100/bx, 500/cs
NON28226	24" U.S. Flag Print	100/bx, 500/cs
NON28228	21" Brush Stroke Print	100/bx, 500/cs
NON28229	24" Brush Stroke Print	100/bx, 500/cs
NON28230R	24" Breast Cancer Ribbon Print	100/bx, 500/cs
NON28233F	24" Butterfly Print	100/bx, 500/cs

**Heavyweight Polypropylene**

NON28235P	21" Multicolor Print	100/bx, 500/cs
NON28234P	24" Multicolor Print	100/bx, 500/cs

**Print Selections**



NON27202XL  
NON27203XL



NON27212M  
NON27213M

**Classic Scrub Wear offers "no show-through" complete coverage**

- Made from breathable, fluid-resistant multi-layer material.
- Anti-static treated.
- Generously cut shirts have one chest pocket and are available in either round neck or V-neck.
- Pants have one rear pocket, reinforced crotch area, and are available with drawstring or elastic waist.
- Available in blue.

When ordering Scrub Wear, add the correct size letter code (S-XXXL for Small through Triple Extra Large) to the product number. For example: NON27212L=Large.

Item No.	Description	Pkg
NON27212S-XXXL	Scrub Shirt, Round Neck	30/cs
NON27202S-XXL	Scrub Shirt, V-Neck	30/cs
NON27203S-XXL	Scrub Pants, Drawstring Waist	30/cs
NON27213S-XXXL	Scrub Pants, Elastic Waist	30/cs

**Classic Breathable Coveralls**

- Made from microporous breathable laminate material.
- Maximum protection - ideal for maternity or the O.R.
- Comfortable for extended periods.
- Available in white.

When ordering, add the correct size letter code (L-XXXL for Large through Triple Extra Large) to the product number. For example: NONCV300L = Large

Item No.	Description	Pkg
NONCV300L-XXL	Straight Sleeve and Ankle	25/cs
NONCV700L-XXXL	Elastic Wrist and Ankle	25/cs

**Classic Multi-Layer Coveralls**

- Made from fluid-resistant material.
- Good combination of comfort and protection.
- Available in white.

When ordering, add the correct size letter code (L-XXL for Large through Double Extra Large) to the product number. For example: NONCV250L = Large

Item No.	Description	Pkg
NONCV250L-XXL	Elastic Wrist and Straight Ankle	25/cs

**Classic Spunbond Coveralls**

- Made from heavyweight polypropylene material
- Breathable fluid-repellent material
- Available in white

When ordering, add the correct size letter code (L-XXXXL for Large through Four Extra Large) to the product number. For example: NONCV150L = Large

Item No.	Description	Pkg
NONCV150L-XXXXL	Elastic Wrist and Straight Ankle	25/cs

**Standard Medium Weight Spunbond Coveralls**

- Made from breathable, medium weight polypropylene material.
- Straight wrists and ankles.
- Available in white.

Item No.	Description	Pkg
CRI3000	Straight Wrist and Straight Ankle, Large	25/cs
CRI3001	Straight Wrist and Straight Ankle, X-Large	25/cs



NONCV300L



NONCV250XL



NONCV150L



CRI3000

# Lab Coats



**NONHW100XL**

## Medline's Lab Coats

Medline's anti-static Lab Coats reach to just below the knee. They feature a four-snap front closure, three reinforced pockets and a back vent for comfortable extended wear. Convenient loop for hanging. Available in white and blue.

### Heavyweight Lab Coats with Traditional Collar and Knit Cuffs

- Made from heavyweight multi-layer material that offers greater protection to fluids than a standard weight lab coat.
- Provides the same comfort as a reusable lab coat.

When ordering, add the correct size letter code (S-XL for Small through Double Extra Large) to the product number. For example: NONHW100L = Large

Item No.	Description	Pkg
NONHW100S-XL	White	30/cs

### Classic Lab Coats with Traditional Collar and Knit Cuffs

- Made from fluid-resistant multi-layer material.
- Available in white or blue.

When ordering, add the correct size letter code (S-XXXL for Small through Triple Extra Large) to the product number. For example: NONSW100L = Large

Item No.	Description	Pkg
NONSW100S-XXXL	White	30/cs
NONSW400S-XXL	Blue	30/cs



**NONSW100L**



**NONSW500XL**

### Classic Lab Coats with Knit Collar and Knit Cuffs

- Made from fluid-resistant multi-layer material.
- Available in white or blue.

When ordering, add the correct size letter code (S-XXXL for Small through Triple Extra Large) to the product number. For example: NONSW500L = Large

Item No.	Description	Pkg
NONSW500S-XXXL	White	30/cs
NONSW600S-XXL	Blue	30/cs



NONRP600L

**Classic Lab Jackets with Knit Collar and Knit Cuffs**

- Made from fluid-resistant multi-layer material.
- Anti-static Lab Jacket reaches just below the waist.
- Five-snap front closure, three reinforced pockets and a back vent for comfortable extended wear.
- Convenient loop for hanging.
- Available in white or blue.

When ordering, add the correct size letter code (S-XXL for Small through Double Extra Large) to the product number.  
For example: NONRP500L = Large

Item No.	Description	Pkg
NONRP500S-XXL	White	30/cs
NONRP600S-XXL	Blue	30/cs

**Heavyweight Apron provides maximum protection against fluids**

- Made from heavyweight 4 mil plastic
- Complete protection for the O.R., E.R., decontamination, morgue or any area where heavy fluid is present.
- Styled to easily pull over the head and tie in back.

Item No.	Description	Pkg
NON24280	41" x 62", Blue	10/bg, 100/cs

**Poly Disposable Aprons effectively protect against liquids and dirt**

- Made from polyethylene material.
- Available in two lengths to protect the torso.
- Convenient pullover styling for quick and easy removal.

Item No.	Description	Pkg
NON24272	Lightweight, 24" x 42", White	100/bx, 1000/cs
NON24274	Medium weight, 28" x 46", White	100/bx, 500/cs
NON24274W	Medium weight, 28" x 46", Individually Wrapped, Dispenser Box, White	50/bx, 250/cs



NON24280

NON24274

# Medline

## Medline = All-Around Protection in Protective Apparel

**Medline offers many options to help you select the right protective apparel for your facility**

To see samples of the materials used in our protective apparel, particularly Isolation Gowns, ask for a copy of our Isolation Gown Fabric brochure.

We can help you organize your protective apparel and other personal protective equipment (PPE) with one of our organizers. For more information, ask for a copy of our Personal Protection Organizers brochure.



**Medline Industries, Inc.**  
One Medline Place  
Mundelein, IL 60060

**Medline United States**  
**1-800-MEDLINE (633-5463)**

[www.medline.com](http://www.medline.com) | [info@medline.com](mailto:info@medline.com)

**Medline Canada**  
**1-800-396-6996**

[www.medline.ca](http://www.medline.ca) | [canada@medline.com](mailto:canada@medline.com)

**Medline México**  
**01-800-831-0898**

[www.medlinemexico.com](http://www.medlinemexico.com) | [mexico@medline.com](mailto:mexico@medline.com)

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure. ©2009 Medline Industries, Inc. Medline, Boundary and Thumbs Up are registered trademarks of Medline Industries, Inc.