



# FALL PREVENTION

ALARMS | WEARABLES | BATH SAFETY | LIFTS  
BED/BEDSIDE | AIDS FOR DAILY LIVING | STRENGTH TRAINING



**Fall Prevention** is a great concern facing health care today. The Centers for Disease Control and Prevention (CDC) has concluded that falls are the leading cause of injury deaths and the most common cause of nonfatal injuries and hospital admissions for trauma among patients 65 years and older.<sup>1</sup> Each year in the United States, nearly one third of older adults experience a fall.<sup>2</sup> Recent compliance regulations, as outlined by the Joint Commission's Patient Safety Goals and the Centers for Medicare & Medicaid Services (CMS) guidelines, have created increased urgency for the prevention of patient/resident falls.<sup>3</sup>

Medline can help health care facilities improve fall management by providing a number of useful products to be used in conjunction with a comprehensive fall management program. One of these products, patient safety alarms, has been highlighted by JCAHO as a tool that can be used to promote improved monitoring of persons at risk for falling.<sup>4</sup> When used properly, patient alarms can alert caregivers when a fall-risk candidate is on the move.

We hope that you will find this product guide to be useful in helping your facility to select products that will help to complete your Fall Prevention Program.

*The following pages contain a sampling of products to address various patient needs. Additional options may be available. Please contact your Medline sales representative or visit [medline.com](http://medline.com) should you require a product not listed in this catalog.*

1 Falls Among Older Adults: An Overview:  
<http://www.cdc.gov/ncipc/factsheets/adultfalls.htm>

2 Costs of Falls Among Older Adults:  
<http://www.cdc.gov/ncipc/factsheets/fallcost.htm>

3 2009 National Patient Safety Goals:  
[http://www.jointcommission.org/GeneralPublic/NPSG/09\\_npsgs.htm](http://www.jointcommission.org/GeneralPublic/NPSG/09_npsgs.htm),  
 Hospital-Acquired Conditions:  
[http://www.cms.hhs.gov/HospitalAcqCond/06\\_Hospital-Acquired\\_Conditions.asp](http://www.cms.hhs.gov/HospitalAcqCond/06_Hospital-Acquired_Conditions.asp)

4 FAQs for the 2007 National Patient Safety Goals:  
[http://www.jointcommission.org/NR/rdonlyres/D4844675-25D7-4B5B-A47D-C549D939F9E5/0/07\\_NPSG\\_FAQs\\_9.pdf](http://www.jointcommission.org/NR/rdonlyres/D4844675-25D7-4B5B-A47D-C549D939F9E5/0/07_NPSG_FAQs_9.pdf)



	Patient Safety Alarms	3-5
	Patient Safety Wear	6
	Bath Safety	7-9
	Lifts/Stand Assists/Slings	10-11
	Beds/Bedside/Mattresses	12-14
	Aids for Daily Living	15
	Strength Training	15

# PREMIUM PRESSURE-SENSING SAFETY ALARMS

The Joint Commission has listed “Reduce the risk of [patient] harm resulting from falls” as one of their National Patient Safety Goals.<sup>4</sup> One component of this goal is the assessment of assistive technologies and protective devices.<sup>4</sup> When used properly, patient alarms can alert caregivers when a fall-risk candidate is on the move.

**Quick Alert® Ultra & Ultra II** Our top sellers. Quick Alert alarms are self-activated, auto-sensing systems that are always working. Each monitor is constructed from heavy-duty plastics that securely house sophisticated electronics. Quick Alert alarms have been our best-selling pressure-sensing alarms for nearly 10 years.



## NEW! Quick Alert Ultra II Item No. MDT8175

- Adjustable delayed alarm setting limits false alarms by patients the flexibility to make brief movements away from their sensor mat
- Volume control button
- Easily connects to a nurse call system with the use of the MDT8000NURSE cable
- Protective silicone cover reduces breakage
- Adjustable strap for use with beds or chairs
- Power indicator light
- Low battery indicator light

## Patient Safety Alarm Selection Guide

For a complete guide to selecting the best alarm for your facility please see the back cover.



## Quick Alert Ultra Item No. MDT8100

Perfect for wheelchair monitoring.

- Contains all the great standard features of the Quick Alert Ultra II except volume control, delayed alarm and nurse call



## Quick Alert Ultra pressure mats & accessories

Item No.	Description	Pkg.
MDT81180C5	Chair sensor mat, 180-day warranty	5/cs
MDT81180B5	Bed sensor mat, 180-day warranty	5/cs
MDT8190C*	Chair sensor mat, 90-day warranty	1/ea 5/cs
MDT8190B*	Bed sensor mat, 90-day warranty	1/ea 5/cs
MDT8130C*	Chair sensor mat, 30-day warranty	1/ea 5/cs
MDT8130B*	Bed sensor mat, 30-day warranty	1/ea 5/cs
MDT8150NURSE	Nurse call cable for Ultra II	Each

\*Items can be ordered by adding 1 for a single mat or 5 for case of 5 mats

# PRESSURE-SENSING SAFETY ALARMS

**MDT8400/8450** All the features you need at a great price. The MDT8400 series patient alarms come packed with some of the most sought-after technological features typically found only on more expensive alarms, including auto-sensing without the need for an on/off switch. The best part is that you can choose from a number of optional accessories as needed.



## Personal Safety Monitor Item No. MDT8400

- Tamper-resistant 2-button design limits patient/resident monitoring interference
- Adjustable delayed alarm setting limits false alarms by patients the flexibility to make brief movements away from their sensor mat
- Volume control
- Easily connects to a nurse call system with the use of the MDT84NURSE cable
- Low battery indicator light
- Extra-durable wheelchair attachment clip included
- Optional: Protective cover Item No. MDT84CP  
Bed bracket Item No. MDT84BRKT

## Patient Safety Alarm Selection Guide

*For a complete guide to selecting the best alarm for your facility please see the back cover.*



## Tamper-Resistant Deluxe Personal Safety Monitor

Item No. MDT8450

*All the great features of the MDT8400 plus more.*

- Built-in protective shock-absorbing rubber bumpers surround monitor
- Wireless remote control pendant for easy resetting
- AC power adapter included — perfect for bedroom monitoring
- Optional: Bed bracket Item No. MDT84BRKT



## Personal Safety Monitor pressure mats & accessories

Item No.	Description	Pkg.
MDT8445TC*	Timed chair sensor mat, 45-day warranty	1/ea 6/cs
MDT8445TB*	Timed bed sensor mat, 45-day warranty	1/ea 6/cs
MDT84180TC*	Timed chair sensor mat, 180-day warranty	1/ea 6/cs
MDT84180TB*	Timed bed sensor mat, 180-day warranty	1/ea 6/cs
MDT84365TC*	Timed chair sensor mat, 1-year warranty	1/ea 6/cs
MDT84365TB*	Timed bed sensor mat, 1-year warranty	1/ea 6/cs
MDT84NURSE	Nurse call cable for MDT8400 & MDT8450	Each

*\*Items can be ordered by adding 1 for a single mat or 6 for case of 6 mats*

## PULL CORD SAFETY ALARMS

Numerous features for varying situations. Medline offers a variety of pull cord alarms - choose the exact feature you want in a proven safety alarm.



### Elite Personal Safety Alarm

Item No. MDT829900-6 & MDT829900-1N

*Pacemaker safe.*

- Magnet-free plastic tab on activation cord makes for safe use near pacemakers
- Volume adjustment
- Interlocking cradle and bracket for easy installation
- 5 alarm tone options
- Battery tester button
- Hook and loop fastener option

### Advantage Magnetic Alarm

Item No. MDT5000 & MDT5000Z

*Value priced - lowest cost per unit in its class.*

- Works from any angle
- Cradle and adjustable strap included
- Battery tester button

### Safety Sentry Alarm

Item No. MDT5100

*Tamper-resistant.*

- No on/off switch prevents tampering
- Volume control
- Tamper-resistant safety locking clip
- Optional: Protective cover Item No. MDT5100CP  
Bed bracket Item No. MDT84BRKT

### Deluxe Personal Safety Alarm

Item No. MDT5200

- Unique magnetic alarm with nurse call compatibility
- No on/off switch prevents tampering
- Volume control
- Tamper-resistant safety locking clip
- Optional: Protective cover Item No. MDT5200CP  
Bed bracket Item No. MDT84BRKT

## WEARABLES (PATIENT SAFETY WEAR)

Our patient safety wear is designed to help you promote safety in a discreet way. Our double-tread slippers and unique Lunar Print Gowns can be used to easily identify at-risk patients without calling unnecessary attention to them, and our transfer belts are a sturdy, simple aid for caregivers who assist patients with walking.



### Medline's Exclusive Slip-Resistant Slippers

For the ideal fall prevention program, choose our standard color-coded double tread slippers for most patients and choose red or yellow slip-resistant slippers for your "at risk" patients. This system will discreetly alert staff to patients that may require additional care during their stay.

Item No.	Description		Pkg.
MDT211218R	One-Size-Fits-Most	Red	48 pr/cs
MDT211218RXL	X-Large	Red	48 pr/cs
MDT211218RY	One-Size-Fits-Most	Yellow	48 pr/cs

### Deluxe Cut Patient Gowns with Tieside Closure for Maximum Coverage

Gowns close completely and comfortably, without embarrassing gaps, for total patient modesty. All gowns feature raglan sleeves for extra room and comfort. Neck and waist ties are made of durable 100% polyester twill tape and reinforced to the body of the gown for added durability.

Item No.	Description	Pkg.
MDT011899	Lunar Print, Demure Cloth	6 dz/cs
MDT021347	IV Gown, Lunar Print, Demure Cloth	6 dz/cs

### Transfer Belts Aid Patients in Walking

These cotton belts are designed to aid caregivers as they assist patients with walking. Fits around patient's waist. Available with metal tooth buckle or plastic quick-release buckle. Washable, 60"L (152cm).



Item No.	Description	Pkg.
MDT821203	White With Blue & Red Pinstripes	Each
MDT828203	White With Blue & Red Pinstripes	6/cs
MDT821204	Pastel Stripes, Metal Buckle	Each
MDT828204C	Pastel Stripes, Metal Buckle	6/cs
MDT821208PAS	Pastel, Cotton, Plastic Buckle	Each



### QualCare HipShield Brief Is Ideal For Those at Risk for Fractures

Here's an economical way to guard against the possibility of hip fracture. This unisex brief has closed-cell foam protection pads that shield the hip and help absorb impact in case of a fall. Pads are perforated so that skin can breathe and won't break out. Brief is made of cotton/spandex. Washable.

Item No.	Description	Pkg.
ALI76223	Small, 30"-32" (76cm-81cm) Waist	Each
ALI76224	Medium, 32"-34" (81cm-86cm) Waist	Each
ALI76225	Large, 36"-38" (91cm-97cm) Waist	Each
ALI76226	Extra Large, 38"-42" (97cm-107cm) Waist	Each

# BATH SAFETY

Medline's bath safety products have two primary goals: to make it easier for patients to get in and out of the bathroom and to help prevent patients from losing their balance during their bathroom routines. In this section, you'll learn about our bathtub bars, grab bars, toilet safety rails, bath and transfer benches and specialty commodes and toilet seats.



## Grab Bars

- Rugged 1¼" (3cm) diameter, steel bars are easy to grip and help reduce the risk of accidents.
- Knurled chrome-plated finish.
- Can be used horizontally or vertically.
- Complete with mounting hardware.

Item No.	Color	Length	Weight Capacity	Pkg.
MDS86012CHR	Knurled Chrome	12" (31cm)	250 lbs (113kg)	3/cs
MDS86016CHR	Knurled Chrome	16" (41cm)	250 lbs (113kg)	3/cs
MDS86018CHR	Knurled Chrome	18" (46cm)	250 lbs (113kg)	3/cs
MDS86024CHR	Knurled Chrome	24" (61cm)	250 lbs (113kg)	3/cs
MDS86032CHR	Knurled Chrome	32" (81cm)	250 lbs (113kg)	3/cs



## Steel Bathtub Bars

- White, enamel-coated, 1" (3cm) steel bar provides reliable patient assistance into and out of the bathtub.
- Crank-style, locking device attaches to tub walls up to 7" (18cm) wide.
- Unit is 6" (15cm) wide and extends 14" (36cm) above the tub wall.
- Equipped with rubber pads to protect tub walls.
- Weighs 12 lbs. (5gm)

Item No.	Description	Height	Length	Bracket Adjustment	Weight Capacity	Pkg.
MDS86321W	Steel Bathtub Bars	14" (36cm)	6" (15cm)	27/8" to 6" (7cm to 15cm)	250 lbs (113kg)	2/cs



## Toilet Safety Rails

- Width adjustable, 2 position arms
- Lightweight, anodized aluminum frame offers safe assistance.
- Waterfall armrests provide added patient comfort.
- Nylon-coated, steel bracket attaches unit securely to toilet.
- Rail legs raise off the floor for cleaning around toilet.
- Legs adjust 10" (25cm).

Item No.	Description	Arm Height From Floor	Depth	Distance Between Arms (Outside)	Weight Capacity	Pkg.
MDS86100RF	Toilet Safety Rails	28" to 32" (71cm x 81cm)	17½" (45cm)	17½" to 20½" (45cm x 52cm)	250 lbs (113kg)	1 pr/cs
MDS86100BRT	Replacement Toilet Safety Rail Bracket			21" to 24" (53cm x 61cm)		Each

# BATH SAFETY



MDS89740RW



MDS89745R

## Bath Benches

- Available in gray (retail package) or white (fully assembled).
- High-density, polyethylene seat and back.
- Durable, aluminum frame is corrosion resistant.
- Angled legs distribute weight over a large area to increase patient security.
- Non-slip tips on all four legs provide extra-stability and safety.
- Height adjusts in 1" (3cm) increments.

Item No.	Description	Seat Height	Seat Depth	Seat Width	Weight Capacity	Pkg.
MDS89740A	Without Back, Assembled, White	14½" to 19½" (37cm to 50cm)	12" (31cm)	20" (51cm)	250 lbs (113kg)	2/cs
MDS89740RW	Without Back*, Gray	14½" to 19½" (37cm to 50cm)	12" (31cm)	20" (51cm)	250 lbs (113kg)	2/cs
MDS89745A	With Back, Assembled, White	14½" to 19½" (37cm to 50cm)	12" (31cm)	20" (51cm)	250 lbs (113kg)	2/cs
MDS89745R	With Back*, Gray	14½" to 19½" (37cm to 50cm)	12" (31cm)	20" (51cm)	250 lbs (113kg)	2/cs

\*Retail packaging requires assembly.



MDS86952XW

## Transfer Benches

- Push-button assembly for back.
- Durable, blow-molded, plastic bench.
- Wide base, and suction-cup tips on two legs in tub offer extra stability and safety.
- Direction of seat back can be reversed easily without tools.
- Anodized, aluminum frame adjusts in 1" increments.
- Fits tubs 11" to 21¼" deep.

Item No.	Description	Seat Height	Seat Width x Depth	Overall Height	Overall Width x Depth At Base	Weight Capacity	Product Weight	Pkg.
MDS86952	Transfer Bench	16½" to 21½" (42cm to 55cm)	17" x 16½" (43cm x 42cm)	34" to 39" (83cm to 99cm)	29" x 21½" (74cm x 55cm)	300 lbs (136kg)	11 lbs (5kg)	Each
MDS86952XW	Bariatric, Plastic	16" to 20" (41cm to 51cm)	17" x 16" (43cm x 41cm)	28" to 32" (71cm to 81cm)	28½" x 20½" (72cm x 52cm)	550 lbs (249kg)	16 lbs (7kg)	Each



# BATH SAFETY



MDS89664



MDS89664XW

## Deluxe, Steel Commodes

- Back rail removes quickly and easily.
- Features a durable, plastic seat with lid.
- Plastic armrests for additional patient comfort.
- Adjustable legs allow seat to be raised or lowered in 1" (3cm) increments to sit from 16" to 22½" (41cm to 51cm) standard, 19½" to 23½" (50cm to 52cm) bariatric.
- Includes 10-quart capacity, commode bucket w/cover & splash shield.
- Non-skid, non-marring, rubber tips provide extra safety.

Commode Accessories	Item No.	Pkg.
Rubber Tips	MDS89664T	1/pk
Seat & Lid for MDS89664	MDS89664SL	4 sets/cs
Seat & Lid for MDS89664XW	MDS89664SLXW	1/set
Set of four splash guards	MDS89664SGC	4/set

Item No.	Description	Base Width	Base Depth Overall	Width Between Arms	Seat Width	Seat Depth Overall	Seat Height Min./Max.	Height Overall Min./Max.	Product Weight	Weight Capacity	Pkg.
MDS89664	Steel, 3-in-1	24" (61cm)	18" (46cm)	19" (48cm)	14" (36cm)	15" (38cm)	15" to 21" (38cm to 53cm)	23½" to 30" (60cm to 76cm)	13.2 lbs (6kg)	300 lbs (136kg)	4/cs
MDS89664XW	Bariatric	28" (71cm)	20½" (52cm)	24" (61cm)	14" (36cm)	15" (38cm)	18¼" to 23½" (46cm to 60cm)	29½" to 32½" (75cm to 83cm)	12.8 lbs (6kg)	650 lbs (295kg)	Each



MDS80314



MDS80316

## Locking, Raised Toilet Seats

- Width adjustable, 2 position arms eliminate the need for toilet safety rails, making it easier for the patient.
- Front clamping mechanism ensures positive locking onto toilet. Helps prevent shifting of toilet seat while sitting or transferring.
- Fits most toilet bowls from 11½" to 15"L (29cm to 38cm).
- Wide, contoured surface provides maximum user comfort.
- Padded support arms offer user a secure hand support and maximum comfort when using toilet seat.
- Raises seat level by 5½" (14cm).

Item No.	Description	Seat Height	Seat Width	Seat Depth	Weight Capacity	Pkg.
MDS80314	Raised Toilet Seat, w/o Arms	5½" (14cm)	16⅞" (44cm)	17⅝" (44cm)	300 lbs (136kg)	3/cs
MDS80316	Raised Toilet Seat, w/ Arms	5½" (14cm)	16⅞" (44cm)	17⅝" (44cm)	300 lbs (136kg)	3/cs

## PATIENT LIFTING/STAND ASSISTS

You can help minimize caregiver injuries by choosing specially designed patient lifts and stand assists. We offer a full line of products that deliver a safe and secure way to lift patients with impaired mobility.



### Battery-Powered Elevating Floor Lift

Attractive white finish, along with a 6 point cradle.

Item No.	Description	Pkg.
MDS400EL	400lb (181 kg) Elevating Floor Lift	Each
MDS600EL	600lb (272 kg) Elevating Floor Lift	Each



### Battery-Powered Stand Assist

Attractive white finish with removable foot platform. Platform can be removed for standing and walking practice. Sling free with purchase of this lift.

Item No.	Description	Pkg.
MDS400SA	400lb (181 kg) Stand Assist	Each
MDS600SA	600lb (272 kg) Stand Assist	Each



MDELSSCALE  
For use on any elevating floor lift  
(MDS88200D, MDS400EL,  
MDS600EL)



### Manual Stand Assist

More affordable alternative to battery-powered stand assists. Patient grabs onto an assist bar and pulls up using their own strength.

Item No.	Description	Pkg.
STA182	Manual Stand Assist	Each

# PATIENT LIFTING/STAND ASSISTS ACCESSORIES



MDSR113



MDSR115



MDSSA3



**FREE!**

## Safe 'n Easy Patient Lifting Program

FREE with any lift purchase. Ask your Medline sales representative or call 1-800-MEDLINE for more information.

### Full Body Solid Fabric Slings

Item No.	Description	Pkg.
MDSR112	450lb (204 kg) Medium	Each
MDSR113	450lb (204 kg) Large	Each
MDSR117	450lb (204 kg) X-Large	Each

### Full Body Mesh Slings

Item No.	Description	Pkg.
MDSR110	450lb (204 kg) Medium	Each
MDSR111	450lb (204 kg) Large	Each
MDSR140	600lb (272 kg) XX-Large	Each

### Full Body Mesh Slings with Commode Opening

Item No.	Description	Pkg.
MDSR114	450lb (204 kg) Medium	Each
MDSR115	450lb (204 kg) Large	Each
MDSR116	450lb (204 kg) X-Large	Each
MDSR141	600lb (272 kg) XX-Large	Each

### Stand Assist Padded Slings

Item No.	Description	Pkg.
MDSSA1	400lb (181 kg) Small	Each
MDSSA2	400lb (181 kg) Medium	Each
MDSSA3	400lb (181 kg) Large	Each
MDSSA4	600lb (272 kg) Bariatric	Each



### Transfer Sheets

Transfer sheets can be used for a variety of patient transfer movements. Repositioning, boosting, and other movements are all accomplished more easily by reducing surface friction through the application of transfer sheets.

Item No.	Description	Pkg.
TS30110	Small, 28.5" x 27.5" (72cm x 70cm)	Each
TS30120	Medium, 42.5" x 31.5" (108cm x 80cm)	Each
TS30130	Large, 53" x 39.5" (136cm x 100cm)	Each

# LOW HEIGHT BED

## Alterra 1232 Hi-Low Bed



MAX  
height of  
26"

LOW  
height of  
7.25"



### Safety and Comfort for the resident

- Four ultra-quiet, dependable DC motors provide power and offer a multitude of comfortable positioning options.
- Unique grid/rib deck for enhanced comfort.

### Ease of use for the caregiver

- Low height of 7.25" (18cm) makes Alterra 1232 bed ideal for Alzheimer's patients and any other resident at risk of falling.
- Battery backup keeps the bed functioning in the event of power outage and also during evacuation.
- High height of 26" (66cm) makes it easy to care for residents while helping to reduce the risk of back injuries.
- Alterra 1232 bed moved from lowest to highest height without any horizontal movement, eliminating wall damage.

Item No.	Description
FCE1232QO	Alterra 1232 Bed Package (includes Alterra 1232 hi-low bed, head section assist rails, head/foot boards in quartered oak finish, mounting hardware, 4 mattress slide stops)
FCE1232B	Alterra 1232 hi-low bed only
FCE1232RSRN	Soft-touch assist rails for FCE1232B (These assist rails can be mounted either on head section or on foot section of the bed deck.)
FCE1232RSR	Head section assist rails for FCE1232B
FCE1232AB	Assist bar for FCE1232B
FCE1232RHFQO	Head/foot boards in quartered oak finish with mounting hardware and mattress slide stops for FCE1232B, (same style as FCE01RHF

### Specifications:

**Height range:**  
7.25" to 26" (18cm to 66cm)

**Bed deck dimensions:**  
80" (203cm) L, 34.5" (88cm) W

**Overall length:** 87.5" (222cm)

**Weight capacity:**  
450 lb. (204kg)

**Movement styles:**  
Head, knee hinge, high/low

**Deck style:**  
Combination grid/rib

**Rolling position:**  
Move at any height

**Motors:** 4 DC

**Casters:** 8 total, 4 locking

**Accessories available:**  
Custom head/footboards, head and foot section half-length assist rails, mattress slide stops, assist bars, staff control



**Warranty information:**

**Electrical components:**  
2 years

**Bed frame:** 10 years

**Welds:** lifetime warranty

### Head & Footboard Options



### Head & Footboard Finish Options



\*Custom finish head/foot boards for FCE1232B and FCE1232HDW must be ordered with this item.

## BEDSIDE SAFETY PRODUCTS

Medline's bedside safety products offer support for independent patients and help cushion any falls. Fall mats help reduce the risk of injury in the event of a fall, and our side rails are designed to offer assistance and additional safety when exiting a bed.



### The Medline Advantage Fall Mat

It helps reduce the possibility of impact-related injuries in the event of an accidental fall from the bed. Made with absorbent, high compression foam, with a durable, easy to clean vinyl cover. Two styles to choose from: folding and non folding

Item No.	Description
<i>Folding</i>	
MDTFM24662BI	24" x 66" x 2" (61cm x 168cm x 5cm)
MDTFM36662BI	36" x 66" x 2" (91cm x 168cm x 5cm)
MDTFM48662BI	48" x 66" x 2" (122cm x 168cm x 5cm)
MDTFM24663BI	24" x 66" x 3" (61cm x 168cm x 8cm)
<i>Non Folding</i>	
MDTFM24662	24" x 66" x 2" (61cm x 168cm x 5cm)
MDTFM36662	36" x 66" x 2" (91cm x 168cm x 5cm)
MDTFM48662	48" x 66" x 2" (122cm x 168cm x 5cm)
MDTFM24663	24" x 66" x 3" (61cm x 168cm x 8cm)



**The Halo Circular Side Rail** design and vertical mounting bar eliminate entrapment concerns and allow exit points for the neck, head, and chest. Also, prevents mattress movement and gaps in the most critical bed zones. Durable steel construction with a 1,000lb (454kg) weight capacity.


Item No.	Description
MDT77101	1 Halo ring w/mounting hardware & stay bar to secure mattress
MDT77102	2 Halo rings w/mounting hardware



**Safety Zone™ Fall Mat** is unusually thin, surprisingly firm and incredibly effective. This innovative design has protected countless nursing home residents from fall-related injuries. With its unique blend of materials and construction, the Safety Zone is Medline's best floor-mat option for maximum resident protection. Heavy duty vinyl cover, bonded to provide better impact distribution. Non-slip base.

Item No.	Description
MDTFLSTRIP	24" x 70" x 1" (61cm x 178cm x 3cm) Beveled edges

# MATTRESSES



## Medline Advantage "Select" PE Mattress

Our most advanced pressure redistribution mattress features our thickest three layered design to completely conform to your resident's body. The soft, sculpted top layer of high resiliency foam redistributes pressure for residents up to 450lbs (204kg) while maintaining its comfort and durability. The visco elastic heel section cradles the bony heels to provide maximum pressure relief. The premium stretch cover is fluid proof, antimicrobial, flame retardant and allows the mattress to conform to your patient while reducing shear and friction. 15 year warranty.


Item No.	Description
MSCADVPE76FR	36" x 76" x 6.75" (91cm x 193cm x 17cm), with Safe-T-Side Raised Edges
MSCADVPE80FR	36" x 80" x 6.75" (91cm x 203cm x 17cm), with Safe-T-Side Raised Edges
MSCADVPE84FR	36" x 84" x 6.75" (91cm x 213cm x 17cm), with Safe-T-Side Raised Edges



## Safe-T-Side Mattress option

► All Medline mattresses can be ordered with the addition of our exclusive Safe-T-Side raised edges to help reduce the risk of patient falls. Both sides have a centered cutout for easy entry, exit and transfer.

► Ask your Medline sales representative or call 1-800-MEDLINE for more information.




## Medline Advantage 1000 Mattress

The Advantage 1000 provides excellent pressure redistribution and will outperform any mattress in its category. It features an Octa-Zoned Surface Sculpted topper to help relieve pressure and shear. The variable depth sculpting provides outstanding comfort while creating horizontal airflow channels to lower surface temperatures and reduce moisture. Nylex top fabric is fluid proof, antimicrobial and flame retardant. 350lb (159kg) weight capacity. 5 year warranty

Item No.	Description
MSCADV1076FR	36" x 76" x 6" (91cm x 193cm x 15cm), with Safe-T-Side Raised Edges
MSCADV1080FR	36" x 80" x 6" (91cm x 203cm x 15cm), with Safe-T-Side Raised Edges
MSCADV1084FR	36" x 84" x 6" (91cm x 213cm x 15cm), with Safe-T-Side Raised Edges

## AIDS FOR DAILY LIVING

From leg lifters to reachers, we have products to assist with the daily living activities that might have become more difficult for patients who have a limited range of motion or who are at risk for falls.



### Sock and Stocking Aid with Ridges

Item No. MDSD1220

- Provides maximum benefits at a moderate price.
- Features the popular three-finger design.
- Sock stays open as the three fingers gently bend to accommodate the contours of the heel and foot.
- Loop strap can be manipulated with only one hand and is length adjustable.

### Stick to Stand safely facilitates sit to stand motion

Item No. MDSSTICK

Developed to assist with the adjustment of beds, chairs, furniture, toileting and other equipment to achieve heights between 100% and 120% of lower leg length.

- Maintains independent transfer and mobility skills.
- Reduces risk of transfer related falls and injuries.
- Reduces the risk of caregiver strain and injury.

### Leg Lift

Item No. MDL1500

- Enables people without lower extremity strength to lift their leg and foot.
- Center section is rigid, with 2 loops on either end to go around the foot and wrist.
- Made of high-strength nylon webbing.

### Raptor® Reacher

Item No. MDL1100

- Low cost, durable, lightweight.
- Provides a tight grip on even the smallest objects.
- 24" (61cm) long.
- Weighs only 6 oz. (170gm)



## STRENGTH TRAINING

The use of strength training products for fall risk patients has been shown to promote increased strength and improved balance.

### Cando Exercise Bands and Tubing


Cando exercise bands and tubing can be used independently or as part of a progressive resistance program (PRE) by simply cutting a strip or knotting it to form a loop. Conveniently color-coded to show its resistance levels using the Thera-Band trademarked color sequence. Low-Powder Cando bands and tubing utilize high-quality latex that feature resistance properties that give its needed uniform expansion and contraction characteristics.

Item No.	Description	Item No.	Description
<i>Cando Bands</i>		<i>Cando Tubing</i>	
MDSP105210	Tan, 6 yd (5.5m)	MDSP105510	Tan, 25 Ft (7.6m)
MDSP105211	Yellow, 6 yd (5.5m)	MDSP105511	Yellow, 25 Ft (7.6m)
MDSP105212	Red, 6 yd (5.5m)	MDSP105512	Red, 25 Ft (7.6m)
MDSP105213	Green, 6 yd (5.5m)	MDSP105513	Green, 25 Ft (7.6m)
MDSP105214	Blue, 6 yd (5.5m)	MDSP105514	Blue, 25 Ft (7.6m)
MDSP105215	Black, 6 yd (5.5m)	MDSP105515	Black, 25 Ft (7.6m)
MDSP105216	Silver, 6 yd (5.5m)	MDSP105516	Silver, 25 Ft (7.6m)
MDSP105217	Gold, 6 yd (5.5m)	MDSP105520	Tan, 100 Ft (30.5m)
MDSP105220	Tan, 50 yd (45.7m)	MDSP105521	Yellow, 100 Ft (30.5m)
MDSP105221	Yellow, 50 yd (45.7m)	MDSP105522	Red, 100 Ft (30.5m)







The following Patient Safety Alarm Selection Guide can be used as a helpful tool to find the Patient Safety Alarm that best meets your needs. See pages 3-5 for detailed product information on all the Patient Safety Alarms listed below.

## Pressure Mat Sensing Systems

	Multiple Volume Control Choices	Multiple Delayed Alarm Choices	Nurse Call compatible	Included Shock-Absorbing Protective Cover	2-Year Warranty
MDT8175 	X	X	X	X	X
MDT8100 				X	X
MDT8450 	(2)	(1)	X		
MDT8400 	(2)	(1)	X		

(Numbers in parenthesis indicate number of adjustments available)

## Pull Cord Systems

	Volume Control Choices	Multiple Alarm Tone Options	Low Battery Indicator Light	Includes Bed and Chair Mounting Hardware	Pacemaker Safe (No Magnet)	Nurse Call Compatible	No On/Off Switch
MDT829900-6 & MDT829900-1N 	X	X	X	X	X		
MDT5200 	X		X	*		X	X
MDT5100 	X		X	*			X
MDT5000 & MDT5000Z 				X			

\*These alarms include a wheelchair mounting clamp only. An optional bed bracket, MDT84BRKT, can also be purchased separately.

1-800-MEDLINE (633-5463)

[www.medline.com](http://www.medline.com)

©2008 Medline Industries, Inc. Raptor is a registered trademark of Maddak, Inc. Safety Zone is a trademark of Avatar Enterprises, Inc. Quick Alert and Medline are registered trademarks of Medline Industries, Inc. One Medline Place, Mundelein, IL 60060  
We reserve the right to correct any errors that may occur in this brochure.  
MKT208354/LIT410R/ZPN6/7.5M

Medline Canada 1-800-396-6996

[www.medline.ca](http://www.medline.ca) [canada@medline.com](mailto:canada@medline.com)

Not all products are currently available for sale in Canada.  
Please contact your Account Representative for additional information.