




Medline Advantage "Select" Mattresses

All Products guaranteed
to ship within 48 hours
or receive a free
pillow for every
mattress ordered²

A New Approach To Pressure Redistribution Mattresses





Medline Advantage “Select” Mattresses

A New Approach To Pressure Redistribution Mattresses

Medline Advantage Mattresses have been setting the standard for value and performance in long term care facilities. With over 35 years of experience manufacturing healthcare mattresses, it is no surprise to our customers why Medline is the industry leader in helping facilities reach their wound prevention and treatment goals.

Today, Medline is proud to introduce our latest improvement - a “Select” group of new Advantage Mattresses. The Advantage “Select” pressure redistribution mattresses are our finest mix of clinical performance and outstanding value. With surface sculpting and premium grades of foam, the new “Select” mattresses are both comfortable and effective in redistributing pressure. Now you can have it all when it comes to choosing a therapeutic mattress.


Advantage “Select” Mattresses – Performance, Value and Patient Comfort!

Medline Advantage “Select” Guide

Medline Item	Risk Level ¹	Weight Capacity	Warranty
Advantage “Select” PE	Single Site Stage III	450 lbs.	15 Year
Advantage “Select” VE	High Risk/Stage I & II	350 lbs.	6 Year Unconditional
Advantage “Select” SE	High Risk/Stage I & II	350 lbs.	4 Year Unconditional

Standard Features on all Medline Advantage “Select”


Three-sided zipper provides easy cover replacement



¹ The resident risk classifications are for guideline purposes only and are not claims for efficacy on individual residents or patients. Mattresses are only one element in proper wound prevention and care. A complete clinical assessment of each resident should be done before selecting any therapeutic support surface.

² All orders of 10 mattresses or less guaranteed to ship within 2 business days. A free healthcare pillow will be provided with each mattress if shipment goes out later than two business days from receipt of order.

Medline Advantage "Select" PE Mattress



Single Site Stage III

Our most advanced pressure redistribution mattress features our thickest three-layered design to completely conform to your resident's body. The soft, sculpted top layer of our finest high resiliency foam redistributes pressure for residents up to 450 lbs while maintaining its comfort and durability.

- The anatomically designed mattress features a zoned core to help cradle bony prominences. The center layer of convoluted foam promotes air flow to help manage heat build up increasing patient comfort.
- The super soft heel section features a convoluted design to provide maximum pressure relief, comfort and support.
- Nylex Super Stretch premium cover allows the mattress to conform to your resident while reducing shear. Fluid proof and antimicrobial.
- Meets federal fire standard 16 CFR 1633
- 450 pound weight capacity

Item	Description	Qty
MSCADVPE76F	36" x 76" x 6.75"	Each
MSCADVPE80F	36" x 80" x 6.75"	Each
MSCADVPE84F	36" x 84" x 6.75"	Each

Medline Advantage "Select" VE

High Risk/Stage I & II

Our high resiliency foam top layer provides excellent pressure redistribution and uncompromised comfort. The sculpted design uses cropped foam to provide air channels to keep patients cool and comfortable.

- The anatomically designed mattress core allows the body to sink into the mattress redistributing pressure away from the bony prominences reducing the risk of skin breakdown.
- The heel section features visco elastic memory foam that cradles the heels allowing for the ultimate in shear reduction and pressure redistribution.
- Nylex top fabric combines durability, performance and reduces shear. Fluid proof and antimicrobial.
- Meets federal fire standard 16 CFR 1633
- 350 pound weight capacity



Item	Description	Qty
MSCADVVE76F	36" x 76" x 6"	Each
MSCADVVE80F	36" x 80" x 6"	Each
MSCADVVE84F	36" x 84" x 6"	Each



Medline Advantage "Select" SE Mattress



High Risk/Stage I & II

When an already stretched budget makes you think you can't upgrade your mattresses, The Advantage "Select" SE is your answer. It will give you excellent pressure and shear reduction at a price that can't be beat.

- Surface sculpted top layer provides pressure redistribution while creating airflow to reduce heat and moisture buildup.
- Special convoluted heel section provides an even softer environment for sensitive heels
- Nylex top fabric combines durability, performance and reduces shear. Fluid proof and antimicrobial.
- Meets federal fire standard 16 CFR 1633
- 350 pound weight capacity

Item	Description	Qty
MSCADVSE76F	36" x 76" x 6.75"	Each
MSCADVSE80F	36" x 80" x 6.75"	Each
MSCADVSE84F	36" x 84" x 6.75"	Each

Medline Advantage Options



Saf-T-Side Raised Edge

Medline's Saf-T-Side raised edges are designed for use when bed side rails are down. Featuring an easy exit-entry cutout which enables residents to sit up and get in and out of bed when they want to, our firm Saf-T-Side raised edges help to prevent resident falls. Available on all Medline Advantage "Select" Mattresses.

When ordering **Raised Edge**, add a "R" to the end of the item number.

All Products guaranteed to ship within 48 hours or receive a free pillow for every mattress ordered²



1-800-MEDLINE (633-5463) www.medline.com
©2007 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. One Medline Place, Mundelein, IL 60060. We reserve the right to correct any errors that may occur within this brochure.

Medline Canada 1-800-396-6996
www.medline.ca canada@medline.com
Not all products are currently available for sale in Canada. Please contact your Account Representative for additional information.