


PROHIBIT[®] Series Face Masks



An Industry First!

Our Exclusive Scented Face Mask
has a Fresh Spearmint Scent.


See page 4 for details



5 Reasons that Make Medline's Prohibit Series® Face Masks Better than the Competition.

1. Incomparable Comfort and Superior Breathability

Masks that restrict air flow aren't comfortable — they're hot. Medline's PROHIBIT Series masks allow air to circulate in both directions while at the same time filtering harmful bacteria, so you stay cool, yet protected.



2. Bacterial Filtration Efficiency (BFE)

All PROHIBIT Series masks have bacterial filtration efficiencies exceeding 95%.* Every mask undergoes testing by Nelson Laboratories, Inc. to ensure that you get the very best protection available.


*Test results available upon request.

3. Particle Filtration Efficiency (PFE)

Surgical and Procedure masks exceed 98% PFE at 0.1 μm . N95 Masks are greater than or equal to 99% PFE at 0.1 μm .

4. Fluid Protection for Greater Safety

Protection from bloodborne pathogens is a serious matter. For the highest level of fluid protection, our Max X series masks pass the ASTM F1862 synthetic blood penetration at 160mm hg. These masks incorporate Drilex media, a revolutionary material prevents fluid impact penetration in either direction, yet remains very breathable.



Drilex media

5. Evaluation Program






Experience what makes Prohibit Series Face Mask better than the rest. Depending on how your account qualifies, you may be eligible to receive FREE product for evaluation. In addition, we will provide the tools and clinical support to work with staff during the evaluation. Conversion incentives may apply. Contact your Medline Sales Representative for details.





We guarantee everything we sell. Should you have any problem with a Medline product, please call us or speak to your Medline Sales Representative. We will take care of the problem right away, either by cheerfully and promptly exchanging the product or by crediting your account.

Color Coding Index

Scented Masks	
Fluid Resistant Masks	
Anti-Fog Masks	
Surgical Masks	
Specialty Surgical Masks	
Isolation & Procedural Masks	
N95 Particulate Respirators	

 All products in this catalog are latex-free.

Table of Contents

Scented Masks	4
Fluid Resistant Masks	4-5
Anti-Fog Masks	6
Surgical Masks	7
Specialty Surgical Masks	8
Isolation Masks	9
Procedural Masks	9
N95 Particulate Respirators	9
Cross Reference	10-11
Mask Accessories	12



An Industry First!

Spearmint scented mask protects against unpleasant odors in the O.R.

Scented Procedure Mask

- Unique scent of spearmint is constant breath of fresh air.
- Fluid-resistant, outer facing provides protection.
- Soft cellulose inner facing increases inner facing comfort.
- Convenient earloops are sonically sealed to the outside of the mask to virtually eliminate possible facial irritation and ensure a secure closure.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27605	Scented Procedure Mask with Earloops, Green	50/bx, 300/cs

BFE > 97%
PFE > 98%



Scented Surgical Mask

- Unique scent of spearmint is constant breath of fresh air.
- Fluid-resistant, outer facing provides protection.
- Soft cellulose inner facing increases inner facing comfort.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27610	Scented Surgical Mask with Ties, Green	50/bx, 300/cs

BFE > 97%
PFE > 98%



Max X Fluid Protection

Our Highest Level of Fluid Protection Ever!


Max X Fluid Protection Masks w/Eyeshield

- Non-glare, fog-free, wraparound eyeshield assures maximum eye and facial protection.
- Unique Drilex filter provides uncompromising fluid resistance.
- Cellulose inner and outer facing offers protection against fluid and inner facing comfort.
- Anti-fog strip of comfortable, soft foam prevents eyeglasses from fogging.

Item #	Description	Pkg
NON27710	Max X w/Eyeshield and Ties, Pattern	25/bx, 100/cs
NON27710EL	Max X w/Eyeshield and Earloops, Pattern	25/bx, 100/cs

BFE > 99%
PFE > 99%

ASTM F1862 160mm hg
Fluid Resistance



Max X Fluid Protection Masks

- Cellulose inner and outer facing offers protection against fluid and inner facing comfort.
- Unique Drilex filter provides uncompromising fluid resistance.
- Anti-fog strip of comfortable, soft foam prevents eyeglasses from fogging, thus ensuring optical clarity.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27712	Max X with Ties, Pattern	50/bx, 300/cs
NON27712EL	Max X with Earloops, Pattern	50/bx, 300/cs

BFE > 99%
PFE > 99%

ASTM F1862 160mm hg
Fluid Resistance



Fluid Resistant Masks

Great fluid protection while filtering out bacteria

BFE > 99%
PFE > 99%

ASTM F1862 120mm hg
Fluid Resistance

Fluid-Resistant Masks with Eyeshield

- Non-glare, fog-free, wraparound eyeshield assures maximum eye and facial protection.
- Polypropylene outer facing for added protection against fluids, plus soft inner cellulose facing for comfort.
- Anti-fog strip of comfortable, soft foam prevents eyeglasses from fogging.

Item #	Description	Pkg
NON27410	Fluid Resistant with Eyeshield and Ties, Purple	25/bx, 100/cs
NON27410EL	Fluid Resistant with Eyeshield and Earloops, Purple	25/bx, 100/cs



BFE > 99%
PFE > 99%

ASTM F1862 120mm hg
Fluid Resistance

Fluid-Resistant Masks

- Polypropylene outer facing for added protection against fluids, plus soft inner cellulose facing for comfort.
- Anti-fog strip of comfortable, soft foam prevents eyeglasses from fogging, thus ensuring optical clarity.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27412	Fluid Resistant with Ties, Purple	50/bx, 300/cs
NON27412EL	Fluid Resistant with Earloops, Purple	50/bx, 300/cs



Prohibit X-Tra®

Offers More Fluid Protection than Traditional Surgical Masks

BFE > 98%
PFE > 98%

ASTM F1862 80mm hg
Fluid Resistance

Prohibit X-Tra® Fluid Protection Masks w/Eyeshield

- Non-glare, fog-free, wraparound eyeshield assures maximum eye and facial protection.
- Cellulose inner and outer facing offers protection against fluid and inner facing comfort.
- Anti-fog strip of comfortable, soft foam prevents eyeglasses from fogging.

Item #	Description	Pkg
NON27405	X-Tra Protection with Eyeshield and Ties, Blue	25/bx, 100/cs



BFE > 98%
PFE > 98%

ASTM F1862 80mm hg
Fluid Resistance

Prohibit X-Tra® Fluid-Protection Masks


- Cellulose inner and outer facing offers protection against fluid and inner facing comfort.
- Anti-fog strip of comfortable, soft foam prevents eyeglasses from fogging, thus ensuring optical clarity.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27408	X-Tra Protection with Ties, Blue	50/bx, 300/cs
NON27408EL	X-Tra Protection with Earloops, Blue	50/bx, 300/cs



Anti-Fog Mask Styles

With four different ways to protect you from "fogging up" Medline gives you more choices.



Anti-Fog Masks with Adhesive Tape

- Two styles to choose from: standard adhesive or aggressive adhesive.
- Non-irritating, adhesive strip forms a tight seal to prevent eyeglass fogging and maintain optical clarity.
- Fluid-resistant facing provides excellent protection.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27371A	Adhesive with Ties, Green	50/bx, 300/cs
NON27361A	Aggressive Adhesive with Ties, Green	50/bx, 300/cs

BFE > 97%
PFE > 98%

Anti-Fog Mask with Film Strip

- A film strip is designed to conform to the face to help seal in air, reducing eyeglass and microscope fogging.
- Cellulose facing for cool comfort.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27379A	Anti-Fog Mask with Film Strip and Ties, Green	50/bx, 300/cs

BFE > 97%
PFE > 98%

Anti-Fog Mask with Foam Strip

- Thick, soft foam strip keeps condensation from escaping out the top of the mask thus preventing eyeglass fogging and ensuring optical clarity.
- Cellulose facing for cool comfort.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27378A	Anti-Fog Mask with Foam Strip and Ties, Blue	50/bx, 300/cs

BFE > 97%
PFE > 98%



Anti-Fog Chamber-Style Mask with Foam Strip

- Large breathing chamber keeps mask away from face for greater comfort and easier breathing; also allows unmuffled speech.
- Thick, foam strip keeps condensation from escaping out the top of the mask, thus preventing eyeglass fogging and ensuring optical clarity.
- Polypropylene outer facing for extra fluid protection along with soft cellulose inner facing for cool comfort.

Item #	Description	Pkg
NON27373A	Chamber Style Mask w/Foam Strip and Ties, Green	50/bx, 300/cs

BFE > 97%
PFE > 98%



Surgical Masks

Our Surgical Masks Give you Peace of Mind in the O.R.

Standard Surgical Mask for your every day needs.

BFE > 97%
PFE > 98%

- NON27376 and NON27600 features a cellulose inner and outer facing keeps the wearer cool and comfortable.
- NON27377 and NON27420 feature a soft, fluid resistant polypropylene inner and outer facing assures comfort on the inside and protection on the outside.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27376	Standard Surgical Mask and Ties, Blue	50/bx, 300/cs
NON27600	Horizontal Tie Surgical Mask, Blue	50/bx, 300/cs
NON27377	Surgical Mask and Ties, Blue	50/bx, 300/cs
NON27420	Surgical Mask with Shield and Ties, Blue	25/bx, 200/cs



Medsoft™ Surgical Mask is our Softest.

BFE > 95%
PFE > 98%

- Our softest mask available.
- Fluid-resistant, outer facing provides protection while the soft inner facing gives the wearer the most comfort against the face.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27400	Soft Surgical Mask and Ties, White	50/bx, 300/cs



MEDLITE™ Surgical Mask feels so light you'll forget your wearing it.

BFE > 97%
PFE > 98%

- Our most comfortable and cool mask, ideal for long procedures.
- Fluid-resistant, outer facing provides protection.
- Soft cellulose inner facing keeps the wearer cool and comfortable.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27402	MEDLITE Surgical Mask and Ties, Blue	50/bx, 300/cs



Cone-Style Mask Enables Easy Donning

BFE > 99%
PFE > N/A

- Fluid-resistant, cup-shaped mask is held away from the nose and mouth for comfort and protection.
- Cone-style enables easy donning.
- Malleable nosepiece allows a custom fit.
- Designed with comfortable, rounded edges.

Item #	Description	Pkg
NON27381	Cone Style Mask with Headband, Blue	50/bx, 300/cs



Chamber-Style Mask Features a Large Breathing Chamber for Greater Comfort

BFE > 97%
PFE > 98%

- Large breathing chamber keeps mask away from face for greater comfort and easier breathing; also allows unmuffled speech.
- An ideal mask for beginning health care professionals as well as veterans.
- Cellulose facing for cool comfort.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27382	Chamber Style Mask and Ties, Green	50/bx, 300/cs



Chamber-Style Mask with Eyeshield Assures Maximum Eye and Facial Protection

- Large breathing chamber keeps mask away from face for greater comfort
- Non-glare, fog-free, wraparound eyeshield assures maximum eye and facial protection.
- Polypropylene outer facing for added protection against fluids, plus soft inner cellulose facing for comfort.
- Anti-fog strip of comfortable foam prevents eyeglasses from fogging.

BFE > 97%
PFE > 98%

Item #	Description	Pkg
NON27411	Chamber Style Mask w/Eyeshield and Ties, Green	25/bx, 100/cs

Hypoallergenic Mask is Ideal for Sensitive Skin

BFE > 97%
PFE > 98%

- Contains no dyes or colorants that could irritate sensitive skin.
- Comfortable cellulose outer and inner facing
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27385	Hypoallergenic Mask with Ties, White	50/bx, 300/cs



Laser Mask Protects Against Laser Plume Particles

BFE > 99%
PFE > 99%

- Special filter blocks laser plume particles as small as 0.1 micron.
- Fluid-resistant outer facing combined with comfortable, soft cellulose inner facing.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27374	Laser Mask with Ties, Gray	50/bx, 300/cs



Did you know?

Medline also offers a full line of disposable protective apparel that complements our fine face masks. Ask your Medline sales representative for details.

Isolation Mask

- Fluid-resistant outer facing provides protection.
- Soft cellulose inner facing increases comfort on the inside.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27120	Isolation Mask with Earloops, Yellow	50/bx, 300/cs

BFE > 97%
PFE > 98%



Basic Isolation Mask

- NON27110 has a tissue inner and outer facing making mask light and comfortable.
- NON27122 features a soft, fluid resistant polypropylene inner and outer facing assures comfort on the inside and protection on the outside.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27110	Basic Isolation Mask with Earloops, Yellow	50/bx, 300/cs
NON27122	Isolation Mask with Earloops, Yellow	50/bx, 300/cs

BFE > 99%
PFE > N/A

Standard Procedure Mask

- NON27375 features a Soft cellulose inner facing increasing comfort on the inside.
- NON27420EL features a soft, fluid resistant polypropylene inner and outer facing assures comfort on the inside and protection on the outside.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27375	Standard Procedure Mask with Earloops, Blue	50/bx, 300/cs
NON27420EL	Procedure Mask with Shield and Earloops, Blue	25/bx, 200/cs

BFE > 95%
PFE > 98%



Basic Procedure Mask

- NON27365 has a tissue inner and outer facing making mask light and comfortable.
- NON27378 features a Soft, fluid resistant polypropylene inner and outer facing assures comfort on the inside and protection on the outside.
- Full-width nosepiece guarantees proper fit.

Item #	Description	Pkg
NON27365	Basic Procedure Mask with Earloops, Blue	50/bx, 300/cs
NON27378	Basic Procedure Mask with Earloops, Blue	50/bx, 300/cs

BFE > 99%
PFE > N/A

N95 Approved Particulate Respirators Meet NIOSH and CDC Guidelines for TB Exposure

Adjustable, Flat Fold and Cone Styles

- N95 approved! Meets NIOSH and CDC guidelines for TB exposure.
- Positive Facial Lock® and malleable nosepiece allow a custom fit for any face, so one size mask fits all. (NON27501 only)
- Micro filter medium protects against micron-size particles.
- Double band holds mask firmly in place for maximum protection.
- All respirators must be selected, fitted and maintained in accordance with OSHA and other applicable regulations.

Item #	Description	Pkg
NON27501	Adjustable Particulate Respirator, Green/White Pattern	35/bx, 210/cs
NON24506	Cone Style Particulate Respirator, Large/Medium, White	20/bx, 240/cs
NON24507	Cone Style Particulate Respirator, Small, White	20/bx, 240/cs
NON24505	Flat Fold Respirator, Dual Headband, White	20/bx, 200/cs



Prohibit Series Competitive Cross Reference

Description	Medline Item #	Medline Packaging	Kimberly-Clark
Isolation mask	NON27110	50/bx, 300/cs	47117
	NON27122	50/bx, 300/cs	
	NON27120	50/bx, 300/cs	
Procedure mask	NON27365	50/bx, 300/cs	47080
	NON27375	50/bx, 300/cs	62356, 62357
	NON27378	50/bx, 300/cs	
	NON27420EL	25/bx, 100/cs	
	NON27605	50/bx, 300/cs	
Standard surgical	NON27376	50/bx, 300/cs	48201, 48210
	NON27600	50/bx, 300/cs	47500
	NON27402	50/bx, 300/cs	48100, 48105
	NON27400	50/bx, 300/cs	48390
	NON27377	50/bx, 300/cs	
	NON27420	25/bx, 100/cs	
	NON27610	50/bx, 300/cs	
Anti-Fog	NON27361A	50/bx, 300/cs	49215
	NON27371A	50/bx, 300/cs	
	NON27379A	50/bx, 300/cs	49235
	NON27378A	50/bx, 300/cs	49310, 49214
Chamber style	NON27382	50/bx, 300/cs	48220, 37525, 48227
	NON27373A	50/bx, 300/cs	49216
	NON27411	25/bx, 100/cs	48228
Hypoallergenic	NON27385	50/bx, 300/cs	59928
Cone	NON27381	50/bx, 300/cs	152
Laser	NON27374	50/bx, 300/cs	47625, 47650
Fluid protection	NON27710	25/bx, 100/cs	48247, 48237, 48248
	NON27710EL	25/bx, 100/cs	47137, 47147
	NON27712	50/bx, 300/cs	48207, 48208, 48297
	NON27712EL	50/bx, 300/cs	47107, 47297
	NON27410	25/bx, 100/cs	62114
	NON27410EL	25/bx, 100/cs	62116, 47298
	NON27412	50/bx, 300/cs	62113
	NON27412EL	50/bx, 300/cs	62115
	NON27405	25/bx, 100/cs	
	NON27408	50/bx, 300/cs	
	NON27408EL	50/bx, 300/cs	
N95 Particulate respirator & accessories	NON24506	20/bx, 240/cs	46727, 46767, 62126 46827, 46867, 62355 47950
	NON24507	20/bx, 240/cs	
	NON24505	50/bx, 300/cs	
	NON27501	35/bx, 210/cs	
	NON24510	each	
Face shields and protective eyewear	NONFS300	24/bx, 96/cs	41204
	NONFS400	24/bx, 96/cs	41205
	NONLENS	250/bx	SV250L
	NONFRAME	100/bx	SV100F
	NON24774	12/bx	
	NON24777A	20/bx	
	NON24777S	20/bx	
	NON24775	36/cs	
	NON24776	36/cs	



Prohibit Series Competitive Cross Reference

Medline Item #	3M	Cardinal	Precept
NON27110 NON27122 NON27120		AT70061,AT70021 AT7004	15100, 15120 15110
NON27365 NON27375 NON27378 NON27420EL NON27605	1826 1820	AT71061,AT71021 AT7511 AT7511-WE	15101 15111
NON27376 NON27600 NON27402 NON27400 NON27377 NON27420 NON27610	1818 1816	AT7505 AT71039 AT71035 AT71135,AT71335	15201, 65-3110 15250 15200, 65-3010 15600, 65-3100
NON27361A NON27371A NON27379A NON27378A	1831 1834 1832	AT72835 AT752005 AT71235,AT72235	15212, 65-3321, 65,3322 15211 15210, 65-3320
NON27382 NON27373A NON27411	1836, 1838 1838FSG	AT51035 AT54535 AT51135,AT54635	15220, 65-3124 65-3346
NON27385	1833	AT73035	15215
NON27381	1800nl	AT7509	65-3333
NON27374		AT75435	65-3310, 65-3311
NON27710 NON27710EL NON27712 NON27712EL NON27410 NON27410EL NON27412 NON27412EL NON27405 NON27408 NON27408EL	1835FS 1840FS 1835 1840 1818FSG	AT74635 AT74631 AT74535 AT74531	15330 15310, 15510 15320 15300, 15301, 15500
NON24506 NON24507 NON24505 NON27501 NON24510	1860 1860S 1870 FT10	N95-ML N95-S	65-3395 65-3395s 65-65000
NONFS300 NONFS400 NONLENSE NONFRAME NON24774 NON24777A NON24777S NON24775 NON24776	1710 1611 1621	F1SHIELD50 H1SHIELD50 ESHIELD50 BFRAME30, WFRAME30	15995 15990



Face Mask Accessories To Complement Our Face Masks

Qualitative Fit Test Kit

- To be used in conjunction with our N95 particulate respirator masks (NON27501, NON24505, NON24506 and NON24507).
- Fit testing is required prior to the initial use of a respirator.
- Saccharin based fit test kit.
- Test is easy to administer and instructions are found within the kit.
- One kit tests approximately 100 staff members.



Item #	Description	Pkg
NON24510	Qualitative Fit Test Kit	each
NON24510A	Fit Test Sensitivity Solution	6 btl/cs
NON24510B	Fit Test Actual Solution	6 btl/cs
NON24510C	Fit Test Bitter Solution	6 btl/cs
NON24510D	Fit Test Hood,	10/cs
NON24510E	Replacement Nebulizers	3/cs

Disposable Face Shields

- Provides the user with maximum eye and facial protection.
- Can be used in conjunction with our Prohibit Face Masks.
- Shields are made of durable 6 ml thick material and are anti-fog to ensure optical clarity.
- Both shields have a latex-free band and comfortable, contoured foam material that rests against the forehead.



Item #	Description	Pkg
NONFS300	Disposable Full Shield, 7.5" long	24/bx, 96/cs
NONFS400	Disposable 3/4 Shield, 6.5" long	24/bx, 96/cs

Fluid Protective Goggles

- Soft, vinyl-frame goggles fit easily over prescription eyewear.
- Four air vents provide additional comfort.
- Available with specially-treated, anti-fog lenses for clear, unobstructed viewing.
- Goggles meet ANSI standards and help your facility meet OSHA requirements for hazardous liquids.



Item #	Description	Pkg
NON24775	Goggles, fog free	36/cs
NON24776	Goggles, standard	36/cs

Safety Glasses

- Safety glasses provide economical and practical protection.
- Built-in brow bar and side shields, plus wide view, optical lenses.
- Comfortable fit.
- Large size glasses can be worn over prescription eyewear.



Item #	Description	Pkg
NON24777A	Safety glasses, visitor, clear, large	20/bx
NON24777S	Safety glasses, visitor, clear, small	20/bx
NON24774	Safety glasses with black strip	12/bx
NONLENS	Disposable lens for use w/NONFRAME, clear	250/bx
NONFRAME	Reusable frame for use NONLENS, assorted colors	100/bx

Face Mask Dispensers

Item #	Description	Pkg
NONFP022	Dispenser for surgical face masks	6/cs
NONFP021	Dispenser for procedure face masks	6/cs



1-800-MEDLINE (633-5463) www.medline.com

©2008 Medline Industries, Inc. Prohibit, Prohibit X-Tra and Medline are registered trademarks Medsoft and Medlite are trademarks of Medline Industries, Inc. One Medline Place, Mundelein, IL 60060. Positive Facial Lock is a registered trademark of Alpha ProTech, Inc. We reserve the right to correct any errors that may occur within this brochure.

Medline Canada 1-800-396-6996

www.medline.ca canada@medline.com

Not all products are currently available for sale in Canada. Please contact your Account Representative for additional information.

MKT211431/LIT168/1.5M/SG6