


Be Independent.™

BEDCADDIE

BECAUSE SITTING UP IN BED CAN BE EASY




PRODUCT #2085



WATCH A DEMO





Ergonomic 5" Rubber Handle

Strap Adjusts Between 48" -84"

Be Independent

FEATURES & BENEFITS

- **3 Hand Grips** – Ladder-like design makes sitting up easy.
- **Versatile** – Attaches to any bed frame.
- **Sturdy Nylon Construction** – Provides long lasting safety.
- **Adjustable**- Customize length to fit your preference.
- **Ergonomic** – Cushioned, no-slip rubber contour grip.
- **Hassle Free Assembly** – Installs and removes in seconds, no tools required.

SPECIFICATIONS

Strap Length: Adjustable 48" – 84"

Length Between Handles: 12"

Handle Grip: 5"

Weight Capacity: 300 lbs

Product Weight: .8 lbs

Package Dimensions: 10" L x 7" W x 1" H; .8 lbs

Case of 10 Dimensions: 8" L x 8" W x 11" H; 8 lbs