Live the Longer Life You Love®

NEW!

5800 STABLE RAIL The Bed Rail With Legs







EASILY SIT UP IN BED

EASILY STAND UP

ORGANIZER KEEPS ITEMS CLOSE BY

Features & Benefits

- Extendable Legs Provides added stability when standing up.
- Ergonomic Cushion Handle Allows for easy transfer in and out of bed.
- 4-Pocket Organizer Provides storage space for handy items.
- Anti-slip grips Secures rail in-between mattress and bed frame.
- Height Adjustment Accomodates any home or hospital bed.
- Easy Installation Installs in seconds with no tools required.

Specifications

Height adjustment from floor to top of handle: 21"- 44" Height adjustment from mattress base to top of handle: 17.5" - 22.5" Rail width: 19"

Height range from floor to horizontal support structure under mattress: 13.5"- 21.5"

Depth of horizontal support structure under mattress: 25"

Weight of Product: 8 lbs

Weight capacity: 300 lbs

Made of: Steel frame with durable powder coat finish Case of 4 Dimensions: 28" L x 20.25 W x 5.75" H; 32 lbs

How Do We Compare	Stable Rail	Other Bed Rails
Tool-Free Assembly		×
Extendable Legs for Added Stability	\checkmark	X
4-Pocket Organizer Pouch Included		X
Provides Assistance In and Out of Bed	S	S
Durable Powder Coat Finish	S	

Easy Assembly (No Additional Tools Required)





Adjust base tubes to accommodate left or right side of bed

Attach base tubes to main handle



Insert base tubes between mattress and box spring



Adjust height of handle and legs as desired