




Surgical Tapes, Adhesive Skin Closures and Self-Adherent Wraps



The most complete line of surgical tapes in the industry



Product Selection Matrix

	Page Number	Latex Free	Non-Occlusive	Waterproof	Sensitive Skin	Single Use	Easy Hand Tear	Adhesion Level
Cloth Tape	2	X	X				X	Mod. ▲
Paper Tape	3	X	X			X	X	Low
Transparent Tape	4	X	X				X	Mod. ▲
Waterproof Tape	5	X			X			High
Elastic Adhesive Bandage	5							High
Elastic Foam Tape	6	X	X					Mod. ▲
Cotton Adhesive	6	X	X				X	High
Ortho-Porous Sports Tape	7			X				High
Megazinc Pink Tape	7	X			X	X		High
Medfix/Medfix EZ Tape	8	X	X			X		High
Co-Flex Med™	9			X				N/A*
Co-Flex NL™	9	X	X			X		N/A*
Medi-Strips Wound Closure	10	X	X					High
Elastic Net Dressing	11	X	X			X		N/A*

* N/A-Product does not stick to skin

▲ Moderate

MEDLINE

Medline Industries, Inc. is the largest privately held national manufacturer and distributor of medical supplies and equipment in the U.S.A. We've been servicing the health care industry for more than 30 years and our growth in wound care has been dramatic.

ORDER WITH CONFIDENCE

With over 700 professional Sales Representatives throughout the U.S., you have an on-site partner to provide solutions. Medline has 23 distribution centers around North America that provide an unmatched 99% service level for Traditional Wound Care products, including all surgical tapes, Medi-Strip™ wound closure strips, and Co-Flex NL™ self-adherent wraps.



Cloth "Silk" Surgical Tape

New Formulation

Cloth "Silk" Surgical Tape



* Compare to 3M™ "Durapore™"

Cloth "Silk" Surgical Tape

Item No.	Description	Pkg
NON260112	½" x 10 yds	24 rl/bx, 12 bx/cs, 288 rl/cs
NON260101	1" x 10 yds	12 rl/bx, 12 bx/cs, 144 rl/cs
NON260102	2" x 10 yds	6 rl/bx, 12 bx/cs, 72 rl/cs
NON260103	3" x 10 yds	4 rl/bx, 12 bx/cs, 48 rl/cs

Single Use Rolls

NON260101S	1" x 1.5 yds	100 rl/bx, 5 bx/cs, 500 rl/cs
NON260102S	2" x 1.5 yds	50 rl/bx, 5 bx/cs, 250 rl/cs



Applications for use:

- Ideal for bulky dressings that require high conformability.
- Post-operative surgical dressings.

- Strong woven rayon/acetate backing provides high tensile strength for confident securement of dressings and devices.
- Easy and clean bi-lateral hand tear quickly customizes sizes.
- Pinked edges reduce fraying when tearing.
- Available in single use rolls that limits cross contamination and waste.
- Hypoallergenic.

Discover Our Reformulated Cloth Tape

Water Repellent	Try it yourself! Water beads on the Medline Cloth Tape. Dressings stay on longer even when exposed to moisture.
Moisture Vapor Transmission	140% greater rate of moisture readily pass through the tape to improve adhesion durability and help prevent skin maceration.
Adhesive Weight	23% reduction in adhesive weight provides a stronger adhesion level with less adhesive mass.
Unwind Tension	75% lower rate of unwind throughout the entire roll makes this easier to use, and helps contribute to a longer shelf life.

Comparisons quoted reflect the improvements made to the Cloth Tape Formulation vs. the previous formula.



www.medline.com
1-800-MEDLINE

Paper Surgical Tape

Now Available in Tan!



latexfree
PRODUCT

* Compare to 3M™ "Micropore™"

Paper Tape: (white)

Item No.	Description	Pkg
NON260012	½" x 10 yds	24 rl/bx, 12 bx/cs, 288 rl/cs
NON260001	1" x 10 yds	12 rl/bx, 12 bx/cs, 144 rl/cs
NON260002	2" x 10 yds	6 rl/bx, 12 bx/cs, 72 rl/cs
NON260003	3" x 10 yds	4 rl/bx, 12 bx/cs, 48 rl/cs

Single Use Rolls:

NON260001S	1" x 1.5 yds	100 rl/bx, 5 bx/cs, 500 rl/cs
NON260002S	2" x 1.5 yds	50 rl/bx, 5 bx/cs, 250 rl/cs

Dispenser Pack:

NON260012D	½" x 10 yds	24 rl/bx, 12 bx/cs, 288 rl/cs
NON260001D	1" x 10 yds	12 rl/bx, 12 bx/cs, 144 rl/cs
NON260002D	2" x 10 yds	6 rl/bx, 12 bx/cs, 72 rl/cs
NON260003D	3" x 10 yds	4 rl/bx, 12 bx/cs, 48 rl/cs

Paper Tape: (tan) **NEW!**

Item No.	Description	Pkg
NON260012T	½" x 10 yds	24 rl/bx, 12 bx/cs, 288 rl/cs
NON260001T	1" x 10 yds	12 rl/bx, 12 bx/cs, 144 rl/cs
NON260002T	2" x 10 yds	6 rl/bx, 12 bx/cs, 72 rl/cs



Applications for use:

- Secures dressings, tubes, lines and I.V. catheters.
- Ideal for light dressings
- Holds well on damp skin
- Tan color helps camouflage dressings on visible areas such as face, neck, arms.
- Gentle adhesion is ideal for general use, pediatric patients, geriatric patients and where frequent re-taping is required.
- Hypoallergenic, breathable tape allows for air and moisture vapor transmission, which helps prevent skin maceration.
- Highly conformable due to exceptional flexibility of backing material.
- Consistently low unwind tension.
- Holds well on damp skin.
- Available in single use rolls that limits cross contamination and waste.



Transparent Surgical Tape



latexfree
PRODUCT

* Compare to 3M™ "Transpore™"



Applications for use:

- Secures dressings, tubes, lines and I.V. catheters.
- Easy to customize widths for NICU patients
- Anchors bulky dressings

- Enhanced perforated pattern provides easy hand tear and rapid transmission of moisture and air.
- Low unwind tension for ease of removal from the roll.
- Easy to write on and easy to clean.
- Durable adhesion that conforms to skin and flexes with the body for comfort.
- A low profile backing material is more comfortable and is less prone to snagging.
- Hypoallergenic.
- Available in single use rolls that limits cross contamination and waste.

Transparent Tape:

Item No.	Description	Pkg
NON260212	½" x 10 yds	24 rl/bx, 12 bx/cs, 288 rl/cs
NON260201	1" x 10 yds	12 rl/bx, 12 bx/cs, 144 rl/cs
NON260202	2" x 10 yds	6 rl/bx, 12 bx/cs, 72 rl/cs
NON260203	3" x 10 yds	4 rl/bx, 12 bx/cs, 48 rl/cs

Single Use Rolls:

NON260201S	1" x 1.5 yds	100 rl/bx, 5 bx/cs, 500 rl/cs
NON260202S	2" x 1.5 yds	50 rl/bx, 5 bx/cs, 250 rl/cs

Excellent Visibility

Exceptionally clear backing material allows for easier viewing of the tubing, monitoring of skin condition and healing progress without removing the tape.



Waterproof Tape



latexfree
PRODUCT

Waterproof Tape:

Item No.	Description	Pkg
NON260512	½" x 10 yds	24 rolls/bx, 12 bx/cs, 288 rl/cs
NON260501	1" x 10 yds	12 rolls/bx, 12 bx/cs, 144 rl/cs
NON260502	2" x 10 yds	6 rolls/bx, 12 bx/cs, 72 rl/cs
NON260503	3" x 10 yds	4 rolls/bx, 12 bx/cs, 48 rl/cs



Applications for use:

- Ideal for general taping needs and outpatient use.
- Bandage wraps, windowed casts, picture frame ostomy flanges, and protect dressings.

- Waterproof, oil resistant and easy to write on.
- Cotton cloth backing with polyethylene coating provides high tensile strength and will not permit moisture to penetrate to skin or dressing
- Strong latex-free adhesive is durable and holds securely.

Elastic Adhesive Bandage



Elastic Adhesive Bandage:

Item No.	Description	Pkg
NON260401	1" x 5 yds	12 rl/bx, 12 bx/cs, 144 rl/cs
NON260402	2" x 5 yds	6 rl/bx, 12 bx/cs, 72 rl/cs
NON260403	3" x 5 yds	4 rl/bx, 12 bx/cs, 48 rl/cs
NON260404	4" x 5 yds	6 rl/bx, 6 bx/cs, 36 rl/cs



Applications for use:

- Controls and reduces post-injury edema by providing firm support and compression.
- Ideal for dynamic/pressure dressings requiring elasticity, support strapping, fractures and muscle injuries.

- 100% cotton elastic cloth tape backing provides high tensile strength.
- Strong rubber-based adhesive permits body movement without loosening of tape.
- Porous construction maintains skin integrity and permits moisture vapor transmission.
- Contains Latex



Elastic Foam Tape



* Compare to 3M™ "Microfoam™"



Foam Tape:

Item No.	Description*	Pkg
NON260601	1" x 5.5 yds	12 rl/bx, 6 bx/cs, 72 rl/cs
NON260602	2" x 5.5 yds	6 rl/bx, 6 bx/cs, 36 rl/cs
NON260603	3" x 5.5 yds	4 rl/bx, 6 bx/cs, 24 rl/cs
NON260604	4" x 5.5 yds	3 rl/bx, 6 bx/cs, 18 rl/cs

* Stretch Length



Applications for use:

- Affix compression dressings to aid hemostasis or control edema (e.g. breast biopsies, orthopedics and abdominal surgical dressings)
- Ideal for securing dressings on difficult to tape areas such as knees, elbows, shoulders, hips and fingers.
- Highly conformable, multi-directional stretch promotes patient comfort and accommodates swelling.
- Moderate adhesion secures dressings over challenging areas and can be used to provide compression applications.
- Smooth, cushioned elastic foam backing protects wound site and dressings.
- Water resistant for greater durability.
- Latex free adhesive and backing leave minimal skin residue.

Cotton Adhesive Tape



* Compare to 3M™ "Cloth Adhesive Tape"



Cotton Adhesive Tape:

Item No.	Description	Pkg
NON260712	½" x 10 yds	24 rl/bx, 12 bx/cs, 288 rl/cs
NON260701	1" x 10 yds	12 rl/bx, 12 bx/cs, 144 rl/cs
NON260702	2" x 10 yds	6 rl/bx, 12 bx/cs, 72 rl/cs
NON260703	3" x 10 yds	4 rl/bx, 12 bx/cs, 48 rl/cs



Applications for use:

- Ideal for positioning body parts, finger splints stabilization, immobilizing digits, or securing tubing.
- Strong, 100% cotton woven backing is durable with high tensile strength.
- High long-term adhesion for confident securement of dressings and tubing.
- Convenient hand tear to customize length and width.
- Latex free adhesive with minimal residue.



Ortho-Porous Sports Tape



Ortho-Porous Sports Tape:

Item No.	Description	Pkg
NON260312	½" x 10 yds	24 rl/bx, 12 bx/cs, 288 rl/cs
NON260301	1" x 10 yds	12 rl/bx, 12 bx/cs, 144 rl/cs
NON260315	1½" x 10 yds	8 rl/bx, 12 bx/cs, 96 rl/cs
NON260302	2" x 10 yds	6 rl/bx, 12 bx/cs, 72 rl/cs
NON260303	3" x 10 yds	4 rl/bx, 12 bx/cs, 48 rl/cs



Applications for use:

- Athletics, immobilization (strains, sprains), rehabilitation, support strapping, patient positioning and restraints.

- Strong woven cotton backing provides high tensile strength and support to highly active patients.
- Powerful rubber-based adhesive is designed for long lasting adhesion with minimal residue.
- Highly porous construction allows the skin to breathe and permits moisture vapor transmission.
- Consistent low unwind tension throughout the roll for ease of use.
- Contains latex.

Megazinc Pink Adhesive Tape



* Compare to 3M™ "Hy-Tape®"

latexfree
PRODUCT

Megazinc Pink Tape

Item No.	Description	Pkg
OMAM55CS	½" x 5 yds	1 rl/bx, 36 bx/cs, 36 rl/cs
OMAM111CS	1" x 5 yds	1 rl/bx, 36 bx/cs, 36 rl/cs
OMAM222CS	2" x 5 yds	1 rl/bx, 36 bx/cs, 36 rl/cs
OMAM55	½" X 5 yds	1 rl/bx
OMAM111	1" X 5 yds	1 rl/bx
OMAM222	2" X 5 yds	1 rl/bx



Applications for use:

- Secures tubes (endotracheal, gastrostomy, nasogastric), epidurals, ostomy, patch testing, IV's and, electrodes.

- Zinc oxide-based adhesive is gentle to sensitive skin.
- Remains intact when wet or washed.
- Easy-tear edges.
- Conformability and flexibility allow for ease of application.
- Easy, painless removal.
- Contraindicated for radiated skin.



Medfix and Medfix EZ Dressing Retention Sheet



* Compare to 3M "Medipore H™ Soft Cloth Surgical Tape"

Medfix Dressing Retention Sheet:

Item No.	Description	Pkg
MSC4002	2" x 11 yds	1 rl/bx
MSC4004	4" x 11 yds	1 rl/bx
MSC4006	6" x 11 yds	1 rl/bx

Medfix EZ Dressing Retention Roll:

Item No.	Description	Pkg
MSC4102	2" x 11 yds	12 rl/bx
MSC4104	4" x 11 yds	12 rl/bx
MSC4106	6" x 11 yds	12 rl/bx
MSC4122	2" x 2 yds	12 rl/bx
MSC4124	4" x 2 yds	12 rl/bx
MSC4126	6" x 2 yds	12 rl/bx



Applications for use:

- Ideal for securing primary dressings, gauze, foams, gel wafers, gastrostomy and other feeding tubes.
- Low sensitivity adhesive for gentle but secure adhesion on fragile skin.
- Medfix Dressing Retention Sheet has printed s-curve release liner makes it easy to measure the desired length, liner removes cleanly and easily.
- Medfix EZ is linerless for faster application and perforated every two inches for easy tearing of desired length.
- Water resistant and acts as an excellent secondary dressing for gauze, foams, gels, etc.
- Hypoallergenic and latex free.



Co-Flex™ Self Adherent Wrap

Easy Hand Tear!



* Compare to 3M™ "Coban™"

Co-Flex NL™ :

Item No.	Description	Pkg
<i>Sterile</i> latexfree PRODUCT		
DYNJ05123	2" x 5 yds stretched, beige	20 rolls/cs
DYNJ05124	3" x 5 yds stretched, beige	20 rolls/cs
DYNJ05129	4" x 5 yds stretched, beige	20 rolls/cs
DYNJ05130	6" x 5 yds stretched, beige	20 rolls/cs

Non-Sterile **latexfree** PRODUCT

* MDS088001	1" x 5 yds stretched, beige	30 rolls/cs
MDS088002	2" x 5 yds stretched, beige	36 rolls/cs
MDS088003	3" x 5 yds stretched, beige	24 rolls/cs
MDS088004	4" x 5 yds stretched, beige	18 rolls/cs
MDS088006	6" x 5 yds stretched, beige	12 rolls/cs

Co-Flex Med™ :

Item No.	Description	Pkg
<i>Sterile, Latex</i>		
DYNJ086002	2" x 5 yds stretched, beige	20 rolls/cs
DYNJ086003	3" x 5 yds stretched, beige	20 rolls/cs
DYNJ086004	4" x 5 yds stretched, beige	20 rolls/cs
DYNJ086006	6" x 5 yds stretched, beige	20 rolls/cs

Non-Sterile, Latex

* MDS086001	1" x 5 yds stretched, beige	30 rolls/cs
MDS086002	2" x 5 yds stretched, beige	36 rolls/cs
MDS086003	3" x 5 yds stretched, beige	24 rolls/cs
MDS086004	4" x 5 yds stretched, beige	18 rolls/cs
MDS086006	6" x 5 yds stretched, beige	12 rolls/cs

* Two 1" rolls per package

Applications for use:



- Secure dressings, splints/armboards, I.V.'s and other tubing.
- Helps control edema and post I.V. hematoma
- Can serve as a compression bandage for hemostasis.

- Self-adherent, cohesive wrap securely sticks to itself and not to skin or other materials and will not slip or loosen over time.
- Easily tears by hand for custom lengths. No need for scissors.
- Porous construction allows ventilation of the wrapped area.
- Retains dressings without tension around digits or limbs.
- Available sterile/non-sterile and with/without latex.
- Each bandage is individually wrapped.

NEW

Kid's Pack Contains Rolls of:

Neon Pink/Purple Hearts
Neon Green/ Dark Blue Dinosaurs
Orange/Black Smiley Face
Light Blue/Dark Blue Trucks and Cars
Yellow/Black Bumblebees
Lilac/Purple Flowers


Colorpack Contains Rolls of:

Pink, Blue, Purple, Light Blue, Neon Green, and Red



Item No.	Description	Pkg
<i>Non-Sterile, Latex Free</i>		
* MDS088001CP	1" x 5 yds stretched, Colorpack	30 rolls/cs
MDS088002CP	2" x 5 yds stretched, Colorpack	36 rolls/cs
MDS088003CP	3" x 5 yds stretched, Colorpack	24 rolls/cs
MDS088004CP	4" x 5 yds stretched, Colorpack	18 rolls/cs
MDS088002KP	2" x 5 yds stretched, Kid's pack	36 rolls/cs
MDS088003KP	3" x 5 yds stretched, Kid's pack	24 rolls/cs

Medi-Strips® Wound Closure Strips



* Compare to 3M™ "Steri-Strips™"

- Constructed of exceptionally porous spun-bond fiber backing to help reduce skin maceration.
- Bi-directional reinforcing nylon filaments provide excellent tensile strength for confident placement and securement especially on large skin surface closures.
- Hypoallergenic and latex-free.
- No skin tackifiers or adhesive enhancers required.
- Perforated card allows for easy delivery of curl resistant strips on clean, dry skin.
- Contraindicated on high tension wounds.

Medi-Strips® Wound Closure Strips:

Item No.	Description	Pkg
<i>Sterile</i>		
NON250318	1/8" x 3"	5/pkg, 50 pk/bx, 4 bx/cs
NON250114	1/4" x 1 1/2"	6/pkg, 50 pk/bx, 4 bx/cs
NON250412	1/2" x 4"	6/pkg, 50 pk/bx, 4 bx/cs
NON250314	1/4" x 3"	3/pkg, 50 pk/bx, 4 bx/cs
NON250414	1/4" x 4"	10/pkg, 50 pk/bx, 4 bx/cs
NON250501	1" x 5"	4/pkg, 25 pk/bx, 4 bx/cs



Applications for use:

- Primary wound/skin closures including lacerations and surgical incisions.
- Procedures where larger incisions need surface skin closure such as cardio-thoracic and orthopedic.
- When early staple or suture removal is desired.

New! Elastic Medi-Strips®



Applications for use:

- Ideal for applications when flexibility and conformability are desired such as on finger, wrist, elbow, shoulder or knee or for plastic surgery.
- Permits stretching with body movement or edema.
- Flexible strips are made of a tan, non woven backing that is soft, highly comfortable and promotes patient comfort.
- Hypoallergenic medical grade adhesive provides durable adhesion while being gentle to the skin.
- Perforated card allows for easy delivery of curl resistant strips on clean, dry skin. Latex-free.

Elastic Medi-Strips®

Item No.	Description	Pkg
NON251314	1/4" x 3"	3/pkg, 50 pk/bx, 4 bx/cs
NON251412	1/4" x 4"	6/pkg, 50 pk/bx, 4 bx/cs

Elastic Net Dressing

The Tape Alternative

* Compare to “Spandage™”



Applications for use:

- Ideal for pediatric, geriatric, fragile and at risk skin that might be damaged by adhesives.
- Areas where tape application or removal would be difficult or cause pain such as cranial procedures, hair replacement, dermabrasion, facial surgery, amputations, and in the sleep lab.
- Secure burn or post-op dressings, IV tubes, stress test or holter monitor fitting to reduce artifacts, secure incontinence pads, decubitus ulcer dressings, hot/cold packs, etc.



Elastic Net Dressing:

Item No.	Size	Application	Max. Stretch Width
NONNET01	1	Small finger, toe	6"
NONNET02	2	Adult finger, toe	7½"
NONNET03	3	Hand, multiple fingers lower arm, wrist, ankle	9"
NONNET04	4	Large hand, foot, upper arm, elbow	11½"
NONNET05	5	Average arm, small leg knee, amputations	12½"
NONNET06	6	Medium/large leg amputations intravenous/armboard	24"
NONNET07	7	Average head, cranium, ears, face, neck	28"
NONNET08	8	Large leg, cranium, ears, face, neck	31"
NONNET09	9	Adult chest, abdomen, breast, shoulder	35"
NONNET10	10	Large chest, abdomen, breast, shoulder	36"
NONNET11	11	XX-large chest, abdomen, breast, shoulder	40"

- Rayon/polyester knitted elastic net tube holds wet and dry dressings firmly in place without damaging skin.
- Soft open weave is low linting, does not constrict, is painless upon removal and allows for maximum airflow.
- Superior elasticity and stretch memory allows frequent inspection and redressing without adhesive.
- Accommodates virtually every part of the body, from finger, head to XX-large torso.
- Easily cut to customize size and applications.
- All sizes are 25 yards stretched length, 1 roll per box.

For additional information and detailed instructions on how to bandage challenging areas, please visit:
www.medline.com/products/elasticnet

Adhesive Bandages

A complete line of bandages for every need.




Sheer-Gard® Adhesive Bandages

Sheer-Gard adhesive bandages are constructed with PVC-free durable polyethylene and perforated for breathability. These bandages feature even adhesion and a non-adherent island pad. Latex-Free.

Plastic, Sterile

Item No.	Description	Pkg
NON25500	¾" x 3"	100/bx, 12 bx/cs, 1200/cs
NON25600	1" x 3"	100/bx, 12 bx/cs, 1200/cs
NON25501	Spot ⅞"	100/bx, 12 bx/cs, 1200/cs
NON25504	X-Lg 2" x 4"	50/bx, 12 bx/cs, 600/cs
NON25509	Jr. ⅝" x 1½"	100/bx, 36 bx/cs, 3600/cs



Comfort Cloth® Adhesive Bandages

Comfort Cloth adhesive bandages are woven to allow them to stretch and conform to the wound. Material naturally breathes. Features even adhesion and a non-adherent island pad. Latex-Free.

Woven Cloth, Sterile

Item No.	Description	Pkg
NON25650	¾" x 3"	100/bx, 12 bx/cs, 1200/cs
NON25660	1" x 3"	100/bx, 12 bx/cs, 1200/cs
NON25510	Knuckle, 3" x 1½"	100/bx, 12 bx/cs, 1200/cs
NON25513	Finger tip, 2½" x 1½"	100/bx, 12 bx/cs, 1200/cs



New! Foam Adhesive Bandage

Designed for the "OR"

Our top of the line adhesive bandage is a favorite of OR nurses everywhere. The wrapper tears cleanly and easily for trouble free product delivery to the sterile field. It has a non-adherent island pad on a dense foam strip that provides greater flexibility and comfort than a traditional adhesive bandage. Excellent for use on joints. Latex-Free.

Foam, Sterile

Item No.	Description	Pkg
NON25950	¾" x 3"	100/bx, 12bx/cs, 1200/cs



Adhesive Bandages

A complete line of bandages for every need.

New! Introducing Medtoons™

Cartoon Adhesive Bandages

Children will adore the smiling bumble bees, soccer playing dogs and snorkeling fish featured on the new Medtoons Cartoon adhesive bandage. Each box includes 50 strips (25 strips each of yellow, and blue) constructed from PVC-free vinyl perforated for increased moisture vapor transmission. An absorbent, non-adherent island pad seals around the wound site to keep out dirt and helps protect granulating tissue. The pressure sensitive, acrylic adhesive is skin friendly and provides long lasting protection. Presented in a convenient dispenser box. Latex-free.



Plastic, Sterile

Item No.	Description	Pkg
NON256130	3/4" x 3"	50/box, 24 bx/cs, 1200/cs

New! Neon Adhesive Bandage

This bright assortment of five neon colors includes 10 plastic strips each of pink, blue, green, purple and yellow in each box. The non-adherent island pad protects granulating tissue, and the hypoallergenic adhesive provides long lasting protection. Latex-free.



Plastic, Sterile

Item No.	Description	Pkg
NON256131	3/4" x 3"	50/box, 24 bx/cs, 1200/cs

New! Pressure Adhesive Bandage

Once in contact with fluid, the absorbent, cellulose island pad expands to over 7 times its normal height helping to provide gentle pressure to stop bleeding, saves nursing time and promotes comfort for dialysis, chemotherapy and cardio cath patients. Latex-free.



Sterile

Item No.	Description	Pkg
NON85100	Large, 1" x 2 3/4"	100/bx, 20 bx/cs, 2000/cs
NON85200	X-Large, 1 1/4" x 2 3/4"	100/bx, 20 bx/cs, 2000/cs

Medline Cross Reference Chart

Cloth Tape	Medline	3M "Durapore"	J&J "Dermicel"	Kendall "Curasilk"
Cloth Tape ½" x 10 yds 24/bx, 12 bx/cs (288/cs)	NON260112	15380 (240/cs)	5143	7137
Cloth Tape 1" x 10 yds 12/bx, 12 bx/cs (144/cs)	NON260101	15381 (120/cs)	5144	7138
Cloth Tape 2" x 10 yds 6/bx, 12 bx/cs (72/cs)	NON260102	15382 (60/cs)	5145	7139
Cloth Tape 3" x 10 yds 4/bx, 12 bx/cs (48/cs)	NON260103	15383 (40/cs)	5146	7140
Cloth Tape Single Use 1" x 1.5 yds 100/bx, 5 bx/cs (500/cs)	NON260101S	1538S-1		
Cloth Tape Single Use 2" x 1.5 yds 100/bx, 5 bx/cs (500/cs)	NON260102S	1538S-2		
Self-Adherent Wrap Non-Sterile with Latex	Medline "Co-Flex"	3M "Coban"		Kendall "Flex-Wrap"
Co-Flex Med, Self-Adherent Wrap, 1" x 5 yds (30/pk)	MDS086001	1581		4581
Co-Flex Med, Self-Adherent Wrap, 2" x 5 yds (36/pk)	MDS086002	1582		4582
Co-Flex Med, Self-Adherent Wrap, 3" x 5 yds (24/pk)	MDS086003	1583		4583
Co-Flex Med, Self-Adherent Wrap, 4" x 5 yds (18/pk)	MDS086004	1584		4584
Co-Flex Med, Self-Adherent Wrap, 6" x 5 yds (12/pk)	MDS086006	1586		4586
Self Adherent Wrap Non-Sterile Latex Free	Medline "Co-Flex"	3M "Coban"		
Co-Flex NL Hand Tear, Latex Free, 1" x 5 yds (30/pk)	MDS088001	2081		
Co-Flex NL Hand Tear, Latex Free, 2" x 5 yds (36/pk)	MDS088002	2082		
Co-Flex NL Hand Tear, Latex Free, 3" x 5 yds (24/pk)	MDS088003	2083		
Co-Flex NL Hand Tear, Latex Free, 4" x 5 yds (18/pk)	MDS088004	2084		
Co-Flex NL Hand Tear, Latex Free, 6" x 5 yds (12/pk)	MDS088006	2086		
Self Adherent Wrap Sterile with Latex	Medline "Co-Flex"	3M "Coban"		Kendall "Flex-Wrap"
Co-Flex Med, Self Adherent Wrap, 2" x 5yds	DYNJ086002 (20/pk)			
Co-Flex Med, Self Adherent Wrap, 3" x 5yds	DYNJ086003 (20/pk)	1583S (24/pk)		4583S (24/pk)
Co-Flex Med, Self Adherent Wrap, 4" x 5yds	DYNJ086004 (20/pk)	1584S (18/pk)		4584S (18/pk)
Co-Flex Med, Self Adherent Wrap, 6" x 5yds	DYNJ086006 (20/pk)	1586S (12/pk)		4586S (12/pk)
Self Adherent Wrap Sterile Latex Free	Medline "Co-Flex"	3M "Coban"		
Co-Flex NL Hand Tear, Latex Free, 2" x 5 yds (20/pk)	DYNJ05123			
Co-Flex NL Hand Tear, Latex Free, 3" x 5 yds (20/pk)	DYNJ05124	2083S (24/pk)		
Co-Flex NL Hand Tear, Latex Free, 4" x 5 yds (20/pk)	DYNJ05129	2084S (18/pk)		
Co-Flex NL Hand Tear, Latex Free, 6" x 5 yds (20/pk)	DYNJ05130	2086S (12/pk)		
Cotton Adhesive Tape Latex Free	Medline	3M "Cloth Adhesive"		Kendall
Cotton Adhesive Tape ½" x 10 yds 24/bx, 12 bx/cs (288/cs)	NON260712	2950-0 (240/cs)		2304 (with Latex)
Cotton Adhesive Tape 1" x 10 yds 12/bx, 12 bx/cs (144/cs)	NON260701	2950-1 (120/cs)		2531 (with Latex)
Cotton Adhesive Tape 2" x 10 yds 6/bx, 12 bx/cs (72/cs)	NON260702	2950-2 (60/cs)		6613 (with Latex)
Cotton Adhesive Tape 3" x 10 yds 4/bx, 12 bx/cs (48/cs)	NON260703	2950-3 (40/cs)		7046 (with Latex)
Elastic Adhesive Bandage	Medline	BSN-JOBST "Elastoplast"	J&J "Elastikon"	Kendall "Expandover"
Elastic Adhesive Bandage 1" x 5 yds 12/bx, 12bx/cs (144/cs)	NON260401	2593 (36/pk)	5172	8884662001
Elastic Adhesive Bandage 2" x 5 yds 6/bx, 12bx/cs (72/cs)	NON260402	2594 (36/pk)	5174	8884662019
Elastic Adhesive Bandage 3" x 5 yds 4/bx, 12bx/cs (48/cs)	NON260403	2595 (36/pk)	5175	8884662027
Elastic Adhesive Bandage 4" x 5 yds 6/bx, 6bx/cs (36/cs)	NON260404	2596 (36/pk)	5177	8884662035



Medline Cross Reference Chart (continued)

Elastic Foam Tape	Medline	3M "Microfoam"		
Elastic Foam Tape 1" x 5.5 yds 12/bx, 6 bx/cs, (72/cs)	NON260601	1528-1		
Elastic Foam Tape 2" x 5.5 yds 6/bx, 6 bx/cs, (36/cs)	NON260602	1528-2		
Elastic Foam Tape 3" x 5.5 yds 4/bx, 6 bx/cs, (24/cs)	NON260603	1528-3		
Elastic Foam Tape 4" x 5.5 yds 3/bx, 6 bx/cs, (18/cs)	NON260604	1528-4		
Elastic Net Dressing	Medline	Glenwood Western	Medi-Tech "Spandage"	Deroyal "Stretch Net" Latex
Elastic Net Dressing Fingers, Toes, Wrist	NONNET01 (Size 1)	GLLF2501 (Size 1)	MT01 (Size 1)	10-7125
Elastic Net Dressing Small Hand, Arm, Leg, Foot	NONNET02 (Size 2)	GLLF2502 (Size 2)	MT02 (Size 2)	10-7126
Elastic Net Dressing Medium Hand, Arm, Leg	NONNET03 (Size 3)	GLLF2503 (Size 3)	MT03 (Size 3)	10-7127
Elastic Net Dressing Large Hand, Arm, Leg, Foot	NONNET04 (Size 4)	GLLF2504 (Size 4)	MT04 (Size 4)	10-7128
Elastic Net Dressing Small Head, Shoulder, Thigh	NONNET05 (Size 5)	GLLF2505 (Size 5)	MT05 (Size 5)	10-7129
Elastic Net Dressing Large Head, Shoulder, Thigh	NONNET06 (Size 6)	GLLF2506 (Size 6)	MT06 (Size 6)	10-7130
Elastic Net Dressing Small Chest, Back, Perineum, Axilla	NONNET07 (Size 7)	GLLF2507 (Size 7)	MT07 (Size 7)	10-7131
Elastic Net Dressing Medium Chest, Back, Perineum, Axilla	NONNET08 (Size 8)	GLLF2508 (Size 8)	MT08 (Size 8)	10-7132
Elastic Net Dressing Large Chest, Back, Axilla, Perineum	NONNET09 (Size 9)	GLLF2509 (Size 9)	MT09 (Size 9)	10-7133
Elastic Net Dressing X-Lg Chest, Back, Axilla, Perineum	NONNET10 (Size 10)	GLLF2510 (Size 10)	MT10 (Size 10)	10-7134
Elastic Net Dressing XX-Lg Chest, Back, Axilla, Perineum	NONNET11 (Size 11)	GLLF2511 (Size 11)	MT11 (Size 11)	10-7135
Medfix Dressing Retention	Medline	3M	Kendall	Smith & Nephew
Medfix 2" x 11 yds 1/bx	MSC4002		9412 (10 yds)	4209 (10 yds 24/cs)
Medfix 4" x 11 yds 1/bx	MSC4004		9414 (10 yds)	4210 (10 yds 24/cs)
Medfix 6" x 11 yds 1/bx	MSC4006		9416 (10 yds)	4216 (10 yds 24/cs)
Medfix EZ 2" x 11 yds 12/bx	MSC4102	2862 (10 yds)		
Medfix EZ 4" x 11 yds 12/bx	MSC4104	2864 (10 yds)		
Medfix EZ 6" x 11 yds 12/bx	MSC4106	2866 (10 yds)		
Medfix EZ 2" x 2 yds 12/bx	MSC4122	2862S (48/cs)		4215 (36/cs)
Medfix EZ 4" x 2 yds 12/bx	MSC4124	2864S (24/cs)		4216 (18/cs)
Medfix EZ 6" x 2 yds 12/bx	MSC4126	2866S (16/cs)		4217 (12/cs)
Megazinc Pink Tape	Medline	Hy-Tape		
Megazinc Pink Tape ½" x 5 yds. (36/cs)	OMAM55CS	105BLF		
Megazinc Pink Tape 1" x 5 yds. (36/cs)	OMAM111CS	110BLF		
Megazinc Pink Tape 2" x 5 yds. (36/cs)	OMAM222CS	120BLF		
Ortho-Porous Sports Tape	Medline	3M	J&J "Zonas"	Kendall "Curity"
Ortho-Porous Tape ½" x10 yds 24/bx, 12 bx/cs (288/cs)	NON260312			
Ortho-Porous Tape 1" x10 yds 12/bx, 12 bx/cs (144/cs)	NON260301		5104	2531
Ortho-Porous Tape 1.5" x10 yds 8/bx, 12 bx/cs (96/cs)	NON260315		5105	
Ortho-Porous Tape 2" x10 yds 6/bx, 12 bx/cs (72/cs)	NON260302		5106	6613
Ortho-Porous Tape 3" x10 yds 4/bx, 12 bx/cs (48/cs)	NON260303		5107	7046
Ortho-Porous Tape 4" x10 yds 3/bx, 12 bx/cs (36/cs)	NON260304			3615

Medline Cross Reference Chart (continued)

Paper Tape	Medline	3M "Micropore"	J&J "Dermiform"	Kendall "Tenderskin"
Paper Tape ½" x10 yds 24/bx, 12 bx/cs (288/cs)	NON260012	15300 (240/cs)	5180	1596
Paper Tape 1" x10 yds 12/bx, 12 bx/cs (144/cs)	NON260001	15301 (120/cs)	5181	1914
Paper Tape 2" x10 yds 6/bx, 12 bx/cs (72/cs)	NON260002	15302 (60/cs)	5182	2419
Paper Tape 3" x10 yds 4/bx, 12 bx/cs (48/cs)	NON260003	15303 (40/cs)	5183	3394
Paper Tape 1" x 1.5 yds 100/bx, 5 bx/cs (500/cs)	NON260001S	1530S-1		
Paper Tape 2" x 1.5 yds 50/bx, 5 bx/cs (250/cs)	NON260002S	1530S-2		
Paper Tape on Dispenser ½" x 10 yds 24/bx, 12 bx/cs (288/cs)	NON260012D	1535-0 (240/cs)		
Paper Tape on Dispenser 1" x 10 yds 12/bx, 12 bx/cs (144/cs)	NON260001D	1535-1 (120/cs)		
Paper Tape on Dispenser 2" x 10 yds 6/bx, 12 bx/cs (72/cs)	NON260002D	1535-2 (60/cs)		
Paper Tape on Dispenser 3" x 10 yds 24/bx, 12 bx/cs (48/cs)	NON260003D	1535-3 (40/cs)		
Tan Paper Tape ½" x 10 yds 24/bx, 12bx/cs (288/cs)	NON260012T	15350 (240/cs)		
Tan Paper Tape 1" x 10 yds 24/bx, 12bx/cs (144/cs)	NON260001T	15351 (120/cs)		
Tan Paper Tape 2" x 10 yds 6/bx, 12bx/cs (72/cs)	NON260002T	15352 (60/cs)		
Transparent Tape	Medline	3M "Transpore"	J&J "Dermiview"	Kendall "Curity"
Transparent Tape ½" x10 yds 24/bx, 12 bx/cs (288/cs)	NON260212	15270(240/cs)		8533
Transparent Tape 1" x10 yds 12/bx, 12 bx/cs (144/cs)	NON260201	15271 (120/cs)	5161 (120/cs)	8534
Transparent Tape 2" x10 yds 6/bx, 12 bx/cs (72/cs)	NON260202	15272 (60/cs)	5162 (60/cs)	8535
Transparent Tape 3" x10 yds 4/bx, 12 bx/cs (48/cs)	NON260203	15273 (40/cs)	5163 (40/cs)	8536
Transparent Tape 1" x1.5 yds 100/bx, 5 bx/cs (500/cs)	NON260201S	1527S-1		
Transparent Tape 1" x10 yds 50/bx, 5 bx/cs (250/cs)	NON260202S	1527-2		
Waterproof Tape	Medline "Latex-free"	3M	J&J "Contains Latex"	Kendall "Contains Latex"
Waterproof Tape ½" x 10 yds 24/bx, 12bx/cs (288/cs)	NON260512		5131	3063
Waterproof Tape 1" x 10 yds 12/bx, 12bx/cs (144/cs)	NON260501		5132	3142
Waterproof Tape 2" x 10 yds 6/bx, 12bx/cs (72/cs)	NON260502		5134	3267
Waterproof Tape 3" x 10 yds 4/bx, 12bx/cs (48/cs)	NON260503		5135	3354
Wound Closure Strips	Medline	3M "Steri-Strip"	J&J "Proxi-Strip"	Kendall "Curi-Strip"
Medi-Strips ½" x 3" 5/pk 50 pk/bx, 4 bx/cs (1000/cs)	NON250318	R1540	100K	9890
Medi-Strips ¼" x 3" 3/pk 50 pk/bx, 4 bx/cs (600/cs)	NON250314	R1541	101K	9892
Medi-Strips ¼" x 1.5" 6/pk 50 pk/bx, 4 bx/cs (1200/cs)	NON250114	R1542	102K	9891
Medi-Strips ¼" x 4" 10/pk 50 pk/bx, 4 bx/cs (2000/cs)	NON250414	R1546	103K	9893
Medi-Strips ½" x 4" 6/pk 50 pk/bx, 4 bx/cs (1200/cs)	NON250412	R1547	104K	9894
Medi-Strips 1" x 5" 4/pk 25 pk/bx, 4 bx/cs (400/cs)	NON250501	R1548		


Designing to Standards

To ensure the highest standards of quality, all Medline tapes are designed and manufactured using quality-controlled professional metrics. These standards have an important role and are incorporated at every stage; scientists and technicians use them in their laboratories; designers use them in their plans; and quality engineers measure them in inspections. All surgical tapes continually undergo rigorous testing, analysis and quality trials to refine and perfect the product. Through this Medline is able to provide improved products that consistently meet your needs for adhesion, strength, conformability, flexibility and maintenance of skin health.

Tests include:

- Test Method of Weight per Unit Area (Unit: g/m²) assesses adhesive mass which creates uniform adhesive distribution
- Test Method of Air Permeability (Unit: cm³/cm²s1) rates air transmission through tape
- Test Method of Water Vapor Permeability (Unit: g/m²/24 hrs.) rates moisture transmission through tape
- Test Method of Ball Tack (Unit:cm) assesses correct level of adhesion
- Test Method of Tensile Strength (Unit: N/2.5 cm) measures tensile strength of tape
- Test Method of Elongation at Break % (Unit: N/2.5 cm) assesses flexibility and stretch that contribute to conformability
- Test Method of Unwinding (Unit: N/width of sample, cm) force required to remove tape from roll
- Test Method of Adhesion on Backing (Unit: N/cm) measures correct level of adhesion on backing
- Test Method of Adhesion on Steel (Unit: N/cm) measures and compares correct level of adhesion on standardized steel plate
- Biocompatibility tests of skin irritation and dermal sensitization assesses potential skin irritation

Test Method of Tensile Strength



Test Method of Water Vapor Permeability



Traditional Wound Care:

Unsurpassed quality in a complete line

Medline offers a comprehensive array of traditional wound care products to meet the full range of our customers' needs. Medline features one of the largest selections of latex-free products, including adhesive bandages that are both latex-free and PVC-free. For our customers' convenience, the latex-free symbol is indicated even on the lowest unit of measure sold. All Medline brand products are manufactured in factories that rigorously follow FDA Quality Systems Regulations. We invite you to consult one of our more than 700 sales representatives or 50 wound care specialists to review all your wound care needs today.

