



# PROTEKT WHEELCHAIR CUSHIONS

Protekt Foam  
Medicare code  
E2601




Protekt Gel  
Medicare code  
E2603



Protekt Ultra  
Bariatric  
Medicare code  
E2608




Protekt Plus  
Medicare code  
E2605




Protekt Plus  
Bariatric  
Medicare code  
E2606



Protekt Ultra  
Medicare code  
E2607




Protekt Gel Bariatric  
Medicare code  
E2604




Protekt Back  
Medicare code  
E2611/ E2612



Protekt Foam Bariatric  
Medicare code E2602






### Protpekt™ Foam Cushion

- Helps in the prevention of pressure sores.
- High density polyurethane foam.
- Stretch cover is breathable and comfortable.
- (2) security ties for added safety.
- 1 year warranty

70001	16"x16"x2"
70002	18"x16"x2"
70003	20"x16"x2"
71001	16"x16"x3"
71002	18"x16"x3"

71003	18"x18"x3"
71004	20"x16"x3"
72001	16"x16"x4"
72002	18"x16"x4"
72003	20"x16"x4"



### Protpekt™ Gel Cushion

- Helps in prevention and treatment of pressure sores.
- Gel bladder designed to eliminate migration
- High density polyurethane foam.
- Stretch cover is breathable and comfortable.
- (2) security ties for added safety.
- 1 year warranty

73001	16"x16"x2"
73002	18"x16"x2"
73003	20"x16"x2"
74001	16"x16"x3"
74002	18"x16"x3"

74003	18"x18"x3"
74004	20"x16"x3"
75001	16"x16"x4"
75002	18"x16"x4"
75003	20"x16"x4"




### Protpekt™ Plus Cushion

- Molded high density pressure sensitive foam.
- Leg troughs promote thigh alignment and postural symmetry.
- Raised front ridge controls sliding.
- Low shear stretch nylon cover.
- Safety strap keeps cushion in place.
- 1 year warranty.

76001	16"x16"x3"
76002	18"x16"x3"

76003	18"x18"x3"
76004	20"x16"x3"



### Protpekt™ Ultra Cushion

- Molded high density pressure sensitive foam.
- Soft foam allows immersion of pelvis for comfort and stability.
- Medial and lateral side supports facilitates positioning.
- Trochanter cut outs facilitates proper pressure redistribution.
- Stretch nylon cover is low shear and breathable
- 1 year warranty.

77001	16"x16"x3"
77002	18"x16"x3"

77003	20"x16"x3"
-------	------------




### Protpekt™ Back Cushion

- Improves seating posture while providing back support and comfort.
- Conforms to provide support to lumbar region.
- Stabilizing board insert prevents hammocking in sling back wheelchair.
- Stretch nylon cover is low shear and breathable.
- Adjustable safety straps for added security.
- 1 year warranty.

79001	16"x17"
79002	18"x17"

79003	20"x19"
79004	22"x19"




### Protpekt™ Foam Bariatric Cushion

- Ultra high density foam.
- Helps in the prevention of pressure sores.
- Stretch cover is comfortable and breathable
- (2) security ties for added safety.
- 1 year warranty

71005	22"x18"x3"
71006	24"x18"x3"
71007	26"x18"x3"

72004	22"x18"x4"
72005	24"x18"x4"
72006	26"x18"x4"



### Protpekt™ Gel Bariatric Cushion

- Ultra high density top foam layer.
- Gel bladder designed to eliminate migration
- Stretch cover is breathable and comfortable.
- (2) security ties for added safety.
- 1 year warranty

74005	22"x18"x3"
74006	24"x18"x3"
74007	26"x18"x3"

75004	22"x18"x4"
75005	24"x18"x4"
75006	26"x18"x4"




### Protpekt™ Plus Bariatric Cushion

- Ultra high density pressure sensitive molded foam.
- Leg troughs promote thigh alignment and postural symmetry.
- Raised front ridge controls sliding.
- Low shear stretch nylon cover.
- Safety strap keeps cushion in place.
- 1 year warranty.

76005	22"x18"x3"
76006	24"x18"x3"

76007	26"x18"x3"
-------	------------



### Protpekt™ Ultra Bariatric Cushion

- Ultra high density pressure sensitive foam.
- Soft foam allows immersion of pelvis for comfort and stability.
- Medial and lateral side supports facilitates positioning.
- Trochanter cut outs facilitates proper pressure redistribution.
- Stretch nylon cover is low shear and breathable
- 1 year warranty.

77004	22"x18"x3"
77005	24"x18"x3"

77006	26"x18"x3"
-------	------------