

INVACARE®  
**MOBILAIRE™**  
50 PSI COMPRESSOR



The durable **Invacare® Mobilaire™ 50 PSI Compressor** is ideal for use with medications that require more flow and pressure for proper nebulization.

**Features**

- Low-noise operation
- Easy-to-read pressure gauge
- Non-marring rubber feet
- Suction cup feet for stable operation
- Durable piston
- Adjustable pressure
- Optional J-tube (IRC5070) available
- UL and CSA listed

**HCPCS Code: E0565**



Pressure gauge  
Suction cup feet



**Invacare Corporation**  
www.invacare.com

**USA**  
One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

**Canada**  
570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2012  
Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols ™ and  
©. All trademarks are  
owned by or licensed  
to Invacare Corporation  
unless otherwise noted.  
Specifications are subject  
to change without  
notification.  
Form No. 11-203  
111146



Model #	Dimensions	Weight	Operating Pressure	Pressure Range	Electrical Requirements	Limited Warranty
IRC607	8" L 10.5" H 8" W	11.3 lb.	7 LPM	0-50 PSI	115 VAC, 60 Hz, 1/8 hp Pole Motor	1 Year

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

