

# INVACARE® CONCENTRATOR ACCESSORIES



Invacare® PreciseRx™ Pediatric Flowmeter



Invacare® Check O2 Plus™ Oxygen Analyzer



The **Invacare® Check O<sub>2</sub> Plus™ Oxygen Analyzer** is a single economical instrument that reports all the operating parameters necessary for performing oxygen concentrator checks. It is an efficient and cost-effective analyzer designed with technicians in mind. It checks oxygen concentration, flow rate and patient outlet pressure on any oxygen concentrator. What's more, it requires no routine calibration, and there's no fuel cell to replace. The Invacare Check O<sub>2</sub> Plus Oxygen Analyzer is three analyzers rolled into one. The **Invacare® PreciseRx™ Pediatric Flowmeter** adapts Invacare oxygen concentrators for low-flow usage. This means you can now use this efficient oxygen modality for pediatric applications.

### Invacare® PreciseR<sub>x</sub>™ Pediatric Flowmeter Features

- Delivers flow rates as low as 1/16 LPM or 50 cc/min. with or without humidifier bottles
- Delivers precise low-flow oxygen delivery: 1/16 - 3/4 LPM or 50 - 750 cc/min.
- Easy to install—no tools required
- Requires no permanent concentrator modification
- Prevents back pressure into concentrator
- Compatible with or without humidifier bottle



PreciseR<sub>x</sub> Pediatric Flowmeter includes: flowmeter, tubing, connectors and humidifier bottle (flow range 1/16 - 3/4 lpm). Invacare PreciseR<sub>x</sub> pediatric flowmeter is designed for use exclusively with Invacare oxygen concentrators. The input pressure is 5 psi (34 kPa).

### Invacare® Check O<sub>2</sub> Plus™ Oxygen Analyzer Features

- Measures three parameters, eliminating the need for multiple instruments: oxygen purity, flow rate and patient outlet pressure
- Accurately measures oxygen purity of the Invacare® XPO2™ Portable Concentrator
- Features simple two-button operation and self-diagnostics
- Two AA batteries provide over 16,000 readings
- Shuts off automatically to extend battery life
- Requires no in-field calibration. Includes connecting tube and protective boot



Check O<sub>2</sub> Plus Oxygen Analyzer includes connecting tube and protective boot

**Invacare Corporation**  
www.invacare.com

**USA**  
One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

**Canada**  
570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2013  
Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols™ and  
®. All trademarks are  
owned by or licensed  
to Invacare Corporation  
unless otherwise noted.  
Specifications are subject  
to change without  
notification.  
Form No. 11-191  
Rev. 04/13  
130346



Model #	Weight	Quantity Per Case	Oxygen Flow	Pressure Range	Oxygen Range	Operating Temperature	Limited Warranty
IRC450	6.4 oz.	1	0-10 LPM ±2%	.5-50 PSI ±.5%	21%-95.6% ±1.5%	59°-104° F (15°-40°C)	3 Years
IRCPF16	2.4 lb.	1	2 LPM	5 ± 0.5 PSI		50-95° F (10°-35°C)	1 Year

