


Protekt™ 500

Gel Infused Foam Pressure Redistribution Mattress

Excess heat is drawn away from the patient providing a constant cool sleeping surface



Gel
infused
foam

 **Proactive**
medical products

Medicare
Code
E0184

270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com

The Protekt™ 500 with (VGPT)™ Visco Gel Polymer Technology is a revolutionary support surface that provides unprecedented pressure redistribution, rapid heat dissipation, and vast increases in breathability compared to traditional foam.

- Heel slope redistributes pressure from heel to the entire leg.
- Precision cut air pockets allows deeper submersion and increased air circulation.
- Gel infused foam.

VGPT™ is open-cell Visco Memory foam infused with gel Polymers. Designed for clinical superiority as well as patient comfort. Some of the many benefits of VGPT™ include:

- **Reduced occurrence of pressure sores** by channeling heat away from the body significantly faster than traditional foam, virtually eliminating perspiration.
- **Superior pressure redistribution** is delivered by providing more support and less cradling at higher levels of compression, redistributing pressure over a large surface area eliminating static fatigue.
- **Protects vulnerable heel area** with the patented **Lure-Quilt™ foam** and a 30 degree slope redistributing pressure to the entire leg and away from the heel.
- **Non-prorated 15 year warranty.**

(VGPT)™ Visco Gel Polymer Technology works by allowing air to pass through open cells, channeling it away from the body instead of trapping it within the closed cells of standard memory foam.

The infusion of gel polymers throughout the foam contributes both cooling properties and the virtual elimination of static fatigue.

As a result of the alliance between open cell foam and gel infusion, the Protekt™ 500 holds its shape better, redistributes pressure over a larger surface area, and successfully channels heat away from the body, keeping patients cool, dry and comfortable.

The heel section of the Protekt™ 500 provides unparalleled pressure mapping results. Designed to guard the vulnerable heel area from incapacitating reoccurring pressure ulcers, the heel section of the Protekt™ 500 utilizes a 30 degree slope to redistribute pressure to the entire lower leg, the top layer is finished with Lure-Quilt™, a patented, highly resilient aerated foam, keeping skin cool and dry.

Optional Raised side rails reduce risk of falls and increase patient safety by providing a firm perimeter for egress and ingress.

Rails are 3 inches high with a 24 inch cutout.

81051	Protekt 500	Pressure Relieving Foam Mattress	36"x76"x6"
81052	Protekt 500	Pressure Relieving Foam Mattress	36"x80"x6"
81053	Protekt 500	Pressure Relieving Foam Mattress	36"x84"x6"
81054	Protekt 500 Plus	Pressure Relieving Foam Mattress w/3" Raised Rails	36"x76"x9"
81055	Protekt 500 Plus	Pressure Relieving Foam Mattress w/3" Raised Rails	36"x80"x9"
81056	Protekt 500 Plus	Pressure Relieving Foam Mattress w/3" Raised Rails	36"x84"x9"



270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com