

Protekt™ Aire Overlays 1000 & 1500

Alternating Pressure Overlays

Protekt™ Aire
1000 Fixed



Protekt™ Aire
1500 Adjustable



 **Proactive**
medical products

Medicare
Code
EO277

270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com

FEATURES AND BENEFITS

PROTEKT™ AIRE 1000 FIXED Alternating Pump & Pad

Aids in the prevention and treatment of pressure ulcers. The design of this bubble pad allows greater skin contact.

- Offering an economical pressure redistribution surface.
- Cycle Time: 10 minutes.
- Lighted On/Off switch.
- Hanging bed hooks.
- Medical grade fire retardant PVC.
- Thickness: 11.5 gauge.
- Size of pad: 35.5"W x 78.5"L.
- Weight capacity: 300 Lbs.
- 18 month warranty.

PROTEKT™ AIRE 1500 ADJUSTABLE Alternating Pressure Pump & Pad

Enables clinician to control interface pressure by adjusting air flow.

- Manual air flow adjustment knob to adjust air flow.
- Cycle Time: 10 minutes.
- Lighted On/Off switch.
- Hanging bed hooks.
- Medical grade fire retardant PVC.
- Thickness: 11.5 gauge.
- Size of pad: 35.5"W x 78.5"L.
- Weight capacity: 300 Lbs.
- 18 month warranty.

80005	Protekt™ Aire 1000	Fixed Alternating Pressure Pump & Pad
80010	Protekt™ Aire 1500	Adjustable Alternating Pressure Pump & Pad (with flaps)
80011	Protekt™ Aire 1000 Bubble Pad	Bubble Pad only for Protekt™ Aire 1000
80012	Protekt™ Aire 1500 Deluxe Bubble Pad	Deluxe Bubble Pad only with flaps for Protekt™ 1500



270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com