

Protekt™ Aire 9000BL

Low Air Loss Mattress System with Pulsation Therapy



 **Proactive**
medical products

Medicare
Code
EO277

270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com

The Protekt™ Aire 9000BL combines “true” low air loss with pulsation therapy, and sets the new standard in the prevention and treatment of pressure ulcers.

The powerful yet quiet user friendly pump produces 1300 LPM (liter per minute) of air flowing directly under the patient. The constant high flow of air keeps the patient cool, dry, and virtually free of perspiration.

FEATURES

- 1275 LPM (liters per minute) true low air loss.
- Provides effective true low air loss relief, or pulse therapy along with continuous low air loss relief.
- In pulse mode every minute the mattress pressure is dropped to 80% of the set pressure for few seconds.
- Max Flow mode with 15 minute timer for rapid inflation of the mattress, inflates in 15 - 50 seconds depending on the size of the mattress.
- In upright (fowler) mode the pump inflates the mattress to 50% higher than the set pressure, to eliminate patient bottoming.
- Lock out key completely locks the control panel to eliminate tampering of the patient settings, including the power switch.
- Audio and visual power failure alarms.
- User-friendly control panel with dome keys and 7 segment LED display for simple user interface, displays current mode, and patient comfort pressure level from Soft 1 to Firm 10 (8 ± 4 to 32 ± 6 mmHg).
- Low friction and low shear force producing, urethane coated breathable nylon top sheet, which is liquid proof and highly vapor permeable.
- 2” convoluted safety foam base to provide additional support to the patient in the event of power failure.
- Weight capacity: 1,000 lbs.
- Mattress dimensions: 36”x80”x10”.
- UL/CSA approved power cord.
- 18 month warranty.

80090	Protekt™ Aire 9000 BL	10” Low Air Loss & Alternating Pressure Mattress System
80091	Protekt™ Aire 9000 BL-Pump	Pump only for Protekt™ Aire 9000 BL
80092	Protekt™ Aire 9000 BL-Mattress	Mattress only for Protekt™ Aire 9000 BL
80093	Protekt™ Aire 9000 BL-Cover	Cover only for Protekt™ Aire 9000 BL
80094	Protekt™ Aire 9000 BL-Bladder	Individual Bladder for Protekt™ Aire 9000 BL



270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com