



Silver Sport 1 Wheelchair


Single Axle



HCPCS
K0001

Item #	Color	Seat Width	Armrest	Product Front Riggings	UOM
SSP118FA-SF	Silver Vein 	18"	Fixed Arm	Swing-away Footrests	1/cs
SSP118FA-ELR	Silver Vein 	18"	Fixed Arm	Swing-away Elevating Leg rests	1/cs

- Powder coated silver vein steel frame. (Figure A)
- Padded armrests.
- Nylon upholstery is durable, lightweight, attractive and easy to clean.
- Urethane tires mounted on composite wheels provide durability and low maintenance.
- Comes with push-to-lock wheel locks.
- Carry pocket on backrest provides additional convenience.
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length front riggings.
- Plastic foot plates.
- Weight Capacity: 300 lbs.



specifications

Overall Width (open)	24"	Overall Length (with riggings)	42"
Overall Width (closed)	12.5"	Weight	41 lbs.
Overall Height	36"	Carton Shipping Weight	48 lbs.
Seat-To-Floor Height	19.5"	Weight Limit	300 lbs.
Back Height	16"	Warranty	Limited 3 Year
Seat Dimensions	18" (W) x 16" (D)		

accessories

Anti-Tippers

STDS833N	Anti-Tippers With Wheels	1/pr
STDS818	-	1/pr

Seat Belts

STDS850	Seat Belts, Auto-Clasp Type Closure	1/ea
STDS851	Seat Belts, Velcro® Type Closure	1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric	1/ea
STDS856	Seat Belts, Velcro® Type Closure, Bariatric	1/ea

Telescoping I.V. Pole Attachment

STDS820	Telescoping IV Pole Attachment	1/ea
---------	--------------------------------	------

Overhead Anti Theft Device

STDS823	Overhead Anti-Theft Device	1/ea
STDS821	Overhead Anti-Theft Device With IV Hooks	1/ea
STDS834	Overhead Anti-Theft Device Single Pole	1/ea

Wheel Lock Extension

STDS801GT	Wheel Lock Extension	1/ea
-----------	----------------------	------

Anti-Fold Bar

STDS806	Anti-Fold Bar	1/ea
---------	---------------	------

Wheelchair Cane/Crutch Holder

STDS1034	Wheelchair Cane/Crutch Holder	1/ea
----------	-------------------------------	------

Heel Loops

STDS831	Heel Loops	1/pr
---------	------------	------

Universal Forward Stabilizer

STDS862	Universal Forward Stabilizer	1/pr
---------	------------------------------	------

Limb Support

WASR	Limb Support, Right	1/ea
WASL	Limb Support, Left	1/ea

Chart Carry Pocket

STDS835	Chart Carry Pocket	1/ea
---------	--------------------	------

Universal Oxygen "E" Cylinder Holder

STDS803	Universal Oxygen "E" Cylinder Holder	1/ea
STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment	1/ea