

Conforming Stretch Gauze

McKesson Conforming Stretch Gauze is designed to secure primary dressings to virtually any body contour by stretching and clinging to itself. The breathable, absorbent and soft cushioning provides additional comfort for patients. This gauze product also features finished edges to prevent unraveling and loose lint. Available in sterile or non-sterile forms.

Features and Benefits

- Rayon and polyester blend provides elasticity while retaining its shape and reducing slippage
- Breathable and absorbent for added comfort
- Secures and stretches to conform to the body's contours
- Finished edges to prevent unraveling
- Needs minimal taping for ease of use
- Low-linting formulation helps minimize the risk of fibers being left behind in the wound bed



ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Conforming Stretch Gauze Bandage, 1" x 1.7 yds., Non-Sterile	16-4101	24/bg, 4bg/cs
Conforming Stretch Gauze Bandage, 2" x 4.1 yds., Non-Sterile	41022000	12/bg, 8bg/cs
Conforming Stretch Gauze Bandage, 3" x 4.1 yds., Non-Sterile	41032000	12/bg, 8bg/cs
Conforming Stretch Gauze Bandage, 4" x 4.1 yds., Non-Sterile	41042000	12/bg, 8bg/cs
Conforming Stretch Gauze Bandage, 1" x 1.7 yds., Sterile	16-4151	1/pk, 24pk/bg, 4bg/cs
Conforming Stretch Gauze Bandage, 2" x 4.1 yds., Sterile	16-4152	1/pk, 12pk/bg, 8bg/cs
Conforming Stretch Gauze Bandage, 3" x 4.1 yds., Sterile	16-4153	1/pk, 12pk/bg, 8bg/cs
Conforming Stretch Gauze Bandage, 4" x 4.1 yds., Sterile	16-4154	1/pk, 12pk/bg, 8bg/cs
Conforming Stretch Gauze Bandage, 6" x 4.1 yds., Sterile	16-4155	1/pk, 6pk/bg, 8bg/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them – making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.