

# GripStep®

## Ultra XL Shoe Cover

- Non-skid in all conditions
- Waterproof
- Water absorbent inner layer
- Safe for medical facilities
- Single use
- Non-linting
- Latex-Free



### Specifications

Construction materials	Outer shell: Low density polyethylene (LDPE) fabric Inner layer: Polypropylene fabric
Material weights: Outer LDPE Coating Inner absorbent fabric Total weight	27 GPM (grams per square meter) 45 GPM 72 GPM
Size	Extra Large—up to size 14 Men’s boots
Color	Navy Blue
Seams	Ultrasonically welded
<b>Item Number</b>	<b>Packaging</b>
GRIP4501	150 Pairs/Case
GRIP4502	25 Pairs/Box; 6 Boxes/case



For more information visit:  
[www.baymedicalcompany.com](http://www.baymedicalcompany.com)  
 1-800-977-7888