

GUARDIAN[®]



Mobility Aids and Bath Safety Products
Quality and Innovation Unified in One Name





GUARDIAN
• SIGNATURE™ •

GUARDIAN®

Choose Guardian for your **Durable Medical**



Table of Contents

Retail Merchandising Program	4-5
Walkers	6-9
Pediatric Walkers	10
Walker Accessories	11-13
Rollators	14-17
Canes	18-21
Crutches	22-25
Knee Walkers	25
Bath Safety	26-28
Bath Benches	28-29
Transfer Benches	30-31
Elevated Toilet Seats/Toilet Safety Rails	32-33
Commodes	34-39
Commode Liners	39

Equipment and a better quality of living.

Guide to symbols used in this catalog.

As you shop through these pages, look for the icons below for additional information on ordering details and product features.



Bariatric Products

The Super Strong symbol highlights products that are able to bear heavier weight capacities — as great as 850lbs (386kgs).



Merchandising Program

An essential point of purchase element for your Retail Merchandising Program. Look for this icon for packaging that is ready to display on your retail shelf.



Latex Free Products

Sourcing latex free products should be quick and easy. Look for this symbol for fast, confident ordering.



Tool-Free

Easy-to-assemble! No tools are needed.

Factory Direct Program

Direct shipping from our factory to your door. By eliminating handling and ordering in bulk, you can cut costs and maximize profits!



Limited Warranty

The Medline Limited Warranty means we stand behind our products, 100%.

We guarantee everything we sell. So rest assured, your Medline product is warranted to be free of defects in materials and workmanship for a lifetime for the original consumer purchaser. All Medline products and devices are built to meet our stringent standards, then carefully inspected prior to shipment. Because we are confident in the material and workmanship of our products, we can assure you of years and years of dependable service.

Some exclusions apply: This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and replacement.

Should you have any problem with a Medline product, please call us. In the event of a defect covered by this warranty, we will use our discretion to take care of the problem right away, either by repairing or replacing the product. If you have a question about your Medline product or this warranty, please call us at 1-800-MEDLINE (1-800-633-5463) or speak to your Medline Sales Representative.

This warranty gives you legal rights, and you may have other rights, which vary from state to state.

Retail Merchandising Program

Our Retail Merchandising Program is designed to help you maximize your profitability. From specially designed planograms, to packaging and signage, discover how this program can be an invaluable way to turn your space into a one-stop shop.



Retail Packaging

Reinforce the quality of the Guardian brand and educate customers.

- Explains the different features and benefits of the product.
- Packaging features a “useful product” section to make customers aware of additional products that are available to enhance their mobility.



Innovative Retail Products

Show unique products that give patients and care givers a superior value compared to traditional products.

- Freedom Transport Chair
- Freedom Rollator
- Excel Translator




Product Layouts and Retail Signage

Maximize the visual impact, emphasis, and appeal of your retail section.

- Suggested items to feature based on analysis of top sellers and profit generators
- Walking Aids – 4' and 8' (1.22m and 2.4m) sections
- Bath Safety – 4' and 8' (1.22m and 2.4m) sections
- Includes valance with hardware and aisle violator

Enhance your showroom floors and boost retail product cash sales.



Tap into our expertise.

We've conducted and compiled research from HME providers across the country. With this valuable information, we developed a Retail Check-Up. The Check-up is a comprehensive analysis to help our customers realize areas of opportunity in retail and maximize their growth.

Capitalize on one of the fastest growing, most profitable parts of the HME business.

Now more than ever, HME providers have great opportunity to increase sales, increase profit margins and get paid faster by growing their retail or cash sale business. The explosive growth of baby boomers, the increased emphasis on home care and the continuation of national competitive bidding makes these cash items all the more important to patients and all the more valuable to HME providers.

Create Immediate Selling Opportunities.

This program offers a variety of point-of-sale pieces to help you generate interest and revenue!

For complete details, contact your Medline sales representative, go online to Medline.com, or call 1-800-MEDLINE.



Walkers

GUARDIAN
SIGNATURE

Easy Care® Walkers

Engineered with quality differences that offer maximum stability and longer life.

- Dual folding mechanism is easy to operate with palms or fingers. Ideal for those with limited dexterity.
- Lightweight aluminum frame with steel cross brace provides durability.
- Folds to only 4" for easy transport and storage.
- Reinforced tube connectors at high-stress areas.
- Feature screws that can be retightened, connect side cross braces.
- Tight tolerance between frame and footpieces for rattle-free use.
- 300lb (136kg) weight capacity.

Item #	Description	Width Inside Hand Grips	Height Adjustment	Pkg.
G30755P	Adult	16½" (42cm)	32½"-39½" (82.5-100cm)	4/cs
G30756P	Junior	16½" (42cm)	27¼"-34¼" (69-87cm)	4/cs
G30757W	Adult, 5" (12.7cm) Wheels	16½" (42cm)	33¼"-39¼" (84-99cm)	4/cs
G30758W	Junior, 5" (12.7cm) Wheels	16½" (42cm)	27½"-34½" (70-88cm)	4/cs



GUARDIAN

Two-Button Folding Walkers

- Easy-to-use, large, push-button mechanisms.
- Each side operates independently to allow easy movement through narrow spaces and greater stability when standing up.
- Front cross braces are located near the top of walker to allow full stride.
- 19½"Wx12½"D (50cmWx32cmD) when open.
- Sturdy 1" (2.5cm) diameter aluminum tubing provides stability while remaining lightweight.
- Contoured, latex-free PVC hand grips for patient comfort.
- Folds up to 4" (10cm) for convenient storage and transportation.
- Height adjusts in 1" (2.5cm) increments.
- 300lb (136kg) weight capacity.

Item #	Description	Width Inside Hand Grips	Height Adjustment	Pkg.
MDS864104	Adult	16" (41cm)	32¼"-39¼" (82-99.5cm)	4/cs
MDS86410W4	Adult, 3" (7.6cm) Wheels	16" (41cm)	32½"-39½" (82.5-100cm)	4/cs
MDS86410W54	Adult, 5" (12.7cm) Wheels	16" (41cm)	32½"-39½" (82.5-100cm)	4/cs
MDS86410J4	Junior	16" (41cm)	24¾"-31¾" (64-81cm)	4/cs
MDS86410JW4	Junior, 3" (7.6cm) Wheels	16" (41cm)	25¼"-32¼" (63-82cm)	4/cs
MDS86410JW54	Junior, 5" (12.7cm) Wheels	16" (41cm)	25¼"-32¼" (64-82cm)	4/cs



GUARDIAN

Paddle Walker

- Features a convenient PUSH or PULL to release folding mechanism that is ideal for patients with manual dexterity limitations.
- Each side operates independently to allow easy movement through narrow spaces and greater stability when standing up.
- Sturdy 1" (2.5cm) diameter aluminum tubing provides stability while remaining lightweight.
- Contoured, latex-free PVC hand grips for patient comfort.
- Folds up to 4" (10cm) for convenient storage and transportation.
- Height adjusts in 1" (2.5cm) increments.
- 300lb (136kg) weight capacity.

Item #	Description	Width Inside Hand Grips	Height Adjustment	Pkg.
MDS86410PAD	Adult	17" (43cm)	32"-39" (81-99cm)	4/cs
MDS86410PAD5	Adult, 5" (12.7cm) Wheels	17" (43cm)	32½"-39½" (82.5-100cm)	4/cs

New!



GUARDIAN

3 in 1 Walker

- Can be used as a stand assist, a walker or a toilet safety rail!
- Dual folding mechanism is easy to operate with palms or fingers. Ideal for those with limited dexterity.
- Feature attractive black accents
- 300lb (136kg) weight capacity.
- Each side operates independently to allow easy movement through narrow spaces and greater stability when standing up.
- Comes standard with a convenient front pouch.

Item #	Description	Width Inside Hand Grips	Height Adjustment	Pkg.
MDS86410UR	Adult	17" (43cm)	32"-39" (81-99cm)	2/cs

New!



MDS86410UR

GUARDIAN SIGNATURE

Red•Dot® Walkers

- Lightweight aluminum frame, with the exclusive Guardian Signature steel "K" cross brace. Chrome plated for added durability.
- Folds to only 4½" (11.4cm) for easy transport and storage.
- Preferred by institutions across the nation.
- Reinforced tube connectors at high-stress areas.
- Feature screws that can be retightened, connect side cross braces.
- Tight tolerance between frame and footpieces for rattle-free use.
- 350lb (159kg) weight capacity.

Item #	Description	Width Inside Hand Grips	Height Adjustment	Pkg.
G07755	Adult	16½" (42cm)	32½"-39½" (82-100cm)	4/cs
G07756	Junior	16½" (42cm)	27¼"-34¼" (69-87cm)	4/cs



G07755

GUARDIAN SIGNATURE

One Button Folding Walker

- One-button folding offers quick, convenient operation and is ideal for those who are weaker on one side.
- Lightweight aluminum frame supports up to 300lbs (136kg) yet weighs only 5lbs (2kg).
- Folds to 4½" for ease of transport and storage.
- Reinforced tube connectors at high-stress areas.
- Feature screws that can be retightened, connect side cross braces.
- Tight tolerance between frame and footpieces for rattle-free use.

Item #	Description	Width Inside Hand Grips	Height Adjustment	Pkg.
G30760P	Adult	16½" (42cm)	31½"-38½" (80-98cm)	4/cs
G30761P	Junior	16½" (42cm)	27¼"-34¼" (69-87cm)	4/cs



G30760P

GUARDIAN

Rigid Walker

- 18" (46cm) between hand grips.
- Height adjusts in one-inch increments from 31"-37½" (79-95cm) from the floor.
- 250lb (113kg) weight capacity.

Item #	Description	Width Inside Hand Grips	Height Adjustment	Pkg.
MDS86204	Rigid Walker	18" (46cm)	31"-37½" (79-95cm)	4/cs



MDS86204

Walkers

GUARDIAN

One-Button Folding Aluminum Walkers

- Rugged, cross bar provides excellent strength without interfering with patient stride.
- Contoured, hand grips maximize patient comfort.
- Deluxe, rubber tips provide excellent traction.
- Folds to a depth of 5" (12.7cm) for convenient storage.
- Height adjusts in 1" (2.5cm) increments.
- 250lb (113kg) weight capacity.



Item #	Description	Width Inside Hand Grips	Width Inside Base	1" (2.5cm) Min./Max. Height Adjustment	Pkg.
MDS86615	Adult	18" (46cm)	21" (53cm)	32"-36" (81-91cm)	2/cs
MDS86617	Junior	18" (46cm)	21" (53cm)	28"-32" (71-81cm)	2/cs

GUARDIAN SIGNATURE

Extra-Duty Walkers

- Lightweight aluminum frame with exclusive Guardian Signature Steel "K" cross brace for extra clearance – weighs only 7½ pounds.
- Chrome plated for added durability
- 19" (48cm) between hand grips.
- Reinforced tube connectors at high-stress areas.
- Feature screws that can be retightened, connect side cross braces.
- Tight tolerance between frame and footpieces for rattle-free use.
- 400lb (181kg) weight capacity.



Item #	Description	Width Inside Hand Grips	Height Adjustment	Pkg.
G07767	Adult	19" (48cm)	32½"-39¼" (82-99cm)	Each
G07768	Junior	19" (48cm)	27½"-34¼" (69-87cm)	Each

GUARDIAN SIGNATURE

Bariatric Folding Walker

- Reinforced aluminum frame with dual steel cross brace provides durability.
- Extra wide frames with up to 22" (56cm) between hand grips accommodate larger users.
- Reinforced tube connectors at high-stress areas.
- Feature screws that can be retightened, connect side cross braces.
- Tight tolerance between frame and footpieces for rattle-free use.
- 650lb (295kg) weight capacity.



Item #	Description	Width Inside Hand Grips	Height Adjustment	Pkg.
G30754B	Adult	22" (56cm)	32"-39" (81-99cm)	2/cs

GUARDIAN

**Bariatric Folding Walkers,
500lb (227kg) Wt. Capacity**



- Easy-to-use, push-button mechanisms.
- Each side operates independently to allow easy movement through narrow spaces and greater stability when standing up.
- Front cross-braces are located near top of walker to allow full stride.
- Sturdy, 1" (2.5cm) diameter, aluminum tubing provides stability, while remaining lightweight.
- Contoured, hand grips for patient comfort.
- Folds to a slim 4½" (11cm) for easy transportation and storage.
- Height adjusts in 1" (2.5cm) increments.



MDS86410XW



MDS86410XWW

Item #	Description	Width Outside Hand Grips	Width Inside Hand Grips	Height Adjustment Min./Max.	Pkg.
MDS86410XW	Bariatric, Two button folding walker	22¾" (58cm)	20¼" (51cm)	31"-41" (79-104cm)	Each
MDS86410XWW	Bariatric, Two button folding walker, with 5" (12.7cm) fixed wheels	22¾" (58cm)	20¼" (51cm)	33½"-43¼" (85-110cm)	Each

GUARDIAN
SIGNATURE

Sidestepper Hemi Walker

- Dual grip hemi-walker offers extra assistance when getting in and out of chairs.
- Provides support and stability to those with limited or no dexterity in one arm or hand.
- Lightweight frame is easy to lift and maneuver.
- Wide base provides more stability than a quad cane.
- 250lb (113kg) weight capacity.

Item #	Description	Width Inside Hand Grips	Height Adjustment	Pkg.
G07770	Adult	19" (48cm)	32"-36" (81-91cm)	2/cs
G07771	Junior	19" (48cm)	28½"-32" (72-81cm)	2/cs



G07770

GUARDIAN

Sidestepper Hemi Walker

- Provides support to patients with use of only one side.
- Folds to just 3" (7.6cm) for easy storage.
- Anodized aluminum frame.
- Deluxe rubber tips provide excellent traction.
- 250lb (113kg) weight capacity.

Item #	Description	Height Adjustment	Pkg.
MDS86620	Hemi Walker/Side Stepper	30"-35" (76-89cm)	Each



MDS86620

Pediatric Walkers

GUARDIAN
SIGNATURE

Strider™ Prescription Walker

- Accommodates most physical requirements with the optional platform attachments.
- Handles adjust without tools to an infinite number of positions to provide stability and posture. 90° bend standard (Tyke models are 45°).
- Extra wide base provides stability, safety and balance.
- Outfitted with 5" front wheels to easily roll over irregular surfaces.

Item #	Description	Width Inside Hand Grips	Height Adjustment	Pkg.
G07782	Tyke-Red	Adjustable	14"-19½" (36-50cm)	Each
G07783	Tweener-Green	Adjustable	24"-36" (61-91cm)	Each



G07782

GUARDIAN
SIGNATURE

Children's Folding Walker

- Lightweight anodized aluminum "A" frame design provides greater balance.
- Height-adjustable "grow with you" design.
- Comfortable vinyl hand grips designed for a child's hand.
- Noise-free silencer on footpieces.
- Folds easily for travel.

Item #	Description	Width Inside Hand Grips	Height Adjustment	Pkg.
G07749	Child-Brite	12¾" (32cm)	24½"-28½" (62-72cm)	Each
G07750	Toddler-Red	10¾" (27cm)	20½"-24½" (52-62cm)	Each



G07749

GUARDIAN
SIGNATURE

Children's Non-Folding Walker

- 3" wheels maneuver over normal floor obstacles: cracks, carpet, or outside play areas.
- Large rear rubber tips provide long braking action.
- Lightweight, sturdy aluminum frame with proportionally large width to length for wide stance and greater stability
- Height-adjustable.
- 80lb (36kg) weight capacity.

Item #	Description	Width Inside Hand Grips	Height Adjustment	Pkg.
G07751	Tyke	13" (33cm)	15½"-20½" (39-52cm)	Each



G07751





Walker Accessories



Front Wheel Attachments

- Front wheel footpieces not only allow you to easily maneuver over surfaces which provide difficulty for walker users, but they also provide a more natural walking stride.
- Available in a choice of 3" (7.6cm), 5" (12.7cm) or swivel casters which allow you to change directions with ease.

Note: Wheel attachments are to be used on the front legs of the walker only.

Item #	Description	Dia.	Pkg.
G07722-8	5" (12.7cm) Wheels, Set of Two	1"	1 pr
G07722G	5" (12.7cm) Wheels, Set of Two	1"	3 pr/cs 
G07723	5" (12.7cm) Wheels, Set of Two	7/8"	1 pr
G07722-8B	5" (12.7cm) Bariatric Wheels, Set of Two	1"	1 pr
G07811	5" (12.7cm) Swivel Caster, Set of Two	1"	1 pr 
G07825-8	3" (7.6cm) Swivel Caster, Set of Two	1"	1 pr 
G07731-8	3" (7.6cm) Wheels, Set of Two	1"	1 pr
G07731G	3" (7.6cm) Wheels, Set of Two	1"	1 pr 
G07729	3" (7.6cm) Wheels, Set of Two	7/8"	1 pr



G07722-8



G07811




G07731-8



Glide Caps

- Included with all wheel footpieces.
- Allows walkers to easily and quietly glide across most non-abrasive surfaces.
- Fits tips with 1½" O.D.

Item #	Description	Pkg.
G07904	Glide Caps	4 caps per bag, 6 bags per case 
G07900CS	Glide Caps, Set of Two	Bulk, 24 pr



G07904



Glide Caps

- Allows walkers to glide easily across most non-abrasive surfaces.
- Rubber tips on rear legs of walker must be removed before adding the MDS86615CG glides.

Item #	Description	Pkg.
MDS86615CG	Glide Cap	6 pr/cs



MDS86615CG

Walker Accessories



Rear Glidebrake Attachments

- Provides safe, gradual rear leg braking action and smooth forward movement when used with front wheel attachments.
- Brake firmly engages the rubber tip on the floor when weight is placed on the walker.
- Choice of easy rolling wheel brakes or Glidebrakes.

Item #	Description	Dia.	Pkg.	
G07732G	Glidebrakes™, Set of Two	1"	3 pr/cs	
G07832-8	Wheel Brakes, Set of Two	1"	1 pr	



G07732G



G07832-8



Footpiece Extensions

- Extend the height of the walker by up to 7".

Item #	Dia.	Pkg.
G07708	1"	1 set (2 extension legs, 2 legs w/ 3" wheels)
G07724	7/8"	1 set (2 extension legs, 2 legs w/ 5" wheels)
G07725	1"	1 set (2 extension legs, 2 legs w/ 5" wheels)
G07712	1"	1 set (4 extension legs)



G07712



Utility Tips

- Heavy-duty rubber tip is reinforced with metal insert for longer wear.
- Fit walkers, commodes and transfer benches.

Item #	Description	Dia.	Pkg.
G01182	Gray	3/4"	8 st/cs
G01202	Gray	1"	8 st/cs
G01212	Gray	1 1/8"	8 st/cs



G01212



Walker Accessories



Folding Flip Tray
G07850



Walker Seat
G07872



Basket Organizer
G07878



Walker Basket
G07715



Front Pouch
G07741



Platform Attachment
G07702 (Adult)
G07703 (Junior)



GUARDIAN
SIGNATURE

Walker Accessory Selection Guide

Model	Description	Size	G07755 G07756 G07767 G07768 Red • Dot	G30755P G30757W G30756P G30758W Easy Care	G30760P/G30761P Single	G07783 G07782 Strider	G07749/G07750 G07751 Pediatric	G30754B Bari
G07731-8	3" (7.62cm) Fixed Wheels - pair	1" (2.54cm)	• • • •	• • • •	•			
G07731G	3" (7.62cm) Fixed Wheels - pair Retail Pkg	1" (2.54cm)	• • • •	• • • •	•			
G07729	3" (7.62cm) Fixed Wheels - pair	7/8" (2.22cm)					•	
G07708	3" (7.62cm) Fixed Wheels Extended - set	1" (2.54cm)	• • • •	• • • •	•			
G07722-8	5" (12.7cm) Fixed Wheels - pair	1" (2.54cm)	• • • •	• Inc • Inc	•			
G07722G	5" (12.7cm) Fixed Wheels - pair Retail Pkg	1" (2.54cm)	• • • •	• Inc • Inc	•			
G07722-8B	5" (12.7cm) Fixed Wheels bariatrics - pair	1" (2.54cm)						•
G07723	5" (12.7cm) Fixed Wheels - pair	7/8" (2.22cm)				Inc Inc		
G07724	5" (12.7cm) Fixed Wheels Extended - set	7/8" (2.22cm)					• •	
G07725	5" (12.7cm) Fixed Wheels Extended - set	1" (2.54cm)	• • • •	• • • •	•			
G07825-8	3" (7.62cm) Swivel Wheels - pair	1" (2.54cm)	• • • •	• • • •	•			
G07820-8	3" (7.62cm) Swivel Wheels w/Glidebrakes	1" (2.54cm)	• • • •	• • • •	•			
G07811	5" (12.7cm) Swivel Wheels - pair	1" (2.54cm)	• •	• •				
G07732G	GlideBrakes - pair	1" (2.54cm)	• • • •	• • • •	•			
G07832-8	Wheel Brakes - pair	1" (2.54cm)	• • • •	• • • •	•			
G07900CS	Glide Caps - Bulk 24 prs/case	1½" (3.81cm)	• • • •	• • • •	•	• •		•
G07904	Glide Caps - 2 pairs/retail pkg	1½" (3.81cm)	• • • •	• • • •	•	• •		•
G07712	Extension Footpieces - set of 4	1" (2.54cm)	• • • •	• • • •	•			
G07796	Platform Attachment	1" (2.54 cm)						
G07797	Platform Attachment	7/8" (2.22 cm)				•		
G07798	Platform Attachment	3/4" (1.9 cm)				•		
G07702	Platform Attachment, Adult	1" (2.54 cm)	• •	• •				
G07703	Platform Attachment, Youth	1" (2.54 cm)	• •	• •				
G07741	Carrying Pouch - Front		• • • •	• • • •	•	• •		
G07850	Fliptray		• •	• • • •				
G07715	Carrying Basket		• •	• • • •		• •		
G07713	Carrying Basket - One Button				•			
G07872	Walker Seat		• •	• • • •				

Rollators

GUARDIAN

Ultralight Rollator

- Weighs only 11lbs (5kgs)!
- Seat height adjustable to accommodate people 4'10"-6'4" (1.5-1.9m).
- Serves as a standard and a junior rollator.
- Nylon bag under the seat with water resistant pouch and shoulder strap.
- Seat measures 12"x12" (30x30cm).
- Loop brakes – pull to stop; push to lock.
- 17" (43cm) distance between handles.
- 250lb (113kg) weight capacity.



Item #	Color	Seat Seat Height	Height Min./Max.	Wheels	Pkg.
MDS86825SLR	Burgundy	18"-23" (46-58cm)	29"-33" (74-84cm)	6" (15cm)	Each
MDS86825SLB	Smoky Blue	18"-23" (46-58cm)	29"-33" (74-84cm)	6" (15cm)	Each



MDS86825SLB

GUARDIAN

Luxe Rollator

- Serves as standard and junior rollator.
- Nylon bag under the seat with water resistant pouch and shoulder strap.
- Seat measures 12"x14" (30.5x35cm).
- Loop brakes-pull to stop; push to lock.
- 19" (48cm) between handles.
- Features a convenient cup holder!
- 300lb (136kg) weight capacity.



Item #	Color	Seat Seat Height	Height Min./Max.	Wheels	Pkg.
MDS86835W	White	18"-23" (46-58cm)	29"-33" (74-84cm)	7" (18cm)	Each

New!



MDS86835W

GUARDIAN

Basic Rollator

- Handles are height adjustable in 1" (2.5cm) increments.
- Includes wire basket underneath seat.
- Comfortable 12"x12" (30x30 cm) padded seat.
- Aluminum frame.
- Weighs only 14lbs (6kgs)!
- Folds easily for compact storage.
- Loop brakes – pull to stop; push to lock.
- 6" (15cm) wheels.
- 250lb (113kg) weight capacity.



MDS86850E



MDS86850EB

Item #	Color	Distance Between Handles	Seat Height	Height Min./Max.	Wheels	Pkg.
MDS86850E	Burgundy	17½" (44cm)	20½" (52cm)	31"-35" (79-89cm)	6" (15cm)	Each
MDS86850EB	Blue	17½" (44cm)	20½" (52cm)	31"-35" (79-89cm)	6" (15cm)	Each



Envoy 480 – Rolling Walker

- Lightweight, aluminum frame has height adjustable handgrips.
- Loop brakes – pull to stop; push to lock.
- Deluxe padded seat and curved padded backrest provide enhanced support and comfort.
- Large 8" (20cm) wheels are great for outdoor terrain.
- Large basket and seat are standard.
- 300lbs (136kg) weight capacity.

Item #	Color	Distance Between Handles	Seat Height	Handle Height Min./Max.	Wheels	Pkg.
G07887B	Blue	18" (46cm)	20" (51cm)	30¾"-36" (78-91cm)	8" (20cm)	Each
G07887R	Rose	18" (46cm)	20" (51cm)	30¾"-36" (78-91cm)	8" (20cm)	Each



G07887R



Retail and Breast Cancer Rollators

- Fully padded seat 14"x13" (35.5x33cm) and padded backrest for total comfort while sitting after a long walk.
- Lightweight aluminum frame.
- Weighs only 17lbs (8kg)!
- Loop brakes – pull to stop; push to lock.
- Large 15"x9" (38x23cm) basket underneath seat.
- Large 8" (20cm) wheels are ideal for outdoor terrain.
- 300lb (136kg) weight capacity.



MDS86825



MDS86825BC



Item #	Color	Distance Between Handles	Seat Height	Height Min./Max.	Wheels	Pkg.
MDS86825	Royal Blue	20" (51cm)	22" (56cm)	32½"-37½" (83-95cm)	8" (20cm)	Each
MDS86825BC	Breast Cancer, Pink Ribbon	20" (51cm)	22" (56cm)	32½"-37½" (83-95cm)	8" (20cm)	Each



Comfort Glide Rollator

- Memory foam seat contours comfortably to the end user.
- Lightweight aluminum frame.
- Loop brakes-pull to stop, push to lock.
- Large 15"x9" (38-23cm) basket underneath seat.
- 8" (20cm) wheels are great for outdoor terrain.
- Comes standard with a convenient cup and cane holder.
- 300lb (136kg) weight capacity.



New!



MDS86826W

Item #	Color	Distance Between Handles	Seat Height	Height Min./Max.	Wheels	Pkg.
MDS86826W	Anodized Blue	20" (51cm)	22" (56cm)	32½"-37½" (83-95cm)	8" (20cm)	Each

Rollators

GUARDIAN
SIGNATURE

Envoy 460 – Economy Rolling Walker

- Weighing 15½lbs (7kg), the lightweight, aluminum frame with 6" (15cm) wheels loads into vehicles with ease.
- Loop brakes – pull to stop; push to lock.
- Padded seat and backrest provide comfort.
- Store personal items securely in basket underneath seat.
- 275lbs (125kg) weight capacity.

Item #	Color	Distance Between Handles	Seat Height	Handle Height Min./Max.	Wheels	Pkg.
G07886G	Green	18½" (47cm)	24" (61cm)	34"-36" (86-96cm)	6" (15cm)	Each
G07886B	Blue	18½" (47cm)	24" (61cm)	34"-36" (86-96cm)	6" (15cm)	Each
G07886RD	Burgundy	18½" (47cm)	24" (61cm)	34"-36" (86-96cm)	6" (15cm)	Each



G07886G

GUARDIAN

Deluxe Rollator

- 36" (91cm) from floor to middle of hand grip in highest position.
- 14"x14" (36x36cm) padded seat.
- Features a convenient handle on the seat for easy folding.
- Folded dimension 14" (36cm).
- Basket dimensions 14"x9"x6" (36x23x15cm).
- Ergonomically designed loop brakes – pull to stop; push to lock.
- Total weight equipped with basket is 16lbs (7.3cm).
- 6" (15cm) non-marking indoor/outdoor wheels.
- 250lb (113kg) weight capacity.

latexfree
PRODUCT

Item #	Color	Distance Between Handles	Seat Height	Height Min./Max.	Wheels	Pkg.
MDS86810	Burgundy	19¾" (50cm)	22" (56cm)	32½"-37" (83-94cm)	6" (15cm)	Each
MDS86810G	Hunter Green	19¾" (50cm)	22" (56cm)	32½"-37" (83-94cm)	6" (15cm)	Each
MDS86810B	Royal Blue	19¾" (50cm)	22" (56cm)	32½"-37" (83-94cm)	6" (15cm)	Each
MDS86810BLK	Black	19¾" (50cm)	22" (56cm)	32½"-37" (83-94cm)	6" (15cm)	Each



MDS86810

Burgundy
MDS86810

Hunter Green
MDS86810G

Royal Blue
MDS86810B

Black
MDS86810BLK

Replacement Parts	Item #	Pkg
Brake System for MDS86810	MDS86810BRK	1 Pair
Basket for MDS86810	MDS86800BSK	Each

GUARDIAN
SIGNATURE

Envoy 380 – 3 Wheel Rolling Walker

- Powder-coated steel frame for long life.
- Tri-wheel design with 8" (20cm) wheels is easy to maneuver.
- Loop brakes – pull to stop; push to lock.
- Folds to only 9" for easy transport and storage.
- 250lbs (113kg) weight capacity.

Item #	Color	Distance Between Handles	Seat Height	Handle Height Min./Max.	Wheels	Pkg.
G07981T	Titanium	14½" (37cm)	N/A	31¼"-36¾" (81-93cm)	8" (20cm)	2/cs

Accessories	Item #	Pkg
Tray and Basket	G07913	2/cs
Tote Bag	G07912	2/cs



G07981T



Envoy 480HD – Heavy Duty Rolling Walker

- Extra wide design, with 21 ½" (55cm) between handles, accommodates larger users.
- Aluminum frame has height adjustable handgrips.
- Loop brakes – pull to stop; push to lock.
- Deluxe padded seat and curved padded backrest provide enhanced support and comfort.
- Large 8" (20cm) wheels are great for outdoor terrain.
- 375lbs (170kg) weight capacity.

Item #	Color	Distance Between Handles	Seat Height	Handle Height Min./Max.	Wheels	Pkg.
G07889T	Titanium	21 ½" (55cm)	20" (51cm)	31"-36 ¾" (79-93cm)	8" (20cm)	Each



G07889T



Envoy 485B – Bariatric Rolling Walker

- Reinforced steel frame offers greater durability and support.
- Extra width between handles (22 ¼") accommodates larger users.
- Handles point forward to promote better posture and walking gait.
- Loop brakes – pull to stop; push to lock.
- Curved padded backrest provides comfort and support.
- Large 8" (20cm) wheels are great for outdoor terrain.
- 500lbs (227kg) weight capacity.

Item #	Color	Distance Between Handles	Seat Height	Handle Height Min./Max.	Wheels	Pkg.
G07890TB	Titanium	22 ¼" (56cm)	22" (56cm)	32 ¼"-34 ¼" (82-87cm)	8" (20cm)	Each



G07890TB



Bariatric Rollator

- Comfortable padded seat – 17 ½"Wx14"D (44cmWx36cmD).
- Comfortable curved backrest.
- Loop brakes – pull to stop; push to lock.
- Reinforced forks to support higher weight capacity.
- Lightweight aluminum makes transporting easy.
- Weighs only 19lbs (8.6kg)!
- Folds easily to place in trunk of car or to store.
- 8" (20cm) wheels.
- 400lb (181kg) weight capacity.



Item #	Color	Distance Between Handles	Seat Height	Height Min./Max.	Wheels	Pkg.
MDS86800XW	Burgundy	23 ½" (60cm)	22" (56cm)	31 ½"-37 ½" (80-95cm)	8" (20cm)	Each



MDS86800XW

GUARDIAN

Bariatric Quad Canes

- Black, all steel construction.
- Double-plated base with extremely-strong welds.
- Solid brass, push button adjustment for extra stability.
- Secured with a locking ring for added safety.
- Comfortable, black, PVC handgrip.
- Dual-Tubing insert for added strength and safety.
- Height adjusts in 1" (2.5cm) increments.
- 500lb (227kg) weight capacity.



MDS86228XWG –
Double footplate and cane



MDS86228XWG

Item #	Base	Adjustment	Footprint	Cane Weight	Pkg.
MDS86222XWG	Small	29"-38" (74-96.5cm)	6"x8" (15x20cm)	3.6lbs (1.6kg)	Each
MDS86228XWG	Large	29"-38" (74-96.5cm)	8"x12" (20x30cm)	4.6lbs (2kg)	Each

GUARDIAN

Bariatric Offset-Handle Cane

- Black, all steel construction.
- Solid brass, double push-button adjustment for extra stability.
- Height adjusts in 1" (2.5cm) increments.
- Secured with a locking ring for added safety.
- Comfortable, PVC handgrip.
- Dual-Tubing insert for added strength and safety.
- 500lb (227kg) weight capacity.



MDS86420XWG

Item #	Description	Color	Adjustment	Cane Weight	Pkg.
MDS86420XWG	Offset Handle	Black	29"-38" (74-96.5cm)	3lbs (1.4kg)	Each

GUARDIAN SIGNATURE

Quad Canes

- Handle inserts and reinforced metal inserts at stress points provide great strength and durability.
- Contoured vinyl handgrip is easier to grasp and offers enhanced durability for long-term use.
- 300lb (136kg) weight capacity.

Item #	Description	Size	Height Adjusts	Pkg.
G05340S	Large Base, Brite	Adult	29"-38" (74-97cm)	2/cs
G05345S	Small Base, Brite	Adult	25"-34" (63-86cm)	2/cs

Quadlite

- Molded resin base is strong yet light-weight, offering a slight flex for self-leveling and stability.
- 300lb (136kg) weight capacity.

G05845S

Item #	Description	Size	Height Adjusts	Pkg.
G05845S	Adult, Brite	Adult	31"-40" (79-100cm)	2/cs
G05846S	Youth, Brite	Youth	31"-40" (79-100cm)	2/cs



GUARDIAN

Aluminum, Quad Canes

- Attractive finishes resist marring and scratching.
- Push-button, adjustment mechanism with locking ring allows quick, secure fitting.
- Height adjusts in 1" (2.5cm) increments.
- Non-skid, rubber tips with metal inserts ensure a secure-grip and long-life.
- No assembly required.
- 250lb (113kg) weight capacity.

Item #	Base	Color	Adjustment	Footprint	Cane Weight	Pkg.
MDS86222G	Small	Black	29"-38" (74-96.5cm)	6"x8" (15x20cm)	2lbs (1kgs)	2/cs
MDS86222BRZ	Small	Bronze	29"-38" (74-96.5cm)	6"x8" (15x20cm)	2lbs (1kgs)	2/cs
MDS86222CHRG	Small	Chrome	29"-38" (74-96.5cm)	6"x8" (15x20cm)	2lbs (1kgs)	2/cs
MDS86228G	Large	Black	30"-39" (76-99cm)	8"x12" (20x30cm)	2.6lbs (1.2kgs)	2/cs
MDS86228BRZ	Large	Bronze	30"-39" (76-99cm)	8"x12" (20x30cm)	2.6lbs (1.2kgs)	2/cs
MDS86228CHRG	Large	Chrome	30"-39" (76-99cm)	8"x12" (20x30cm)	2.6lbs (1.2kgs)	2/cs

**RETAIL
PACKAGING**



MDS86228 series

GUARDIAN SIGNATURE

Offset-Handle Canes

- Offers an anatomically correct design that relieves palm fatigue and provides extra comfort for the long term user.
- Contoured vinyl handgrip is easier to grasp and offers enhanced durability for long-term use.
- 300lb (136kg) weight capacity.

Item #	Description	Color	Adjustment	Cane Weight	Pkg.
G05365	Offset Handle	Bronze Slim	31"-40" (79-100cm)	.75lbs (0.3kg)	6/cs
G05355	Offset Handle	Brite	31"-40" (79-100cm)	.75lbs (0.3kg)	6/cs
G05356	Offset Handle	Black	31"-40" (79-100cm)	.75lbs (0.3kg)	6/cs
G05555	Offset Handle with Sure Grip	Brite	31"-40" (79-100cm)	.75lbs (0.3kg)	6/cs



G05356

G05555

GUARDIAN

Offset Handle Aluminum Canes

- Contoured handles with foam hand grips provide a comfortable, secure-grip, complete with convenient wrist strap.
- Attractive finishes resist marring and scratching.
- Handle height adjusts in 1" (2.5cm) increments.
- Secured with a locking ring for added safety.
- 250lb (113kg) weight capacity.

**RETAIL
PACKAGING**

Item #	Color	Adjustment	Cane Weight	Pkg.
MDS86420GU	Black	29"-38" (74-96.5cm)	0.8lbs (0.4kg)	6/cs
MDS86420BRZG	Bronze	29"-38" (74-96.5cm)	0.8lbs (0.4kg)	6/cs
MDS86420CHRG	Chrome	29"-38" (74-96.5cm)	0.8lbs (0.4kg)	6/cs
MDS86420BG	Blue Ice	29"-38" (74-96.5cm)	0.8lbs (0.4kg)	6/cs
MDS86420GG	Green Ice	29"-38" (74-96.5cm)	0.8lbs (0.4kg)	6/cs



MDS86420 Series

GUARDIAN SIGNATURE

Standard Handle Cane

- Standard handle allows user to retain cane over an arm while opening doors or performing other tasks.

Item #	Description	Color	Adjustment	Cane Weight	Pkg.
G30372	Adult	Bronze	31"-40" (79-100cm)	.75lbs (0.3kg)	6/cs



GUARDIAN SIGNATURE

Standard Handle Cane

- Handle inserts and reinforced metal inserts at the stress points provide greater strength and durability.
- 300lb (136kg) weight capacity.
- Contoured vinyl handgrip is easier to grasp and offers enhanced durability for long-term use.

Item #	Description	Color	Adjustment	Cane Weight	Pkg.
G05370	Adult	Chrome	31"-40" (79-100cm)	.75lbs (0.3kg)	6/cs



**RETAIL
PACKAGING**

GUARDIAN

Standard Handle Canes

- Standard handle allows user to retain cane over an arm while opening doors or performing other tasks.
- Contoured handles with foam hand grip provide a comfortable secure-grip.
- 250lb (113kg) weight capacity.

Item #	Description	Color	Adjustment	Cane Weight	Pkg.
MDS86415G	Adult	Black	30"-38" (76-96.5cm)	0.8lbs (0.4kg)	6/cs
MDS86415CHR	Adult	Chrome	30"-38" (76-96.5cm)	0.8lbs (0.4kg)	6/cs

MDS86415CHR MDS86415G

GUARDIAN

Folding Aluminum Cane

- Height adjusts in 1" (2.5cm) increments.
- Comfortable, contour grip.
- Anodized, aluminum shaft.
- Folds to 11" (28cm).
- Locks into place.
- 250lb (113kg) weight capacity.

Item #	Color	Adjustment	Cane Weight	Pkg.
MDS86430WG	Black	32"-36" (81-91cm)	1lb (0.5kg)	6/Cs



MDS86430WG


**RETAIL
PACKAGING**



Fashion Canes

- Available in seven stylish patterns
- Adjusts in 1" (2.5cm) increments.
- 250lb (113kg) weight capacity.

Item #	Description	Adjustment	Cane Weight	Pkg.
MDS86435LP	Leopard, T handle	29"-38" (74-96.5cm)	0.8lbs (0.4kg)	6/cs
MDS86435FL	Floral, T handle	29"-38" (74-96.5cm)	0.8lbs (0.4kg)	6/cs
MDS86420RS	Rose, Offset Handle	29"-38" (74-96.5cm)	0.8lbs (0.4kg)	6/cs
MDS86420CM	Camo, Offset Handle	29"-38" (74-96.5cm)	0.8lbs (0.4kg)	6/cs
MDS86420PR	Purple, Offset Handle	29"-38" (74-96.5cm)	0.8lbs (0.4kg)	6/cs
MDS86420LB	Light Blue, Offset Handle	29"-38" (74-96.5cm)	0.8lbs (0.4kg)	6/cs
MDS86435B	Black, T handle	29"-38" (74-96.5cm)	0.8lbs (0.4kg)	6/cs



**RETAIL
PACKAGING**

From Left to Right:
MDS86435LP
MDS86435FL
MDS86420RS
MDS86420CM
MDS86420PR
MDS86420LB



Wooden Canes

- Attractive, classic style.
- Available in black or mahogany finishes.
- 300lb (136kg) weight capacity.

Item #	Base Dia.	Color	Adjustment	Cane Weight	Pkg.
G05390B	7/8" (2.2cm)	Black	36" (91cm)	0.8lbs (0.4kg)	6/cs
G05390M	7/8" (2.2cm)	Mahogany	36" (91cm)	0.8lbs (0.4kg)	6/cs



Wood Canes

- Available in natural or walnut finish.
- Rubber tip provides good traction.
- 300lb (136kg) weight capacity.

Item #	Base	Color	Adjustment	Cane Weight	Pkg.
MDS80212	7/8" (2cm)	Natural	37" (94cm) long	0.8lbs (0.4kg)	6/cs
MDS80217	1" (2.5cm)	Walnut	37" (94cm) long	0.8lbs (0.4kg)	6/cs



Cane Tips

- Heavy-duty rubber tip is reinforced with metal insert for longer wear.

Item #	Description	Dia.	Pkg.
G30561-PC	Black	5/8" (1.6cm)	1 pr
G30572-PC	Gray	3/4" (2cm)	1 pr
G30581-PC	Black	7/8" (2cm)	1 pr



Quad Cane Tips

- Heavy-duty rubber tip is reinforced with metal insert for longer wear.
- Special tread for skid resistance and secure traction.

Item #	Description	Dia.	Pkg.
G02162	Gray	1/2" (1.3cm)	4/pkg
G02172	Gray	5/8" (1.6cm)	4/pkg



G02162



Standard Replacement Tips

Replacement Part	Description	Size	Color	Packaging
MDS80212T	Cane Tips for Wood Canes	7/8" (2cm) ID	Black	1 pr
MDS80217T	Cane Tips for Wood Canes	1" (2.5cm) ID	Black	1 pr
MDS86424G	Tips for Small Base Canes	1/2" (1.3cm)	Black	2/bx, 6 bx/cs
MDS86425G	Tips for Large Base Canes	5/8" (1.6cm)	Black	2/bx, 6 bx/cs
MDS86426G	Replacement Tips for Offset and Standard Canes	3/4" (2cm)	Black	2/bx, 6 bx/cs
MDS86427G	Replacement Hand Grips for Offset-Handle Canes	4" (10cm)	Black	2/bx, 6 bx/c



Crutches



Forearm Crutches

- Telescoping parts are manufactured with internal bushings and external lock nuts to provide totally silent usage.
- Vinyl-coated, tapered, contoured arm cuffs
- Heavy-wall, high strength aluminum tubing.
- Black components hide signs of wear and scratches

Item #	Description	Patient Height	Pkg.	Weight Capacity
G05160	Tall Adult	5'10"-6'6" (178-198cm)	1 pr	300lbs (136kg)
G05161	Adult	5'-6'2" (152-183cm)	1 pr	300lbs (136kg)




G05161



Aluminum, Forearm Crutches

- Vinyl-coated, tapered, contoured arm, cuffs mold to user for extra comfort.
- Pivoting cuff bracket is secured to a composite post for added durability.
- Comfortable, contoured, vinyl hand grips provide long wear.
- Leg and forearm sections adjust independently for a custom fit.
- All telescoping parts are manufactured with internal bushings and external lock nuts to provide silent usage.
- Recommended for long-term disability.
- All wear items (arm pads, hand grips and tips) are warranted for 6 months.

Item #	Description	Patient Height	Pkg.	Weight Capacity
MDS805160	Tall Adult	5'10"-6'6" (1.8-2m)	1 pr	250lbs (113kg)
MDS805161	Adult	5'0"-6'2" (1.5-1.9m)	1 pr	250lbs (113kg)
MDS805162	Youth	4'2"-5'2" (1.3-1.6m)	1 pr	250lbs (113kg)
CNI600A*	Bariatric (special order)* 		1 pr	600lbs (272kg)
MDS80516CF	Universal Replacement Forearm Cuff		1 pr	
MDS80266RW	Crutch Tips for Bariatric Forearm Crutch		6pr/cs	

*Special Measurement Required



MDS80516CF

MDS805161



Children's Forearm Crutches

- Telescoping parts are manufactured with internal bushings and external lock nuts to provide totally silent usage.
- Vinyl-coated, tapered, contoured arm cuffs mold to user for extra comfort
- Heavy-wall, high strength aluminum tubing.
- Black components hide signs of wear and scratches
- Available to fit lil tykes from 2 feet tall.

Item #	Description	Patient Height	Pkg.	Weight Capacity
G05162	Youth-Brite	4'2"-5'2" (127-158cm)	1 pr	300lbs (136kg)
G05163	Child-Brite	3'2"-4'6" (97-137cm)	1 pr	200lbs (91kg)



G05163



Red•Dot® Crutches

- Clip-lock adjusted handle and push-button adjusted height allows for 10 second adjustments.
- 1" adjustment increments on both handle and height provided a precise fit which reduces muscle fatigue.
- High density latex free rubber cushions and hand grips provide long wear and greater comfort.
- How-to-Fit and Use guides included as part of risk management program.
- 300lbs (136kg) weight capacity.



Item #	Description	Patient Height	Adj. Range	Pkg.
G90-214-8	Tall	5'10"-6'6" (178-198cm)	52"-60" (132-150cm)	8 pr/cs
G91-214-8	Adult	5'1"-5'9" (155-175cm)	44"-52" (112-132cm)	8 pr/cs
G92-214-8	Youth	4'6"-5'2" (132-157cm)	37"-45" (94-114cm)	8 pr/cs



G91-214-8



Bariatric Crutches

- Heavy-duty steel construction provides durability and supports up to 550lbs (249kg).
- Push-button adjusted footpieces with fitting scale allow for accurate, fast fitting.
- 1" (2.54cm) adjustment increments on both handle and height provide a precise fit which reduces muscle fatigue.



Item #	Description	Patient Height	Adj. Range	Pkg.
G60314B	Tall Adult	5'10"-6'6" (1.8-2m)	52"-60" (132-152cm)	1 pr/cs
G61314B	Adult	5'1"-5'9" (155-175cm)	44"-52" (112-132cm)	1 pr/cs



G61314B



Bariatric Steel Crutches

- 650lb (295kg) weight capacity.

Item #	Description	Patient Height	Adj. Range	Pkg.
MDS80334XW	Tall Adult	5'10"-6'6" (1.8-2m)	54"-62" (1.4-1.6m)	1 pr/cs
MDS80335XW	Adult	5'1"-5'9" (1.5-1.8m)	46"-54" (1.2-1.4m)	1 pr/cs



I-beam post



MDS80334XW

Crutches

GUARDIAN

Aluminum Crutches

- Our push-button crutches are lighter than wood crutches thanks to their aluminum construction.
- Double-extruded, center tube provides load-bearing strength where it is needed the most.
- Special footpiece utilizes durable rivets to secure the aluminum tubing from the inside, providing more stability than plastic connectors.
- Patented, I-beam, center post limits footpiece rotation to a 1/4 turn, always keeping the push-button in view and facilitating easy adjustment.
- Long, virtually unbreakable stair deflector helps protect against the crutch slipping on stairs. It helps increase stress tolerance at all angles without cracking.
- Non-skid, rubber tips provide good traction. Built-in, metal ring helps prevent excessive wear and ensures maximum durability.
- Packaged completely accessorized, no assembly required.
- 300lb (136kg) weight capacity, 250lbs (113kg) for child.



Item #	Description	Patient Height	Adjustment Range	Pkg.
MDSV80534	Tall Adult	5'10"-6'6" (1.8-2m)	52½"-60½" (1.3-1.5m)	8 pr/cs
MDSV80535	Med Adult	5'2"-5'10" (1.6-1.8m)	44½"-52½" (113-1.3m)	8 pr/cs
MDSV80536	Youth	4'6"-5'2" (1.4-1.6m)	36½"-44½" (93-113cm)	8 pr/cs
MDS80337Z	Child	4'2"-4'6" (1.3-1.4m)	31"-37" (79-94cm)	2 pr/cs



MDS80535



Push-Button Aluminum Crutch Accessories

Item #	Description	Pkg.
MDS80265W	Crutch Tips, Gray – Super (Extra-Large)	2/bx, 6 bx/cs
MDS80266RW	Crutch Tips, Gray Universal	2/bx, 6 bx/cs
MDS80267W	Crutch Underarm Pad, Foam, Gray* *Does not fit MDS80334XW or MDS80335XW	2/bx, 6 bx/cs
MDS80268W	Crutch Hand Grip, Gray	2/bx, 6 bx/cs
MDS80269	Crutch Kit, Gray	6 kits/ea
MDSPPC107	Crutch Pillows, Charcoal. Helps avoid the negative side effects associated with crutch therapy, chaffing and friction under arms, auxiliary nerve damage, bruised and painful underarms and hands. Includes 1 pair underarm pads and 1 pair palm pads.	2/bx



MDS80265W

GUARDIAN

Children's Crutches

- Push-button adjusted footpiece with fitting guide on crutch footpiece provides accurate fitting.
- 1" adjustment increments on both handle and height provide a precise fit, reducing fatigue.
- Latex free rubber components
- 200lbs (91kg) weight capacity.

Item #	Description	Approximate Patient Height	Height Adjusts	Pkg.
G53314-8	Children's Crutches	4'0"-4'6" (1.2-1.4m)	31½"-37½" (80-95cm)	8 pr/cs



G53314-8



Wood Crutches

- 250lb (113kg) weight capacity.

Item #	Description	Patient Height	Pkg.
MDS80889	Adult, Xlarge, w/ tip, pad, and closed grip	5'8"-7'0" (1.7-2.1m)	10 pr/cs
MDS80891	Adult, Large, w/ tip, pad, and closed grip	5'5"-6'3" (1.6-2m)	10 pr/cs
MDS80895	Medium, w/ tip, pad, and closed grip	4'10"-5'6" (1.5-1.7m)	10 pr/cs
MDS80899	Youth, w/ tip, pad, and closed grip	4'1"-4'10" (1.2-1.5m)	10 pr/cs
MDS80900	Child, w/ tip, pad, and closed grip	3'4"-4'0" (1-1.2m)	10 pr/cs
MDS2633C	Child	26"-33" (66-84cm)	10 pr/cs
MDS4252C	Adult, Medium	42"-52" (1-1.3m)	10 pr/cs
MDS4860C	Adult, Wood w/ Access	48"-60" (1.2-1.5m)	10 pr/cs



MDS80891



Weil Knee Walker

- Provides comfortable alternative to crutches.
- Large basket standard.
- Turns for easy use at home or outside.
- Turning front wheels.
- Attractive burgundy color using Triple Coat painting process for better durability.
- Locking brakes.
- 300lb (136kg) weight capacity.

Knee Walker Measurements

Transportable Position	12½" (32cm)
Lowest Position	17½" (44cm)
Highest Position	24" (44.5m)
Increments	2" (5cm)



MDS86000



Item #	Stability (front wheel separation)	Height Adjustment	Dimensions	Knee Rest Dimensions	Wheels	Pkg.
MDS86000	22" (56cm)	4'10"-6'4" (1.5-2m)	36"Lx22½"W (91cmLx57cmW)	20½"Lx8½" (52cmLx22cmW)	8" (20cm)	Each



Basic Knee Walker

- A great alternative for customers looking for a cost effective knee walker.
- Can be used universally, do not need additional hardware to change from left or right leg.
- Eye catching green color.
- 300lb (136kg) weight capacity.

Item #	Stability (front wheel separation)	Height Adjustment	Dimensions	Knee Rest Dimensions	Wheels	Pkg.
MDS81000	20½" (52cm)	2'8"-3'4" (0.8-1m)	33"Lx20½"W (84cmLx52cmW)	15"Lx9"W (38cmLx23cmW)	6" Frt, 8" Bk (15cm, 20cm)	Each



MDS81000

Bath Safety



GUARDIAN
SIGNATURE

Ergonomic Hand Held Shower

- Easily adapts to existing shower arm, delivering standard or pulsating spray.
- Extra long 86" (218 cm) hose for use while seated.
- On/Off control and massage control are conveniently positioned at your fingertips.
- Diverter valve allows you to use existing shower head.

Item #	Description	Pkg.
G98035	Hand held shower with diverter valve	4/cs



RETAIL
PACKAGING

G98035

GUARDIAN
SIGNATURE

Deluxe Tub Grab Bar

- Clamp design, including stainless steel components.
- Patented step through design provides an unobstructed entry/exit path allowing better hand positioning.
- Bi-level handgrip for maximum flexibility and safety.
- Tool free installation and large rubber-lined pads won't mar tub surface.
- Adjusts to fit tubs 2¾" - 6¼" (7-16cm) thick.

Item #	Description	Pkg.
G98006	Deluxe Plastic Tub Grab Bar	3/cs



RETAIL
PACKAGING

G98006

GUARDIAN
SIGNATURE

Bathtub Safety Rails

- For a safe grip around the tub.
- White powder-coated steel with foam-lined clamps that won't mar tub surface.
- Mounting bracket adjusts 3½" to 6" (9-16cm) to fit most standard tubs.

Item #	Description	Pkg.
G98080	8"Hx11"W (20x27cm)	3/cs



RETAIL
PACKAGING


G98080



Knurled Chrome Wall Grab Bars

- Adjustable oval flange rotates to best position for secure mounting.
- 1" (2.5 cm) diameter, 18 gauge chrome plated steel.
- Knurled finish provides a steady grasp even when wet.
- 300lb (136kg) weight capacity.

Item #	Length	Pkg.
G81020KE	12" (30cm)	3/cs
G81030KE	16" (41cm)	3/cs
G81040KE	18" (46cm)	3/cs
G81050KE	24" (61cm)	3/cs
G81060KE	32" (81cm)	3/cs
G81070KE	36" (91cm)	3/cs



G81060KE



White Finish Wall Grab Bar

- White painted steel matches most decors.
- 1" (2.5cm) diameter, 18 gauge steel with 2½" (6cm) flange.
- White screw covers included.
- 300lb (136kg) weight capacity.

Item #	Length	Pkg.
G87520	12" (30cm)	3/cs
G87530	16" (41cm)	3/cs
G87540	18" (46cm)	3/cs
G87550	24" (61cm)	3/cs
G87560	32" (81cm)	3/cs
G87570	36" (91cm)	3/cs



G87560

GUARDIAN

Grab Bars

- Rugged 1¼" (3cm) diameter, steel bars are easy to grip and help reduce the risk of accidents.
- Can be used horizontally or vertically.
- Complete with mounting hardware.
- 250lb (113kg) weight capacity.

Item #	Color	Length	Diameter	Pkg.
MDS86012CHRG	Knurled Chrome	12" (30.5cm)	1¼" (3cm)	3/cs
MDS86016CHRG	Knurled Chrome	16" (40.6cm)	1¼" (3cm)	3/cs
MDS86018CHRG	Knurled Chrome	18" (45.7cm)	1¼" (3cm)	3/cs
MDS86024CHRG	Knurled Chrome	24" (61cm)	1¼" (3cm)	3/cs
MDS86032CHRG	Knurled Chrome	32" (81.3cm)	1¼" (3cm)	3/cs




Knurled Chrome Grab Bar Family



Safe-T-Pole

- Floor-to-ceiling pole that safely assists standing, sitting, climbing stairs, or transferring.
- No tools required for installation. Screw mechanism with padded bracket fits securely without marring surface.
- Made of white powder-coated, heavy gauge steel for added durability.

Item #	Description	Diameter	Pkg.
G98120B	Base	1½" (4cm)	Each
G98120P	Pole	1½" (4cm)	Each



G98120B
G98120P

CURAD®

Cast Protectors

- Ideal for keeping casts, bandages, or wounds dry while showering.
- Easy one-handed application.
- No need for tape or straps.
- Reusable.



Item #	Description	Pkg.
CUR200AAA	Adult Arm Cast Protector	2/bx, 6 bx/cs
CUR200ALL	Adult Leg Cast Protector	2/bx, 6 bx/cs
CUR100KAA	Kids Arm Cast Protector	2/bx, 6 bx/cs
CUR100KLL	Kids Leg Cast Protector	2/bx, 6 bx/cs



GUARDIAN

Bariatric Bath Benches, 550lb (249kg) Weight Capacity

- Comfortable, blow-molded seat with integrated handles for additional support.
- Non-slip feet secure bench in bathtub.
- Dual cross brace for extra stability with extremely durable welds.
- High-strength, aluminum frame to prevent rusting.
- Aircraft rivets will not loosen.
- Completely assembled.
- 550lb (249kg) weight capacity.



Item #	Description	Seat Height	Seat Depth	Seat Width	Pkg.
MDS89740AXW	Bariatric Bath Bench w/o Back	14"-17" (35.5-43cm)	12" (30.5cm)	20" (51cm)	Each
MDS89745AXW	Bariatric Bath Bench w/ Back	16½"-20½" (42-52cm)	21½" (54.6cm)	18" (46cm)	Each

GUARDIAN SIGNATURE


Shower Chair/Stool

For those who require a basic chair with height adjustability and an extra wide seating surface.

- Sturdy, lightweight, anodized aluminum frame resists rusting.
- Legs adjust in 1" (2.54cm) increments with just a push of a button to ideal height between 13" (33cm) and 20" (51cm).
- Available with a back for extra support and leverage, or backless for easier access in bathing.
- Textured, blow-molded polyethylene wide track seat and backrest with built-in hand-holds provide extra safety.
- Rubber feet grip even wet surfaces.
- 250lb (113kg) weight capacity.



G30402



G30403

Item #	Description	Color	Seat Height Range	Seat Depth	Seat Width	Pkg.
G30400-4	Chair with back, assembled	Gray	13"-20" (33-51cm)	12" (30.5cm)	19" (48cm)	4/cs
G30401-4	Stool without back, assembled	Gray	13"-20" (33-51cm)	12" (30.5cm)	19" (48cm)	4/cs
G30402	Chair with back, unassembled	White	13"-20" (33-51cm)	12" (30.5cm)	19" (48cm)	3/cs
G30403	Stool without back, unassembled	White	13"-20" (33-51cm)	12" (30.5cm)	19" (48cm)	3/cs





Padded Shower Chair

- Extra thick padding on seat and backrest provides additional comfort for long-term use.
- Extra deep 16½" seat for greater stability and comfort.
- Legs adjust in 1" increments to ideal height between 18" (46cm) and 22" (56cm).
- Rubber tips grip bathtub surface.
- 300lb (136kg) weight capacity.

Item #	Description	Seat Height Range	Seat Depth	Seat Width	Pkg.
G98027	Padded Shower Chair	18"-22" (46-59cm)	16¾" (42.5cm)	17" (43cm)	Each



G98027



Bath Benches

- Available in gray (retail package) or white (fully assembled).
- High-density, polyethylene seat and back.
- Durable, aluminum frame is corrosion resistant.
- Angled legs distribute weight over a large area to increase patient security.
- Non-slip tips on all four legs provide extra-stability and safety.
- Height adjusts in 1" (2.5cm) increments.
- 250lb (113kg) weight capacity.



MDS89740RWG



MDS89745R

Item #	Description	Seat Height	Seat Depth	Seat Width	Pkg.
MDS89740A	Without Back, Assembled, White	14½"-19½" (37-49.5cm)	12" (30.5cm)	20" (51cm)	2/cs
MDS89740RWG	Without Back, Gray	14½"-19½" (37-49.5cm)	12" (30.5cm)	20" (51cm)	2/cs
MDS89745A	With Back, Assembled, White	14½"-19½" (37-49.5cm)	12" (30.5cm)	20" (51cm)	2/cs
MDS89745RG	With Back, Gray	14½"-19½" (37-49.5cm)	12" (30.5cm)	20" (51cm)	2/cs



Composite Plastic Bath Benches

- Easy assembly, no tools required.
- Precise, height-adjustable legs adjust in one-inch increments.
- Built-in, handheld, shower holder on each side of chair.
- Optional, removable arms to help aid in patient balance.
- Durable and rustproof all plastic construction.
- Non-slip finish on seat.
- Non-slip feet.
- 300lb (136kg) weight capacity.



MDS89750R



MDS89755R



Replacement Tips	Item #	Pkg
Replacement Non-Slip Feet Tips	MDS89755T	4/set

Item #	Description	Seat Height	Seat Depth	Seat Width	Product Weight	Pkg.
MDS89750R	Without Back with Arms, White	17"-21½" (43-54.6cm)	16¾" (42.5cm)	22" (55.9cm)	6lbs (2.7kg)	2/cs
MDS89755R	With Back with Arms, White	17"-21½" (43-54.6cm)	16¾" (42.5cm)	22" (55.9cm)	7lbs (3.2kg)	2/cs

Transfer Benches

GUARDIAN
SIGNATURE

Padded Transfer Bench

- Vinyl covered foam backrest and seating surface offer a highly resilient cushioned surface for comfort.
- Lightweight, folding aluminum frame is height adjustable in 1/2" (1.27cm) increments.
- Suction cup footpieces secure the bench to the tub floor.
- Backrest adjusts for right or left facing tubs.
- Side arm and seat handle provide extra stability and leverage when transferring.
- 300lb (136kg) weight capacity.



G98338F

Item #	Description	Color	Seat Height	Seat Depth	Seat Width	Pkg.
G98338F	With Back, Assembled	Beige	17"-22" (43-56cm)	16½" (42cm)	26½" (67cm)	2/cs

GUARDIAN
SIGNATURE

Padded Transfer Bench with Commode Opening

- A commode opening affords the user greater comfort and better access for cleaning.
- Thick, high density foam cushions and backrest for comfortable seating that won't bottom out.
- Suction cup footpieces secure the bench to tub floor.
- Lightweight, folding aluminum frame is height adjustable in 1/2" (1.27cm) increments to accommodate most tubs.
- 300lb (136kg) weight capacity.



G98013A

Item #	Description	Seat Height	Seat Depth	Seat Width	Pkg.
G98013A	With Back, Assembled, Beige	17"-22" (43-56cm)	16½" (42cm)	26½" (67cm)	Each

GUARDIAN

Deluxe-Padded Transfer Bench

- Comfortable, 2" (5cm) cushioned seat and back are covered in stain-resistant, light gray vinyl.
- Wide-base, and suction-cup tips on two legs in tub offer extra stability and safety.
- Direction of seat back can be reversed easily without tools.
- Anodized, aluminum frame.
- Fits tubs 11"-21¼" (28-54cm) deep.
- 300lb (136kg) weight capacity.



MDS86945P

Replacement Parts	Item #	Pkg
Suction-cups w/ Leg Extensions	MDS86952LSC	1 pr
Rubber Tips	MDS89664T	1/pk

Item #	Seat Height	Seat Width x Depth	Overall Height	Overall Width x Depth At Base	Product Weight	Pkg.
MDS86945P	20½"-25½" (52-65cm)	16"x17¼" (41x44cm)	34"-39" (86-99cm)	31"x21½" (79x55cm)	14½lbs (6.6kg)	2/cs



Non-Padded Transfer Bench

- Textured, molded backrest and seating surface offer a non-slip surface with holes for water drainage.
- Lightweight, anodized aluminum frame is height adjustable in 1/2" (1.27cm) increments.
- Backrest attaches for right or left facing tubs easily without tools.
- Side arm and seat back handle provide extra stability and leverage when transferring.
- 300lb (136kg) weight capacity.



G98308A

Item #	Description	Seat Height	Seat Depth	Seat Width	Pkg.
G98308A	With Back, Assembled, Gray	17"-22" (43-56cm)	16¼" (41cm)	26½" (67cm)	2/cs



Non-Padded Transfer Bench

- Suction cup footpieces secure the bench to tub floor.
- Lightweight, anodized aluminum frame is height adjustable in 1/2" (1.27cm) increments to accommodate most tubs.
- Backrest adjusts without tools for right or left facing tubs.
- Built-in soap dish and hand-held shower slot for user convenience.
- 300lb (136kg) weight capacity.



G98309G

Item #	Description	Color	Seat Height	Seat Depth	Seat Width	Pkg.
G98309G	With Back, Some Assembly Required	White	17"-21½" (43-55cm)	21" (53cm)	31½" (80cm)	3/cs



Transfer Bench

- Push-button assembly for back.
- Durable, blow-molded, plastic bench.
- Wide base, and suction-cup tips on regular transfer bench. Two legs in tub offer extra stability and safety.
- Direction of seat back can be reversed easily without tools.
- Anodized, aluminum frame adjusts in 1" (2.5cm) increments.
- Fits tubs 11"-21¼" (28-54cm) deep.



MDS86952XW

Item #	Description	Seat Height	Seat Width x Depth	Overall Height	Overall Width x Depth At Base	Weight Capacity	Product Weight	Pkg.
MDS86952	Transfer Bench	16½"-21½" (42-55cm)	17"x16½" (43x42cm)	34"-39" (86-99cm)	29"x21½" (74x55cm)	300lbs (136kg)	11lbs (5kg)	2/cs
MDS86952XW	Bariatric	16"-20" (41-51cm)	17"x16" (43x41cm)	28"-32" (71-81cm)	28½"x20½" (72x52cm)	600lbs (272kg)	16lbs (7.3kg)	Each

Elevated Toilet Seats/Toilet Safety Rails



G30270A



Locking Raised Toilet Seats

- Raises seat height 5" (13cm) above the toilet, arms 11" (28cm).
- Easy-to-use locking system requires no tools and fits most toilets.
- Ergonomically correct seating surface for extra comfort.
- Available with wrap around foam padded armrests for a sure grip, comfort and support.
- Lightweight polypropylene resin is easy to clean.

Item #	Weight Capacity	Package
G30270A (with padded arms)	300lbs (136kg)	3/cs
G30260 (without arms)	250lbs (113kg)	3/cs



G30300



Toilet Safety Frame

- Handles are adjustable and rotate back to allow a wide range of comfortable and secure positions.
- Closed cell foam armrests provide a comfortable, sure grip.
- Adjustable height 26"-31" (66-74cm) accommodates standard or elevated toilet seats.
- Easy-to-clean aluminum frame mounts securely onto bowl with adjustable bracket 18"-24" (46-61cm).
- Overall Depth: at arm: 16³/₄" (43cm), at leg 12" (30cm)

Item #	Ht. Adjustment	Width Adjustment	Wt. Capacity	Pkg.
G30300	26"-31" (66-74cm)	8"-24" (46-61cm)	250lbs (113kg)	2 sets/cs
G30300-1	26"-31" (66-74cm)	8"-24" (46-61cm)	250lbs (113kg)	2 sets/cs



G30250



Economy Raised Toilet Seat

- Universal design is made of easy-to-clean polypropylene resin.
- Non Slip pads offer additional stability.
- Orthopedically designed to provide user with stability and comfort.
- Shaped to fit most toilets.

Item #	Weight Capacity	Package
G30250	250lbs (113kg)	3/cs



MDS80315



Universal Raised Toilet Seat

- Seat fits securely over any toilet bowl.
- Vinyl-covered brackets allow height adjustments from 4³/₄"-6³/₄" (12-17cm) above standard seat and prevent chipping or scratching of porcelain.
- Plastic seat and plastic splash guard make cleaning easy.
- 250lb (113kg) weight capacity.

Item #	Weight Capacity	Package
MDS80315 (Universal)	250lbs (113kg)	Each



Lock-On Elevated Toilet Seat with Arms

- Adds 5" (12.7cm)
- Arms assist in getting on and off the seat
- Locking device secures the seat to the toilet
- No tools required

Item #	Weight Capacity	Pkg.
MDS80320S	300lbs (136kg)	Each



Bolt-down Elevated Toilet Seat with Arms

- Elevates seat by 4" (10cm)
- 4½" (11.4cm) molded-in arms assist in getting on and off the seat
- Hardware included

Item #	Description	Wt. Capacity	Pkg.
MDS80322S	Standard	300lbs (136kg)	Each
MDS80322E	Elongated	300lbs (136kg)	Each



Hinged Elevated Toilet Seat

- 4" (10cm) elevation
- Easy to clean
- Hinged design can be lifted and lowered like a regular toilet
- Uses existing seat and lid so as to not affect the decor of the room

Item #	Description	Wt. Capacity	Pkg.
MDS80324S	Standard	300lbs (136kg)	Each
MDS80324E	Elongated	300lbs (136kg)	Each



Extra Wide Tall-Ette with Arms and Legs

- Steel frame
- 4" (10cm) elevation
- Width between arms is 18" (46cm)
- Soft foam armrests assist in sitting and rising

Item #	Weight Capacity	Pkg.
MDS80326	600lbs (272kgs)	Each



Arthro Tall-Ette Elevated Toilet Seats

- 4" (10cm) elevation allows for easy sitting and rising

Item #	Description	Wt. Capacity	Pkg.
Right Side Sloped			
MDS80328RS	Standard	300lbs (136kg)	Each
MDS80328RE	Elongated	300lbs (136kg)	Each
Left Side Sloped			
MDS80328LS	Standard	300lbs (136kg)	Each
MDS80328LE	Elongated	300lbs (136kg)	Each



Total Hip Replacement Tall-Ette Elev. Toilet Seat

- 4" (10cm) elevation permits for easy sitting and rising
- Specifically designed to meet the needs of people who have undergone unilateral or bilateral total hip replacement

Item #	Description	Wt. Capacity	Pkg.
MDS80328TS	Standard	300lbs (136kg)	Each
MDS80328TE	Elongated	300lbs (136kg)	Each

Commodes

GUARDIAN
SIGNATURE

Non-Folding 3 in 1 Steel Commode

- 3 in 1 design allows for use as a bedside commode, or raised toilet seat and safety frame over a standard toilet.
- Roomy 19" (48cm) between arms.
- Wide base depth provides more stability.
- Clip-on seat can be easily removed for cleaning.
- Steel frame is rust resistant.
- G30211-4 chrome plated finish offers greater durability.
- Includes pail, lid and splash guard.
- Height adjusts 16½"-22½" (42-57cm).



G30213-4

Item #	Frame Material	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Max. Wt. Capacity	Packaging
G30213-4	Steel	16½"-22½" (42-57cm)	16" (41cm)	19" (48cm)	21¼" (54cm)	350lbs (159kg)	Each
G30211-4	Chromed Steel	16½"-22½" (42-57cm)	16" (41cm)	19" (48cm)	21¼" (54cm)	350lbs (159kg)	Each

GUARDIAN

Non-Folding Steel Commodes

- Back rail removes quickly and easily.
- Features a durable, plastic seat with lid.
- Plastic armrests for additional patient comfort.
- Adjustable legs allow seat to be raised or lowered in 1" (2.5cm) increments.
- Includes 12-quart capacity, commode bucket with cover & splash shield.
- Non-skid, non-marring, rubber tips provide extra safety.



MDS89664

Commode Accessories	Item #	Pkg
Rubber Tips	MDS89664T	1/pk
Seat & Lid for MDS89664	MDS89664SL	4 sets/cs
Set of four splash guards	MDS89664SGC	4/set

Item #	Description	Base Width	Base Depth Overall	Width Between Arms	Seat Width	Seat Depth Overall	Seat Height Min./Max.	Height Overall Min./Max.	Product Weight	Weight Capacity	Pkg.
MDS89664	Steel, 3-in-1	24" (61cm)	18" (46cm)	19" (48cm)	14" (36cm)	15" (38cm)	15"-21" (38-53cm)	23½"-30" (59.7-76cm)	13.2lbs (6kgs)	350lbs (159kgs)	4/cs



GUARDIAN

Non-Folding Aluminum Commode

- High-strength, aluminum frame prevents rusting.
- Durable, plastic seat with lid.
- Plastic armrests for additional, patient comfort.
- Adjustable legs allow the seat to be raised or lowered in 1" (2.5cm) increments to sit from 16"-20" (41-51cm).
- Includes 10-quart capacity, commode bucket with cover and splash shield.
- Non-skid, non-marring, rubber tips provide extra safety.
- 300lb (136kg) weight capacity.



MDS89664A

Item #	Base Width	Overall Base Depth	Width Between Arms	Seat Width	Overall Seat Depth	Seat Height Min./Max.	Overall Height Min./Max.	Pkg.
MDS89664A	25½" (65cm)	16½" (42cm)	18½" (47cm)	14" (35.5cm)	15" (38cm)	16"-20" (41-51cm)	23"-27" (58.5-68.5cm)	4/cs

GUARDIAN

3 in 1 Steel Commode

- 3 in 1 design allows for use as a bedside commode, or raised toilet seat and safety frame over a standard toilet.
- Folding frame design assembles in seconds without tools and reduces storage space by 35% while making delivery easy. Allows user to refold for transport or storage.
- Roomy 19" between arms.
- Wide base depth provides more stability.
- Clip-on seat can be easily removed for cleaning.
- Steel frame is rust resistant.
- Includes pail, lid and splash guard.
- Height adjusts 16½"-22½".



G30213-4F

Item #	Frame Material	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	max. wt. Capacity	Packaging
G30213-4F	Steel, Folding	16½"-22½" (42-57cm)	16" (41cm)	19" (48cm)	21¼" (54cm)	350lbs (159kg)	Each

Commodes

GUARDIAN

Folding, Elongated Seat

- Features an elongated seat for added maneuverability and comfort.
- 3 in 1 design allows for use as a bedside commode, or raised toilet seat and safety frame over a standard toilet.
- Roomy 19" (48cm) between arms.
- Clip-on seat can be easily removed for cleaning.
- Steel frame is rust resistant.
- Includes pail, lid and splash guard.
- 400lb (181kg) weight capacity.



MDS89664EL

Item #	Base Width	Overall Base Depth	Width Between Arms	Seat Width	Overall Seat Depth	Seat Height Min./Max.	Overall Height Min./Max.	Pkg.
MDS89664EL	22" (56cm)	19" (48cm)	19" (48cm)	14" (35.5cm)	16½" (42cm)	17"-23" (43-58cm)	26"-32" (66-81cm)	4/cs

GUARDIAN SIGNATURE

3 in 1 Aluminum Commode

- Aluminum frame is lightweight and rust-resistant.
- Folding frame design assembles in seconds without tools and reduces storage space by 35% while making delivery easy. Allows user to refold for transport or storage.
- 3 in 1 design allows for use as a bedside commode, or raised toilet seat and safety frame over a standard toilet.
- Full 19" (48cm) between arms.
- Wide base depth provides more stability.
- Clip-on seat can be easily removed for cleaning.
- Includes pail, lid and splash guard
- Height adjusts 17½"-21½" (44-55cm).



G30210-4F

Item #	Frame Material	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Max. Wt. Capacity	Packaging
G30210-4F	Aluminum	17½"-21½" (44-54cm)	16" (41cm)	19" (48cm)	21¼" (54cm)	250lbs (113kg)	4/cs

GUARDIAN SIGNATURE

Non-Padded Drop-Arm Commode

- Designed to accommodate lateral transfers.
- Hard, chrome plated steel frame is rust resistant and offers greater durability.
- Foam handgrips provide a more secure grip.
- Backrest can be removed for use as an over-the-toilet safety frame.
- Pail, lid and splash guard included.
- Height adjusts 18¼"-23¼" (46-59cm).



G98202

Item #	Frame Material	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Max. Wt. Capacity	Packaging
G98202	Chromed Steel	18¼"-23¼" (46-59cm)	18"-20" (46-51cm)	18½" (47cm)	22" (56cm)	300lbs (136kg)	Each



GUARDIAN
SIGNATURE

Padded Drop-Arm Commode

- Features thick, resilient foam padded seat and back for increased comfort.
- Designed to accommodate lateral transfers.
- Hard chrome plated steel frame is rust resistant and offers greater durability.
- Foam handgrips provide a more secure grip.
- Height adjust 18¼"-23¼" (46-58.5cm).
- Backrest can be removed for use as an over-the-toilet safety frame.



G98204

Item #	Frame Material	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Max. Wt. Capacity	Packaging
G98204	Chromed Steel	20¼"-25¼" (51-64cm)	18"-20" (46-51cm)	18lbs (47cm)	22" (56cm)	350lbs (159kg)	Each

GUARDIAN

Drop-Arm Steel Commodes

- Single-handed, arm release.
- 12-quart pail w/cover included.
- Chrome plated steel for added durability.

Replacement Accessories	Item #	Pkg
Rubber Tips	MDS89664T	4/cs
Commode Bucket	MDS89664SL	4/cs



MDS89668

Item #	Description	Base Width	Base Depth Overall	Width Between Arms	Seat Width	Seat Depth Overall	Seat Height Min./Max.	Height Overall Min./Max.	Weight Capacity	Pkg.
MDS89668	Drop-Arm Steel	20¼" (51.4cm)	20½" (52cm)	17½" (44cm)	14" (35.5cm)	15" (38cm)	17"-21" (43-53cm)	28¾"-33¼" (73-84cm)	250lbs (113kg)	Each

Commodes



GUARDIAN SIGNATURE

Heavy Duty Commode

- Extra wide frame has a full 21.12" between arms.
- Large base depth provides more stability.
- Clip-on seat can be easily removed for cleaning.
- Height adjusts 18"-22" (46-56cm)



G30214-2

Item #	Frame Material	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Max. Wt. Capacity	Packaging
G30214-2	Painted Steel	17½"-21½" (46-56cm)	16" (41cm)	21 1/8" (54cm)	23¾" (60cm)	400lbs (182kg)	Each

GUARDIAN

Bariatric Deluxe, Steel Commodes

- Back rail removes quickly and easily.
- Features a durable, plastic seat with lid.
- Plastic armrests for additional patient comfort.
- Adjustable legs allow seat to be raised or lowered in 1" (2.5cm) increments.
- Includes 12-quart capacity, commode bucket with cover & splash shield.
- Non-skid, non-marring, rubber tips provide extra safety.

Commode Accessories	Item #	Pkg
Rubber Tips	MDS89664T	1/pk
Seat & Lid for MDS89664XW	MDS89664SLXW	1/set
Set of four splash guards	MDS89664SGC	4/set



MDS89668XW



MDS89664XW

Item #	Description	Base Width	Base Depth Overall	Width Between Arms	Seat Width	Seat Depth Overall	Seat Height Min./Max.	Height Overall Min./Max.	Weight Capacity	Pkg.
MDS89664XW	Bariatric	28" (71cm)	20½" (52cm)	24" (61cm)	14" (36cm)	15" (38cm)	18¼"-23½" (46-59.7cm)	29½"-32½" (75-82.6cm)	650lbs (295kgs)	Each
MDS89668XW	Bariatric Platform	29" (74cm)	19½" (49.5cm)	28" (71cm)	28" (71cm)	19" (48cm)	19¼" (49cm Only)	25½" (64cm Only)	850lbs (386kg)	Each



Bariatric Commode

- Extra wide frame with 23½" between arms to accommodate larger users.
- Reinforced steel frame offers greater durability.
- Height adjusts 16½"-23½".
- 2 per case.

G30216B



Item #	Frame Material	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Max. Wt. Capacity	Packaging
G30216B	Painted Steel	16½"-23½" (42-62cm)	18½" (47cm)	23½" (60cm)	27" (69cm)	500lbs (227kg)	2/cs



Replacement Commode Bucket

- Constructed of high-impact polyethylene.
- Fits all non-platform Medline commodes.
- Complete with handle and lid.
- Diameter of bucket is 11" (28cm).
- 12-quart/2-gallon capacity.

MDS80306B



Item #	Description	Pkg.
MDS80306B	Commode Bucket	6/cs




Commode Liners

- Make commode clean-up easy and hygienic.
- Gelling agent solidifies the waste for convenient disposal
- Great add-on sale to any commode



Item #	Description	Pkg.
MDS89664LINER	Commode Liners	6 Boxes/cs



MDS89664LINER

Medline Makes Mobility More Mobile



Medline Industries, Inc.

One Medline Place
Mundelein, IL 60060

Medline United States
1-800-MEDLINE (633-5463)

www.medline.com | info@medline.com

Medline Canada
1-800-396-6996

www.medline.ca | canada@medline.com

Medline México
01-800-831-0898

www.medlinemexico.com | mexico@medline.com

Some products may not be available for sale in Mexico or Canada. Please contact your Account Representative for additional information.

We reserve the right to correct any errors that may occur within this brochure. ©2010 Medline Industries, Inc.

Medline, Guardian and Guardian Signature are registered trademarks of Medline Industries, Inc.

MKT210072/LIT124R/10M/DP6