


CURAD[®]

WE HELP HEAL

Detectable Bandages For the Food Industry

Built-in safety features:


- Bright blue color for high visibility
- Hidden metal marker for metal detection



CURAD® Detectable Adhesive Bandages



Detectable Metal, Highly Visible Color



Curad blue bandages meet the needs of the food service industry with their high visibility and innovative, long-lasting adhesive that's gentle on the skin.

Detectable metal hidden behind the non-adherent wound pad easily triggers metal detectors used by food packagers. This routine process scans for dangerous metal debris such as screws or knife tips that can accidentally fall into food during processing. In addition, the highly identifiable blue color stands out in all different foods.

Choose from two types of bandages:



Revolutionary woven Flex-Fabric™
Moves with the skin for greater comfort and durability.



Standard plastic backing is effective and economical.

➔ All of our bandages seal around all 4 sides to protect from germs and dirt.

Item #	Description	Dimensions	Qty
Flex-Fabric			
NON25650BL	Woven blue detectable bandage	¾" x 3" (19 x 76 mm)	100pcs/bx, 12bx/cs
NON25660BL	Woven blue detectable bandage	1" x 3" (25 x 76 mm)	100pcs/bx, 12bx/cs
NON25513BL	Fingertip, Woven blue detectable bandage	1¾" x 2" (44 x 50 mm)	100pcs/bx, 12bx/cs
NON25510BL	Knuckle, Woven blue detectable bandage	1½" x 3" (38 x 76 mm)	100pcs/bx, 12bx/cs
Plastic			
NON25500BL	Plastic blue detectable bandage	¾" x 3" (19 x 76 mm)	100pcs/bx, 12bx/cs
NON25600BL	Plastic blue detectable bandage	1" x 3" (25 x 76 mm)	100pcs/bx, 12bx/cs



Medline Industries, Inc.

One Medline Place
Mundelein, IL 60060

Medline United States
1-800-MEDLINE (633-5463)

www.medline.com | info@medline.com

Medline Canada
1-800-396-6996

www.medline.ca | canada@medline.com

Medline México
01-800-831-0898

www.medlinemexico.com | mexico@medline.com

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.
©2010 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. Curad is a registered trademark of Medline Industries, Inc.