

Everest & Jennings®
DURA-gel™
SP III

The Dura-Gel™ SP III (Skin Protection) wheelchair cushion reduces surface pressure through the use of non-aqueous Dura-Gel™ placed in three pressure zones. The cushion uses a high-density base foam in combination with a visco-elastic topper for added comfort and also positions the user with a moderate slope from front to back.



- ❑ For use with moderate risk patients
- ❑ Reduces tissue interface pressures at critical points
- ❑ Features a specialized Dura-Gel™ center layer surrounded by comfortable foam layers on the top and bottom of the cushion
- ❑ The soft, stretchable lightweight cover reduces shear/friction and is fluid resistant and has a strap at the back to secure the cushion firmly to the chair
- ❑ Accommodates a wide range of body shapes, sizes and weights
- ❑ 350 lb maximum weight capacity
- ❑ Lightweight
- ❑ 18 month limited warranty



Latex Free

Designed to meet the requirements of HCPCS Code: E2603

- 8930166** 16" x 16" x 3"
- 8930186** 18" x 16" x 3"
- 8930188** 18" x 18" x 3"
- 8930206** 20" x 16" x 3"
- 8930208** 20" x 18" x 3"

Designed to meet the requirements of HCPCS Code: E2604

- 8930228** 22" x 18" x 3"
- 8930220** 22" x 20" x 3"
- 8930248** 24" x 18" x 3"

