


BASELINE[®]

EVALUATION INSTRUMENTS

12-0290 / 12-0291

12-0292 / 12-0293

PNEUMATIC BULB DYNAMOMETER



12-0291

Squeeze bulb dynamometer measures PSI.

12-0290: 0-30 PSI, without max pointer, 0.5 PSI gradations

12-0291: 0-30 PSI, with max pointer, 0.5 PSI gradations

12-0292: 0-15 PSI, without max pointer, 0.5 PSI gradations

12-0293: 0-15 PSI, with max pointer, 0.5 PSI gradations

For units with max pointer, max pointer is re-set to zero by turning center knob.

CAUTION: this product contains natural rubber latex, which may cause allergic reactions.



**Manufacturer and Master Distributor of
Physical Therapy and Rehabilitation Products**

Baseline[®] is a registered trademark of Goldberg.

Fabrication Enterprises, Inc.

PO Box 1500, White Plains, NY 10602 USA

tel: 800-431-2830 914-345-9300

fax: 800-634-5370 914-345-9800

info@FabricationEnterprises.com

Authorized CE Representative:
RMS UK Ltd.
28 Trinity Road
Nailsea, Somerset BS48 4NU
United Kingdom

