

BASELINE[®]

EVALUATION INSTRUMENTS

-30° to 180° plastic finger goniometer 12-1012

- permits observation of joint's axis of motion, and ROM measurements of MCPs, PIPs and DIPs
- measures a range of 30° hyperextension through 120° flexion
- primary scale reads in 1° increments
- goniometer arm has a linear scale in inches and centimeters
- goniometer arm remains at measurement until reset

other sizes and style plastic goniometers available

12-inch 360°	12-1000*
8-inch 360°	12-1001*
6-inch 360°	12-1002*
6-inch 180° pocket.....	12-1005*
6-inch 360° rulongmeter....	12-1006*
3-inch -30° to 180° finger...	12-1012*

*available standard or with white
HiRes[™] strip behind numbers
(add HR to end of product code
for HiRes[™] model)



Fabrication Enterprises Inc.
P.O. Box 1500
White Plains, NY 10602 (USA)
tel: 1-914-345-9300
fax: 1-914-345-9800
www.Fab-Ent.com



Authorized CE representative:
RMS UK Ltd.
28 Trinity Road
Nailsea, Somerset BS48 4NU (UK)

ver 10/11