


BASILINE®

EVALUATION INSTRUMENTS

180° digit goniometer



Readings can be taken by the therapist using only one hand; this frees up the therapist's other hand to stabilize the patient during measurement or to record the measurement.

Measures 110° flexion through 40° hyper-extension in 2° increments.

Arm has a 10cm/4" scale.

Made of sturdy ABS plastic with white markings on black background.

Authorized CE representative:
RMS UK Ltd.

28 Trinity Road, Nailsea, Somerset BS48 4NU
United Kingdom



Fabrication Enterprises Inc.
P.O. Box 1500
Little Plains, NY 10602 (USA)
tel: 1-914-345-9300
fax: 1-914-345-9800
www.Fab-Ent.com

Baseline® is a trademark of Goldberg
©2011 FEI, all rights reserved

ver 3/11