



Contents

Stock kit # _____

- 1.52cm² flat circular probe
- 2 pound (.90 kg) gauge #12-1440
- 5 pound gauge (2.27 kg) #12-1441
- 25 pound gauge (11.34 kg) #12-1442
- 60 pound gauge (27.22 kg) #12-1443

Baseline[®] Dolorimeter

Evaluate pain threshold. Standardized (1.52cm²) flat circular probe is pushed against subject's skin until pain threshold is reached. Gauge reads in pounds and kilograms.

Distributed By:

Manufactured By:

Fabrication Enterprises, Inc.
PO Box 1500
White Plains, NY 10602

tel: 800-431-2830 914-345-9300
fax: 800-634-5370 914-345-9800
info@FabricationEnterprises.com