




SQUID Active Cold Compression Therapy

The Squid active cold compression therapy system provides intermittent sequential cold compression by combining an air compression wrap and cold gel pack.

Patients can use Squid to improve circulation and reduce pain and swelling on affected areas. This system is an effective alternative to traditional cryotherapy tools in a clinical setting.



shown with Squid pump, ankle wrap and gel pack

Key Features:

- Use for temporary relief of aches and pains
- Portable, air pressure pump is small and lightweight (under 0.45 kg)
- 4 preprogrammed levels of compression from low (30 mmHg) to high (85 mmHg)
- Intermittent sequential compression
- Rechargeable battery
- Compression wraps are ergonomically designed to fit specific body areas
- Soft gel packs provide maximum cold and comfort without the mess of ice water
- Wraps and gel packs work with both Squid Go and Squid clinical units



Squid Go home unit

- Rechargeable battery allows for **15 uses** on one charge
- Squid Go set includes pump, wrap, and gel pack



Upper Extremity Wrap

Upper Extremity Wrap	Set
wrist	13-2480
elbow	13-2473
shoulder, left, regular*	13-2474
shoulder, left, large*	13-2475
shoulder, right, regular*	13-2476
shoulder, right, large*	13-2477

*only offers intermittent compression

Mid Body Wrap

back	13-2472
------	---------

Lower Extremity Wrap

leg, regular	13-2478
leg, large	13-2479
ankle, regular	13-2470
ankle, large	13-2471

Squid clinical unit

- Rechargeable battery allows for **50 uses** on one charge
- Squid clinical set includes pump, wrap and gel pack



Upper Extremity Warp

Upper Extremity Warp	Set
wrist	13-2465
elbow	13-2425
shoulder, left, regular*	13-2450
shoulder, left, large*	13-2457
shoulder, right, regular*	13-2453
shoulder, right, large*	13-2460

*only offers intermittent compression

Mid Body Wrap

back	13-2420
------	---------

Lower Extremity Wrap

knee	13-2430
leg, regular	13-2435
leg, large	13-2441
leg circumferential, regular	13-2439
leg circumferential, large	13-2445
ankle, regular	13-2410
ankle, large	13-2414

Wrap and Gel Pack

13-2468
13-2428
13-2452
13-2459
13-2455
13-2462
13-2422
13-2433
13-2438
13-2442
13-2447
13-2448
13-2413
13-2417

Carry Bag



- Use to store and transport Squid Go and clinical pump and wraps

13-2469 carry bag

Squid Sizing Chart

	Chest Circumference	Men's US Shoe Size	Women's US Shoe Size
Shoulder Regular	81 cm - 114 cm	7-13	6-14
Shoulder Large	114 cm - 144 cm	13-20	14-20
	Wrap Length	Knee Circumference	Circumference on Top of Wrap
Leg Regular	35 cm	12 cm - 45 cm	12 cm - 50 cm
Leg Large	42 cm	45 cm - 58 cm	50 cm - 66 cm
		Circumference on Bottom of Wrap	
		12 cm - 40 cm	
		40 cm - 53 cm	