

CANDo® **Massage Gun**

The CanDo® Massage Gun is a cordless state-of-the-art massage device that helps to reduce muscle soreness and tension, promote blood circulation, optimize mobility, and improve the general condition of the body's musculoskeletal system.

Key Features

- Powerfully delivers customized, handheld percussion massages
- Use to relieve muscle cramps, trigger points, knots and sore muscles
- Effectively promotes better circulation, optimizes mobility and reduces stress and tension
- Easily swap out massage heads to target different muscle groups
- Battery provides 6-8 hours of use per charge
- Carrying case and six interchangeable massage heads (bullet, fork, wedge, flat, air and round) are included
- Handheld device weighs less than 5 lbs

Three operational modes:

Fitness Mode

Delivers constant pressure and speed at a preset speed level ranging from 1,200 RPMs to 3,200 RPMs. Ideal for fitness enthusiasts and athletes.

Medical Mode


Designed for point to point local treatment. Speed and pressure will increase from 1,000 to 1,900 RPMs then decrease back to 1,000 RPMs within 6 seconds. This pattern will continue to cycle automatically.

Health Care Mode

Designed for those looking to massage away daily stress and tension. Offers a lower speed and torque than other settings. Preset speed levels range from 1,000 RPMs to 2,000 RPMs.



Six interchangeable massage heads included



14-1770 CanDo® massage gun 350.00



FABRICATION ENTERPRISES INC
PO Box 1500, White Plains, NY 10602 (USA)
tel: 800-431-2830 / +1-914-345-9300
fax: 800-634-5370 / +1-914-345-9800
FabEnt.com sales@fab-ent.com