

CANDo® Treatment Furniture

CanDo® Platform Mounted Parallel Bars

- parallel bars are mounted to a durable 1.5" hardwood base
- height is easily adjustable and numbered with height indicator strips
- width adjustable
- 400 lb weight capacity
- intended for physical therapy use only
- width adjusts from 18" to 28"; height adjusts 26" to 39"



15-4239	7 foot length
15-4240	10 foot length
15-4241	12 foot length

CanDo® Floor Mounted Parallel Bars

- parallel bars are designed to mount directly to the floor
- height is easily adjustable and numbered with height indicator strips
- width adjustable
- 400 lb weight capacity
- intended for physical therapy use only
- height adjusts from 26" to 39"



15-4242	7 foot length
15-4243	10 foot length
15-4244	12 foot length

CanDo® Staircases

- steel fasteners and metal braces provide support and stability
- molded, solid hardwood handrails offer patient comfort
- safety treads on each step
- two-sided staircase patients can ascend and descend without turning
- one-sided staircase saves space
- 500 lb weight capacity



two-sided

one-sided

15-4230	two-sided staircase	54"L x 34"W x 48"H
15-4238	one-sided staircase	56"L x 34"W x 60"H

CanDo® Adjustable Balance Platform

- use for balance and stability training support
- adjustable handrails are mounted to a solid platform
- platform features rugged indoor / outdoor carpet
- 36"L x 36"W x 26-39"H



15-4254	adjustable balance platform
---------	-----------------------------

CanDo® Weight and Mirror Rack

- holds 6 band or tubing rolls up to 50 yds (hollow core center only)
- dumbbell holding rods are angled for extra stability; holds up to 22 dumbbells
- white pegboard can hold up to 48 cuff weights
- mirror features ANSI safety backing
- easily move rack across surfaces with 4 swivel wheels (2 are able to lock)
- dimensions: 31.5"L x 18"W x 72"H
- accessories include:
 - 16 cuff weights (2 ea: 1, 1.5, 2, 2.5, 3, 4 & 5 lbs; 1 ea: 7.5, 10 lbs)
 - 20 dumbbells (2 ea: 1 - 10 lbs)
 - 5 - 50 yard CanDo® Exercise Band Rolls (yellow - black)

15-4255	with accessories
15-4234	without accessories (Rack Only)



upholstery color options

black

slate blue

royal blue

warm gray

desert tan

finish options

natural

dark cherry

CanDo® Treatment Table with Adjustable Back and Shelf

- Features hardwood legs, laminate shelf and adjustable backrest
- Scratch resistant topcoat finish
- Seamless round corner top
- 2" firm foam padding
- 400 lb weight capacity
- Choose upholstery color and finish from table above



15-4231 72"L x 27"W x 31"H
 15-4232 72"L x 30"W x 31"H
 15-4233 78"L x 30"W x 31"H

CanDo® Treatment Table with Adjustable Back

- Features hardwood legs and adjustable backrest
- Scratch resistant topcoat finish
- Seamless round corner top
- 2" firm foam padding
- Heavyweight, knit-backed upholstery
- 400 lb weight capacity
- Choose upholstery color and finish from table above



15-4245 72"L x 27"W x 31"H
 15-4246 72"L x 30"W x 31"H

CanDo® Treatment Table with Flat Top and Shelf

- Features solid hardwood legs and full laminated shelf
- H-brace construction
- Seamless round corner top
- 2" firm foam padding
- 400 lb weight capacity
- Choose upholstery color and finish from table above



15-4249 72"L x 27"W x 31"H
 15-4250 72"L x 30"W x 31"H

CanDo® H-Brace Treatment Table with Flat Top

- Features solid hardwood legs
- H-brace construction
- Seamless round corner top
- 2" firm foam padding
- 400 lb weight capacity
- Choose upholstery color and finish from table above



15-4247 72"L x 27"W x 31"H
 15-4248 72"L x 30"W x 31"H

CanDo® Upholstered Mat Table



- Features solid hardwood legs and plywood frame
- Heavy duty, upholstered mat platform
- 6 leg design with center support adds stability
- 2" ultra-firm, foam padding
- Gray edge bumper guard helps prevent tears
- 700 lb weight capacity
- Choose upholstery color from table above

15-4235 84"L x 48"W x 18"H
 15-4236 84"L x 60"W x 18"H
 15-4237 96"L x 72"W x 18"H

CanDo® Hi-Lo Mat Platform Table with Upholstered Top

- Electric powered foot-controlled height adjusts from 22" to 34"
- Ultra-strong uniframe base
- Powder-coated steel frame
- Custom bumper strip protects top's edges from tears
- 2" foam padding
- 600 lb weight capacity
- Choose upholstery color from table above



15-4251 84"L x 48"W
 15-4252 84"L x 60"W
 15-4253 96"L x 72"W

CanDo and CANDo are registered trademarks of Goldberg