


Orfit® Classic Sheets

- Non-coated material sticks easily to itself and to outrigger attachments
- Temporarily bonds during fabrication; permanent bonding possible with dry heat
- Turns transparent when heated for visibility
- Unlimited potential for fabrication of very large to very small orthoses
- Available in soft (high stretch) and stiff (low stretch)
- Soft version offers increased, accurate molding around anatomic corners; stiff version offers increases support for larger orthoses



| perf | stretch | sheet | case qty |
|-------|---------|-------------------|--------------------------|
| solid | soft | 1/16" x 18" x 24" | 24-5620-1 24-5620-4 4 ea |
| solid | soft | 1/12" x 18" x 24" | 24-5622-1 24-5622-4 4 ea |
| solid | soft | 1/8" x 18" x 24" | 24-5628-1 24-5628-4 4 ea |
| solid | soft | 3/16" x 24" x 36" | 24-5634-1 24-5634-4 4 ea |
| solid | stiff | 1/8" x 18" x 24" | 24-5636-1 24-5636-4 4 ea |
| solid | stiff | 1/8" x 24" x 36" | 24-5638-1 24-5638-4 4 ea |
| solid | stiff | 3/16" x 24" x 36" | 24-5640-1 24-5640-4 4 ea |
| 25% | soft | 1/12" x 18" x 24" | 24-5625-1 24-5625-4 4 ea |
| 25% | soft | 1/8" x 18" x 24" | 24-5629-1 24-5629-4 4 ea |
| 3.5% | soft | 1/8" x 18" x 24" | 24-5630-1 24-5630-4 4 ea |
| 3.5% | soft | 3/16" x 24" x 36" | 24-5635-1 24-5635-4 4 ea |
| 3.5% | stiff | 1/8" x 24" x 36" | 24-5639-1 24-5639-4 4 ea |
| 13% | soft | 1/16" x 18" x 24" | 24-5621-1 24-5621-4 4 ea |
| 13% | soft | 1/12" x 18" x 24" | 24-5623-1 24-5623-4 4 ea |

Orfit® Classic Precuts

- Make splint creation quick and easy
- Precut pieces are anatomically shaped
- Easily remoldable
- Available in 3 sizes: small, medium and large



| perf | thick | gauntlet thumb post | 1 piece | case qty |
|------|-------|-------------------------------------|------------|----------------|
| | | 13% 1/16" small | 24-5650-1 | 24-5650-2 2 ea |
| | | 13% 1/16" medium | 24-5651-1 | 24-5651-2 2 ea |
| | | 13% 1/16" large | 24-5652-1 | 24-5652-2 2 ea |
| | | complete thumb | | |
| | | 13% 1/12" small | 24-5653-1 | 24-5653-2 2 ea |
| | | 13% 1/12" medium | 24-5654-1 | 24-5654-2 2 ea |
| | | 13% 1/12" large | 24-5655-1 | 24-5655-2 2 ea |
| | | wrist + thumb | | |
| | | 13% 1/12" small | 24-5656-1 | 24-5656-2 2 ea |
| | | 13% 1/12" medium | 24-5657-1 | 24-5657-2 2 ea |
| | | 13% 1/12" large | 24-5658-1 | 24-5658-2 2 ea |
| | | gauntlet immobilization | | |
| | | 13% 1/12" small | 24-5659-1 | 24-5659-2 2 ea |
| | | 13% 1/12" medium | 24-5660-1 | 24-5660-2 2 ea |
| | | 13% 1/12" large | 24-5661-1 | 24-5661-2 2 ea |
| | | intrinsic resting hand | | |
| | | solid 1/8" small | 24-5666-1 | 24-5666-2 2 ea |
| | | solid 1/8" medium | 24-5667-1 | 24-5667-2 2 ea |
| | | solid 1/8" large | 24-5668-1 | 24-5668-2 2 ea |
| | | anti-spastic + thumb | | |
| | | solid 1/8" small | 24-5669-1 | 24-5669-2 2 ea |
| | | solid 1/8" medium | 24-5670-1 | 24-5670-2 2 ea |
| | | solid 1/8" large | 24-5671-1 | 24-5671-2 2 ea |
| | | intrinsic anti-spastic hand | | |
| | | 3.5% 1/8" small | 24-5672-1 | 24-5672-2 2 ea |
| | | 3.5% 1/8" medium | 24-5673-1 | 24-5673-2 2 ea |
| | | 3.5% 1/8" large | 24-5674-1 | 24-5674-2 2 ea |
| | | rheumatoid arthritis resting | | |
| | | 3.5% 1/8" small | 24-5675-1 | 24-5675-2 2 ea |
| | | 3.5% 1/8" medium | 24-5676-1 | 24-5676-2 2 ea |
| | | 3.5% 1/8" large | 24-5677-1 | 24-5677-2 2 ea |
| | | ulnar drift | | |
| | | solid 1/8" left | 24-5678L-1 | |
| | | solid 1/8" right | 24-5678R-1 | |

Orfit® NS Soft Sheets

- Coated material is non-stick, designed to prevent accidental bonding
- Allows for temporary bonding during splint creation
- Permanent bonding possible when coating removed and heat is applied
- Soft version is high stretch and allows for increased accurate molding around anatomic corners

| perf | sheet | case qty |
|-------|-------------------|--------------------------|
| solid | 1/16" x 18" x 24" | 24-5680-1 24-5680-4 4 ea |
| solid | 1/12" x 18" x 24" | 24-5682-1 24-5682-4 4 ea |
| solid | 1/8" x 18" x 24" | 24-5688-1 24-5688-4 4 ea |
| 25% | 1/12" x 18" x 24" | 24-5685-1 24-5685-4 4 ea |
| 25% | 1/8" x 18" x 24" | 24-5689-1 24-5689-4 4 ea |
| 3.5% | 1/12" x 18" x 24" | 24-5686-1 24-5686-4 4 ea |
| 3.5% | 1/8" x 18" x 24" | 24-5690-1 24-5690-4 4 ea |
| 13% | 1/16" x 18" x 24" | 24-5681-1 24-5681-4 4 ea |
| 13% | 1/12" x 18" x 24" | 24-5683-1 24-5683-4 4 ea |
| 13% | 3/32" x 18" x 24" | 24-5687-1 24-5687-4 4 ea |



Orfit® Colors NS Sheets

- Non-stick coating allows for temporary bonding during application
- Choice of unique colors
- Designed to appeal to patients of all ages
- Material has a high gloss and a shiny finish
- Available precuts make splint creation quick and easy



| perf | sheet | case qty |
|--|-----------|----------------|
| solid sonic silver 1/12" x 18" x 24" | 24-5770-1 | 24-5770-4 4 ea |
| solid sonic silver 1/8" x 18" x 24" | 24-5772-1 | 24-5772-4 4 ea |
| solid atomic blue 1/12" x 18" x 24" | 24-5774-1 | 24-5774-4 4 ea |
| solid atomic blue 1/8" x 18" x 24" | 24-5776-1 | 24-5776-4 4 ea |
| solid gold 1/12" x 18" x 24" | 24-5778-1 | 24-5778-4 4 ea |
| solid gold 1/8" x 18" x 24" | 24-5780-1 | 24-5780-4 4 ea |
| solid violet 1/12" x 18" x 24" | 24-5782-1 | 24-5782-4 4 ea |
| solid hot green 1/12" x 18" x 24" | 24-5784-1 | 24-5784-4 4 ea |
| solid dominant black 1/12" x 18" x 24" | 24-5786-1 | 24-5786-4 4 ea |
| solid dynamic red 1/12" x 18" x 24" | 24-5788-1 | 24-5788-4 4 ea |
| solid ocean blue 1/12" x 18" x 24" | 24-5790-1 | 24-5790-4 4 ea |
| 3.5% sonic silver 1/8" x 18" x 24" | 24-5773-1 | 24-5773-4 4 ea |
| 3.5% atomic blue 1/8" x 18" x 24" | 24-5777-1 | 24-5777-4 4 ea |
| 3.5% gold 1/8" x 18" x 24" | 24-5781-1 | 24-5781-4 4 ea |
| 13% sonic silver 1/12" x 18" x 24" | 24-5771-1 | 24-5771-4 4 ea |
| 13% atomic blue 1/12" x 18" x 24" | 24-5775-1 | 24-5775-4 4 ea |
| 13% gold 1/12" x 18" x 24" | 24-5779-1 | 24-5779-4 4 ea |
| 13% violet 1/12" x 18" x 24" | 24-5783-1 | 24-5783-4 4 ea |
| 13% hot green 1/12" x 18" x 24" | 24-5785-1 | 24-5785-4 4 ea |
| 13% dominant black 1/12" x 18" x 24" | 24-5787-1 | 24-5787-4 4 ea |
| 13% dynamic red 1/12" x 18" x 24" | 24-5789-1 | 24-5789-4 4 ea |
| 13% ocean blue 1/12" x 18" x 24" | 24-5791-1 | 24-5791-4 4 ea |

Orfit® Colors NS Precuts

perf thick

Precut: wrist + thumb



wrist + thumb

| | small | medium | large |
|--------------------------|-----------|-----------|-----------|
| 13% 1/12" atomic blue | 24-5818-1 | 24-5819-1 | 24-5820-1 |
| 13% 1/12" ocean blue | 24-5806-1 | 24-5807-1 | 24-5808-1 |
| 13% 1/12" sonic silver | 24-5830-1 | 24-5831-1 | 24-5832-1 |
| 13% 1/12" dominant black | 24-5842-1 | 24-5843-1 | 24-5844-1 |

Precut: gauntlet immobilization



gauntlet immobilization

| | | | |
|--------------------------|-----------|-----------|-----------|
| 13% 1/12" atomic blue | 24-5821-1 | 24-5822-1 | 24-5823-1 |
| 13% 1/12" ocean blue | 24-5809-1 | 24-5810-1 | 24-5811-1 |
| 13% 1/12" sonic silver | 24-5833-1 | 24-5834-1 | 24-5835-1 |
| 13% 1/12" dominant black | 24-5845-1 | 24-5846-1 | 24-5847-1 |

Orflight® Sheets

- Non-coated material is ideal for creating lightweight, supportive splints
- Easy to trim edges and soft foamy feel for patient comfort



| perf | sheet | case qty |
|------------------------|-----------|----------------|
| solid 1/8" x 18" x 24" | 24-5722-1 | 24-5722-4 4 ea |
| 3.5% 1/8" x 18" x 24" | 24-5723-1 | 24-5723-4 4 ea |
| 13% 1/16" x 18" x 24" | 24-5720-1 | 24-5720-4 4 ea |
| 13% 3/32" x 18" x 24" | 24-5721-1 | 24-5721-4 4 ea |

Orflight® Precuts



| perf | thick | 1 piece | case qty |
|---|-------------|-----------|----------------|
| Precut: rheumatoid arthritis resting | | | |
| 3.5% | 1/8" small | 24-5940-1 | 24-5940-2 2 ea |
| 3.5% | 1/8" medium | 24-5941-1 | 24-5941-2 2 ea |
| 3.5% | 1/8" large | 24-5942-1 | 24-5942-2 2 ea |

Orflight® NS Sheets

- Non-stick coated material is ideal for creating lightweight, supportive splints
- Easy to trim edges and soft foamy feel for patient comfort
- Available in black and atomic blue colors



| perf | sheet | case qty |
|------------------------------------|-----------|----------------|
| solid black 1/8" x 18" x 24" | 24-5742-1 | 24-5742-4 4 ea |
| 3.5% black 1/8" x 18" x 24" | 24-5743-1 | 24-5743-4 4 ea |
| 13% black 1/16" x 18" x 24" | 24-5740-1 | 24-5740-4 4 ea |
| 13% black 3/32" x 18" x 24" | 24-5741-1 | 24-5741-4 4 ea |
| solid atomic blue 1/8" x 18" x 24" | 24-5762-1 | 24-5762-4 4 ea |
| 3.5% atomic blue 1/8" x 18" x 24" | 24-5763-1 | 24-5763-4 4 ea |
| 13% atomic blue 1/16" x 18" x 24" | 24-5760-1 | 24-5760-4 4 ea |
| 13% atomic blue 3/32" x 18" x 24" | 24-5761-1 | 24-5761-4 4 ea |

Orfit® Eco Sheets

- Non-coated, off-white material is ideal for large and rigid splints
- Made of high quality recycled material
- Appealing smooth finish resists fingerprints



| perf | sheet | case qty | |
|-------------------------|-----------|-----------|------|
| solid 3/32" x 18" x 24" | 24-5860-1 | 24-5860-4 | 4 ea |
| solid 1/8" x 18" x 24" | 24-5862-1 | 24-5862-4 | 4 ea |
| solid 1/8" x 24" x 36" | 24-5864-1 | 24-5864-4 | 4 ea |
| 3.5% 3/32" x 18" x 24" | 24-5861-1 | 24-5861-4 | 4 ea |
| 3.5% 1/8" x 18" x 24" | 24-5863-1 | 24-5863-4 | 4 ea |
| 3.5% 1/8" x 24" x 36" | 24-5865-1 | 24-5865-4 | 4 ea |

Orfit® Eco Precuts

- make splint creation quick and easy



| perf | 1 piece | case qty | |
|-----------------------|-----------|-----------|------|
| Dorsal Cock-Up | | | |
| solid 1/8" x-small | 24-5662-1 | 24-5662-2 | 2 ea |
| solid 1/8" small | 24-5663-1 | 24-5663-2 | 2 ea |
| solid 1/8" medium | 24-5664-1 | 24-5664-2 | 2 ea |
| solid 1/8" large | 24-5665-1 | 24-5665-2 | 2 ea |

Orfit® Eco Black NS Sheets

- Non-stick, coated black material is designed for large and rigid splints
- Made of high quality recycled materials
- Appealing smooth finish resists fingerprints
- Does not stick to bandages or dressings



| perf | sheet | case qty | |
|-------------------------|-----------|-----------|------|
| solid 3/32" x 18" x 24" | 24-5870-1 | 24-5870-4 | 4 ea |
| solid 1/8" x 18" x 24" | 24-5872-1 | 24-5872-4 | 4 ea |
| solid 1/8" x 24" x 36" | 24-5874-1 | 24-5874-4 | 4 ea |
| 3.5% 3/32" x 18" x 24" | 24-5871-1 | 24-5871-4 | 4 ea |
| 3.5% 1/8" x 18" x 24" | 24-5873-1 | 24-5873-4 | 4 ea |
| 3.5% 1/8" x 24" x 36" | 24-5875-1 | 24-5875-4 | 4 ea |

Orfit® Eco Black NS Precuts

- make splint creation quick and easy



| perf | 1 piece | case qty | |
|------------------------------------|-----------|-----------|------|
| Intrinsic Anti-Spastic Hand | | | |
| 3.5% 1/8" small | 24-5945-1 | 24-5945-2 | 2 ea |
| 3.5% 1/8" medium | 24-5946-1 | 24-5946-2 | 2 ea |
| 3.5% 1/8" large | 24-5947-1 | 24-5947-2 | 2 ea |

Orfit® Natural NS Soft Sheets

- Non-stick, anti-bacterial coating
- Ideal for hand splints
- Features a smooth white finish



| perf | sheet | case qty | |
|-------------------------|-----------|-----------|------|
| solid 1/16" x 18" x 24" | 24-5710-1 | 24-5710-4 | 4 ea |
| solid 3/32" x 18" x 24" | 24-5712-1 | 24-5712-4 | 4 ea |
| solid 1/8" x 18" x 24" | 24-5714-1 | 24-5714-4 | 4 ea |
| 25% 1/8" x 18" x 24" | 24-5716-1 | 24-5716-4 | 4 ea |
| 1% 1/8" x 18" x 24" | 24-5715-1 | 24-5715-4 | 4 ea |
| 13% 1/16" x 18" x 24" | 24-5711-1 | 24-5711-4 | 4 ea |
| 13% 3/32" x 18" x 24" | 24-5713-1 | 24-5713-4 | 4 ea |

Orfit® Flex NS Sheets

- Non-stick, coated material
- Features high drapability and low resistance to stretch
- Anti-bacterial coating



| perf | sheet | case qty | |
|-------------------------|-----------|-----------|------|
| solid 1/8" x 18" x 24" | 24-5902-1 | 24-5902-4 | 4 ea |
| solid 1/8" x 24" x 36" | 24-5904-1 | 24-5904-4 | 4 ea |
| solid 3/32" x 18" x 24" | 24-5900-1 | 24-5900-4 | 4 ea |
| 1% 1/8" x 18" x 24" | 24-5903-1 | 24-5903-4 | 4 ea |
| 1% 1/8" x 24" x 36" | 24-5905-1 | 24-5905-4 | 4 ea |
| 1% 3/32" x 18" x 24" | 24-5901-1 | 24-5901-4 | 4 ea |

Orfibrace® NS Sheets

- Non-stick, coated material maintains perfect shape form and size for increased client comfort
- Anti-bacterial coating
- Cost effective when splints need to be remade



| perf | sheet | case qty | |
|------------------------|-----------|-----------|------|
| solid 1/8" x 24" x 36" | 24-5880-1 | 24-5880-4 | 4 ea |
| 1% 1/8" x 24" x 36" | 24-5881-1 | 24-5881-4 | 4 ea |

Orfit® Allfit

- Small pellets of thermoplastic material
- Easily mold into any desired shape or pattern



| | bag | box of 6 bags | |
|-------------------------|-----------|---------------|------|
| 1 bag, 1.1 lbs (500 gm) | 24-5910-1 | 24-5910-6 | 6 ea |

Orficast™ Thermoplastic Tape

- Ideal for creating quick and easy immobilization splints
- Tape features thermoplastic threads knitted with cotton weave
- Feels soft, yet comfortable on client
- Offers full support and immobilization
- Product will adhere to self when warmed, even more so with dry heat
- **Orficast™ More** has more thermoplastic threads and is more rigid; allows for more orthotic fabrication options

Orficast™ (0.05" thick)

| | roll | case qty |
|-----------------------------|-----------|------------------|
| 1/19" x 1" x 9.8 ft blue | 24-5600-1 | 24-5600-12 12 ea |
| 1/19" x 1" x 9.8 ft black | 24-5602-1 | 24-5602-12 12 ea |
| 1/19" x 2.5" x 9.8 ft blue | 24-5601-1 | 24-5601-6 6 ea |
| 1/19" x 2.5" x 9.8 ft black | 24-5603-1 | 24-5603-6 6 ea |

Orficast™ More (0.10" thick)

| | | |
|-----------------------------|-----------|----------------|
| 7/64" x 2.5" x 9.8 ft blue | 24-5610-1 | 24-5610-6 6 ea |
| 7/64" x 2.5" x 9.8 ft black | 24-5613-1 | 24-5613-6 6 ea |
| 7/64" x 5" x 9.8 ft blue | 24-5611-1 | 24-5611-6 6 ea |
| 7/64" x 5" x 9.8 ft black | 24-5614-1 | 24-5614-6 6 ea |
| 7/64" x 6" x 9.8 ft blue | 24-5612-1 | 24-5612-6 6 ea |
| 7/64" x 6" x 9.8 ft black | 24-5615-1 | 24-5615-6 6 ea |



Orfit® NS Strips

- Strips are made of non-stick, coated material to prevent accidental bonding
- Use for corrective finger and hand orthoses
- Adhere to other orthoses to increase rigidity

| | |
|---|-----------|
| 1/12" x 18" x 1/5" sonic silver (10 pc) | 24-5850-1 |
| 1/12" x 18" x 1/5" gold (10 pc) | 24-5851-1 |
| 1/12" x 18" x 1/5" sonic silver (5 pc), gold (5 pc) | 24-5852-1 |
| 1/8" x 18" x 4/5" gold (10 pc) | 24-5853-1 |
| 1/8" x 18" x 4/5" atomic blue (10 pc) | 24-5854-1 |
| 1/8" x 18" x 4/5" gold (5 pc), atomic blue (5 pc) | 24-5855-1 |



Orfit® Padding and Liner Luxofoam

- Self-adhesive, soft PVC foam used for padding small surfaces
- Measures: 19.3" x 38.2" x 0.11"



| | |
|---------|--------------------------|
| 24-5929 | luxofoam, non perforated |
| 24-5930 | luxofoam, perforated |

Orfilastic

- Terryloth lining with self-adhesive open cell foam backing
- Measures: 20" x 39.4" x 0.11"



| | |
|---------|------------------|
| 24-5931 | orfilastic |
| 24-5932 | orfilastic black |

Thermoplastic and Elastic Edging Strip

- Use as an edging material on splints, external immobilization devices and rehabilitation aids
- Measures: 0.04" x 0.75" x 118"



| | |
|---------|------------------------------|
| 24-5933 | thermoplastic, natural color |
| 24-5934 | elastic soft velvet, black |

Orfitube™ Outrigger Attachments

- Attach outriggers to splinting material to create mobilization orthoses (not shown)

| | |
|---------|---|
| 24-5911 | orfitube kit |
| 24-5925 | radial wrist extension splint kit, small |
| 24-5926 | radial wrist extension splint kit, medium |
| 24-5927 | radial wrist extension splint kit, large |
| 24-5912 | orfitubes (10 pc), bending wires (2 pcs) |
| 24-5914 | elastic thread, max resistance 7 oz, 30 ft |
| 24-5915 | elastic thread, max resistance 14 oz, 30 ft |
| 24-5916 | elastic thread, max resistance 21 oz, 30ft |
| 24-5917 | orfitube ends (20 pc) |
| 24-5918 | single orfitube adjuster |
| 24-5919 | single orfitube adjuster (5 pc), allen wrench |
| 24-5920 | orfitube bending tool |
| 24-5921 | finger coil spring (10 pc) |
| 24-5922 | knuckle bender coil spring (10 pc) |
| 24-5923 | finger extension assist (10 pc) |
| 24-5924 | wrist extension coil spring (10 pc, 5 left & 5 right) |
| 24-5928 | stainless steel spring wire, 45 ft |

- Secure to any type of base splint
- Precoiled wiring saves time in energy in the creation of more difficult splints



- **Orfitube™ Kit** includes: 10 orfitubes that are 330mm (13") long; 3 different types of elastic thread on a spool with different degrees of elastic strength; 1 spool of polyester monofilament for an optimal gliding effect; 2 bending wires
- **Radial Wrist Extension Splint Kit** includes: 2 Orfit® Classic Precuts (3.2 mm (1/8") non perforated, stiff); and 2 lateral coil springs

