



Bantam Extra Small Packages

Bantam Extra Small User Specifications

Approximate Height Range	28"-40"
Weight Limit	50 lbs.
Seat Depth Range (seat to back surface)	7"-12"
Seat to Foot Plate Range (from seat pivot)	4.5"-16.5"

Patterns and Colors

Hygienic Cover Pattern (select one):

- Space Race (standard)
- Rainbow Unicorn
- Sea Creatures
- Fancy Pants Elephants
- All Smiles Animals

Frame Color (select one):

- Green (standard)
- Purple
- White
- Orange

Minimum Support Package

The minimum support package includes: (no substitutions)

Supine, Easy-Adjust Seat Depth, Black Molded Tray, Velcro® Positioning Belt, Head Support, Planar Seat & Back, Gas Spring Lift with Foot Pedal, Knee Pads, 5" Front Wheels, 3" Rear Locking Swivel Casters, Chest Strap, Foot Plates.

Foot Holders (select one):

- Foot Holders 7.75"Lx3.25"W (standard)
- Foot Holders 6"Lx3"W
- Foot Holders 9.75"Lx4"W

Knee Pads (select one):

- Knee Pads 3.25" (inside width) (standard)
- Knee Pads 2.5"
- Knee Pads 4.25"

Back (select one):

- Planar Back 11"-13"H Range (height from seat) (standard)
- Planar Back 15"-17"H Range (height from seat)

Positioning Belt (select one):

- Velcro® Positioning Belt (fits hip circumference up to 23") (standard)
- Velcro® Positioning Belt (fits hip circumference up to 33")

Head Support (select one):

- Head Support-Form to Fit 5"Hx10"W (standard)
- Head Support 5"Hx8"W
- Head Support 6"Hx10"W

Moderate Support Package

The moderate support package includes: (no substitutions)

Supine, Easy-Adjust Seat Depth, Front Swivel Casters, Black Molded Swing-Away Tray, Forearm Wings, Multi-Adjustable Foot Plates, Secure Foot Straps, Hygienic Covers for (Seat, Back, Head Support), Swing-Away Knee Pads, Form to Fit Seat, Form to Fit Back, Velcro® Positioning Belt, X-Style Chest Vest, High Mount Chest Vest Bracket, Head Support, Gas Spring Lift with Foot Pedal, 3" Rear Locking Swivel Casters.

Foot Holders (select one):

- Foot Holders 7.75"Lx3.25"W (standard)
- Foot Holders 6"Lx3"W
- Foot Holders 9.75"Lx4"W

Knee Pads (select one):

- Knee Pads 3.25" (inside width) (standard)
- Knee Pads 2.5" (inside width)
- Knee Pads 4.25" (inside width)

Foot Straps (select one):

- Secure Foot Straps 10"L (length over top of foot) (standard)
- Secure Foot Straps 15"L (length over top of foot)

Back (select one):

- Form to Fit Back 11"-13"H Range (standard)
- Form to Fit Back 15"-17"H Range

Positioning Belt (select one):

- Velcro® Positioning Belt (fits hip circumference up to 23") (standard)
- Velcro® Positioning Belt (fits hip circumference up to 33")

Chest Vest (select one):

- X-Style Chest Vest 9.5"Lx9"W (standard)
- X-Style Chest Vest 11"Lx9.5"W

Head Support (select one):

- Head Support-Form to Fit 5"Hx10"W (standard)
- Head Support 5"Hx8"W
- Head Support 6"Hx10"W

Maximum Support Package

The maximum support package includes: (no substitutions)

Supine, Easy-Adjust Seat Depth, Dual Control, Gas Spring Lockout, Front Swivel Casters, Clear Swing-Away Tray, Forearm Wings, Multi-Adjustable Foot Plates, Secure Foot Straps, Hygienic Covers for (Seat, Back, Head Support), Swing-Away Knee Pads, Hip Supports, Calf Pad, Lateral Supports, Velcro® Positioning Belt, X-Style Chest Vest, High Mount Chest Vest Bracket, Head Support, Push Handle, Gas Spring Lift with Foot Pedal, Planar Seat & Back, 3" Rear Locking Swivel Casters.

Foot Holders (select one):

- Foot Holders 7.75"Lx3.25"W (standard)
- Foot Holders 6"Lx3"W
- Foot Holders 9.75"Lx4"W

Knee Pads (select one):

- Knee Pads 3.25" (inside width) (standard)
- Knee Pads 2.5" (inside width)
- Knee Pads 4.25" (inside width)

Back (select one):

- Planar Back 11"-13"H Range (height from seat) (standard)
- Planar Back 15"-17"H Range (height from seat)

Foot Straps (select one):

- Secure Foot Straps 10"L (length over top of foot) (standard)
- Secure Foot Straps 15"L (length over top of foot)

Lateral Supports (select one):

- Lateral Supports 5"-9" Range-Flat Pad (standard)
- Lateral Supports 8"-12" Range-Curved Pad
- Lateral Supports 8"-12" Range-Flat Pad
- Lateral Supports 5"-9" Range-Curved Pad

Positioning Belt (select one):

- Velcro® Positioning Belt (fits hip circumference up to 23") (standard)
- Velcro® Positioning Belt (fits hip circumference up to 33")

Chest Vest (select one):

- X-Style Chest Vest 9.5"Lx9"W (standard)
- X-Style Chest Vest 11"Lx9.5"W

Head Support (select one):

- Head Support-Form to Fit 5"Hx10"W (standard)
- Head Support 5"Hx8"W
- Head Support 6"Hx10"W

Moderate Support Package (Mobile)

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Tall Mobile, Mobile Tray, Multi-Adjustable Foot Plates, Secure Foot Straps, Hygienic Covers for (Seat, Back), Swing-Away Knee Pads, Hip Supports, Velcro® Positioning Belt, X-Style Chest Vest, High Mount Chest Vest Bracket, Gas Spring Lift with Foot Pedal, Planar Seat & Back, 3" Rear Locking Casters, Chest Strap.

Foot Holders (select one):

- Foot Holders 7.75"Lx3.25"W (standard)
- Foot Holders 6"Lx3"W
- Foot Holders 9.75"Lx4"W

Knee Pads (select one):

- 3.25" (inside width) (standard)
- 2.5" (inside width)
- 4.25" (inside width)

Back (select one):

- Planar Back 11"-13"H Range (height from seat) (standard)
- Planar Back 15"-17"H Range (height from seat)

Foot Straps (select one):

- Secure Foot Straps 10"L (length over top of foot) (standard)
- Secure Foot Straps 15"L (length over top of foot)

Positioning Belt (select one):

- Velcro® Positioning Belt (fits hip circumference up to 23") (standard)
- Velcro® Positioning Belt (fits hip circumference up to 33")

Chest Vest (select one):

- X-Style Chest Vest 9.5"Lx9"W (standard)
- X-Style Chest Vest 11"Lx9.5"W