



Bantam Small Packages

Bantam Small User Specifications

Approximate Height Range	36"-54"
Weight Limit	100 lbs.
Seat Depth Range (seat to back surface)	11"-16"
Seat to Foot Plate Range (from seat pivot)	4.5"-16.5"

Patterns and Colors

Hygienic Cover Pattern (select one):

- ☐ Space Race (standard)
- ☐ Rainbow Unicorn
- ☐ Sea Creatures
- ☐ Fancy Pants Elephants
- ☐ All Smiles Animals

Frame Color (select one):

- ☐ Green (standard)
- ☐ Purple
- ☐ White
- ☐ Orange

Minimum Support Package

The minimum support package includes: (no substitutions)

Supine, Easy-Adjust Seat Depth, Black Molded Tray, Velcro® Positioning Belt, Head Support, Planar Seat & Back, Gas Spring Lift with Foot Pedal, Knee Pads, 5" Front Wheels, 3" Rear Locking Swivel Casters, Chest Strap, Foot Plates.

Foot Holders (select one):

- ☐ Foot Holders 9.75"Lx4"W (standard)
- ☐ Foot Holders 7.75"Lx3.25"W
- ☐ Foot Holders 11.75"Lx5"W

Knee Pads (select one):

- ☐ Knee Pads 4.25" (inside width) (standard)
- ☐ Knee Pads 3.25" (inside width)
- ☐ Knee Pads 5" (inside width)

Back (select one):

- ☐ Planar Back 14"-16"H Range (height from seat) (standard)
- ☐ Planar Back 18"-20"H Range (height from seat)

Positioning Belt (select one):

- ☐ Velcro® Positioning Belt (fits hip circumference up to 29") (standard)
- ☐ Velcro® Positioning Belt (fits hip circumference up to 39")

Head Support (select one):

- ☐ Head Support-Form to Fit 5"Hx10"W (standard)
- ☐ Head Support 5"Hx8"W
- ☐ Head Support 6"Hx10"W

Moderate Support Package

The moderate support package includes: (no substitutions)

Supine, Easy-Adjust Seat Depth, Front Swivel Casters, Black Molded Swing-Away Tray, Forearm Wings, Multi-Adjustable Foot Plates, Secure Foot Straps, Hygienic Covers for (Seat, Back, Head Support), Swing-Away Knee Pads, Form to Fit Seat, Form to Fit Back, Positioning Belt, X-Style Chest Vest, High Mount Chest Vest Bracket, Head Support, Gas Spring Lift with Foot Pedal, 3" Rear Locking Swivel Casters.

Foot Holders (select one):

- ☐ Foot Holders 9.75"Lx4"W (standard)
- ☐ Foot Holders 7.75"Lx3.25"W
- ☐ Foot Holders 11.75"Lx5"W

Knee Pads (select one):

- ☐ Knee Pads 4.25" (inside width) (standard)
- ☐ Knee Pads 3.25" (inside width)
- ☐ Knee Pads 5" (inside width)

Foot Straps (select one):

- ☐ Secure Foot Straps 10"L (length over top of foot) (standard)
- ☐ Secure Foot Straps 15"L (length over top of foot)

Back (select one):

- ☐ Form to Fit Back 12.5"-14.5"H (range) (standard)
- ☐ Form to Fit Back 16.5"-18.5"H (range) Range

Positioning Belt (select one):

- ☐ Velcro® Positioning Belt (fits hip circumference up to 29") (standard)
- ☐ Velcro® Positioning Belt (fits hip circumference up to 39")

Chest Vest (select one):

- ☐ X-Style Chest Vest 9.5"Lx9"W (standard)
- ☐ X-Style Chest Vest 11"Lx9.5"W

Head Support (select one):

- ☐ Head Support-Form to Fit 5"Hx10"W (standard)
- ☐ Head Support 5"Hx8"W
- ☐ Head Support 6"Hx10"W

Maximum Support Package

The maximum support package includes: (no substitutions)

Supine, Easy-Adjust Seat Depth, Dual Control, Gas Spring Lockout, Front Swivel Casters, Clear Swing-Away Tray, Forearm Wings, Multi-Adjustable Foot Plates, Secure Foot Straps, Hygienic Covers for (Seat, Back, Head Support), Swing-Away Knee Pads, Hip Supports, Calf Pad, Lateral Supports, Positioning Belt, X-Style Chest Vest, High Mount Chest Vest Bracket, Head Support, Push Handle, Gas Spring Lift with Foot Pedal, Planar Seat & Back, 3" Rear Locking Swivel Casters.

Foot Holders (select one):

- ☐ Foot Holders 9.75"Lx4"W (standard)
- ☐ Foot Holders 7.75"Lx3.25"W
- ☐ Foot Holders 11.75"Lx5"W

Knee Pads (select one):

- ☐ Knee Pads 4.25" (inside width) (standard)
- ☐ Knee Pads 3.25" (inside width)
- ☐ Knee Pads 5" (inside width)

Back (select one):

- ☐ Planar Back 14"-16"H Range (height from seat) (standard)
- ☐ Planar Back 18"-20"H Range (height from seat)

Foot Straps (select one):

- ☐ Secure Foot Straps 10"L (length over top of foot) (standard)
- ☐ Secure Foot Straps 15"L (length over top of foot)

Lateral Supports (select one):

- ☐ Lateral Supports 8"-12" Range-Flat Pad (standard)
- ☐ Lateral Supports 8"-12" Range-Curved Pad
- ☐ Lateral Supports 5"-9" Range-Curved Pad
- ☐ Lateral Supports 5"-9" Range-Flat Pad

Positioning Belt (select one):

- ☐ Velcro® Positioning Belt (fits hip circumference up to 29") (standard)
- ☐ Velcro® Positioning Belt (fits hip circumference up to 39")

Chest Vest (select one):

- ☐ X-Style Chest Vest 9.5"Lx9"W (standard)
- ☐ X-Style Chest Vest 11"Lx9.5"W

Head Support (select one):

- ☐ Head Support-Form to Fit 5"Hx10"W (standard)
- ☐ Head Support 5"Hx8"W
- ☐ Head Support 6"Hx10"W

Moderate Support Package (Mobile)

The moderate support package includes: (no substitutions)

Easy-Adjust Seat Depth, Tall Mobile, Mobile Tray, Multi-Adjustable Foot Plates, Secure Foot Straps, Hygienic Covers for (Seat, Back), Swing-Away Knee Pads, Hip Supports, Positioning Belt, X-Style Chest Vest, High Mount Chest Vest Bracket, Gas Spring Lift with Foot Pedal, Planar Seat & Back, 3" Rear Locking Casters, Chest Strap.

Foot Holders (select one):

- ☐ Foot Holders 9.75"Lx4"W (standard)
- ☐ Foot Holders 7.75"Lx3.25"W
- ☐ Foot Holders 11.75"Lx5"W

Knee Pads (select one):

- ☐ Knee Pads 4.25" (inside width) (standard)
- ☐ Knee Pads 3.25" (inside width) (inside width)
- ☐ Knee Pads 5" (inside width) (inside width)

Back (select one):

- ☐ Planar Back 14"-16"H Range (height from seat) (standard)
- ☐ Planar Back 18"-20"H Range (height from seat)

Foot Straps (select one):

- ☐ Secure Foot Straps 10"L (length over top of foot, pair) (standard)
- ☐ Secure Foot Straps 15"L (length over top of foot, pair)

Positioning Belt (select one):

- ☐ Velcro® Positioning Belt (fits hip circumference up to 29") (standard)
- ☐ Velcro® Positioning Belt (fits hip circumference up to 39")

Chest Vest (select one):

- ☐ X-Style Chest Vest 9.5"Lx9"W (standard)
- ☐ X-Style Chest Vest 11"Lx9.5"W