

Designed to provide the best absorbency to protect a mattress or chair from urinary incontinence and other spills. These are washable and reusable for long periods of use.

Features:

- Washable and quick drying
- Triple layer protection
- Quilted cotton comfort top layer
- Extra absorbant soaker layer in middle
- Protective vinyl barrier on bottom with knitted backing
- Rounded corners with reinforced stitching for durability
- Absorbency is approximately: 1,500 ml when it's laying flat on a smooth surface



Item	Description	Size	UOM
D0095-1724R-1	Chair Pad	17" x 24"	1 Ea
D0095-4235R-1	Bed Pad for full	42" x 35"	1 Ea
D0095-5435R-1	Bed Pad for Queen	54" x 35"	1 Ea
D0095-7835R-1	Bed Pad for King	78" x 35"	1 Ea
D0095	Bed Pad for partial	29"x35"	1 Ea
D0095-3	Bed Pad Pack	29"x35"	3/Pck

Care Instructions:

- Machine wash hot water
- Rinse
- Tumble dry low heat
- No softener or bleach

Note: patient's weight and positioning on top can alter the absorbency amount along with eventual wear and tear which can also have an effect. Also, it's not recommended that this quantity of liquid or even a portion of it is left under a patient, as the patient's comfort and well being will be effected when positioned on a wet surface for an extended period of time.