

Old thermometer got your temperature rising?



Four easy steps.
One accurate reading.

Medline's Tympanic Thermometer with Quick Probe Release



Tympanic Thermometer with Quick Probe Release

Get fast, accurate and easy-to-read results with Medline's Tympanic Thermometer with Quick Probe Release. This ear thermometer is designed as a comfortable and convenient method for taking a person's body temperature. No embarrassing or invasive procedures; you simply place the thermometer in the ear to get a reading within seconds. It's perfectly suited for both home and clinical use on patients of all ages. Plus, its unique mechanical probe release helps reduce cross contamination, allowing you to easily change probe covers with every use!

Why an ear thermometer?

It turns out that the eardrum is an extremely accurate point to measure body temperature from because it is recessed inside the head (just like your tongue). Instead of touching the eardrum with a thermometer, the tympanic thermometer uses an infrared sensor to measure the temperature of energy radiating from the eardrum.

Innovative, quick-response ear thermometer features:

- One-touch operation
- Mechanical probe cover quick-release allows easy probe cover change
- Includes 100 probe covers
- Safe to use on patients of any age
- Readings in Fahrenheit and Celsius
- Takes an accurate, non-invasive temperature from the tympanic membrane



Fast. Simple. Noninvasive. It's easy to see why Medline's Tympanic Thermometer with Quick Probe Release is ideal for use on patients of all ages at home, in hospitals and in long-term care facilities.

Item Number	Description
MDS9700	Thermometer, each
MDS9701	Probe covers, 100/bx


To order your Tympanic Thermometer with Quick Probe Release, please contact your local Medline representative, call 1.800.MEDLINE or visit us online at www.medline.com.



Apply disposable cover onto the tympanic thermometer.



Turn device "on" by pressing button on front.



Place thermometer in ear and depress top button to take reading.



After reading, flip the blue tab to release the disposable probe cover without contact.

Medline Industries, Inc.
One Medline Place
Mundelein, Illinois 60060
www.medline.com
1.800.MEDLINE

**Simple Solutions.
Significant Results.™**

