



TRADITIONAL **WOUND CARE**

Unsurpassed quality in a complete line of products



Why Medline?

Medline is a one-stop shop with the products and services you need to achieve immediate results and sustain improvements for the long haul. How do we do it?

- We're a manufacturer and distributor, unlike other companies. By eliminating the middleman, we save you time and money.
- We're knowledgeable about our products. If you call us with a question, you get an answer. Simple as that.
- We have sales reps, inventory and warehouses all over the world. So we're there when and where our customers need us.
- We provide ongoing education and competency reviews that ensure your staff has the skill levels necessary to do their jobs.
- We understand that ongoing nursing education helps to improve the quality of patient care and employee morale while empowering staff to deliver care consistent with current standards of practice and government guidelines.
- That's why we offer a variety of courses and programs designed to help nurses learn and continue to grow professionally.
- Medline offers a comprehensive array of traditional wound care products to meet the full range of our customers' needs. We are one of the largest Providers of wound care products in the country.
- To ensure the highest standards of quality, all Medline brand products are manufactured in factories that follow FDA Quality System Regulations.
- All of our sterile gauze products are also produced to meet USP Type VII regulations.
- We are proud to offer one of the largest selections of latex-free products.

We manufacture 70 percent of the products we sell and have 30 distribution sites strategically located countrywide.

As a result, we can serve 98.2 percent of the marketplace within two days.

Table of Contents



Abdominal (ABD) Pads
See page 8



**Sof-Form®
Conforming Bandages**
See page 12



Bulkee II® Bandages
See page 10

Gauze Sponges

Sterile Woven Gauze Sponges	4
Non-Sterile Woven Gauze Sponges	4
X-Ray Detectable Woven Sponges	5
Non-Woven Avant Deluxe® Gauze Sponges	6
Drain Sponges	7
Non-Woven Avant Gauze® Sponges	7
Sterile Post-Op Gauze Sponges	8

Bandages/Rolls

Abdominal (ABD) Pads	8
Eye Pads	8
Cotton Balls	9
Cotton Rolls	9
Non-Adherent Pads	9
Bulkee II® Wound Dressings	10-11
Sof-Form® Conforming Bandages	12
Elastic Net	13
Elastic Bandages	14-15
Co-Flex Self-Adherent Wrap	20
Wytex Undercast Padding	23

Impregnated Gauze

Petrolatum	16
Xeroform	16
Oil Emulsion Gauze	17
Packing Strips	17

Tapes & Adhesives

Tapes	18-19
Medi-Strips® Wound Closure Strips	19
Co-Flex Cohesive Bandage	20-21
Adhesive Bandages	24-28

Miscellaneous

Cotton and Rayon Balls	9
Undercast Padding	23
First Aid Kits	29
Cross-Reference Guide	30-39

Woven Gauze Sponges

The best for traditional wound care

These 100% cotton sponges are ideal for wound dressings, wound packing and general wound care. Our strictly controlled manufacturing process assures clean, debris-free packing with folded edges to prevent unraveling. Sterile 2's are packed in coated paper envelopes to reduce fiber debris. This easy-open envelope is also carefully sealed to prevent the intrusion of dust and contaminants. Sterile 10's feature strong rigid trays. Also available in non-sterile bulk packs. Latex-Free.



Gauze Sponges, Sterile, 2's (USP Type VII)

Item#	Size	Quantity
NON21420	2" x 2", 8-ply (5.08 x 5.08 cm)	50 pk/bx, 30bx/cs, 1,500 pks/cs
NON21422	3" x 3", 12-ply (7.62 x 7.62 cm)	40 pk/bx, 30bx/cs, 1,200 pks/cs
NON21448	4" x 4", 8-ply (10.16 x 10.16 cm)	50 pk/bx, 12bx/cs, 600 pks/cs
NON21424	4" x 4", 12-ply 10.16 x 10.16 cm	25 pk/bx, 24bx/cs, 600 pks/cs
NON21842	8" x 4", 12-ply 20.3 x 10.16 cm	25 pk/bx, 20bx/cs, 500 pks/cs

Gauze Sponges, Sterile 10's (Trays) (USP Type VII)

Item#	Size	Quantity
NON21426	4" x 4", 12-ply (10.16 x 10.16 cm)	10/tray, 128 trays/cs, 1280/cs
NON21428	4" x 4", 16-ply (10.16 x 10.16 cm)	10/tray, 128 trays/cs, 1280/cs
NON21425	4" x 4", 8-ply (10.16 x 10.16 cm)	10/tray, 128 trays/cs, 1280/cs

Gauze Sponges, Non-Sterile (Bulk)

Item#	Size	Quantity
NON25208	2" x 2", 8-ply (5.08 x 5.08 cm)	200/sleeve, 25 sleeves/cs, 5000/cs
NON25212	2" x 2", 12-ply (5.08 x 5.08 cm)	200/sleeve, 40 sleeves/cs, 8000/cs
NON25312	3" x 3", 12-ply (7.62 x 7.62 cm)	200/sleeve, 20 sleeves/cs, 4000/cs
NON25408	4" x 4", 8-ply (10.16 x 10.16 cm)	200/sleeve, 20 sleeves/cs, 4000/cs
NON25412	4" x 4", 12-ply (10.16 x 10.16 cm)	200/sleeve, 10 sleeves/cs, 2000/cs
NON25416	4" x 4", 16-ply (10.16 x 10.16 cm)	200/sleeve, 10 sleeves/cs, 2000/cs
NON25812	4" x 8", 12-ply (10.16 x 10.16 cm)	200/sleeve, 10 sleeves/cs, 2000/cs

X-Ray Detectable Woven Gauze Sponges

Low linting, fluffy and highly absorbent, our x-ray detectable sponges are made from 100% cotton. Available in multiple sizes and plies. This product features dual x-ray detectable barium sulfate strands in a unique figure 8 design (for easy detection) layered into the sponges. Latex-Free.

Gauze Sponges, Sterile (Trays)

(USP Type VII)

Item#	Size	Quantity
NON21430LF	4" x 4", 16-ply (10.16 x 10.16 cm)	10 sponges/tray, 128 trays/cs, 1280/cs
NON21432LF	8" x 4", 12-ply (20.3 x 10.16 cm)	10 sponges/tray, 80 trays/cs, 800/cs
NON21433	4" x 4", 32-ply (10.16 x 10.16 cm)	10 sponges/tray, 48 trays/cs, 480/cs
NON21438	8" x 4", 24-ply (20.3 x 10.16 cm)	10 sponges/tray, 20 trays/cs, 200/cs

Gauze Sponges, Non-Sterile (Bulk)

Item#	Size	Quantity
NON25430	4" x 4", 16-ply (10.16 x 10.16 cm)	10 sponges/band, 100/sleeve, 2000/cs
NON25432	8" x 4", 12-ply (20.3 x 10.16 cm)	10 sponges/band, 100/sleeve, 2000/cs

100% Cotton Drain Sponges

If you're looking for natural breathability in a drain sponge, cotton is the right choice. Made of 100% cotton, they are pre-cut for your convenience.

100% Cotton Drain Sponges, Sterile

Item#	Size	Quantity
NON255000	4" x 4", 16-ply (10.16 x 10.16 cm)	1/pk 25/bx, 12 bx/cs

NON21430LF 4" x 4" Sterile Woven Gauze Sponges

X-Ray Detectable



NON255000 4" x 4" Sterile Cotton Drain Sponges



Our deluxe line of non-woven sponges

Avant Gauze®

Medline has a complete line of non-woven gauze sponges. Constructed from a blend of rayon and polyester, these sponges offer unique benefits compared to traditional cotton gauze:

- low linting
- increased absorbency
- fast wicking
- economical
- less traumatic to granulating tissue

Avant Deluxe™ Gauze

This is our deluxe line of non-woven sponges that are a unique rayon and polyester blend to provide strength and absorbency. The increased material weight helps to provide a more efficient sponge than our traditional Avant line by reducing the number of sponges that have to be used for each procedure. These non-woven sponges provide increased wicking ability as well as less lint versus the traditional woven sponges. Latex-Free.

Gauze Sponges, Sterile 2's

Item#	Size	Quantity
NON264442	4" x 4", 4-ply, (10.16 x 10.16 cm)	50 pk/bx, 12 bx/cs, 600 pks/cs

Gauze Sponges, Non-Sterile (Bulk)

Item#	Size	Quantity
NON26223	2" x 2", 3-ply (5.08 x 5.08 cm)	200/sleeve, 40 sleeves/cs, 8000/cs
NON26224	2" x 2", 4-ply (5.08 x 5.08 cm)	200/sleeve, 40 sleeves/cs, 8000/cs
NON26334	3" x 3", 4-ply (7.62 x 7.62 cm)	200/sleeve, 20 sleeves/cs, 4000/cs
NON26443	4" x 4", 3-ply (10.16 x 10.16 cm)	200/sleeve, 10 sleeves/cs, 2000/cs
NON26444	4" x 4", 4-ply (10.16 x 10.16 cm)	200/sleeve, 10 sleeves/cs, 2000/cs

NON26223	2" x 2" Non-Sterile Avant Deluxe Gauze
NON21444	4" x 4" Sterile Avant Non-woven Gauze



latexfree
PRODUCT

Avant Gauze® Non-Woven Sponge

Non-woven dressing

Avant Gauze is Medline's standard non-woven dressing. Compared to traditional cotton gauze, Avant Gauze produces less lint and is less traumatic to fragile granulating wound tissue because it's virtually non-adherent. It also has significantly faster wicking ability which promotes a healthy healing environment. Discover the perfect balance between strength and absorbency in this rayon/polyester blend gauze. Latex-Free.



Gauze Sponges, Sterile 2's

Item#	Size	Quantity
NON21224	2" x 2", 4-ply (5.08 x 5.08 cm)	50 pk/bx, 30 bx/cs, 1500 pks/cs
NON21334	3" x 3", 4-ply (7.62 x 7.62 cm)	40 pk/bx, 30 bx/cs, 1200 pks/cs
NON21444	4" x 4", 4-ply (10.16 x 10.16 cm)	50 pk/bx, 12 bx/cs, 600 pks/cs
NON21446	4" x 4", 6-ply (10.16 x 10.16 cm)	25 pk/bx, 12 bx/cs, 300 pks/cs

Gauze Sponges, Sterile 10's (Trays)

Item#	Size	Quantity
NON21449	4" x 4", 4-ply, (10.16 x 10.16 cm)	10/tray, 128 trays/cs, 1280/cs

Gauze Sponges, Non-Sterile (Bulk)

Item#	Size	Quantity
NON25223	2" x 2", 3-ply (5.08 x 5.08 cm)	200/sleeve, 25 sleeves/cs, 5000/cs
NON25224	2" x 2", 4-ply (5.08 x 5.08 cm)	200/sleeve, 40 sleeves/cs, 8000/cs
NON25334	3" x 3", 4-ply (7.62 x 7.62 cm)	200/sleeve, 20 sleeves/cs, 4000/cs
NON25443	4" x 4", 3-ply (10.16 x 10.16 cm)	200/sleeve, 20 sleeves/cs, 4000/cs
NON25444	4" x 4", 4-ply (10.16 x 10.16 cm)	200/sleeve, 10 sleeves/cs, 2000/cs

Avant Gauze® Drain Sponge

Special fenestration for easy application

This sponge features an extra fenestration which allows it to lay flat when wrapped around a drainage tube. Without it, the sponge will pucker, decreasing effectiveness by allowing fluids to pool around the drain tube. It's a simple feature that most other drain sponges do not offer. Includes all the benefits of non-woven Avant Gauze. Each 6-ply rayon/polyester sponge is highly absorbent, low-linting and fast wicking to wound beds. Latex-Free.

Drain Sponges, Sterile 2's

Item#	Size	Quantity
NON256000	4" x 4", 6-ply (10.16 x 10.16 cm)	25 pk/bx, 300 pk/cs, 12 bx/cs
NON256001	2" x 2", 6-ply (5.08 x 5.08 cm)	35 pk/bx, 700 pk/cs, 20 bx/cs

NON256000 4" x 4" Sterile Avant Gauze Drain Sponge



NON21442 4" x 4" Sterile Post-Op Gauze



Post-Op Gauze

More absorbent than cotton gauze sponges, these post-op sponges are filled with extra-absorbent cellulose fibers for those situations that require increased absorbency. The outer layer is constructed of our absorbent non-woven Avant Gauze fibers. Packed in sterile 2's. Latex-Free.

Sterile 2's

Item#	Size	Quantity
NON21442	4" x 4" (10.16 x 10.16 cm)	25 pk/bx, 24 bx/cs, 600 pks/cs
NON21441	4" x 3" (10.16 x 7.62 cm)	25 pk/bx, 24 bx/cs, 600 pks/cs

Non-Sterile (Bulk)

Item#	Size	Quantity
NON25442	4" x 4" (10.16 x 10.16 cm)	100/sleeve, 20 sleeves/cs, 2000/cs
NON25441	4" x 3" (10.16 x 7.62 cm)	100/sleeve, 20 sleeves/cs, 2000/cs

Eye Pads

These sterile Eye Pads are made with a soft non-woven cover, over an absorbent cotton filler. Compressed edges help reduce linting. Latex-Free.

Sterile 1's

Item#	Size	Quantity
NON21600	1 ⁵ / ₈ " x 2 ⁵ / ₈ " (4.1 x 6.7 cm)	1/pk, 50 pk/bx, 600/cs
NON21601	2 ¹ / ₈ " x 2 ⁵ / ₈ " (5.4 x 6.7 cm)	1/pk, 50 pk/bx, 600/cs

NON21450 5" x 9" Sterile ABD Pads



Abdominal (ABD) Pads

Features extra absorbency and wicking

Medline's super absorbent abdominal pads feature a soft non-woven outer layer that quickly wicks fluid to a cellulose center. This thick layer of cellulose quickly absorbs and disperses fluids laterally to prevent pooling. Hydrophobic back impedes strike-through and all edges are sealed to prevent linting. Latex-Free.

Sterile 1's

Item#	Size	Quantity
NON21450	5" x 9" (12.7 x 22.9 cm)	25/bx, 16 bx/cs, 400/cs
NON21453	8" x 7 ¹ / ₂ " (20.3 x 19.1 cm)	20/bx, 12 bx/cs, 240/cs
NON21454	8" x 10" (20.3 x 25.4 cm)	18/bx, 20 bx/cs, 360/cs


Non-Sterile (Bulk)

Item#	Size	Quantity
NON21451	5" x 9" (12.7 x 22.9 cm)	576 pads/cs
NON21452	8" x 7 ¹ / ₂ " (20.3 x 19.1 cm)	576 pads/cs
NON21456	8" x 10" (20.3 x 25.4 cm)	432 pads/cs
NON21457	12" x 16" (30.5 x 40.6 cm)	144 pads/cs

Multi-Trauma Dressing (Sterile)

NON21459	10" x 30" (25.4 x 76.2 cm)	50/cs
DYNDRAINPACK	Heavy drainage pack/sterile	48 ea/cs



NON25710 3" x 4" Sterile Non-Adherent Pads*Non-Adherent Pads*

This sterile, absorbent rayon/polyester pad will not interfere with granulating tissue by sticking to the lightly draining wound. The soft perforated film allows air to circulate and the absorbent pad allows the fluid to collect. It can easily be cut to fit any shape without separating. Latex-Free.

Sterile 1's

Item#	Size	Quantity
NON25700	2" x 3", Sterile, 1's (5.08 x 7.62 cm)	100/bx, 12 bx/cs, 1200/cs
NON25710	3" x 4", Sterile, 1's (7.62 x 10.16 cm)	100/bx, 12 bx/cs, 1200/cs
NON25720	3" x 8", Sterile, 1's (7.62 x 20.3 cm)	50/bx, 12 bx/cs, 600/cs

Cotton Roll


A one pound of 100% cotton wrapped and neatly rolled in paper. Highly absorbent surgical grade cotton is ideal for wound cleansing, application of medications, and padding of wound area. Latex-Free.

Sterile

Item#	Size	Quantity
NON6026	1 lb. (.5 kg)	25 roll/cs

Non-Sterile

Item#	Size	Quantity
NON6027	1 lb. (.5 kg)	25 roll/cs

MDS21460 1" Cotton Balls*Cotton Balls and Rayon Balls*

Soft and absorbent. Designed for many uses around any facility.

Sterile

Item#	Size	Quantity
NON6101	Cotton, med, 1" diameter (2.54 cm)	65 ea/bx, 48 bx/cs, 3120 ea/cs

Non-Sterile

Item#	Size	Quantity
MDS21460	Cotton, med, 1" diameter (2.54 cm)	2000 ea/bag, 2 bags/cs, 4000 ea/cs
MDS21461	Cotton, med, 1" diameter (2.54 cm)	200 ea/bag, 20 bags/cs, 4000 ea/cs
MDS21462	Cotton, lg, 1 1/4" diameter (3.18 cm)	1000 ea/bag, 2 bags/cs, 2000 ea/cs
MDS21463	Cotton, lg, 1 1/4" diameter (3.18 cm)	50 ea/bag, 80 bags/cs, 4000 ea/cs
MDS21470	Rayon, med, 1" diameter (2.54 cm)	2000 ea/bag, 2 bags/cs, 4000 ea/cs
MDS21472	Rayon, lg, 1 1/4" diameter (3.18 cm)	1000 ea/bag, 2 bags/cs, 2000 ea/cs

NON6026 Cotton Roll

Bulkee II® Wound Dressings

Our finest and most popular cotton bandages



Bulkee II® 100% Cotton Gauze Bandages

Compare to Kerlix®*

Ideal for both primary and secondary wound care dressings. When used as a primary dressing, the open weave design provides fast wicking, absorbency and aeration. Bulkee II has a tightly finished edge to reduce linting; and crimps for increased cushioning and patient comfort. The 6-ply Bulkee II bandage stretches and conforms to a variety of patient needs. It performs exceptionally well in burn care applications. Latex-Free.

Sterile (Singles)

Item#	Size	Quantity
NON25865	4.5" x 4.1 yds, stretched (11.4 cm x 3.75 m)	100 rolls/cs
NON25861	3.4" x 3.6 yds, stretched (8.6 cm x 3.30 m)	96 rolls/cs
NON25850	2.25" x 3 yds, stretched (5.7 cm x 2.74 m)	96 rolls/cs

Non-Sterile

Item#	Size	Quantity
NON25855	4.5" x 4.1 yds, stretched (11.4 cm x 3.75 m)	100 rolls/cs

Bulkee® Lite Conforming Bandage

Compare to Kerlix® Lite*

This 100% crimped cotton 3-ply woven bandage is ideal for a secondary dressing to secure IV's, dressings, splints, or for providing mild compression. The bandage provides cushioning and a light stretch, which helps it to conform to difficult body contours. Latex-Free.

NON27499 6" x 4.5 yds Sterile Bulkee Lite Bandage



Sterile (Singles)

Item#	Size	Quantity
NON27496	2" x 3.5 yds, stretched (5.08 cm x 3.20 m)	12/bx, 8 bx/cs, 96/cs
NON27497	3" x 4.1 yds, stretched (7.62 cm x 3.75 m)	12/bx, 8 bx/cs, 96/cs
NON27498	4" x 4.1 yds, stretched (10.16 cm x 3.75 m)	12/bx, 8 bx/cs, 96/cs
NON27499	6" x 4.5 yds, stretched (15.24 cm x 4.11 m)	6/bx, 8 bx/cs, 48/cs

Non-Sterile (Bulk)

Item#	Size	Quantity
NON27492	2" x 3.5 yds, stretched (5.08 cm x 3.2 m)	12/bag, 8 bx/cs, 96/cs
NON27493	3" x 4.1 yds, stretched (7.62 cm x 3.75 m)	12/bx, 8 bx/cs, 96/cs
NON27494	4" x 4.1 yds, stretched (10.16 cm x 3.75 m)	12/bx, 8 bx/cs, 96/cs
NON27495	6" x 4.5 yds, stretched (15.24 cm x 4.11 m)	6/bx, 8 bx/cs, 48/cs

Bulkee® Super Fluff Sponges

All the great features of Bulkee II® in a sponge

Compare to Kerlix® Fluff Sponges*

Ideal loft, bulk and absorbency make the Bulkee sponge an excellent primary dressing. It effectively manages exudate in and around the wound. Because of its strength and absorbency, this sponge is ideal for pre-op prepping and cleansing, debriding and packing. Latex-Free.

NON25853 6" x 6.75" Fluff Sponges



Sterile

Item#	Size	Quantity
NON25852	6" x 6.75", (15.24 x 17.1 cm) diagonally	2 sponges/pk, 20 pk/bx, 12 bx/cs, 480 sponges/cs
NON25853	6" x 6.75", (15.24 x 17.1 cm) diagonally	5 sponges/tray, 120 trays/cs, 600 sponges/cs
NON25854	6" x 6.75", (15.24 x 17.1 cm) diagonally	10 sponges/tray, 48 trays/cs, 480 sponges/cs
NON25856	7.75" x 8.75", (19.7 x 22.2 cm) diagonally	2 sponges/pk, 25 pk/bx, 16 bx/cs, 400 sponges/cs
NON25844	4" x 4", 12-ply sponges	2 sponges/pk, 25 pk/bx, 24 bx/cs, 1200 sponges/cs
NON25884	4" x 4", 12-ply sponges	10 sponges/tray, 128 trays/cs, 1280 sponges/cs

Non-Sterile

Item#	Size	Quantity
NON1167B	6" x 6.75" (15.24 x 17.1 cm)	100/bag, 1200/cs

*This product is not manufactured by Tyco Healthcare Group, makers of Kerlix dressings, bandages and sponges.

Sof-Form[®] Conforming Bandages

Compare to Kendall Conform^{®*}

This premier quality bandage gently secures to any body contour without restriction of movement and accommodates an increase or decrease in swelling. The uniquely constructed weave and finished edges help eliminate unraveling and lint. The combination of rayon and polyester help the bandage to provide stretch memory when applied. It's a patient-safe dressing that renders a neat application each time. Ideal for secondary dressing applications. Latex-Free.



Sterile (Singles)

Item#	Size	Quantity
NON254955	1" x 75", relaxed (2.54 x 190.52 cm)	12/bx, 8 bx/cs, 96/cs
NON25496	2" x 75", relaxed (5.08 x 190.52 cm)	12/bx, 8 bx/cs, 96/cs
NON25497	3" x 75", relaxed (7.62 x 190.52 cm)	12/bx, 8 bx/cs, 96/cs
NON25498	4" x 75", relaxed (10.16 x 190.52 cm)	12/bx, 8 bx/cs, 96/cs
NON25499	6" x 75", relaxed (15.2 x 190.52 x cm)	6/bx, 8 bx/cs, 48/cs

Non-Sterile (Bulk)

Item#	Size	Quantity
NON25491	1" x 75", relaxed (2.54 x 190.52 cm)	12/bag, 8 bags/cs, 96/cs
NON25492	2" x 75", relaxed (5.08 x 190.52 cm)	12/bag, 8 bags/cs, 96/cs
NON25493	3" x 75", relaxed (7.62 x 190.52 cm)	12/bag, 8 bags/cs, 96/cs
NON25494	4" x 75", relaxed (10.16 x 190.52 cm)	12/bag, 8 bags/cs, 96/cs
NON25495	6" x 75", relaxed (15.2 x 190.52 cm)	6/bag, 8 bags/cs, 48/cs

*This product is not manufactured or distributed by Tyco Healthcare Group, makers of Conform gauze bandages.

Tubular Gauze Bandages

Made of knitted 100% cotton, these white tubular bandages are available in a range of sizes to facilitate bandaging contoured parts of the body. They also can be used for dressing retention, for simple fractures, and for splinting. For use with applicators. Supplied in 50-yd roll.

NONTUB058 Tubular Gauze Bandages

Tubular Gauze Bandages

Item#	Description	Size	Quantity
NONTUB058	Small finger, toe	50 yd, 5/8" (1.59 cm)	36/cs
NONTUB100	Large finger, toe	50 yd, 1" (2.54 cm)	36/cs
NONTUB112	Small arm, hand, foot	50 yd, 1 1/2" (3.81 cm)	24/cs
NONTUB258	Arm, lower leg, small thigh	50 yd, 2 5/8" (7.94 cm)	12/cs
NONTUB358	Leg, thigh, head, shoulder	50 yd, 3 5/8" (11.11 cm)	6/cs



Elastic Net

Pre-cut Lengths for Individual use now available!

Compare to Spandage™

NONNET04 Elastic Net (Size 4)



Elastic Net

25 yards stretched length, 1 roll per box

Item#	Size	Max Circ. Stretch Width	Description
NONNET01	1	6" (15.2 cm)	Fingers, toes, wrists
NONNET02	2	7 ² / ₃ " (19.5 cm)	Small hand, arm, leg, foot
NONNET03	3	9 ¹ / ₃ " (23.7 cm)	Med hand, arm, leg, foot
NONNET04	4	11 ² / ₃ " (29.6 cm)	Lrg hand, arm, leg, foot
NONNET05	5	12 ³ / ₈ " (31.4 cm)	Small head, shoulder, thigh
NONNET06	6	24" (61 cm)	Lrg head, shoulder, thigh
NONNET07	7	28 ² / ₃ " (72.8 cm)	Small chest, back, perineum, axilla
NONNET08	8	31 ¹ / ₃ " (79.6 cm)	Med chest, back, perineum, axilla
NONNET09	9	35" (88.9 cm)	Lrg chest, back, perineum, axilla
NONNET10	10	36 ² / ₃ " (93.1 cm)	XL chest, back, perineum, axilla
NONNET11	11	40" (101.6 cm)	XXL chest, back, perineum, axilla
NONNET12	12	56" (142.2 cm)	XXXL chest, back, perineum, axilla
NONNET22	22	80" (203.2 cm)	Bariatric Size

TIP! If you lose your applicator, simply take a 20cc syringe (without the needle), remove the plunger and place the netting over the syringe. Place the digit inside the syringe and slowly move netting on to the digit.

TIP! Measure circumference of appendage to select correct product size.

Favored by wound care nurses everywhere. This tubular retainer net holds dressings firmly in place without compression. It is an essential product for use with fragile, geriatric, pediatric, and at-risk skin, or those allergic to tape.

- Rayon/polyester knitted elastic net tube holds dressings firmly in place and reduces dressing time
- Soft, open weave does not constrict
- Easily cut to customize length
- Accommodates virtually any part of the body from finger to double extra large torso
- Latex-free and neoprene free

NONNET924P Pre-Cut Elastic Net



Pre-Cut Elastic Net

Individually wrapped pre-cut sizes of Elastic Net provide ease of use, ideal for home health settings.

Individual 24" Lengths (Relaxed)

Item#	Size	Maximum Circumference Stretch Width	Number per Case
NONNET124P	1	6" (15.2 cm)	25
NONNET224P	2	7 ² / ₃ " (19.5 cm)	25
NONNET324P	3	9 ¹ / ₃ " (23.7 cm)	25
NONNET424P	4	11 ² / ₃ " (29.6 cm)	25
NONNET524P	5	12 ³ / ₈ " (31.4 cm)	25
NONNET624P	6	17" (43.2 cm)	25
NONNET724P	7	28 ² / ₃ " (72.8 cm)	25
NONNET824P	8	31 ¹ / ₃ " (79.6 cm)	25
NONNET924P	9	35" (88.9 cm)	25
NONNET1024P	10	41" (104.1 cm)	25
NONNET1124P	11	46" (116.8 cm)	25
NONNET1224P	12	56" (142.2 cm)	25
NONNET2224P	22	80" (203.2 cm)	25

Individual 36" Lengths (Relaxed)

NONNET636P	6	17" (43.2 cm)	50
NONNET736P	7	28 ² / ₃ " (72.8 cm)	50

Elastic Bandages

A complete selection for all your needs



Matrix Elastic Bandage

Matrix is our finest elastic bandage. Ideal for low compression situations, the polyester/cotton weave provides excellent elasticity and recoverability. Features an easy-to-use hook and loop closure at both ends and is color-coded for easy size identification on both ends of roll. Latex-Free.

Sterile (Singles)

Item#	Size	Quantity
DYNJ05152LF	2" x 5 yds, stretched (5.08 cm x 4.57 m)	20/cs
DYNJ05153LF	3" x 5 yds, stretched (7.62 cm x 4.57 m)	20/cs
DYNJ05154LF	4" x 5 yds, stretched (10.16 cm x 4.57 m)	20/cs
DYNJ05156LF	6" x 5 yds, stretched (15.2 cm x 4.57 m)	20/cs
DYNJ05157LF	6" x 10 yds, stretched (15.2 cm x 9.144 m)	20/cs
DYNJ05158LF	6" x 15 yds, stretched (15.2 cm x 13.716 m)	20/cs

Non-Sterile

Item#	Size	Quantity
MDS087002LF	2" x 5 yds, stretched (5.08 x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS087003LF	3" x 5 yds, stretched (7.62 x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS087004LF	4" x 5 yds, stretched (10.16 x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS087006LF	6" x 5 yds, stretched (15.2 x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS087104LF	4" x 10 yds, stretched (10.16 cm x 9.144 m)	4/bx, 5 bx/cs, 20/cs
MDS087106LF	6" x 10 yds, stretched (15.2 cm x 9.144 cm)	4/bx, 5 bx/cs, 20/cs
MDS087156LF	6" x 15 yds, stretched (15.2 cm x 13.716 m)	4/bx, 5 bx/cs, 20/cs

Soft-Wrap® Elastic Bandage

Knitted, breathable, suited for general use support and compression. This 90% cotton and 10% polyester bandage serves as an ideal secondary bandage. Clips included. Latex-Free.

Sterile (Singles)

Item#	Size	Quantity
DYNJ05125	2" x 5 yds, stretched, beige (5.08 cm x 4.57 m)	20/cs
DYNJ05126	3" x 5 yds, stretched, beige (7.62 cm x 4.57 m)	20/cs
DYNJ05131	4" x 5 yds, stretched, beige (10.16 cm x 4.57 m)	20/cs
DYNJ05132	6" x 5 yds, stretched, beige (15.2 cm x 4.57 m)	20/cs

Non-Sterile

Item#	Size	Quantity
MDS046002	2" x 5 yds, stretched, beige (5.08 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS046003	3" x 5 yds, stretched, beige (7.62 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS046004	4" x 5 yds, stretched, beige (10.16 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS046006	6" x 5 yds, stretched, beige (15.2 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs

MDS077003 Swift-Wrap Non-Sterile Elastic Bandage*Swift-Wrap® Elastic Bandage*

High quality elastic bandage with convenient hook and loop closure, eliminating the need for clips or tape. Cotton/ polyester knit construction, ideal for high compression application on a variety of conditions. Latex-Free.

Sterile (Singles)

Item#	Size	Quantity
DYNJ05144	2" x 5 yds, stretched, white (5.08 cm x 4.57 m)	20/cs
DYNJ05145	3" x 5 yds, stretched, white (7.62 cm x 4.57 m)	20/cs
DYNJ05146	4" x 5 yds, stretched, white (10.16 cm x 4.57 m)	20/cs
DYNJ05147	6" x 5 yds, stretched, white (15.2 cm x 4.57 m)	20/cs

Non-Sterile

Item#	Size	Quantity
MDS077002	2" x 5 yds, stretched, white (5.08 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS077003	3" x 5 yds, stretched, white (7.62 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS077004	4" x 5 yds, stretched, white (10.16 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS077006	6" x 5 yds, stretched, white (15.2 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs

MDS055004 Sure-Wrap Non-Sterile Elastic Bandage*Sure-Wrap® Elastic Bandage*

Sure-Wrap is our most cost-effective elastic bandage. It delivers consistent pressure and is appropriate for applications requiring higher compressions. Clips included. 10% Cotton/90% Polyester. Latex-Free.


Non-Sterile (White)

Item#	Size	Quantity
MDS057002	2" x 5 yds, stretched, white (5.08 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS057003	3" x 5 yds, stretched, white (7.62 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS057004	4" x 5 yds, stretched, white (10.16 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS057006	6" x 5 yds, stretched, white (15.20 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs

Non-Sterile (Beige)

Item#	Size	Quantity
MDS055002	2" x 5 yds, stretched, beige (5.08 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS055003	3" x 5 yds, stretched, beige (7.62 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS055004	4" x 5 yds, stretched, beige (10.16 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs
MDS055006	6" x 5 yds, stretched, beige (15.2 cm x 4.57 m)	10/bx, 5 bx/cs, 50/cs

Impregnated Gauze



Petrolatum Gauze Dressing

An effective, non-adherent primary dressing that maintains a moist wound bed for autolytic healing. U.S.P. cotton gauze is evenly impregnated with U.S.P. white petrolatum to create a soothing, conforming dressing that won't stick to the wound site. Ideal for use on non-draining wounds including tube sites, stomas, skin grafts, circumcision, burns, umbilical bandage, non-infected wounds, pressure ulcers (stage 2-4), fistulas or abrasions. Latex-Free.

Sterile

Item#	Size	Quantity
NON251136	1" x 36"	12ea/bx, 12bx/cs, 144/cs
NON251136Z	1" x 36" (2.54 x 4.57 cm)	12ea/bx
NON251180	1" x 8"	12ea/bx, 12bx/cs, 144/cs
NON251180Z	1" x 8" (2.54 x 20.3 cm)	12ea/bx
NON251272	1/2" x 72"	12ea/bx, 12bx/cs, 144/cs
NON251272Z	1/2" x 72" (1.3 x 182.9 cm)	12ea/bx
NON251318	3" x 18"	12ea/bx, 12bx/cs, 144/cs
NON251318Z	3" x 18" (7.62 x 45.7 cm)	12ea/bx
NON251336	3" x 36"	12ea/bx, 12bx/cs, 144/cs
NON251336Z	3" x 36" (7.62 x 91.4 cm)	12ea/bx
NON251390	3" x 9"	12ea/bx, 12bx/cs, 144/cs
NON251390Z	3" x 9" (7.62 x 22.9 cm)	12ea/bx
NON251636	6" x 36"	12ea/bx, 12bx/cs, 144/cs
NON251636Z	6" x 36" (15.2 x 91.4 cm)	12ea/bx

Xeroform Petrolatum Gauze Dressing

Petrolatum gauze with the addition of 3% Xeroform (Bismuth Tribromophenate), a deodorizing agent. Ideal for use on non-draining to lightly draining wounds including donor sites, colostomies, cancer resections, burns abrasions, skin graft sites, and all surgical incisions. Contraindicated for heavily draining wounds. Latex-Free.

Sterile

Item#	Size	Quantity
NON253180	1" x 8"	12ea/bx, 12bx/cs, 144/cs
NON253180Z	1" x 8" (2.54 x 20.3 cm)	12ea/bx
NON253220	2" x 2"	12ea/bx, 12bx/cs, 144/cs
NON253220Z	2" x 2" (5.08 x 5.08 cm)	12ea/bx
NON253440	4" x 4"	12ea/bx, 12bx/cs, 144/cs
NON253440Z	4" x 4" (10.16 x 10.16 cm)	12ea/bx
NON253590	5" x 9"	12ea/bx, 12bx/cs, 144/cs
NON253590Z	5" x 9" (12.7 x 22.9 cm)	12ea/bx

NON250330 3" x 3" Oil Emulsion Dressing

Oil Emulsion Dressing, Non-Occlusive

Compare to Adaptic®* Non-Adhering Dressing

Knitted, high porosity cellulose acetate is evenly impregnated with U.S.P. white petrolatum in an emulsion blend that permits the flow of exudates without adhering to granulating tissue and minimizes skin maceration. Highly conformable dressing is ideal for most draining wounds including burns, suture lines, lacerations, skin ulcers, episiotomy incisions and abrasions. To be used directly on the wound, beneath an absorbent secondary dressing. Latex-Free.

Sterile

Item#	Size	Quantity
NON250330	3" x 3"	50/bx, 12 bx/cs, 600/cs
NON250330Z	3" x 3" (7.62 x 7.62 cm)	50/bx, 12 bx/cs, 600/cs
NON250381	3" x 8"	24/bx, 6 bx/cs, 144/cs
NON250381Z	3" x 8" (7.62 x 20.3 cm)	24/bx, 6 bx/cs, 144/cs
NON250317	3" x 16"	36/bx, 6 bx/cs, 216/cs
NON250317Z	3" x 16" (7.62 x 40.6 cm)	36 pk/bx
NON250383	3" x 8"	36 packs of 3 dressings/bx, (7.62 x 20.3 cm) (3/pkg) 6 bx/cs, 216 pks/cs

*This product is not manufactured or distributed by Johnson & Johnson Corporation, distributors of Adaptic non-adherent dressings.

NON255015 1" x 5 yd Plain Packing Strips**NON256025** 1" x 5 yd Iodoform Packing Strips

Packing Strips

These tightly woven, 100% cotton gauze packing strips are ravel resistant, and won't fray or lint. Used for packing or as drainage conduits in wound care management applications, such as nasal or sinus packing, boils, abscesses, fistulas and other draining or tunneling wounds. Strips are available plain or impregnated with iodoform, an antiseptic agent. Packaged as one 5 yard strip per bottle. Cut to desired length, won't fray. Latex-Free.

Sterile, Plain

Item#	Size	Quantity
NON255015	1" x 5 yd (2.54 cm x 4.57 m)	12/cs
NON255025	2" x 5 yd (5.08 cm x 4.57 m)	12/cs
NON255125	1/2" x 5 yd (1.3 cm x 4.57 m)	12/cs
NON255145	1/4" x 5 yd (0.6 cm x 4.57 m)	12/cs

Sterile, Iodoform

Item#	Size	Quantity
NON256015	1" x 5 yd (2.54 cm x 4.57 m)	12/cs
NON256025	2" x 5 yd (5.08 cm x 4.57 m)	12/cs
NON256125	1/2" x 5 yd (1.3 cm x 4.57 m)	12/cs
NON256145	1/4" x 5 yd (0.6 cm x 4.57 m)	12/cs

Adhesive Tapes

Re-engineered for just the right adhesion

NON260001 1" x 10 yds Medfix Paper Tape



NON260101 1" x 10 yds Medfix Cloth Tape



Paper Tape Compare to Micropore®

A gentle, breathable tape for sensitive skin, recommended especially for pediatric and geriatric patients. Flexible and conformable. Adheres well to damp skin. Also available in tan. Latex-Free.

Item#	Size	Quantity
NON260012	1/2" x 10 yds (1.3 cm x 9.14 m)	24 rl/bx, 12 bx/cs, 288 rl/cs
NON260001	1" x 10 yds (2.54 cm x 9.14 m)	12 rl/bx, 12 bx/cs, 144 rl/cs
NON260002	2" x 10 yds (5.08 cm x 9.14 m)	6 rl/bx, 12 bx/cs, 72 rl/cs
NON260003	3" x 10 yds (7.62 cm x 9.14 m)	4 rl/bx, 12 bx/cs, 48 rl/cs

Single Use Rolls

Item#	Size	Quantity
NON260001S	1" x 1.5 yds (2.54 cm x 1.37 m)	100 rl/bx, 5 bx/cs, 500 rl/cs
NON260002S	2" x 1.5 yds (5.08 cm x 1.37 m)	50 rl/bx, 5 bx/cs, 250 rl/cs

Cloth "Silk" Tape Compare to Durapore®

With its strong woven construction, this cloth "silk-like" tape is perfect for affixing dressings. It's water repellent and easy to tear. Latex-Free.

Item#	Size	Quantity
NON260112	1/2" x 10 yds (1.3 cm x 9.14 m)	24 rl/bx, 12 bx/cs, 288 rl/cs
NON260101	1" x 10 yds (2.54 cm x 9.14 m)	12 rl/bx, 12 bx/cs, 144 rl/cs
NON260102	2" x 10 yds (5.08 cm x 9.14 m)	6 rl/bx, 12 bx/cs, 72 rl/cs
NON260103	3" x 10 yds (7.62 x 9.14 m)	4 rl/bx, 12 bx/cs, 48 rl/cs

Also Available:

Elastic Foam Tape (compare to Microfoam)	NON260601/2/3/4
Waterproof Tape (compare to Wetpruf)	NON260501/2/3
Elastic Adhesive Bandage (compare to Elastoplast Elastic Adhesive Bandage)	NON260401/2/3



Please consult our **Surgical Tape Brochure** for information on our comprehensive line of tapes.

NON260201 1" x 10 yds Medfix Transparent Tape



Transparent Tape Compare to Transpore™

A breathable perforated plastic tape that permits skin examination without tape removal. Provides a clean and easy bilateral tear. Excellent tape to hold facial dressings or for IV sites and tubing retention. Latex-Free.

Item#	Size	Quantity
NON260212	1/2" x 10 yds (1.3 cm x 9.14 m)	24 rl/bx, 12 bx/cs, 288 rl/cs
NON260201	1" x 10 yds (2.54 cm x 9.14 m)	12 rl/bx, 12 bx/cs, 144 rl/cs
NON260202	2" x 10 yds (5.08 cm x 9.14 m)	6 rl/bx, 12 bx/cs, 72 rl/cs
NON260203	3" x 10 yds (7.62 cm x 9.14 m)	4 rl/bx, 12 bx/cs, 48 rl/cs

Cotton Adhesive Tape NON260701/2/3
(compare to Cloth Adhesive Tape)

Ortho Porous Sports Tape NON260301/2/3
(compare to Zonas®)

Megazinc Pink™ Adhesive OMAM111
(compare to Hy-Tape®)

Medfix EZ MSC4102
(compare to Medipore H Soft Cloth)

Medi-Strips® Reinforced Wound Closure Strips

Adhesive wound closures with an exceptionally porous spun-bond fiber backing, includes reinforced filaments for high tensile strength. The medical grade adhesive affords excellent adhesion while being gentle to the skin. Medi-Strips are used in addition to, or in place of, sutures. Each size is in a color-coded package for quick and easy size identification. Latex-Free.

NON250314 1/4" x 3"

Medi-Strips Reinforced Wound Closure Strips



Reinforced, Sterile

Item#	Size	Quantity
NON250318	1/8" x 3" (0.32 x 7.62 cm)	5/pk, 50 pk/bx, 4 bx/cs
NON250114	1/4" x 1 1/2" (0.6 x 3.8 cm)	6/pk, 50 pk/bx, 4 bx/cs
NON250412	1/2" x 4" (1.3 x 10.16 cm)	6/pk, 50 pk/bx, 4 bx/cs
NON250314	1/4" x 3" (0.6 x 7.62 cm)	3/pk, 50 pk/bx, 4 bx/cs
NON250414	1/4" x 4" (0.6 x 10.16 cm)	5 strips/card, 2 cards/pk, 50 pk/bx, 4 bx/cs
NON250501	1" x 5" (2.54 x 12.7 cm)	4/pk, 25 pk/bx, 4 bx/cs

Medfix Adhesive Moleskin



A roll of tan colored flannel with a napped cotton backing on a zinc oxide-based adhesive. Used to cushion the skin in an effort to help manage blisters, corns, calluses and chafing. Cut to desired size. Contains natural rubber latex.

Item#	Size	Quantity
NON5147	9" x 4 yds (22.9 cm x 3.66 m)	1 roll/tube, 12 tubes/cs

Self-Adherent Wraps

- Self-adherent, cohesive wrap securely sticks to itself and not to skin or other materials and will not slip or loosen over time.
- Easily tears by hand for custom lengths. No need for scissors.
- Porous construction allows ventilation of the wrapped area.
- Retains dressings without tension around digits or limbs.
- Available sterile/non-sterile and with/without latex.
- Each bandage is individually wrapped.

Co-Flex® NL

This is a latex-free self-adherent, hand tear compression wrap.

Sterile

Item#	Size	Quantity
DYNJ05123	2" x 5 yds, stretched, tan (5.08 cm x 4.57 m)	20 rolls/cs
DYNJ05124	3" x 5 yds, stretched, tan (7.62 cm x 4.57 m)	20 rolls/cs
DYNJ05129	4" x 5 yds, stretched, tan (10.16 cm x 4.57 m)	20 rolls/cs
DYNJ05130	6" x 5 yds, stretched, tan (15.2 cm x 4.57 m)	20 rolls/cs

Non-Sterile

Item#	Size	Quantity
MDS088001	1" x 5 yds, stretched, tan (2.54 cm x 4.57 m)	30 rolls/cs
MDS088002	2" x 5 yds, stretched, tan (5.08 cm x 4.57 m)	36 rolls/cs
MDS088003	3" x 5 yds, stretched, tan (7.62 cm x 4.57 m)	24 rolls/cs
MDS088004	4" x 5 yds, stretched, tan (10.16 cm x 4.57 m)	18 rolls/cs
MDS088006	6" x 5 yds, stretched, tan (15.2 cm x 4.57 m)	12 rolls/cs

Co-Flex® Med

This is a latex based self-adherent, hand-tear compression wrap.

Sterile

Item#	Size	Quantity
DYNJ086002	2" x 5 yds, stretched, tan (5.08 cm x 4.57 m)	20 rolls/cs
DYNJ086003	3" x 5 yds, stretched, tan (7.62 cm x 4.57 m)	20 rolls/cs
DYNJ086004	4" x 5 yds, stretched, tan (10.16 x 4.57 m)	20 rolls/cs
DYNJ086006	6" x 5 yds, stretched, tan (15.2 cm x 4.57 m)	20 rolls/cs

Non-Sterile

Item#	Size	Quantity
MDS086001	1" x 5 yds, stretched, tan (2.54 cm x 4.57 m)	30 rolls/cs
MDS086002	2" x 5 yds, stretched, tan (5.08 cm x 4.57 m)	36 rolls/cs
MDS086003	3" x 5 yds, stretched, tan (7.62 cm x 4.57 m)	24 rolls/cs
MDS086004	4" x 5 yds, stretched, tan (10.16 cm x 4.57 m)	18 rolls/cs
MDS086006	6" x 5 yds, stretched, tan (15.2 cm x 4.57 m)	12 rolls/cs



Kid's Pack Contains Rolls of:

- Neon Pink/Purple Hearts
- Neon Green/Dark Green Dinosaurs
- Orange/Black Smiley Face
- Light Blue/Dark Blue Trucks and Cars
- Yellow/Black Bumblebees
- Lilac/Purple Flowers


Colorpack Contains

Rolls of: Pink, Blue, Purple, Light Blue, Neon Green, and Red



Non-Sterile, Latex Free

Item#	Description	Quantity
MDS088001CP	1" x 5 yds, stretched, Colorpack (2.54 cm x 4.57 m)	30 rolls/cs
MDS088002CP	2" x 5 yds, stretched, Colorpack (5.08 cm x 4.57 m)	36 rolls/cs
MDS088003CP	3" x 5 yds, stretched, Colorpack (7.62 cm x 4.57 m)	24 rolls/cs
MDS088004CP	4" x 5 yds, stretched, Colorpack (10.16 cm x 4.57 m)	18 rolls/cs
MDS088002KP	2" x 5 yds, stretched, Kid's pack (5.08 cm x 4.57 m)	36 rolls/cs
MDS088003KP	3" x 5 yds, stretched, Kid's pack (7.62 cm x 4.57 m)	24 rolls/cs



Co-Flex[®] LF² latexfree PRODUCT

A unique quick stick self-adhering bandage with a soft foam layer. Latex-Free but performs better than latex. Hand tear.


Sterile

Item#	Size	Quantity
DYNJ089002	2" x 5 yds (5.08 cm x 4.57m), stretched, tan	20 rolls/cs
DYNJ089003	3" x 5 yds (7.62 cm x 4.57 m), stretched, tan	20 rolls/cs
DYNJ089004	4" x 5 yds (10.16 cm x 4.57 m), stretched, tan	20 rolls/cs
DYNJ089006	6" x 5 yds (15.2 cm x 4.57 m), stretched, tan	20 rolls/cs

Non-Sterile

Item#	Size	Quantity
MDS089001	1" x 5 yds (2.54 cm x 4.57 m), stretched, tan	36 rolls/cs
MDS089002	2" x 5 yds (5.08 cm x 4.57 m), stretched, tan	30 rolls/cs
MDS089003	3" x 5 yds (7.62 cm x 4.57 m), stretched, tan	24 rolls/cs
MDS089004	4" x 5 yds (10.16 cm x 4.57 m), stretched, tan	18 rolls/cs
MDS089006	6" x 5 yds (15.2 cm x 4.57 m), stretched, tan	12 rolls/cs
MDS089001CP	1" x 5 yds, stretched, Colorpack (2.54 cm x 4.57 m)	36 rolls/cs
MDS089002CP	2" x 5 yds, stretched, Colorpack (5.08 cm x 4.57 m)	30 rolls/cs
MDS089003CP	3" x 5 yds, stretched, Colorpack (7.62 cm x 4.57 m)	24 rolls/cs
MDS089004CP	4" x 5 yds, stretched, Colorpack (10.16 cm x 4.57 m)	18 rolls/cs

NON5129 Montgomery Straps



! TIP! Laces are not included. Plain packing strips or gauze can be used and purchased separately.

Montgomery Straps

Medfix Montgomery Straps are adhesive straps used to facilitate frequent dressing changes without having to remove and reapply tape. Breathable, strong and comfortable. Helps reduce the risk of patient skin trauma associated with frequent taping. Can cut to size. Latex-Free.

Item#	Size	Quantity
NON5129	7 ¹ / ₄ " x 11 ¹ / ₈ "	24 straps per box (12 pairs), (18.4 x 28.3 cm) 6 bx/cs, 144 straps per case

Adhesive Surgical Dressing

This all-in-one dressing helps reduce nursing time by assembling in one dressing a sterile sponge, non-adhesive pad and cloth "silk" tape. A non-adherent contact layer that permits free flow of drainage, protects the wound site, and will not stick to granulating tissue. A soft, highly absorbent rayon/cellulose pad for light to medium draining wounds. The tape backing holds the pad in place on two sides. Easily cut to customized sizes. Ideal for long surgical incisions, lacerations, abrasions and invasive lines. Latex-Free.

Sterile 1's

Item#	Size	Quantity
NON4311	4" x 6" (4" x 3" pad) (10.16 x 15.2 cm)	25/bx, 8 bx/cs, 200 ea/cs
NON4313	8" x 6" (8" x 3" pad) (20.3 x 15.2 cm)	25/bx, 4 bx/cs, 100 ea/cs

NON4311 Adhesive Surgical Dressing



NON216645 6" x 6" Sterile Burn Dressing



Burn Dressing

45-Ply Cotton Dressing Provides Loft and Extra Cushioning of the Wound

100% Cotton wide mesh dressings absorb wound drainage and cushion the site. Specifically designed for use on burns. Available in a variety of sizes, ply, sewn or unsewn edges. Sterile and Latex-Free.

Sterile

Item#	Size	Quantity
NON7911	18" x 18" 10-ply Sewn (45.7 x 45.7 cm)	1 ea/pk, 80 pk/cs
NON7912	18" x 18" 2-ply Sewn (45.7 x 45.7 cm)	2 ea/pk, 100 pk/cs
NON7913	18" x 18" 4-ply Sewn (45.7 x 45.7 cm)	5 ea/pk, 40 pk/cs
NON7915	18" x 18" 50-ply Unsewn (45.7 x 45.7 cm)	1 ea/pk, 20 pk/cs
NON7916	36" x 36" 18-ply Unsewn (91.4 x 91.4 cm)	1 ea/pk, 20 pk/cs
NON216645	6" x 6" 45-ply Sewn (15.2 x 15.2 cm)	6 ea/tray, 20/tray/cs

NON7911 Fine Mesh Layered Burn Dressing



NONTRI200 36" x 36" Triangular Bandage



Triangular Bandage

- Non-sterile individually packaged, Latex-Free
- Includes 2 safety pins
- Conveniently unfolds for arm sling
- Ideal for EMS and First Aid Kits

Item#	Size	Quantity
NONTRI200	36" x 36" x 51" (91.4 x 91.4 x 129.5 cm)	12/bg, 240/cs
NONTRI300	40" x 40" x 56" (101.6 x 101.6 x 142.2 cm)	12/bg, 240/cs

Nursing Pad

This disposable, nursing pad is soft and highly absorbent. Helps wick fluids away from the skin. Latex-Free.

Item#	Size	Quantity
HSPBP2412	Disposable nursing pad	12/bx, 288/cs

Sterile Disposable Drape

Create a sterile field with a sterile paper drape. Repellent plastic liner helps prevent strike through. Latex-Free.

Item#	Size	Quantity
NON21001	18" x 26" Paper drape (45.7 x 66 cm)	50/bx, 300/cs
NON21002	18" x 26" Paper drape (45.7 x 66 cm) 2.75" Fenestration	50/bx, 300/cs

MDS066004 WyTex Undercast Padding*WyTex Undercast Padding*

This product is constructed of 100% cotton fibers to provide a soft environment under a hard cast. WyTex is a low linting product with more loft to provide a better feel for the patient. Latex-Free.

Sterile (Singles)

Item#	Size	Quantity
DYNJ066002	2" x 4 yds (5.08 cm x 3.66 m)	20 rolls/cs
DYNJ066003	3" x 4 yds (7.62 cm x 3.66 m)	20 rolls/cs
DYNJ066004	4" x 4 yds (10.16 cm x 3.66 m)	20 rolls/cs
DYNJ066006	6" x 4 yds 15.2 cm x 3.66 m	20 rolls/cs

Non-Sterile

Item#	Size	Quantity
MDS066002	2" x 4 yds (5.08 cm x 3.66 m)	24/bg, 3 bg/cs, 72/cs
MDS066003	3" x 4 yds (7.62 cm x 3.66 m)	12/bg, 6 bg/cs, 72/cs
MDS066004	4" x 4 yds (10.16 cm x 3.66 m)	12/bg, 6 bg/cs, 72/cs
MDS066006	6" x 4 yds (15.2 cm x 3.66 m)	6/bg, 6/cs, 36/cs

MDS067002 Syn-Tex Undercast Padding*Syn-Tex Undercast Padding*

These synthetic cast padding rolls are ideal for use with all forms of synthetic hard casting material. 100% polyester needle punched construction absorbs excess moisture and dries quickly to retain a comfortable fit. Latex-Free.

Non-Sterile

Item#	Size	Quantity
MDS067002	2" x 4 yds (5.08 cm x 3.66 m)	12/bg, 6bg/cs, 72/cs
MDS067003	3" x 4 yds (7.62 cm x 3.66 m)	12/bg, 6bg/cs, 72/cs
MDS067004	4" x 4 yds (10.16 cm x 3.66 m)	12/bg, 6bg/cs, 72/cs
MDS067006	6" x 4 yds (15.2 cm x 3.66 m)	6/bg, 6/cs, 36/cs

Arm Sleeve

This protective arm sleeve is designed to provide an extra layer of protection to at-risk skin. Available by the pair. Latex-Free.

Item#	Size	Quantity
NONSLEEVE	14" long with thumb loop (35.6 cm)	1 pair/pk

Adhesive Bandages

A complete line of bandages for every need.

Jr., Plastic
NON25509



3/4" Flex-Fabric
NON25650



3/4" Plastic
NON25500



Finger Tip
Flex Fabric
NON25513



Spot, Plastic
NON25501

X-Large, Plastic
NON25504



Knuckle,
Flex-Fabric
NON25510



Sheer-Gard[®] Adhesive Bandages

Sheer-Gard adhesive bandages are constructed with PVC-free durable polyethylene and perforated for breathability. These bandages feature even adhesion and a non-adherent island pad. Latex-Free.



Plastic, Sterile

Item#	Size	Quantity
NON25500	3/4" x 3" (1.9 x 7.62 cm)	100/bx, 12 bx/cs, 1200/cs
NON25600	1" x 3" (2.54 x 7.62 cm)	100/bx, 12 bx/cs, 1200/cs
NON25501	Spot 7/8" (2.2 cm)	100/bx, 12 bx/cs, 1200/cs
NON25504	X-Lg 2" x 4" (5.08 x 10.16 cm)	50/bx, 12 bx/cs, 600/cs
NON25509	Jr. 3/8" x 1 1/2" (1 x 3.8 cm)	100/bx, 36 bx/cs, 3600/cs

Comfort Cloth[®] Adhesive Bandages

Comfort Cloth adhesive bandages are woven to allow them to stretch and conform to the wound. Material naturally breathes. Features even adhesion and a non-adherent island pad. Latex-Free.



Woven Cloth, Sterile

Item#	Size	Quantity
NON25650	3/4" x 3" (1.9 x 7.62 cm)	100/bx, 12 bx/cs, 1200/cs
NON25660	1" x 3" (2.54 x 7.62 cm)	100/bx, 12 bx/cs, 1200/cs
NON25510	Knuckle, 3" x 1 1/2" (7.62 x 3.8 cm)	100/bx, 12 bx/cs, 1200/cs
NON25513	Finger tip, 2 1/8" x 1 1/2" (5.4 x 3.8 cm)	100/bx, 12 bx/cs, 1200/cs
NON25650B	3/4" x 3" (1.9 x 7.62 cm)	2025/bx, 4 bx/cs, 8100/cs
NON25660B	1" x 3" (2.54 x 7.62 cm)	1800/bx, 4 bx/cs, 7200/cs

Introducing Medtoons®

Children will adore the smiling bumble bees, soccer playing dogs and snorkeling fish featured on the new Medtoons Cartoon adhesive bandage. Each box includes 50 strips (25 strips each of yellow, and blue) constructed from PVC-free vinyl perforated for increased moisture vapor transmission. An absorbent, non-adherent island pad seals around the wound site to keep out dirt and helps protect granulating tissue. The pressure sensitive, acrylic adhesive is skin friendly and provides long lasting protection. Presented in a convenient dispenser box. Latex-Free.

Sterile

Item#	Size	Quantity
NON256130	¾" x 3" (1.9 x 7.62 cm)	50/box, 24 bx/cs, 1200/cs

NON256130 Cartoon Adhesive Bandage



Pressure Adhesive Bandage

Once in contact with fluid, the absorbent, cellulose island pad expands to over 7 times its normal height helping to provide gentle pressure to stop bleeding. Waterproof PE backing prevents strike through. Saves nursing time and promotes comfort for dialysis, chemotherapy and cardio cath patients. Latex-Free.

Sterile

Item#	Size	Quantity
NON85100	Large, 1" x 2¾" (2.54 x 7 cm)	100/bx, 20 bx/cs, 2000/cs
NON85200	X-Large, 1¼" x 2¾" (3.2 x 7 cm)	100/bx, 20 bx/cs, 2000/cs

NON85100 Large Pressure Adhesive Bandage

unactivated



activated



latexfree
PRODUCT

Neon Adhesive Bandage

This bright assortment of five neon colors includes 10 plastic strips each of pink, blue, green, purple and yellow in each box. The non-adherent island pad protects granulating tissue, and the adhesive provides long lasting protection. Latex-Free.


Sterile

Item#	Size	Quantity
NON256131	¾" x 3" (1.9 x 7.62 cm)	50/bx, 24 bx/cs, 1200/cs

NON256131 Neon Adhesive Bandage



latexfree
PRODUCT



CURAD® Plastic Strips Are Ideal for Everyday Use

CURAD® Plastic Strips are ideal multi-purpose bandages for the whole family. Different from other plastic strips, our strips are made from flexible plastic reinforced with a layer of non-woven material for softness and strength. A skin-friendly adhesive doesn't leave a sticky, itchy residue on the skin, so you are ensured maximum comfort as you heal. Sterile.

Item#	Description	Pkg
CUR02278	80 - 3/4" x 2 7/8" (1.9 x 7.3 cm)	80 ea/bx, 24 bx/cs
CU02278	80 - 3/4" x 2 7/8" (1.9 x 7.3 cm)	80 ea/bx, 24 bx/cs
CUR45157	20 - 7/8" Spot (2.2 cm) 20 - 3/4" x 2 1/8" (1.9 x 5.4 cm) 30 - 3/4" x 2 7/8" (1.9 x 7.3 cm) 10 - 1 1/8" x 2 7/8" (2.9 x 7.3 cm)	80 ea/bx, 24 bx/cs
CUR43021	20 - 1" x 3 1/4" (2.54 x 8.2 cm)	20 ea/bx, 24 bx/cs
CUR44010	30 - 3/4" x 3" (1.9 x 8.2 cm)	30 ea/bx, 24 bx/cs

CURAD® Extreme Lengths and Extreme Hold Bandages Are Designed to Protect Hardworking Hands

Curad® Extreme Lengths bandages are extra long for extra hold. Unlike other bandages, they wrap around the finger twice to stay on strong! Made with flex-fabric and sure stick adhesive, this finger bandage stays in place! Curad® Extreme – Stay on strong to protect your hardworking hands! Sterile.

Curad® Extreme Hold are bandages designed to protect hardworking hands. They are made with a heavy-duty fabric for maximum protection, and the sure stick adhesive means the bandage will stay on until you want it off! Curad® Extreme Hold has a bandage for every size cut or scrape. There are special shaped bandages for knuckles and fingertips, as well as regular and extra large bandages. It is the ideal assortment for your hardworking life. Sterile.

Item#	Description	Pkg
CUR01101	3/4" x 4 3/4" (1.9 x 12.1 cm)	20 ea/bx, 24 bx/cs
CU01101	3/4" x 4 3/4" (1.9 x 12.1 cm)	20 ea/bx, 24 bx/cs
CUR14924	4 - 2" x 3 3/4" (5.08 x 9.5 cm), 4 Knuckle, 4 Fingertip, 18 - 7/8" x 3" (2.2 x 7.62 cm)	30 ea/bx, 24 bx/cs
CUR14926	10 - 3 3/4" x 2" (9.5 x 5.08 cm)	10 ea/bx, 24 bx/cs
CUR14925	20 - 1" x 3" (2.5 x 7.62 cm)	20 ea/bx, 24 bx/cs

CURAD Hydro Heal® is great for advanced first aid.

These hydrocolloid strips are constructed with a smooth satin backing and a tapered edge for longer wear time. Sterile.

Item#	Description	Pkg
CUR076123	6 - 0.82" x 1.6" (2.1 x 4.2 cm)	6 ea/bx, 24 bx/cs

CURAD® Sheer Strips Are the Ideal Bandage for Less Visible Protection

Our sheer strips are designed to blend with your skin for a more appealing look and help you hide the bandage while you heal. They are also made with a skin-friendly adhesive that doesn't leave a sticky, itchy residue on the skin, so you are ensured maximum comfort as you heal. Sterile.

Item#	Description	Pkg
CUR02279	80 - 3/4" x 27/8" (1.9 x 7.3 cm)	80 ea/bx, 24 bx/cs
CU02279	80 - 3/4" x 27/8" (1.9 x 7.3 cm)	80 ea/bx, 24 bx/cs
CUR45243	20 - 7/8" round (2.2 cm)	
	20 - 3/4" x 21/8" (1.9 x 5.4 cm)	
	30 - 3/4" x 27/8" (1.9 x 7.3 cm)	
	10 - 11/8" x 27/8" (2.9 x 7.3 cm)	80 ea/bx, 24 bx/cs
CUR02276	10 - 1" x 3" (2.5 x 8.2 cm)	
	20 - 3/4" x 3" (1.9 x 7.62 cm)	
	20 - 5/8" x 21/4" (1.5 x 5.7 cm)	80 ea/bx, 24 bx/cs
CUR02279	10 - 2" x 33/4" (5.08 x 9.5 cm)	80 ea/bx, 24 bx/cs



CURAD® Flex-Fabric™ Strips Are Designed For People With Active Lifestyles Who Require a Comfortable Bandage to “Stay Put” Even When Wet

Our Flex-Fabric™ strips are made of woven fabric designed to stretch, bend, breathe and stay on comfortably, even when wet. The water-repellent coating keeps the bandages dry and secure. CURAD® uses a skin-friendly adhesive that doesn't leave a sticky, itchy residue on the skin, so you are ensured maximum comfort as you heal. Sterile.

Item#	Description	Pkg
CUR45245	3/4" x 21/2" (1.9 x 6.4 cm)	40 ea/bx, 24 bx/cs
CUR40659	3/4" x 3" (1.9 x 7.62 cm)	10 ea/bx, 24 pk/cs
CUR0700	Fabric Assorted	100 ea/bx, 24 bx/cs
CUR45246	10 - Fingertip 13/4" x 2" (44 mm x 51 mm)	
	10 - Knuckle 11/2" x 3" (38 mm x 76 mm)	20 ea/bx, 24 bx/cs

CURAD Sensitive Skin® Strips Are the Ideal Bandage for People With Sensitive Skin

Our sensitive skin adhesive bandages have been specially designed to be gentle to your skin. Our ultralite fabric is breathable, comfortable and hypoallergenic. The absorbent non-stick pad is ouchless to remove. Innovative adhesive stays on longer. Sterile.

Item#	Description	Pkg
CUR45230	50 - 7/8" Spots (2.2 cm)	50 ea/bx, 24 bx/cs
CUR45168	12 - 7/8" Spots (2.2 cm)	
	12 - 3/4" x 21/2" (1.9 x 6.4 cm)	
	6 - 11/8" x 21/2" (2.9 x 6.4 cm)	30 ea/bx, 24 bx/cs
CU45168	12 - 7/8" Spots (2.2 cm)	
	12 - 3/4" x 21/2" (1.9 x 6.4 cm)	
	6 - 11/8" x 21/2" (2.9 x 6.4 cm)	30 ea/bx, 24 bx/cs
CUR45227	30 - 3/4" x 21/2" (1.9 x 6.4 cm)	30 ea/bx, 24 bx/cs
CUR45167	10 - 2" x 27/8" (5.08 x 7.3 cm)	10 ea/bx, 24 bx/cs





CURAD® Ouchless® Non-Stick Pads

CURAD® Ouchless® Non-Stick Pads are wound covers for post surgical wounds, lightly bleeding and draining wounds, large surface cuts, scrapes, and burns. Our non-stick pads are made of a soft, perforated film bonded to a highly absorbent and breathable pad. The film is ventilated to allow the wound to breathe and absorb any fluids. The pad can be cut or trimmed to fit any wound. Our sterile pads are available with or without adhesive tabs.

Item#	Description	Pkg
Sterile pads without adhesive tabs		
CUR47395	2" x 3" (5.08 x 7.62 cm)	10 ea/bx, 12 bx/cs
CUR47396	2" x 3" (5.08 x 7.62 cm)	20 ea/bx, 12 bx/cs
CUR47397	3" x 4" (7.62 x 10.16 cm)	10 ea/bx, 12 bx/cs
CUR47398	3" x 4" (7.62 x 10.16 cm)	20 ea/bx, 12 bx/cs
CUR47399	3" x 8" (7.62 x 20.3 cm)	8 ea/bx, 24 bx/cs
Sterile pads with adhesive tabs		
CUR47146	2" x 3" (5.08 x 7.62 cm)	10 ea/bx, 12 bx/cs
CUR47147	2" x 3" (5.08 x 7.62 cm)	20 ea/bx, 12 bx/cs
CU47147	2" x 3" (5.08 x 7.62 cm)	20 ea/bx, 12 bx/cs
CUR47148	3" x 4" (7.62 x 10.16 cm)	10 ea/bx, 12 bx/cs

Mediplast® is effective against Corn, Calluses and Warts

Our adhesive cut-to-fit pads are soft, convenient and easy to use. The active ingredient used in Mediplast® is Salicylic Acid (40%), the highly effective ingredient widely recommended by physicians and pharmacists for the removal of corns, calluses and common warts. Adhesive contains Latex.

Item#	Description	Pkg
CUR01496	2" x 3" (5.08 x 7.62 cm) pads	25 ea/bx, 6 bx/cs

All Purpose CURAD® Waterproof Tape



Our tape is waterproof, oil resistant and has a cotton cloth backing with polyethylene coating which provides high tensile strength that will not permit moisture to penetrate to the skin or dressing. The strong adhesive is durable and holds securely. Ideal for bandage wraps, picture frame ostomy flanges, protection of dressings for general taping needs and outpatient use.

Item#	Description	Pkg
CUR47440	1/2 x 5 yards (1.3 x 4.57 cm)	1 ea/rl, 24 rl/cs
CUR47441	1/2 x 10 yards (1.3 x 9.14 cm)	1 ea/rl, 24 rl/cs

STERILE GAUZE SPONGES	MEDLINE	J&J	COVIDIEN (Kendall)	CARDINAL HEALTH (Allegiance)
Gauze 2x2, 8-ply Sterile 2's 50 pk/bx, 30 bx/cs (3000/cs)	NON21420	2318	1806	GZ2208-2
Gauze 3x3, 12-ply Sterile 2's 40 pk/bx, 30 bx/cs (2400/cs)	NON21422		1903	N/A
Gauze 4x4, 8-ply Sterile 2's 50 pk/bx, 12 bx/cs (1200/cs)	NON21448	2315	2187	GZ4408-2
Gauze 4x4, 12-ply Sterile 2's 25 pk/bx, 24 bx/cs (1200/cs)	NON21424	2317	3033	GZ4412-2
Gauze 4x4, 12-ply Sterile 10's 10/tray, 128 pk/cs (1280/cs)	NON21426	2322 960/cs	6939	GZ4412-10
Gauze 4x4, 16-ply Sterile 10's 10/tray, 128 pk/cs (1280/cs)	NON21428	2329 720/cs	7605	GZ4416-10
Gauze 4x4, 8-ply Sterile 10's 10/tray, 128 pk/cs (1280/cs)	NON21425	2321	6318	
Gauze 4x8, 12-ply Sterile 2's 25 pk/bx, 20 bx/cs (1000/cs)	NON21842		2259	
NON-STERILE GAUZE SPONGES	MEDLINE	J&J	COVIDIEN	CARDINAL HEALTH
Gauze 2x2, 8-ply N/S 200/pk, 25 pk/cs (5000/cs)	NON25208	7612	2146	GZ2208
Gauze 2x2, 12-ply N/S 200/pk, 40 pk/cs (8000/cs)	NON25212		2252	GZ2212
Gauze 3x3, 12-ply N/S 200/pk, 20 pk/cs (4000/cs)	NON25312		2346	GZ3312
Gauze 4x4, 8-ply N/S 200/pk, 20 pk/cs (4000/cs)	NON25408	7610	2556	GZ4408
Gauze 4x4, 12-ply N/S 200/pk, 10 pk/cs (2000/cs)	NON25412	7623	2634	GZ4412
Gauze 4x4, 16-ply N/S 200/pk, 10 pk/cs (2000/cs)	NON25416		2733	GZ4416
Gauze 8x4, 12-ply N/S 200/pk, 10 pk/cs (2000/cs)	NON25812		2835	GZ8412
X-RAY DETECTABLE GAUZE SPONGES	MEDLINE	J&J	COVIDIEN	CARDINAL HEALTH
Gauze X-Ray Detect, 4x4, 16-ply, Sterile 10's 10/tray, 128 trays/cs (1280/cs)	NON21430LF	2515 720/cs	7317 1350/cs	GT4416-10
Gauze X-Ray Detect, 8x4, 12-ply, Sterile 10's 10/tray, 80 trays/cs (800/cs)	NON21432LF	2520 480/cs	7318	GT8412-10 48 tr/cs, 480/cs
Gauze X-Ray Detect, 4x4, 32-ply, Sterile 10's 10/tray, 48 trays/cs (480/cs)	NON21433		7320	
Gauze X-Ray Detect, 8x4, 24-ply, Sterile 10's 10/tray, 20 trays/cs (200/cs)	NON21438		7321	
COTTON DRAIN SPONGES	MEDLINE			
100% Cotton, 4x4, 16-ply, Sterile 1/pk 25/bx, 12 bx/cs	NON255000			

NON-WOVEN SPONGES - AVANT DELUXE	MEDLINE	J&J "Nu-Gauze"		
Avant Deluxe Gauze Non-Woven 4x4, 4-ply Sterile 2's, 50 pk/bx, 12 bx/cs (1200/cs)	NON264442	2337		
Avant Deluxe Gauze Non-Woven 2x2, 3-ply N/S 200/pk, 40 pk/cs (8000/cs)	NON26223	7635		
Avant Deluxe Gauze Non-Woven 2x2, 4-ply N/S 200/pk, 40 pk/cs (8000/cs)	NON26224	7632		
Avant Deluxe Gauze Non-Woven 3x3, 4-ply N/S 200/pk, 20 pk/cs (4000/cs)	NON26334	7633		
Avant Deluxe Gauze Non-Woven 4x4, 3-ply N/S 200/pk, 10 pk/cs (2000/cs)	NON26443	7636		
Avant Deluxe Gauze Non-Woven 4x4, 4-ply N/S 200/pk, 10 pk/cs (2000/cs)	NON26444	7634		
NON-WOVEN GAUZE SPONGES - AVANT	MEDLINE		COVIDIEN "Versalon"	CARDINAL HEALTH
Avant Gauze Non-Woven, 2x2, 4-ply, Sterile 2's 50 pk/bx, 30 bx/cs (3000/cs)	NON21224		8042	NW2204-2
Avant Gauze Non-Woven, 3x3, 4-ply, Sterile 2's 40 pk/bx, 30 bx/cs (2400/cs)	NON21334		8043	NW3304-2 1200/cs
Avant Gauze Non-Woven, 4x4, 4-ply, Sterile 2's 50 pk/bx, 12 bx/cs (1200/cs)	NON21444		8044	NW4404-2
Avant Gauze Non-Woven, 4x4, 6-ply, Sterile 2's 25 pk/bx, 12 bx/cs (600/cs)	NON21446		7084	NW4406-2
Avant Gauze Non-woven, 4x4, 4 ply Sterile 10's 10/tray 128 trays per case	NON21449		8047	N/A
Avant Gauze Non-Woven, 2x2, 4-ply, Non-Sterile 200/pk, 40 pk/cs (8000/cs)	NON25224		9022	NW2204
Avant Gauze Non-Woven, 3x3, 4-ply, Non-Sterile 200/pk, 20 pk/cs (4000/cs)	NON25334		9023	NW3304 8000/cs
Avant Gauze Non-Woven, 4x4, 4-ply, Non-Sterile 200/pk, 10 pk/cs (2000/cs)	NON25444		9024	NW4404
DRAIN SPONGES - AVANT	MEDLINE	J&J	COVIDIEN	CARDINAL HEALTH
Avant Gauze Non-Woven Drain Sponges, 4x4, 6-ply, Sterile 2's, 25 pk/bx, 12 bx/cs (600/cs)	NON256000	2391	7086	DR4406-2
Avant Gauze Non-Woven Drain Sponges, 2x2, 6-ply, Sterile 2's, 35 pk/bx, 20 bx/cs (1400/cs)	NON256001	2392	7087	VS2206-2
POST-OP SPONGES	MEDLINE	J&J "Topper"	COVIDIEN "Cover"	CARDINAL HEALTH
Post-Op Sponges, 4x4, Sterile 2's 25 pk/bx, 24 bx/cs (1200/cs)	NON21442	2436	2913	CV0404-2
Post-Op Sponges, 4x3, Sterile 2's 25 pk/bx, 24 bx/cs (1200/cs)	NON21441	2435/2345	3157	CV0403-2
ABD PADS	MEDLINE	J&J "Surgipad"	COVIDIEN "Tendersorb/Curity"	CARDINAL HEALTH
ABD Pads 5x9 Sterile 1's, 25/bx, 16 bx/cs (400/cs)	NON21450	2145	9190/7196 256/cs/400/cs	AB0509-1
ABD Pads 8x7.5 Sterile 1's, 20/bx, 12 bx/cs (240/cs)	NON21453	2144	9192/7197 200/cs/240/cs	AB0875-1 500/cs
ABD Pads 8x10 Sterile 1's, 18/bx, 20 bx/cs (360/cs)	NON21454	2148 320/cs	9194/7198 144/cs/360/cs	AB0810-1
ABD Pads 5x9 Non-Sterile Bulk (576/cs)	NON21451		8190/6196	AB0509
ABD Pads 8x7.5 Non-Sterile Bulk (576/cs)	NON21452		8192/6197	AB0875
ABD Pads 8x10 Non-Sterile Bulk (432/cs)	NON21456		8194/6198	AB0810



BULKEE II GAUZE BANDAGES	MEDLINE	J&J "Kling Fluff"	COVIDIEN (Kendall) "Kerlix"	CARDINAL HEALTH (Allegiance) "Hi-Bulk"
Bulkee II Gauze Rolls 4.5" x 4.1 yds Sterile 1's, (100/cs)	NON25865	6930	6715	FR4541-1
Bulkee II Gauze Rolls 3.4" x 3.6 yds Sterile 1's, (96/cs)	NON25861	6964	6725	FR3436-1
Bulkee II Gauze Rolls 4.5" x 4.1 yds Non-Sterile, (100/cs)	NON25855	6937	1892	FR4541
BULKEE FLUFF SPONGES	MEDLINE	J&J "Kling Fluff"	COVIDIEN "Kerlix Super"	CARDINAL HEALTH
Bulkee Super Fluff Sponge, 6" x 6.75" Diag. Sterile 2's, 20/bx, 12 bx/cs (480/cs)	NON25852	2335 20/bx, 9 bx/cs	2585	FLFM-2
Bulkee Super Fluff Sponge, 6" x 6.75" Diag. Sterile 5's, 5/tray, 120 trays/cs (600/cs)	NON25853	2336 5/pk, 90 pk/cs	3085	FLFM-5
Bulkee Super Fluff Sponge, 6" x 6.75" Diag. Sterile 10's, 10/tray, 48 trays/cs (480/cs)	NON25854	2334	7310	FLFM-10
Bulkee Super Fluff Sponge, Lg 7.75" x 8.75" Diag. Sterile 2's, 25 pks/bx, 16 bx/cs (400/cs)	NON25856		4308	
Bulkee 4x4, 12 ply Sponge, Sterile 2's 25 pks/bx, 24 bx/cs (1200/cs)	NON25844		5072	
Bulkee 4x4, 12 ply Sponge, Sterile 10's 10 each/tray, 128 trays/cs (1280/cs)	NON25884		6120	
BULKEE LITE	MEDLINE	J&J "Kling"	COVIDIEN "Kerlix Lite"	CARDINAL HEALTH
Bulkee Lite Conforming bandages 2" x 3.5 yds non-sterile, 12/bx, 8 bx/cs (96/cs)	NON27492	6902	8052	PR0402
Bulkee Lite Conforming bandages 3" x 4.1 yds non-sterile, 12/bx, 8 bx/cs (96/cs)	NON27493	6903	8053	PR0403
Bulkee Lite Conforming bandages 4" x 4.1 yds non-sterile, 12/bx, 8 bx/cs (96/cs)	NON27494	6904	8054	PR0404
Bulkee Lite Conforming bandages 6" x 4.5 yds non-sterile, 6/bx, 8 bx/cs (48/cs)	NON27495	6906	8056	PR0406
Bulkee Lite Conforming bandages 2" x 3.5 yds Sterile 1's, 12/bx, 8 bx/cs (96/cs)	NON27496	6922	8072	PR0402-1
Bulkee Lite Conforming bandages 3" x 4.1 yds Sterile 1's, 12/bx, 8 bx/cs (96/cs)	NON27497	6923	8073	PR0403-1
Bulkee Lite Conforming bandages 4" x 4.1 yds Sterile 1's, 12/bx, 8 bx/cs (96/cs)	NON27498	6924	8074	PR0404-1
Bulkee Lite Conforming bandages 6" x 4.5 yds Sterile 1's, 6/bx, 8 bx/cs (48/cs)	NON27499	6926	8076	PR0406-1

STERILE SOF-FORM	MEDLINE	J&J	COVIDIEN (Kendall)	CARDINAL HEALTH (Allegiance)	
Sof-Form Bandages, 1" x 75" Sterile 1's, 12/bx, 8 bx/cs (96/cs)	NON254955				
Sof-Form Bandages, 2" x 75" Sterile 1's, 12/bx, 8 bx/cs (96/cs)	NON25496	6992	2231	GR0402-1	
Sof-Form Bandages, 3" x 75" Sterile 1's, 12/bx, 8 bx/cs (96/cs)	NON25497	6993 72/cs	2232	GR0403-1	
Sof-Form Bandages, 4" x 75" Sterile 1's, 12/bx, 8 bx/cs (96/cs)	NON25498	6994 72/cs	2236	GR0404-1	
Sof-Form Bandages, 6" x 75" Sterile 1's, 6/bx, 8 bx/cs (48/cs)	NON25499	6996 48/cs	2238	GR0406-1	
NON-STERILE SOF-FORM	MEDLINE	J&J "Sof-Kling"	COVIDIEN "Conform"	CARDINAL HEALTH	
Sof-Form Bandages, 1" x 75" Non-Sterile, 12/bx, 8 bx/cs (96/cs)	NON25491	6981	2239	GR0401	
Sof-Form Bandage, 2" x 75" Non-Sterile, 12/bx, 8 bx/cs (96/cs)	NON25492	6982	2242	GR0402	
Sof-Form Bandage, 3" x 75" Non-Sterile, 12/bx, 8 bx/cs (96/cs)	NON25493	6983	2244	GR0403	
Sof-Form Bandage, 4" x 75" Non-Sterile, 12/bx, 8 bx/cs (96/cs)	NON25494	6984	2247	GR0404	
Sof-Form Bandage, 6" x 75" Non-Sterile, 6/bx, 8 bx/cs (48/cs)	NON25495	6986	2249	GR0406	
ELASTIC NET, 25 yards LATEX-FREE	MEDLINE	GLENWOOD/ WESTERN MEDICAL "Surgilast"	MEDI-TECH "Spandage"	DEROYAL 1% LATEX "Stretch Net"	MEDICAL ACTION "SePro Net"
Fingers, toes, wrist, Size 1	NONNET01	GLLF2501	MT01	10-7421	53110LF
Small hand, arm, leg, foot, Size 2	NONNET02	GLLF2502	MT02	10-7422	53120LF
Medium hand, arm, leg, Size 3	NONNET03	GLLF2503	MT03	10-7433	53130LF
Large hand, arm, leg, foot, Size 4	NONNET04	GLLF2504	MT04	10-7424	53140LF
Small head, shoulder, thigh, Size 5	NONNET05	GLLF2505	MT05	10-7425	53150LF
Large head, shoulder, thigh, Size 6	NONNET06	GLLF2506	MT06	10-7427	53160LF
Small chest, back, peineum. axilla, Size 7	NONNET07	GLLF2507	MT07	10-7428	53170LF
Medium chest, back, perineum, axilla, Size 8	NONNET08	GLLF2508	MT08	10-7429	53180LF
Large chest, back, axilla, perineum, Size 9	NONNET09	GLLF2509	MT09	10-7430	53190LF
X-lg chest, back, axilla, perineum, Size 10	NONNET10	GLLF25010	MT010	10-7431	53210LF
XX-lg chest, back, axilla, perineum, Size 11	NONNET11	GLLF25011	MT011	10-7432	53220LF
3X-lg chest, back, axilla, perineum, Size 12	NONNET12	GLLF25012	MT012	10-7433	
Custom size	NONNET22		MT022		
ELASTIC NET, PRE-CUT INDIVIDUAL SIZES					
Finger, toe & wrist, Size 1, 24" Long	NONNET124P		MT1X24		
Small hand, arm, leg, foot, Size 2, 24" Long	NONNET224P		MT2X24		
Medium hand, arm, leg, foot, Size 3, 24" Long	NONNET324P		MT3X24		
Large hand, arm, leg, foot, Size 4, 24" Long	NONNET424P		MT4X24		
Small head, shoulder, thigh, Size 5, 24" Long	NONNET524P		MT5X24		
Large head, shoulder, thigh, Size 6, 24" Long	NONNET624P		MT6X24		
Large head, shoulder, thigh, Size 6, 36" Long	NONNET636P		MT6X36		
Small chest, back, perineum, axilla, Size 7, 24" Long	NONNET724P		MT7X24		
Small chest, back, perineum, axilla, Size 7, 36" Long	NONNET736P		MT7X36		
Medium chest, back, perineum, axilla, Size 8, 24" Long	NONNET824P		MT8X24		



ELASTIC NET, PRE-CUT INDIVIDUAL SIZES	MEDLINE		MEDI-TECH "Spandage"	
Large chest, back, perineum, axilla, Size 9, 24" Long	NONNET924P		MT9X24	
X-large chest, back, perineum, axilla, Size 10, 24" Long	NONNET1024P		MT10X24	
XX-large chest, back, perineum, axilla, Size 11, 24" Long	NONNET1124P		MT11X24	
3X-large chest, back, perineum, axilla, Size 12, 24" Long	NONNET1124P		MT12X24	
Custom Size, Size 22, 24" Long	NONNET1124P		MT22X24	
MATRIX ELASTIC BANDAGE, STERILE	MEDLINE	CONCO/HARTMANN		CARDINAL HEALTH (Allegiance)
Matrix Elastic Bandage, 2" x 5 yd, latex free, 20/cs	DYNJ05152LF	59720000		23593-12LF
Matrix Elastic Bandage, 3" x 5 yd, latex free, 20/cs	DYNJ05153LF	59730000		23593-13LF
Matrix Elastic Bandage, 4" x 5 yd, latex-free, 20/cs	DYNJ05154LF	59740000		23593-14LF
Matrix Elastic Bandage, 6" x 5 yd, latex-free, 20/cs	DYNJ05156LF	59760000		23593-16LF
Matrix Elastic Bandage, 6" x 10 yd, latex-free, 20/cs	DYNJ05157LF	59790000		
Matrix Elastic Bandage, 6" x 15 yd, latex-free, 20/cs	DYNJ05158LF			23593-S960
MATRIX ELASTIC BANDAGE, NON-STERILE	MEDLINE	CONCO/HARTMANN		CARDINAL HEALTH
Matrix Elastic Bandage, 2" x 5 yd, latex-free 10/bx, 5 bx/cs, 50/cs	MDS087002LF	59120000		23593-02LF
Matrix Elastic Bandage, 3" x 5 yd, latex-free 10/bx, 5 bx/cs, 50/cs	MDS087003LF	59130000		23593-03LF
Matrix Elastic Bandage, 4" x 5 yd, latex-free 10/bx, 5 bx/cs, 50/cs	MDS087004LF	59140000		23593-04LF
Matrix Elastic Bandage, 6" x 5 yd, latex-free 10/bx, 5 bx/cs, 50/cs	MDS087006LF	59160000		23593-06LF
Matrix Elastic Bandage, 4" x 10 yd, latex-free 4/bx, 5 bx/cs, 20/cs	MDS087104LF	59180000		
Matrix Elastic Bandage, 6" x 10 yd, latex-free 4/bx, 5 bx/cs, 20/cs	MDS087106LF	59190000		
Matrix Elastic Bandage, 6" x 15 yd, latex-free 4/bx, 5 bx/cs, 20/cs	MDS087156LF			23593-960
SOFT WRAP, STERILE	MEDLINE	CONCO/HARTMANN		
Soft-Wrap, 2" x 5 yd, latex-free, 20/cs	DYNJ05125	27200000		
Soft-Wrap, 3" x 5 yd, latex-free, 20/cs	DYNJ05126	27300000		
Soft-Wrap, 4" x 5 yd, latex-free, 20/cs	DYNJ05131	27400000		
Soft-Wrap, 6" x 5 yd, latex-free, 20/cs	DYNJ05132	27600000		
SOFT WRAP, NON-STERILE	MEDLINE	CONCO/HARTMANN	COVIDIEN	CARDINAL HEALTH
Soft-Wrap, 2" x 5 yd, latex-free 10/bx, 5 bx/cs, 50/cs	MDS046002	38200000	3302	23606-020
Soft-Wrap, 3" x 5 yd, latex-free 10/bx, 5 bx/cs, 50/cs	MDS046003	39300000	3550	23606-030
Soft-Wrap, 4" x 5 yd, latex-free 10/bx, 5 bx/cs, 50/cs	MDS046004	38400000	3616	23606-040
Soft-Wrap, 6" x 5 yd, latex-free 10/bx, 5 bx/cs, 50/cs	MDS046006	39600000	3930	23606-060

SWIFT WRAP, STERILE	MEDLINE	CONCO/HARTMANN		
Swift Wrap, 2" x 5 yd, 20/cs	DYNJ05144	59920000		
Swift Wrap, 3" x 5 yd, 20/cs	DYNJ05145	59930000		
Swift Wrap, 4" x 5 yd, 20/cs	DYNJ05146	59940000		
Swift Wrap, 6" x 5 yd, 20/cs	DYNJ05147	59960000		
SWIFT WRAP, NON-STERILE	MEDLINE	CONCO/HARTMANN	COVIDIEN	
2" x 5 yd, 10/bx, 5 bx/cs, 50/cs	MDS077002	59520000	4202	
3" x 5 yd, 10/bx, 5 bx/cs, 50/cs	MDS077003	59530000	4203	
4" x 5 yd, 10/bx, 5 bx/cs, 50/cs	MDS077004	59540000	4204	
6" x 5 yd, 10/bx, 5 bx/cs, 50/cs	MDS077006	59560000	4206	
SURE WRAP, NON-STERILE	MEDLINE	CONCO/HARTMANN	COVIDIEN	CARDINAL HEALTH
2" x 5 yd, beige, latex-free, 10/bx, 5 bx/cs, 50/cs	MDS055002	54200000	4422	23602-02LF
3" x 5 yd, beige, latex-free, 10/bx, 5 bx/cs, 50/cs	MDS055003	54300000	4423	23602-03LF
4" x 5 yd, beige, latex-free, 10/bx, 5 bx/cs, 50/cs	MDS055004	54400000	4424	23602-04LF
6" x 5 yd, beige, latex-free, 10/bx, 5 bx/cs, 50/cs	MDS055006	54600000	4426	23602-06LF
2" x 5 yd, white, 10/bx, 5 bx/cs, 50/cs	MDS057002	16200000		23602-021
3" x 5 yd, white, 10/bx, 5 bx/cs, 50/cs	MDS057003	16300000		23602-031
4" x 5 yd, white, 10/bx, 5 bx/cs, 50/cs	MDS057004	16400000		23602-041
6" x 5 yd, white, 10/bx, 5 bx/cs, 50/cs	MDS057006	16600000		23602-061
PETROLATUM GAUZE DRESSINGS	MEDLINE	J&J	COVIDIEN	CARDINAL HEALTH
Petrolatum Gauze Dressing, 1" x 36" 12/bx, 12 bx/cs (144/cs)	NON251136	2046 72/cs	888441260 72/cs	PG136
Petrolatum Gauze Dressing, 1" x 8" 12/bx, 12 bx/cs (144/cs)	NON251180	2045 72/cs	8884417601 200/cs	PG18
Petrolatum Gauze Dressing, 1/2" x 72" 12/bx, 12 bx/cs (144/cs)	NON251272	2053 72/cs	8884411600 72/cs	PG572
Petrolatum Gauze Dressing – overwrap 1/2" x 72" 12/bx, 12 bx/cs (144/cs)	NON252272			
Petrolatum Gauze Dressing, 3" x 18" 12/bx, 12 bx/cs (144/cs)	NON251318	2048 72/cs	8884414600 72/cs	PG318
Petrolatum Gauze Dressing – overwrap 3" x 18" 12/bx, 12 bx/cs (144/cs)	NON252318			
Petrolatum Gauze Dressing – 3" x 36" 12/bx, 12 bx/cs (144/cs)	NON251336	2049 72/cs	8884415600 72/cs	PG336
Petrolatum Gauze Dressing – 3" x 9" 12/bx, 12 bx/cs (144/cs)	NON251390	2047 72/cs	8884413605 0200/cs	PG39
Petrolatum Gauze Dressing – 6" x 36" 12/bx, 12 bx/cs (144/cs)	NON251636	2051 72/cs	8884416600 72/cs	PG636
XEROFORM	MEDLINE	J&J	COVIDIEN	CARDINAL HEALTH
Xeroform Petrolatum Gauze Dressing 1" x 8" 12/bx, 12 bx/cs (144/cs)	NON253180	2006 200/cs	8884433301	PGX18
Xeroform Petrolatum Gauze Dressing overwrap 1" x 8" 12/bx, 12 bx/cs (144/cs)	NON255180			
Xeroform Petrolatum Gauze Dressing 2" x 2" 12/bx, 12 bx/cs (144/cs)	NON253220		8884433400	
Xeroform Petrolatum Gauze Dressing 4" x 4" 12/bx, 12 bx/cs (144/cs)	NON253440		8884433400	
Xeroform Petrolatum Gauze Dressing 5" x 9" 12/bx, 12 bx/cs (144/cs)	NON253590	2007 200/cs	8884433605	PGX59

OIL EMULSION DRESSING	MEDLINE	J&J "Adaptic"	COVIDIEN (Kendall)	CARDINAL HEALTH (Allegiance)	
Oil Emulsion Gauze, 3" x 16", 36/bx, 6 bx/cs (216/cs)	NON250317	2014	6114		
Oil Emulsion Gauze, 3" x 3", 36/bx, 6 bx/cs (216/cs)	NON250330	2012	6112	OED33	
Oil Emulsion Gauze, 3" x 8", 36/bx, 6 bx/cs (216/cs)	NON250381	2015	6113	OED38	
Oil Emulsion Gauze, 3" x 8", 36/bx, 6 bx/cs (216/cs) (3 dressings/pkg), 648 dressings/cs	NON250383	2013	6115		
PAPER TAPE	MEDLINE	J&J "Dermiform"	COVIDIEN "Tenderskin"	CARDINAL HEALTH	3M "Micropore"
Paper Tape 1/2" x 10 yds, 24/bx, 12 bx/cs (288/cs)	NON260012	5180	1596		1530-0
Paper Tape 1" x 10 yds, 12/bx, 12 bx/cs (144/cs)	NON260001	5181	1914	32261	1530-2 120/cs
Paper Tape 2" x 10 yds, 6/bx, 12 bx/cs (72/cs)	NON260002	5182	2419	32262	1530-2 60/cs
Paper Tape 3" x 10 yds, 4/bx, 12 bx/cs (48/cs)	NON260003	5183	3394		1530-3
CLOTH TAPE	MEDLINE	J&J "Dermicel"	COVIDIEN "Curasilk"	CARDINAL HEALTH	3M "Durapore"
Cloth Tape 1/2" x 10 yds, 24/bx, 12 bx/cs (288/cs)	NON260112	5143	7137		1538-0 240/cs
Cloth Tape 1" x 10 yds, 12/bx, 12 bx/cs (144/cs)	NON260101	5144	7138	32266 120/cs	1538-1 120/cs
Cloth Tape 2" x 10 yds, 6/bx, 12 bx/cs (72/cs)	NON260102	5145	7139	32267 60/cs	1538-2 60/cs
Cloth Tape 3" x 10 yds, 4/bx, 12 bx/cs (48/cs)	NON260103	5146 40/cs	7140	32268 40/cs	1538-3
TRANSPARENT TAPE	MEDLINE	J&J "Dermiview"	COVIDIEN "Curity"	CARDINAL HEALTH	3M "Transpore"
Transparent Tape 1/2" x 10 yds 24/bx, 12 bx/cs (288/cs)	NON260212		8533		15270-0 240/cs
Transparent Tape 1" x 10 yds 12/bx, 12 bx/cs (144/cs)	NON260201	5161 120/cs	8534	32271 120/cs	1527-1 120/cs
Transparent Tape 2" x 10 yds 6/bx, 12 bx/cs (72/cs)	NON260202	5162 60/cs	8535	32272 60/cs	1527-2 60/cs
Transparent Tape 3" x 10 yds 4/bx, 12 bx/cs (48/cs)	NON260203	5163 40/cs	8536		1527-3 40/cs
WOUND CLOSURE STRIPS	MEDLINE	J&J "Proxi-Strip"	COVIDIEN "Curi-Strip"		3M "Steri-Strip"
Medi-Strips 1/8" x 3" 5/pk, 50 pk/bx, 4 bx/cs (200 pk/cs)	NON250318	100K	9890		R1540
Medi-Strips 1/4" x 3" 3/pk, 50 pk/bx, 4 bx/cs (200 pk/cs)	NON250314	101K	9892		R1541
Medi-Strips 1/4" x 1.5" 6/pk, 50 pk/bx, 4 bx/cs (200 pk/cs)	NON250114	102K	9891		R1542
Medi-Strips 1/4" x 4" 10/pk, 50 pk/bx, 4 bx/cs (200 pk/cs)	NON250414	103K	9893		R1546
Medi-Strips 1/2" x 4" 6/pk, 50 pk/bx, 4 bx/cs (200 pk/cs)	NON250412	104K	9894		R1547
Medi-Strips 1" x 5" 4/pk, 25 pk/bx, 4 bx/cs (100 pk/cs)	NON250501				R1548

ADHESIVE BANDAGES - PLASTIC	MEDLINE	J&J "Band-Aid"	COVIDIEN (Kendall) "Curity"	CARDINAL HEALTH (Allegiance)
Adhesive Bandage 3/4" x 3" - Sheer-Gard 100/bx, 12 bx/cs (1200/cs)	NON25500	4634	3065 3600/cs	BS3403 9600/cs
Adhesive Bandage 1" x 3" - Sheer-Gard 100/bx, 12 bx/cs (1200/cs)	NON25600	4644 1500/cs	3270 3600/cs	BS0103 9600/cs
Adhesive Bandage X-Lg 2" x 4" - Sheer-Gard 50/bx, 12 bx/cs (600/cs)	NON25504	5716 2x4.5" 900/cs	7539	BS0245 1200/cs
Adhesive Bandage 3/8" x 1-1/2" - Sheer-Gard "Junior" 100/bx, 36 bx/cs (3600/cs)	NON25509	4633		
Adhesive Bandage 7/8" Spot - Sheer-Gard 100/bx, 12 bx/cs (1200/cs)	NON25501	4930 3600/cs	6403 3600/cs	BS0078 20,000/cs
Adhesive Bandage 5" x 3" - Medtoons 50/bx, 24 bx/cs (1200/cs)	NON256130		9519	
ADHESIVE BANDAGES - WOVEN CLOTH	MEDLINE	J&J "Band-Aid"	COVIDIEN	CARDINAL HEALTH
Adhesive Bandage 3/4" x 3" - Comfort Cloth 100/bx, 12 bx/cs (1200/cs)	NON25650	4434	5024	BF3403 7680/cs
Adhesive Bandages 1" x 3" - Comfort Cloth 100/bx, 12 bx/cs (1200/cs)	NON25660	4444	6135	BF0103 7680/cs
Adhesive Bandage Knuckle - Comfort Cloth 100/bx, 12 bx/cs (1200/cs)	NON25510	4438	6411	BF0001 4400/cs
Adhesive Bandage Fingertip - Comfort Cloth 100/bx, 12 bx/cs (1200/cs)	NON25513	4436	6409	BF0002 4200/cs
ADHESIVE BANDAGE – MEDTOON	MEDLINE	J&J "Band-Aid"	COVIDIEN	CARDINAL HEALTH
Adhesive Bandage 3/4" x 3" – Medtoon 50/bx, 24 bx/cs (1200/cs)	NON256130	AB67	9519 91300	1075737 1073737
PACKING STRIPS	MEDLINE	COVIDIEN "Curity"	DYNAREX	CARDINAL HEALTH
Packing Strip, Plain, 1" x 5 yd	NON255015	7633	3403	59321
Packing Strip, Plain, 2" x 5 yd	NON255025	7634	3404	59421
Packing Strip, Plain, 1/2" x 5 yd	NON255125	7632	3402	59221
Packing Strip, Plain, 1/4" x 5 yd	NON255145	7631	3401	59121
Packing Strip, Iodoform, 1" x 5 yd	NON256015	7833	3413	59346
Packing Strip, Iodoform, 2" x 5 yd	NON256025	7834	3414	59446
Packing Strip, Iodoform, 1/2" x 5 yd	NON256125	7831	3412	59246
Packing Strip, Iodoform, 1/4" x 5 yd	NON256145	7832	3411	59146
EYE PADS	MEDLINE	J&J	COVIDIEN	CARDINAL HEALTH
1-5/8" x 2-5/8" Pad, 1's, 50/bx, 12 bx/cs, 600/cs	NON21600	8771	2841	EP2581
2-1/8" x 2-5/8" Pad, 1's, 50/bx, 12 bx/cs, 600/cs	NON21601	8773		EP2582



NON-ADHERENT PADS, STERILE	MEDLINE	J&J "Release"	COVIDIEN (Kendall) "Telfa"	CARDINAL HEALTH (Allegiance)
2" x 3", 1's 100/box, 12 bx/cs, 1200/cs	NON25700	2052	1961 2400/cs	NAO302-1
3" x 4", 1's 100/box, 12 bx/cs, 1200/cs	NON25710	2054/4731	2132/1050 2400/cs	NAO304-1
3" x 8", 1's 50/box, 12 bx/cs, 600/cs	NON25720	2058	1238	NAO308-1
ADHESIVE SURGICAL DRESSING, STERILE	MEDLINE	J&J	COVIDIEN	
4" x 6" (4" x 3" pad), 1's, latex-free 25/box, 8 bx/cs, 200/cs	NON4311	4311	5112	
8" x 6" (8" x 3" pad), 1's, latex-free 25/box, 4 bx/cs, 100/cs	NON4313	4313	5113	
SELF-ADHERENT BANDAGES, STERILE	MEDLINE	ANDOVER "Co-Flex NL"		3M "Coban"
Co-Flex NL™, 2" x 5 yds stretched, tan, latex-free 20 rolls/cs	DYNJ05123			
Co-Flex NL™, 3" x 5 yds stretched, tan, latex-free 20 rolls/cs	DYNJ05124	5300S		2083S 24/pk
Co-Flex NL™, 4" x 5 yds stretched, tan, latex-free 20 rolls/cs	DYNJ05129	5400S		2084S 18/pk
Co-Flex NL™, 6" x 5 yds stretched, tan, latex-free 20 rolls/cs	DYNJ05130	5600S		2086S 12/pk
SELF-ADHERENT BANDAGES, NON-STERILE	MEDLINE	ANDOVER "Co-Flex NL"		3M "Coban"
Co-Flex NL™, 1" x 5 yds stretched, tan, latex-free 30 rolls/cs	MDS088001	5100TN		2081
Co-Flex NL™, 2" x 5 yds stretched, tan, latex-free 36 rolls/cs	MDS088002	5200TN		2082
Co-Flex NL™, 3" x 5 yds stretched, tan, latex-free 24 rolls/cs	MDS088003	5300TN		2083
Co-Flex NL™, 4" x 5 yds stretched, tan, latex-free 18 rolls/cs	MDS088004	5400TN		2084
Co-Flex NL™, 6" x 5 yds stretched, tan, latex-free 12 rolls/cs	MDS088006	5600TN		2086

NEED SOMETHING YOU CAN'T FIND IN THIS BROCHURE?

This brochure contains hundreds of Medline manufactured products. But if you can't find something you need, please ask your Medline Sales Representative for assistance, call **1-800-MEDLINE (1-800-633-5463)** or consult our online catalog at **www.medline.com**.

SELF-ADHERENT BANDAGES, STERILE	MEDLINE		COVIDIEN (Kendall) "Flex-Wrap"	3M "Coban"
Co-Flex Med, 2" x 5 yds stretched, tan, 20 rolls/cs	DYNJ086002			
Co-Flex Med, 3" x 5 yds stretched, tan, 20 rolls/cs	DYNJ086003		4583S	1583S
Co-Flex Med, 4" x 5 yds stretched, tan, 20 rolls/cs	DYNJ086004		4584S, 18/pk	1584S, 18/pk
Co-Flex Med, 6" x 5 yds stretched, tan, 20 rolls/cs	DYNJ086006		4586S, 12/pk	1586S, 12/pk
SELF-ADHERENT BANDAGES, NON-STERILE	MEDLINE	ANDOVER "Co-Flex Med"	COVIDIEN "Flex-Wrap"	3M "Coban"
Co-Flex Med, 1" x 5 yds stretched, tan, 36 rolls/cs	MDS086001	7100TN	4581	1581
Co-Flex Med, 2" x 5 yds stretched, tan, 36 rolls/cs	MDS086002	7200TN	4582	1582
Co-Flex Med, 3" x 5 yds stretched, tan, 24 rolls/cs	MDS086003	7300TN	4583	1583
Co-Flex Med, 4" x 5 yds stretched, tan, 18 rolls/cs	MDS086004	7400TN	4584	1584
Co-Flex Med, 6" x 5 yds stretched, tan, 12 rolls/cs	MDS086006	7600TN	4586	1586
SELF-ADHERENT BANDAGE, NON-STERILE, LATEX-FREE	MEDLINE	ANDOVER "Co-Flex LF2"		
Co-Flex LF2, 1" x 5 yds stretched, tan, 36 rolls/cs	MDS089001	9100		
Co-Flex LF2, 2" x 5 yds stretched, tan, 30 rolls/cs	MDS089002	9200		
Co-Flex LF2, 3" x 5 yds stretched, tan, 24 rolls/cs	MDS089003	9300		
Co-Flex LF2, 4" x 5 yds stretched, tan, 18 rolls/cs	MDS089004	9400		
Co-Flex LF2, 6" x 5 yds stretched, tan, 12 rolls/cs	MDS089006	9600		
SELF-ADHERENT BANDAGE, STERILE, LATEX-FREE	MEDLINE	ANDOVER "Co-Flex LF2"		
Co-Flex LF2, 2" x 5 yds stretched, tan, 20 rolls/cs	DYNJ089002	9200S		
Co-Flex LF2, 3" x 5 yds stretched, tan, 20 rolls/cs	DYNJ089003	9300S		
Co-Flex LF2, 4" x 5 yds stretched, tan, 20 rolls/cs	DYNJ089004	9400S		
Co-Flex LF2, 6" x 5 yds stretched, tan, 18 rolls/cs	DYNJ089006	9600S		
WYTEX, STERILE	MEDLINE		COVIDIEN	
Undercast padding, 2" x 4 yds, 20 rolls/cs	DYNJ066002		2283	
Undercast padding, 3" x 4 yds, 20 rolls/cs	DYNJ066003		2394	
Undercast padding, 4" x 4 yds, 20 rolls/cs	DYNJ066004		2502	
Undercast padding, 6" x 4 yds, 20 rolls/cs	DYNJ066006		2554	
WYTEX, NON-STERILE	MEDLINE		COVIDIEN	
Undercast padding, 2" x 4 yds, 24/bg, 3 bg/cs, 72/cs	MDS066002		1418	
Undercast padding, 3" x 4 yds, 12/bg, 6 bg/cs, 72/cs	MDS066003		2059	
Undercast padding, 4" x 4 yds, 12/bg, 6 bg/cs, 72/cs	MDS066004		3175	
Undercast padding, 6" x 4 yds, 6/bg, 6/cs, 36/cs	MDS066006		3489	
SYN-TEX, NON-STERILE	MEDLINE	J&J	COVIDIEN	3M
Undercast padding, Synthetic, 2" x 4 yds, 12/bg, 6 bg/cs, 72/cs	MDS067002	9052	6242	CMW02
Undercast padding, Synthetic, 3" x 4 yds, 12/bg, 6 bg/cs, 72/cs	MDS067003	9033	6243	CMW03
Undercast padding, Synthetic 4" x 4 yds, 12/bg, 6 bg/cs, 72/cs	MDS067004	9034	6244	CMW04
Undercast padding, Synthetic 6" x 4 yds, 6/bg, 6/cs, 36/cs	MDS067006	9036	6246	CMW06



CURAD® Sesame Street Plastic Strips Are the Ideal Bandages for Everyday Use for Kids

Kids can collect all 30 unique and colorful designs. Our strips are made from soft, flexible plastic and a skin-friendly adhesive that doesn't leave a sticky, itchy residue on the skin, so you are ensured maximum comfort as you heal. 20 assorted sizes per box. Sterile.

Sesame Street® and associated characters, trademarks and design elements are owned and licensed by Sesame Workshop. ©2008 Sesame Workshop. All Rights Reserved.

Item#	Description	Pkg
CUR47069	6 - 1 1/8" x 2 1/8", 14 - 3/4" x 2 1/8" (2.9 x 5.4 cm, 1.9 x 5.4 cm)	20 ea/bx, 24 bx/cs

CURAD® Curious George™ Adhesive Bandage

Relive the adventures of our favorite monkey with Curad Curious George adhesive bandages. Each box includes 20 strips of assorted sizes (constructed from water resistant vinyl). An absorbent, non-stick island pad seals out dirt and germs and is "ouchless" to remove. A must-have bandage for any pediatrics unit.



Item#	Description	Pkg
CUR47067	6 - 1 1/8" x 2 1/8" (2.9 x 5.4 cm), 20 ea/bx, 24 bx/cs 14 - 3/4" x 2 1/8" (1.9 x 5.4 cm)	

CURAD® My Little Pony™ Adhesive Bandage

Little patients' tears turn to smiles when My Little Pony™ brand bandages are applied. Each box includes 20 strips of assorted sizes (constructed from water resistant vinyl). High quality color printing brings the ponies alive on each adhesive bandage. An absorbent, non-stick island pad seals out dirt and germs and is "ouchless" to remove. A must-have bandage for any pediatrics unit.



Item #	Description	Pkg
CUR14927	6 - 1 1/8" x 2 1/8" (2.9 x 5.4 cm), 20 ea/bx, 24 bx/cs 14 - 3/4" x 2 1/8" (1.9 x 5.4 cm)	



CURAD® Camp Camo Adhesive Bandages

Our authentic green, pink and blue camo patterns are printed on woven fabric that stretches with movement. Strong adhesive offers durable protection. The absorbent non-stick island pad is ouchless to remove. Ideal for knees, elbows and looking cool!

Fabric, Sterile

Item #	Description	Pkg
CUR45701	3/4" x 3", Green (1.9 x 7.62 cm)	25/bx, 24 bx/cs
CUR45702	3/4" x 3", Assorted (1.9 x 7.62 cm) Pink (13 each) Blue (12 each)	25 bx/cs, 600/cs
CUR45704	3 3/4" x 2" (9.5 x 5.08 cm)	12 ea/bx, 24 bx/cs

Camo Tape

Item #	Description	Pkg
CUR45705	1 1/2" x 10 yds (3.8 cm x 9.1 m)	1 rl/bx, 24 rls/cs



Medline Industries, Inc.
One Medline Place
Mundelein, IL 60060

Medline United States
1-800-MEDLINE (633-5463)
www.medline.com | info@medline.com

Medline Canada
1-800-396-6996
www.medline.ca | canada@medline.com

Medline México
01-800-831-0898
www.medline.mexico.com | mexico@medline.com

We reserve the right to correct any errors that may occur within this brochure. Some products may not be available for sale in Mexico or Canada. Please contact your Account Representative for additional information. ©2009 Medline Industries, Inc. Medline is registered trademarks of Medline Industries, Inc. Adaptic is a registered trademark of Johnson & Johnson Corporation. Avant Gauze is a Registered Trademark of Medline Industries, Inc. Band-Aid is a registered trademark of Johnson & Johnson Corporation. Bulkee II is a Registered Trademark of Medline Industries, Inc. Bulkee Lite is a Registered Trademark of Medline Industries, Inc. Coban is a registered trademark of 3M. Co-Flex is a Registered Trademark of Andover Coated Products, Inc. Comfort Cloth is a Registered Trademark of Medline Industries, Inc. Conform is a Registered Trademark of Tyco Healthcare Group Curad is a Registered Trademark of Medline Industries, Inc. Curasilk is a registered trademark of Beiersdorf Aktiengesellschaft Corporation Curi-strip is a registered trademark of Tyco Healthcare Group Curity is a registered trademark of Tyco Healthcare Group. Dermiform is a registered trademark of Johnson & Johnson Corporation. Durapore is a registered trademark of 3M. Dynarex is a registered trademark of Dynarex Corporation. Elastoplast is a registered trademark of Beiersdorf Aktiengesellschaft Corporation. Flex-Fabric is a Trademark of Medline Industries, Inc. Flex-Wrap is a registered trademark of Tyco Healthcare Group. Hy-Tape is a registered trademark of Hy-Tape International. Kerlix is a registered trademark of Tyco Healthcare Group. Kerlix Lite is a registered trademark of Tyco Healthcare Group. Kerlix Super is a registered trademark of Tyco Healthcare Group. Kling Fluff is a registered trademark of Johnson & Johnson Corporation. Mediplast is a Registered Trademark of Medline Industries, Inc. Medi-strip is a registered trademark of Medline Industries, Inc. Medtoons is a Registered Trademark of Medline Industries, Inc. Megazinc Pink is a Trademark of Omega Medical Products Corporation. Micropore is a registered trademark of 3M. Nu-Gauze is a registered trademark of Johnson & Johnson Corporation. Ouchless is a registered trademark of Johnson & Johnson Corporation. Proxi-strip is a registered trademark of Johnson & Johnson Corporation. Sheer-Gard is a registered trademark of Medline Industries, Inc. Sof-Form is a registered trademark of Medline Industries, Inc. Soft-Wrap is a registered trademark of Medline Industries, Inc. Spandage is a trademark of Medi-Tech International Corporation Sten-strip is a registered trademark of 3M Sure-Wrap is a registered trademark of Medline Industries, Inc. Surgilast is a registered trademark of Western Medical Inc. Swift-Wrap is a registered trademark of Medline Industries, Inc. Telfa is a registered trademark of Tyco Healthcare Group. Tenderskin is a registered trademark of Tyco Healthcare Group. Tendersorb is a registered trademark of Tyco Healthcare Group. Zonas is a registered trademark of Johnson & Johnson Corp. Curad Hydro Heal and Curad Sensitive Skin are registered trademarks of Medline Industries, Inc. Transpore is a registered trademark of 3M. Versalon is a registered trademark of Tyco Healthcare Group. MKT209080 / LIT175 / 5M / TKZ15