



# OR FLUID MANAGEMENT

Fluid Control and Disposal Solutions

FLUID MANAGEMENT  
PRODUCTS FOR THE O.R.  
**LOOK NO FURTHER  
THAN MEDLINE**

Medline is your one-stop provider  
for all your fluid management needs.





# TABLE OF CONTENTS

Suction Canisters	4
Suction Canister Carriers	5
Suction Canister Holders	5
Jumbo Jug™	5
Solidifiers	6
Medline Absorbent Pads	7
QuickSuite® Floor Suction Pad	8
QuickSuite Floor Mats	9
QuickSuite Floor Suction	9
Medline Yankauers	10
Frazier Suction Instrument	10
Medline Lipo-Suction Clear, Vinyl Tubing	11
Sterile, Non Conductive Connecting Tubing	11



## Medline Suction Canisters

We've heard you loud and clear. You told us your #1 issue with suction cannisters is finding, connecting and misplacing spare parts like elbows. This new elbow-less Medline cannister makes that issue a thing of the past.

- Horizontal ports eliminate the need for elbows \*
- Patent pending lid design simplifies use, reduces O.R. turnaround time
- Available in two popular sizes, 1500 cc and 2400 cc



DYNDSC2400



DYNDSC1500

### ORDERING INFORMATION

	Item #	Description	Pkg
Suction Canisters	DYNDSC2400	2400 cc With Turret Lid	27/cs
	DYNDSC1500	1500 cc With Turret Lid	45/cs

## Suction Canister Carriers

- Carries up to four disposable suction canisters of either 1500 cc or 2400 cc size
- Coating provides easy to clean surface
- Adjustable in height using an easy two button push method

Ask about our **FREE ACCESSORY** program!  
Eligible customers can receive the suction canister carriers and holders for **FREE!**



DYNDSCC

## Canister Holders

**Holds either 1500 cc or 2400 cc suction canisters. Insert into slide bracket on the carrier.**

- Coating provides easy to clean surface
- Strong metal construction



DYNDSCH

## JumboJug™ Large Volume Suction Canisters

- Integrated ergonomic handle reduces the risk of fluid contact and makes transporting simple and efficient
- Jumbo collection capacity: 13,000 ml (13 Liter)
- Fits standard kick bucket stands
- Saves surgical time by eliminating frequent canister changes
- Works great with Medline's New QuickSuite® Floor Suction products



Multi-Connection

JumboJug™  
DYDJJ13JK

### ORDERING INFORMATION

	Item #	Description	Pkg
<b>Suction Canister Carrier</b>	DYNDSCC	Suction Canister Carrier	1 ea/cs
<b>Suction Canister Holders</b>	DYNDSCH	Suction Canister Holder	20 ea/cs
<b>Suction Canister Walmount</b>	DYNDSCW	Suction Canister Walmount	20 ea/cs
<b>JumboJug™ Large Volume Suction Canisters</b>	DYNDJJ13K	JumboJug™, 13,000 ml	6/cs

## Liqui-Loc™

Liqui-Loc solidifiers improve safety, speed up O.R. turnover times and even contribute to a greener environment.

- Two different delivery methods—Dissolvable PVA packs and bottles
- Virtually eliminates employee exposure to aerosolized bloodborne pathogens
- Reduces OR turnover time—Canister contents are solidified when the procedure is complete
- You can add Liqui-Loc to your canister at the end of the procedure, if this is your protocol
- Earth-friendly—PVA packs eliminate bottle disposal forever



Liqui-Loc can be added to the suction canister prior to the procedure for safety and convenience.

Color-coded bottles for easy identification.

### ORDERING INFORMATION

	Item #	Description	Pkg
<b>Solidifier/ Non Disinfectant</b>	DYNS500B	Solidifier, 500 cc	96 ea/cs
	DYNS1500B	Solidifier, 1500 cc	96 ea/cs
	DYNS3000B	Solidifier, 3000 cc	96 ea/cs
	DYNS16000B	Solidifier, 16000 cc	20 ea/cs
<b>Solidifier/Disinfectant</b>	DYNS500DB	Solidifier/Disinfect, 500 cc	96 ea/cs
	DYNS3000DB	Solidifier/Disinfect, 3000 cc	96 ea/cs
	DYNS2000DB	Solidifier/Disinfect, 3000 cc	96 ea/cs
	DYNS1500DB	Solidifier/Disinfect, 1500 cc	96 ea/cs
	DYNS1200DB	Solidifier/Disinfect, 1200 cc	96 ea/cs

**Eliminates bottle waste and O.R. turnover time.**

Choose the new Liqui-Loc™ disolvable PVA Packs from Medline.



DYND3000P



DYND1500P



DYND1000P



DYND500P



**ORDERING INFORMATION**

	Item #	Description	Pkg
<b>Solidifier/ Non Disinfectant</b>	DYND500P	Solidifier, 500 cc	96 ea/cs
	DYND3000P	Solidifier, 3000 cc	96 ea/cs
	DYND1500P	Solidifier, 1500 cc	96 ea/cs
	DYND1000P	Solidifier, 1000 cc	96 ea/cs
<b>Solidifier/Disinfectant</b>	DYND500DP	Solidifier/Disinfect, 500 cc	96 ea/cs
	DYND3000DP	Solidifier/Disinfect, 3000 cc	96 ea/cs
	DYND1500DP	Solidifier/Disinfect, 1500 cc	96 ea/cs



# QuickSuite® Floor Management Products

## The New QuickSuite Floor Suction Pad

- Combines the functionality of absorption and suction to efficiently remove large amounts of fluid from the O.R. floor
- Quiet and easy-to-use
- Works with standard wall suction
- Can be rolled over without loss of suction
- Red Bag friendly.

Absorbs and suctions at the same time.




### ORDERING INFORMATION

	Item #	Description	Pkg
<b>Suction Pad</b>	DYNDORP3648AS	Yellow, 36" x 48" (91 cm x 122 cm)	5/cs



## QuickSuite Floor Mats

- Helps keep floors clean, dry, and slip-free during fluid-intense procedures
- Quickly absorb fluids while conforming to uneven floor surfaces without bunching
- Mats are offered in 4 varieties depending on need and preference
- Choose from one or two-sided absorbency
- Many different sizes to mix and match based on your needs.



Absorbent,  
Impervious Pad  
DYNDORP4640A



Absorbent Pad  
DYNDORP3640A



Absorbent, Impervious  
DYNDORP3036A

## QuickSuite Floor Suction

- Disc-shaped devices that quietly suction fluid from the operating room floor
- Deluxe and basic models available
- Maneuvered with the push of a foot and include flexible tubing for easy repositioning
- Work with standard wall suction devices



Basic Floor Suction  
DYNDQSF1000



Deluxe Floor Suction  
DYNDQSF2000

### ORDERING INFORMATION

	Item #	Description	Pkg
<b>Absorbent, Impervious Pad</b>	DYNORP4640A	Blue, 46" x 40" (117 cm x 102 cm)	10/cs
<b>Absorbent Pad</b>	DYNDORP3640A	Yellow, 36" x 40" (91 cm x 102 cm)	30/cs
	DYNDORP3672A	Yellow, 36" x 72" (91 cm x 183 cm)	12/cs
<b>Absorbent, Impervious</b>	DYNDORP3036A	White w/ Blue Backing, 30" x 36" (76 cm x 91 cm)	20/cs
<b>Basic Floor Suction</b>	DYNDQSF1000	Blue	20/cs
<b>Deluxe Floor Suction</b>	DYNDQSF2000	Yellow	20/cs

## Medline Yankauers

- Slip-resistant handle for secure grip
- Smooth inner surface for maximum, non-occluding flow
- Transparent construction for easy assessment of aspirate
- Shatterproof construction aids patient safety
- Ribbed connector for various sizes of connecting tubing
- Optional control vent for suctioning control



DYND50138



DYND50130



DYND52130



DYND50140



DYND50143



DYND50144

## Introducing our Frazier Suction Instrument

- Sterile obturator to help clear obstructions
- Malleable aluminum shaft can be bent to the desired position without occlusion
- Control vent for continuous or intermittent suctioning
- Ergonomic handle for improved grip and ease of use
- Universal connector easily connects to many different sizes of suction tubing



DYNDFR12S

### ORDERING INFORMATION

	Item #	Description	Pkg
<b>Bulb Tip Yankauers, Sterile</b>	DYND50130	Without Control Vent, Rigid	50/cs
	DYND52130	Without Control Vent, Flexible	50/cs
	DYND50132	With Control Vent, Rigid	50/cs
	DYND52132	With Control Vent, Flexible	50/cs
<b>Flange Tip Yankauers, Sterile</b>	DYND50140	Without Control Vent, Rigid	50/cs
	DYND50142	With Control Vent, Rigid	50/cs
	DYND50143	Standard Capacity, No Vent, Flexible	50/cs
	DYND50144	Fine Capacity, No Vent, Flexible	50/cs
<b>Bulb Tip Yankauers with Tubing, Sterile</b>	DYND50135	6' (1.8 m) Tubing, Without Control Vent	20/cs
	DYND50137	6' (1.8 m) Tubing, With Control Vent	20/cs
	DYND50138	10' (3.0 m) Tubing, Without Control Vent	20/cs
<b>Frazier Suction Instrument</b>	DYNDFR8S	8Fr Frazier Suction, with Control Vent, Sterile	40/cs
	DYNDFR10S	10Fr Frazier Suction, with Control Vent, Sterile	40/cs
	DYNDFR12S	12Fr Frazier Suction, with Control Vent, Sterile	40/cs
<b>Poole Suctioning Instruments, Sterile</b>	DYND50180	Poole Tip	50/cs
	DYND50172	Poole Tip with 72" (1.8 m) Non-Conductive Vinyl Tubing	20/cs


## Medline Liposuction Clear, Vinyl Tubing

Heavy-duty liposuction tubing with 1/4" (0.6 cm) wall thickness for strength and flexibility.

- Universal fit
- Highly durable
- Easy-to-use flexibility
- Utilize with pressures up to 29.5 in Hg

"We trialed Medline liposuction tubing in three different cases today with three different plastic surgeons. . . The physicians and the nurses liked it better than the special order tubing we were using because it is longer, more flexible and easier to use. We will now use Medline's liposuction tubing."


HealthSouth Outpatient Surgery Center  
Birmingham, Alabama



DYND50110

## Sterile, Non-Conductive Connecting Tubing

- Available in 5 mm (3/16") or 6 mm (1/4") inner diameter tubing
- Female connectors
- For pressures up to 22 in Hg (Mercury)



DYND50211



DYND50221

### ORDERING INFORMATION

	Item #	Description	Pkg
<b>Medline Liposuction Clear, Vinyl Tubing</b>	DYND50110	3/8" (0.95 cm) x 5/8" (1.59 cm) x 10' (3.0 m)	10/cs
	DYND50112	3/8" (0.95 cm) x 5/8" (1.59 cm) x 12' (3.7 m)	10/cs
<b>Sterile, Non-Conductive Connecting Tubing</b>	DYND50211	3/16" (0.48 cm) x 20" (50.8 cm)	50/cs
	DYND50216	3/16" (0.48 cm) x 6' (1.8 m), With Male Connector	50/cs
	DYND50221	3/16" (0.48 cm) x 10' (3.0 m), With Male Connector	50/cs
	DYND50223	3/16" (0.48 cm) x 12' (3.7 m), With Male Connector	20/cs
	DYND50226	3/16" (0.48 cm) x 20' (6.1 m)	20/cs
	DYND50240	3/16" (0.48 cm) x 20" (50.8 cm) & 6' (1.8 m), Tandem Tube Kit With Male Connector	50/cs
	DYND50246	1/4" (0.64 cm) x 6' (1.8 m), With Male Connector	50/cs
	DYND50251	1/4" (0.64 cm) x 10' (3.0 m), With Male Connector	50/cs
	DYND50252	1/4" (0.64 cm) x 12' (3.7 m), With Male Connector	20/cs
	DYND50253	1/4" (0.64 cm) x 20' (6.1 m), With Male Connector	20/cs
	DYND50230	1/4" (0.64 cm) x 100' (30.5 m), Non-Sterile, Bubble Connector Every 3' (0.9 m) Bulk (No Female Connector)	1/bx





**Medline Industries, Inc.**

One Medline Place  
Mundelein, IL 60060

**Medline United States**  
**1-800-MEDLINE (633-5463)**

[www.medline.com](http://www.medline.com) | [info@medline.com](mailto:info@medline.com)

**Medline Canada**  
**1-800-396-6996**

[www.medline.ca](http://www.medline.ca) | [canada@medline.com](mailto:canada@medline.com)

**Medline México**  
**01-800-831-0898**

[www.medlinemexico.com](http://www.medlinemexico.com) | [mexico@medline.com](mailto:mexico@medline.com)

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure. ©2010 Medline Industries, Inc. Medline and QuickSuite are registered trademarks of Medline Industries, Inc. Liqui-Loc and JumboJug are trademarks of Medline Industries, Inc.